BD RPR Card Antigen Suspension

For Use in Performing the **Macro-Vue™** RPR Card Tests
(Qualitative and Quantitative) for the Serologic Detection of Syphilis

English: pages 1 – 2 Italiano: pagine 4 – 5
Français : pages 2 – 3 Português: páginas 5 – 6
Deutsch: Seiten 3 – 4 Español: páginas 6 – 7

0210761JAA(03)
2015-08

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкциите. / Pokyny vám poskytné místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyttä lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нусқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie ziskate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakt närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinize temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

The **BD Macro-Vue™** RPR Card Tests (Qualitative and Quantitative) are for the serologic detection of syphilis. The RPR Card Antigen Suspension is a component of the test.

SUMMARY AND EXPLANATION

RPR Card Antigen Suspension is a reagent component part used in performing the **BD Macro-Vue** RPR (Rapid Plasma Reagin) 18 mm Circle Test which is a nontreponemal testing procedure for the serologic detection of syphilis and is described in *A Manual of Tests for Syphilis*.¹ The test techniques described in this manual are the most widely used in the United States and have been thoroughly evaluated and standardized.

PRINCIPLES OF THE PROCEDURE

The RPR Card Antigen Suspension is a carbon-particle cardioliipin antigen^{2,3,4} which detects “reagin,” an antibody-like substance present in sera from syphilitic persons, and occasionally from persons with other acute or chronic conditions.

When a specimen contains reagin, flocculation occurs with a coagglutination of the carbon particles of the RPR Card antigen, which appear as black clumps against the white background of the plastic-coated card. This coagglutination can be read macroscopically. By contrast, nonreactive specimens appear to have an even light-gray color.

REAGENTS

The ingredients* of the RPR Card Antigen Suspension are¹ 0.003% cardioliipin, 0.020 – 0.022% lecithin, 0.09% cholesterol, 0.0125 M EDTA, 0.01 M Na₂HPO₄, 0.01 M KH₂PO₄, 0.2% sodium azide (preservative), 0.02% charcoal (specially prepared, Becton Dickinson), 10% choline chloride, w/v, and deionized/distilled water.

*Adjusted and/or supplemented as required to meet performance criteria.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

RPR Card Antigen Suspension is recommended only for use in the **BD Macro-Vue** RPR Card Test procedures. Prior to using the antigen, it is essential that the direction sheet accompanying each **BD Macro-Vue** RPR Card Test Kit be reviewed.

WARNING

H302 Harmful if swallowed. **H315** Causes skin irritation. **H319** Causes serious eye irritation. **P101** If medical advice is needed, have product container or label at hand. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P264** Wash thoroughly after handling. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage: Refrigerate (2 to 8 °C) the RPR Card Antigen Suspension. Once the antigen has been placed in the *dispensing bottle* (Cat. No. 270509) and refrigerated (2 to 8 °C), the reactivity remains satisfactory for approximately three months, or until the expiration date, if it occurs sooner. Label the dispensing bottle with the antigen lot number, expiration date and date antigen was placed in the bottle.

Storage of the antigen in bright sunlight or temperatures above 30 °C (86°F) should be avoided; such conditions may cause a rough appearance of the antigen when used with nonreactive sera. If the ampule of antigen is frozen during shipment, it can be reconstituted once by warming to room temperature; repeated freezing and thawing should be avoided. Immediate use of a refrigerated antigen may result in decreased sensitivity of the test. Therefore, upon removal from the refrigerator, allow the antigen to warm to room temperature (23 to 29 °C) before use.

The antigen suspension should not be used beyond the expiration date.

PROCEDURE

Materials Provided: RPR Card Antigen Suspension, with pre-installed ampule opener, a reagent component part of the **BD Macro-Vue** RPR Card Test Kits.

Materials Required But Not Provided:

1. Antigen plastic dispensing bottles (Cat. No. 270509).
2. Dispensing square cut 20 gauge needles (Cat. No. 273509).
3. **BD Macro-Vue** RPR Card Test cards (Cat. Nos. 271849, 271949 or 272001).
4. Control sera with established patterns of graded reactivity which should be included in each day's testing to confirm optimal reactivity of the antigen suspension.
5. A rotator, 100 ± 2 rpm, circumscribing a circle of 2 cm in diameter, with automatic timer, friction drive and a cover containing a moistened blotter or sponge. See section “Availability” in **BD Macro-Vue** RPR Card Test Kit direction sheets for RPR Card Test Control Cards and **BD Macro-Vue** Card Test Rotator with Humidifying Cover.

6. Saline (0.9%) for use in quantitative testing. Prepare by adding 900 mg dry NaCl, ACS to 100 mL with deionized/distilled water.
7. Serum nonreactive to syphilis in 0.9% Saline; required for diluting test specimens giving a Reactive result at the 1:16 dilution.

Also required are the necessary equipment and labware used in serology testing laboratory for preparation, storage and handling of serologic specimens.

TEST PROCEDURE, QUALITY CONTROL, RESULTS, LIMITATIONS OF THE PROCEDURE, EXPECTED VALUES AND PERFORMANCE CHARACTERISTICS

This package is a reagent component part of the **BD Macro-Vue™** RPR Card Test Kits. Refer to the direction sheet included in all kits.

AVAILABILITY

Cat. No. Description

	RPR Card Antigen Suspension (dated and refrigerated product):
270333	3 mL, ampule, for use with the 18 mm Circle Card Test, Box of three.
270309	3 mL, ampule, for use with the 18 mm Circle Card Test, Box of ten.

REFERENCES

1. Creighton, E.T. 1990. Rapid plasma reagin (RPR) 18-mm circle card test, p. 99-108. In S.A. Larsen, E. F. Hunter, and S. J. Kraus, (ed.), A manual of tests for syphilis, 8th ed. American Public Health Association, Washington, D.C.
2. Portnoy, J., J. H. Brewer, and A. Harris. 1962. Rapid plasma reagin card test for syphilis and other treponematoses. Pub. Health Rep. 77:645-652.
3. Portnoy, J. 1963. Modifications of the rapid plasma reagin (RPR) card test for syphilis, for use in large scale testing. Am. J. Clin. Pathol. 40:473-479.
4. Falcone, V.H., G.W. Stout, and M.B. Moore, Jr. 1964. Evaluation of rapid plasma reagin (circle) card test. Public Health Rep. 79:491-495.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Benex Limited