
Revisions

Rev from	Rev to	ECO #
F	G	5508-10

Notes:


- BD Cat. Number 421571, 421570
- Blank (Sheet) Size: Length: TBD Width: TBD
 Number of Pages: TBD Number of Sheets: 1
 Page Size: Length 8" Width 5.25" Final Folded Size: 4" x 5.25"
- Style (see illustrations below): #4


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 1570000011		Category and Description Package Insert, Clay Adams Brand Sedi Stain	Sheet: 1 of 7 Scale: N/A	A

BD Clay Adams™ Sedi Stain Concentrated Stain

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INTENDED USE

Sedi-Stain Concentrated Stain is a stabilized modification of the Sternheimer-Malbin stain for use in the microscopic examination of urine sediment.

SUMMARY AND EXPLANATION

The microscopic examination of the urinary sediment is generally recognized to be a most valuable diagnostic technique. Next to an actual tissue biopsy of the kidney, the microscopic findings in the urinary sediment are the clearest indicators of intrinsic renal disease.¹

Consistent and reliable interpretation of unstained specimens is often difficult due to the refractive nature of many of the formed elements, the requirement for reduced illumination, and the need for experience on the part of the observer.^{2,3}

In 1951, Sternheimer and Malbin⁴ reported on the use of a stain for urinary sediment which was more convenient to use than previous staining methods, and which provided clear delineation of all the formed elements in the sediment, as well as valuable diagnostic information.

The original formulation, however, required filtration and often formed precipitates if stored beyond three months.

Sedi-Stain Concentrated Stain is a stabilized modification of the Sternheimer-Malbin urinary stain. This highly selective formula stains blood cells, casts, and other formed elements in urinary sediment in a distinctive fashion which permits rapid and accurate identification.

PRINCIPLES OF THE PROCEDURE

The chemical and physical properties of the various formed elements of urinary sediment cause the crystal violet and safranin to be taken up in varying proportions. The resultant distinctive staining permits rapid and accurate identification.

Of particular significance are the internal staining characteristics such as the presence or absence of nuclei or cytoplasmic granules which would allow for differential staining and physical properties of cellular elements in urine.

REAGENTS

Sedi-Stain Concentrated Stain for urinary sediment is a stabilized modification of Sternheimer-Malbin stain.

Sedi-Stain Concentrated Stain does not form a sediment such as that formed in ordinary S-M stain and does not require filtration as required by ordinary S-M stain.

Sedi-Stain Concentrated Stain contains:

Crystal Violet.....	0.10%	Ethyl Alcohol (SD-3A).....	10.00%
Safranin O.....	0.25%	Water and Stabilizers.....	89.62%
Ammonium Oxalate.....	0.03%		

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Storage: Normal room temperature storage conditions are adequate to protect the stability of **Sedi-Stain Concentrated Stain**. Exposure to light and normal variations in ambient temperature will not affect the stain.

Product Deterioration: If a precipitate develops in the stain, it should be discarded.

SPECIMEN COLLECTION AND PREPARATION

WARNING: Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"^{5,6} and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

The urine specimen should be a freshly voided sample collected in a clean container which is closed immediately and processed as soon as possible. There is no need for elaborate cleansing procedures for routine microscopic examination unless the specimen is contaminated by vaginal discharge or hemorrhage.⁹

Specimens which cannot be examined immediately should be refrigerated, but must not be frozen. If examination is to be delayed beyond 4 hours, it is suggested that a drop of 40% formalin be added to prevent growth of microorganisms and preserve urinary sediments.¹⁰

Chloroform should not be used as a preservative since it settles to the bottom of the container and may interfere with microscopic examination.

Caution should be exercised in selecting preservatives which will not interfere with other tests to which the specimen may be subjected.

PROCEDURE

1. Use a freshly voided specimen collected in a clean container.
2. Mix the specimen thoroughly and pour into a centrifuge tube.
3. Centrifuge for 5 minutes at about 400 g's (e.g., 1500 RPM, 6" radius).
4. Decant the supernatant without disturbing the sediment.
5. Add 1-2 drops of stain to the sediment in the tube. (Occasionally, a technician may prefer to use 3 drops of stain.)
6. Thoroughly mix the contents of the tube by flicking the bottom of the tube sharply with the index finger several times, or by mixing with a disposable pipette.
7. Transfer 1 drop of the stained sediment onto a microscope slide. A coverslip may be placed over the drop to facilitate handling and to provide a uniform layer.
8. As an alternative to steps 5, 6, and 7, transfer 1 drop of the unstained sediment from the centrifuge tube onto a microscope slide and add up to one drop of stain directly to the sediment on the slide. Small quantities of stain may be transferred to the sediment by using a glass rod instead of the dropper.
9. Examine microscopically. Report average number of red blood cells and white blood cells per high power field. Report average number of casts and other formed elements per low power field, as examined under dimmed light.

When using **BD Sedi-Cal™ Urine Centrifuge Tubes**:

Follow directions which accompany tubes. Be SURE TO: Add only 1 drop of stain to the sediment. Use a coverslip when reading microscopically on the plastic cap to ensure a constant depth of specimen.

User Quality Control

Quality control requirements must be followed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

LIMITATIONS OF THE PROCEDURE

Microscopic examination of urinary sediment is a semi-quantitative procedure. In cases where exact counts of leukocytes, bacteria, casts, etc., are required, alternative methods may be used.

EXPECTED VALUES

Some red blood cells, white blood cells, and casts are excreted by normal individuals, but they are seen only occasionally in urinary sediments examined microscopically.

From 0-3 red blood cells per high powered field (males), 0-5 red blood cells per high powered field (females), 0-5 white blood cells per high powered field and 0-4 hyaline casts per low powered field are accepted as normal.¹¹

PERFORMANCE CHARACTERISTICS

Sternheimer and Malbin,⁴ using a stain containing both crystal violet and safranin, demonstrated that valuable diagnostic information could be obtained by staining the sediment of centrifuged urine. Pyuria and hematuria were identified based on the differential staining characteristics of white blood cells, epithelial elements and casts that took up the stain well as compared to red blood cells that did not; also the ratio of white blood cells to red blood cells was quickly estimated. White blood cells stained either deep red to violet or pale blue. The violet staining cells were uniform in size, contained a dark red or deep purple nucleus and violet granules.

They occurred commonly in lower urinary tract infections without renal involvement. Dead bacteria stained dark purple, while living bacteria remained either unstained or stained pink; mycelia and spores of fungi appeared light and purple. *Trichomonas* parasites were either colorless or pale blue. Swollen, pale blue-staining cells with granules showing Brownian movement (identified as leukocytes by a positive peroxidase reaction) were present in inflammatory renal diseases, such as acute cystitis, kidney abscesses and particularly in advanced cases of pyelonephritis.

AVAILABILITY

Cat. No.	Description
421570	Sedi-Stain™ Concentrated Stain , 12.5 mL
421571	Sedi-Stain™ Concentrated Stain , 62 mL
420597	Sedi-Cal™ Centrifuge Tubes with caps, 4 mL, packaged 100 tubes/shelf pack.

REFERENCES

1. Fuller, J.B., (1963) "The Role of the Urine in Diagnosis" in Audiovisual Seminar Program No. 100 of Am. Soc. Clin. Path., p.2.
2. Lynch, M.J. et al., (1969) Medical Laboratory Technology and Clinical Pathology, W.B. Saunders Co., Philadelphia, pp. 101 & 102.
3. An Atlas of Use of Sternheimer-Malbin staining technique in examination of urinary sediments, Abbot Laboratories, North Chicago, Ill., (1960, 1961).
4. Sternheimer, R. and Malbin, B., (1951) "Clinical Recognition of Pyelonephritis with a New Stain for Urinary Sediments," Am. J. Med., 11, 312
5. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
6. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
7. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
8. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
9. Kark, R.M., et al., (1963) A Primer of Urinalysis, Harper & Row, New York, p. 5.
10. Newell, J.E. and Duke, E., (1962) "The Routine Examination of Urine in the Laboratory," Am. Soc. Clin. Path., p.8.
11. Jacobs, D.S., et al., (1990) Laboratory Test handbook, 2nd Edition, Lexi-Comp, Inc. Stow, Ohio, p.937

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Benex Limited