
**Revisions**

Rev from	Rev to	ECO #
1003	0104	2493-04

**Notes:**

- BD Cat. Number 240905, 240906
- Blank (Sheet) Size : Length: N/A      Width: N/A  
 Number of Pages: N/A      Number of Sheets: N/A  
 Page Size: Length N/A      Width N/A      Final Folded Size: N/A
- Style (see illustrations below): N/A

- See Specification Control Number 240905 for Material Information
- Ink Colors: Printed two sides  Yes     No  
 No. of Colors: N/A
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 <b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 201359F		Category and Description Package Insert Link2 H. Pylori	Sheet: 1 of 27 <hr/> Scale: 1:1	<b>A</b>

# BD LINK2™ *H. Pylori* Rapid Test

For the qualitative detection of IgG antibodies to *Helicobacter pylori* (*H. pylori*) in SERUM

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Italiano: pagine 15 - 19  
 Español: páginas 20 - 24

201359F  
 2004/01

CLIA Moderately Complex Procedure

See symbol glossary at end of insert. / Se symbolglossaret i slutningen af indlægssedlen. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Consulte o glossário de símbolos no fim do folheto informativo. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

## INTENDED USE

LINK2™ *H. pylori* Rapid Test is a qualitative immunochemical membrane assay for the detection of IgG antibodies to *H. pylori* in serum by visual interpretation. The test is for professional use as an aid in the diagnosis of *H. pylori* infection in adult symptomatic patients. The LINK2 *H. pylori* Rapid Test Quality Control Kit, available separately, is designed to be used in conjunction with the LINK2 *H. pylori* Rapid Test to verify test performance.

## SUMMARY AND EXPLANATION

*Helicobacter pylori* (formerly called *Campylobacter pylori*) was first isolated from human gastric biopsy material in 1982.<sup>1</sup> Epidemiological studies indicate that *H. pylori* is present in approximately 30% of adults in the industrialized nations. Prevalence increases with age (more than 50% of individuals may be infected above the age of 50) and correlates with low socio-economic status in childhood. Transmission is most probably by the fecal-oral or oral-oral route.

Colonization by this gram-negative, microaerophilic bacterium is characterized by acute inflammatory reaction, with infiltration of the lamina propria by mononuclear (and frequently polymorphonuclear) cells. Infection is self-limiting in a few individuals, but the majority develop an active, chronic (predominantly antral) gastritis of varying severity. *H. pylori* infection is found in more than 90% of duodenal ulcer patients and in around 75% of all peptic ulcer sufferers.

*H. pylori* infection is also more common in gastric cancer patients. The risk of gastric cancer has been estimated to be six-fold higher in *H. pylori* infected populations than in uninfected populations.<sup>2</sup>

Although the actual mechanism by which *H. pylori* infection promotes ulcer formation is unknown, several studies have concluded that eradication of *H. pylori* in ulcer patients leads to an eight-fold reduction in the rate of recurrence found with conventional anti-ulcer therapies.<sup>3-5</sup>

An accurate diagnosis of *H. pylori* infection is of value in targeting antibiotic-based eradication regimens in symptomatic patients. Since *H. pylori* is rarely found in patients with Type A (autoimmune) or Type C (reactive) gastritis, or in reflux esophagitis, the detection of infection can aid in the differential diagnosis of peptic ulcer disease.

The LINK2 *H. pylori* Rapid Test is an immunoassay for the detection of antibodies against *Helicobacter pylori*, the causative agent of Type B chronic gastritis.<sup>1,2</sup> The test has been designed to give a positive result in patients who are or have been infected with *H. pylori*. The test is easy to perform and requires no specialized equipment or instrumentation. Visual interpretation provides an accurate qualitative result in the physician's office, hospital or laboratory. Diagnosis of *H. pylori* infection by antibody immunoassay can reduce the number of patients requiring endoscopy.<sup>6</sup>

The LINK2 *H. pylori* Rapid Test Quality Control Kit allows the user to implement a program of quality control testing to ensure the validity of test results obtained with the LINK2 *H. pylori* Rapid Test.

## PRINCIPLES OF THE PROCEDURE

LINK2 *H. pylori* Rapid Test is a dry chemistry immunological test which detects antibodies to *H. pylori*.

The test device contains a strip of plastic-backed membrane held in a result "window." In contact with this strip is a porous pad containing blue latex particles which have been coated with highly purified *H. pylori* antigen. The membrane has been coated with two lines in the window area (i.e., highly purified *H. pylori* antigen on the test line and a control line which captures the blue latex particles).

When a serum sample is introduced into the test device, it solubilizes the colored latex dried into the porous pad. Any antibodies to *H. pylori* in the serum bind to the antigen-coated latex particles. The serum then moves by capillary action across the membrane. Any anti-*H. pylori* present in the serum binds to the antigen on the test line. The presence of bound antibodies is revealed by the colored latex reagent. Latex particles then migrate along the test strip where they are captured by the control line. The presence of two lines in the result window indicates a positive test while only a single control line is a negative test.

The LINK2 *H. pylori* Rapid Test Quality Control Kit, available separately, comprises one bottle each of positive and negative control. Each contains human serum. The serum used to prepare the positive control has been found to contain

a moderate level of antibodies to *H. pylori* while that used for the negative control has no antibodies to *H. pylori* and little activity. The positive control produces a positive test result due to the presence of antibodies to *H. pylori* and helps validate both the LINK2 *H. pylori* Rapid Test and operator performance. The negative control does not contain antibodies to *H. pylori* and produces a negative test result, which again assists in verifying test performance.

#### REAGENTS

The LINK2 *H. pylori* Rapid Test Kit includes:

<b>Test Devices</b>	15	Each foil pouch contains one test device and a orange desiccant package.
<b>Sample Collection Combs</b>	15	Each plastic bag contains five sample collection combs.
<b>Test Result Interpretation Card</b>	1	Card

**Storage Instructions:** LINK2 *H. pylori* Rapid Test kits are stable until the expiration date when stored between 15 - 30°C. Do not freeze.

The LINK2 *H. pylori* Rapid Test Quality Control Kit includes:

Control +	1 mL	Ready to use positive control serum with 0.1% sodium azide (preservative).
Control -	1 mL	Ready to use negative control serum with 0.1% sodium azide (preservative).

**Storage and Stability:** Store the LINK2 *H. pylori* Rapid Test Quality Control Kit refrigerated at 2 - 8°C. Do not freeze. The unopened kit is stable until the labeled expiration dating printed on the outer box. After opening, the controls are stable for 3 months when stored at 2 - 8°C or the original expiration date, whichever comes first.

#### Warnings and Precautions:

For *in vitro* Diagnostic Use.

Do not use test components or controls beyond the expiration date. Do not open foil pouch until just before use.

Do not use test if orange desiccant package is completely green.

Always wear gloves when taking blood samples, using the controls and performing the test. Proper attention should also be given to potential infectivity of the samples used in the test. All blood products and materials that come in contact with them should be handled as potentially infectious.

**WARNING:** Because no test method can offer complete assurance that HIV, hepatitis B virus, or other infectious agents are absent, SPECIMENS AND CONTROL SOLUTIONS SHOULD BE HANDLED AS THOUGH CAPABLE OF TRANSMITTING AN INFECTIOUS DISEASE. "Standard Precautions"<sup>7-10</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

After completion of the test, all components should be disposed of as biohazardous waste.

Positive and Negative Controls contain sodium azide, which is very toxic by inhalation, in contact with skin and if swallowed. Contact with acids liberates very toxic gas. After contact with skin wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposing, flush with a large volume of water to prevent azide buildup.

#### PROCEDURE

**Materials Provided:** see "Reagents".

**Materials Required But Not Provided:** Timing device.


##### Sample Collection and Test Procedures

Fresh serum should be used as a sample type for the test.


Approximately 50 µL of test sample is needed.

**To use a serum sample:**

1. Dip the collection comb into serum which has been separated from red blood cells after centrifugation.
2. Make certain that the collection comb is completely filled with serum. Do not overfill the comb.

- Carefully push the collection comb into the test device until the comb locks into place. Place the test device on a flat surface. Begin timing the test. Results are read in 5 min.

**To use the LINK2 *H. pylori* Rapid Test Quality Control Kit:**

- Remove positive and negative control serum bottles from the refrigerator and allow to warm to room temperature (18 - 30°C) for 15 - 30 min prior to use.
- Mix a control bottle gently by repeated inversion. Remove the cap.
- Holding a sample collection comb in one hand, carefully invert the control bottle until the nozzle is touching the teeth of the comb. Gently squeeze the bottle so that serum fills the gaps between the teeth of the comb; if necessary, move the nozzle over the surface of the comb to facilitate filling.
- Insert the collection comb into the LINK2 *H. pylori* Rapid Test device until the comb locks into place. Place the test device on a flat surface. Begin timing the test. Results are read in 5 min.

**Quality Control**

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

A procedural control is built into the LINK2 *H. pylori* Rapid Test: the presence of a control line in the test result window provides assurance that sufficient sample amount was collected and that the latex particles in the strip were properly rehydrated and flowed through the test and control areas, thereby indicating that the test is performing properly and the procedure has been run correctly. If the control line does not appear, the test result is invalid.

A built-in negative background control is provided by the clearing of background color in the test window. This verifies the test has been performed correctly. This area should be straw colored to light pink within 5 min and not interfere with reading the test result. A clear background in the test window is considered an internal negative background control. If background color remains in the test window which interferes with your ability to read the test, the test result is invalid.


The LINK2 *H. pylori* Rapid Test Quality Control Kit is designed to be used in conjunction with the LINK2 *H. pylori* Rapid Test to verify test performance. The positive control produces TWO distinct blue lines in the test result window indicating that the LINK2 *H. pylori* Rapid Test is performing satisfactorily. The negative control produces ONE distinct blue line above the "C" in the test result window indicating that the LINK2 *H. pylori* Rapid Test is performing satisfactorily. Neither the shade nor the intensity of the control line should be used as a reference for the appearance of the quality control test line.

A positive and negative external control must be tested with each test kit, with each change in operator within a test kit and whenever recommended storage temperature ranges have been exceeded. If controls do not perform as expected, contact BD Technical Services in the U.S.A. at (800) 638-8663 or your local BD office.

**INTERPRETATION OF TEST RESULTS**

A positive result may appear within a few minutes when the liquid reaches the far end of the test window, otherwise read the test at 5 min. Results should be interpreted at a normal reading distance.


**Positive:** Two distinct blue lines in the result window, one above the letter C and one above the letter T, indicate the presence of IgG antibodies to *H. pylori*.

**Negative:** One distinct blue line in the result window above the letter C indicates the absence of detectable IgG antibodies.

Any other result is invalid and the test should be repeated. If the sample fails to flow across the test window, the test is invalid and should be repeated.

*Positive test results which appear after 5 min are not valid.*

**LIMITATIONS OF THE PROCEDURE**

The LINK2 *H. pylori* Rapid Test should be used only to evaluate patients with clinical signs and symptoms suggestive of gastrointestinal disease. The test is not intended for use with asymptomatic patients. Performance characteristics for persons under the age of 18 have not been established with this test.

The results obtained from this test are intended to be an aid in diagnosis only. Each physician must interpret the results in conjunction with the patient's history, physical findings and other diagnostic procedures.

A negative test result indicates that antibodies to *H. pylori* are either not present or at levels undetectable by the test. A positive test indicates the presence of antibodies to *H. pylori*, but determination of an active or inactive infection cannot be made.

The **LINK2 *H. pylori* Rapid Test** is qualitative. The intensity of the sample test line does not necessarily correlate to the concentration of antibody in the blood.

The effect of hematocrit on the **LINK2 *H. pylori* Rapid Test** was assessed using venous blood from four volunteers. Hematocrits were altered by removing and adding plasma. Hematocrits ranged from 22.5% to 64.8%. Altering the hematocrit did not affect the accuracy of the test.

Test results are not affected by elevated levels of creatinine or bilirubin or by hemolysis.

The controls provided in the **LINK2 *H. pylori* Rapid Test Quality Control Kit** are designed to be used only with the **LINK2 *H. pylori* Rapid Test**. Their validity with other similar tests cannot be guaranteed.

Both controls are prepared from human serum rather than whole blood. The positive control has been manufactured to equate to a sample in the middle of the range of *H. pylori* positive subjects and yields a test line that is moderately intense.

**EXPECTED RESULTS**

Epidemiological studies indicate that *H. pylori* is present in around 20 to 30% of adults in the industrialized nations. Prevalence increases with age (more than 50% of individuals may be infected above the age 50) and correlates with low socio-economic status in childhood. Close person to person contact in childhood is an important determinant of sero-prevalence of *H. pylori* in adulthood, suggesting that the infection is transmitted directly from one person to another and may be commonly acquired early in life.<sup>11</sup>

A study to determine the presence of *H. pylori* in asymptomatic adults was conducted at a U.S. hospital and at a medical device manufacturer. Ninety-eight adults of varying ages participated. Each participant was evaluated using the **LINK2 *H. pylori* Rapid Test** and an ELISA serum test. Of these participants, 7.1% (7 out of 98) were found to be positive by both methods.

**PERFORMANCE CHARACTERISTICS**

Two studies were performed to evaluate the performance of the **LINK2 *H. pylori* Rapid Test**.

One study was performed at four U.S. sites and one UK site with 229 patients to evaluate the performance of the **LINK2 *H. pylori* Rapid Test** to the true condition of the patient. Each patient involved in this study was evaluated with the following tests for *H. pylori*: **LINK2 *H. pylori* Rapid Test**, ELISA *H. pylori* antibody detection assay, CLO-test for active *H. pylori* infection, and histology analysis for detection.

Because *H. pylori* bacteria are sometimes difficult to collect during biopsy, the actual bacteria may not be sampled by the physician during the biopsy collection. Antibody detection tests are more likely to determine if an infection is present or has recently occurred, provided that the patient is not immunosuppressed and is actually producing antibodies to *H. pylori* bacteria.

For this study, a "true positive" is defined as having a positive ELISA test, as well as either a positive CLO-test or a positive histology analysis.

A "true negative" is defined as having a negative ELISA test, and both the CLO-test and histology analysis must also be negative.

Any results not meeting the criteria for "true negative" or "true positive" are ruled as indeterminate.

Combining data from the four U.S. sites and the one UK site, 229 capillary samples were collected in order to compare the **LINK2 *H. pylori* Rapid Test** results to the true condition of the patient. Twenty-nine samples were classified as indeterminate. The combined results for the remaining 200 samples are shown in Table 1.

**Table 1**  
**LINK2 *H. pylori* Capillary Result**

	Positive	Negative
True Positive	70	9
True Negative	18	103

Clinical Sensitivity (95% CI): 88.6% (79.5%, 94.7%)

Clinical Specificity (95% CI): 85.1% (77.5%, 90.9%)

The second study was performed with the same U.S. sites used in the first study combined with one additional U.S. site. The purpose was to compare the performance of serum to the performance of capillary blood with the **LINK2 *H. pylori* Rapid Test**.

Combining data from the sites, 274 serum **LINK2 *H. pylori* Rapid Test** results were compared to capillary **LINK2 *H. pylori* Rapid Test** results. The combined results are shown in Table 2.

**Table 2**  
**LINK2 *H. pylori* Result**

		Serum	
		Positive	Negative
Capillary	Positive	102	8
	Negative	16	148

Comparative Sensitivity (95% CI): 92.7% (86.2%, 96.8%)

Comparative Specificity (95% CI): 90.2% (84.7%, 94.3%)

Accuracy 91.2%

**Reproducibility:** Reproducibility testing was performed at three external sites. Each day for 3 days, 10 blinded samples were tested at each of the three sites. These blinded samples consisted of negatives, weak positives and positives. Additionally, three patients traveled to each site for blinded capillary blood testing.

During reproducibility testing, at least 90% of negatives were found to be negative and 100% of positives were found to be positive.

**Cross Reactivity:** Studies were conducted to investigate the extent to which the antibodies detected by the LINK2 *H. pylori* Rapid Test cross-react with proteins from other related bacteria.

Crude bacterial sonicates at concentrations ranging from 0 to 2500 µg/mL, were added to a *H. pylori* antigen detection assay and their inhibition effects assessed. LINK2 *H. pylori* Rapid Test antigens were used as a comparison, since they compete in the assay and hence inhibit binding.

The results showed that proteins from the following bacteria do not inhibit the binding of LINK2 *H. pylori* antigen and, therefore, do not cross-react with the LINK2 *H. pylori* Rapid Test.

*Borrelia burgdorferi*    *Campylobacter jejuni*

*Campylobacter coli*    *Escherichia coli*

*Campylobacter fetus*

#### AVAILABILITY

Cat. No.	Description
240905	LINK2™ <i>H. pylori</i> Rapid Test Kit, 15 tests
240906	LINK2™ <i>H. pylori</i> Rapid Test Quality Control Kit

#### REFERENCES

- Warren, J.R. and B.J. Marshall. 1983. Unidentified curved bacillus on gastric epithelium in active gastritis. *Lancet* 1:1273-1275.
- The EUROGAST Study Group: 1993. An international association between *Helicobacter pylori* infection and gastric cancer. *Lancet* 341: 8847.
- Fiocca, R., et al. 1991. Duodenal ulcer relapse after eradication of *Helicobacter pylori*. *Lancet* 337: 1614.
- Graham, D.Y., et al. 1992. Effect of treatment of *Helicobacter pylori* infection on the long-term recurrence of gastric or duodenal ulcer. *Ann. Intern. Med.* 116: 705.
- Borody, T., et al. 1992. *Helicobacter pylori* reinfection 4 years post eradication. *Lancet* 339: 1295.
- Sobola, G. et al. 1991. Screening dyspepsia by serology to *Helicobacter pylori*. *Lancet* 338: 94.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
- Webb, P.M., et al. 1994. Relation between infection with *Helicobacter pylori* and living conditions in childhood; evidence for person to person transmission in early life. *Br. Med. J.* 308: 750.

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BD, BD Logo and LINK2 are trademarks of Becton, Dickinson and Company. ©2004 BD.
BENEX Limited