**Revisions**

Rev from	Rev to	ECO #
0104	0504	2873-04

**Notes:**


- BD Cat. Number 440454
- Blank (Sheet) Size : Length: 11"      Width: 8.5"  
 Number of Pages: 12      Number of Sheets: 3  
 Page Size: Length 8.5"      Width 5.5"      Final Folded Size: 8.5" x 5.5"
- Style (see illustrations below): # 5


#1

#2


#3


#4

#5

#6

#7

- See Specification Control Number 3300809 for Material Information
- Ink Colors: Printed two sides  Yes       No  
 No. of Colors: 1      PMS# Standard Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 3300809JAA		Category and Description Package Insert, BD ProbeTec ET Urine Processing Kit	Sheet: 1 of 13 Scale: N/A	A

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## INTENDED USE

The **BD ProbeTec™ ET Urine Processing Kit** contains urine processing pouches (UPPs) which are designed to remove amplification inhibitors for testing male and female urine specimens with the **BD ProbeTec ET Chlamydia trachomatis** and *Neisseria gonorrhoeae* amplified DNA assays on the **BD ProbeTec ET System**. The UPP may be added at the collection site or laboratory test site.

## SUMMARY AND EXPLANATION

When in contact with a urine specimen for at least 2 h, the UPP removes inhibitors of Strand Displacement Amplification, thereby increasing the likelihood of detecting *C. trachomatis* and *N. gonorrhoeae* with the **BD ProbeTec ET System**. When placed in the urine specimen at the collection site, the UPP also stabilizes the specimen for room temperature transport.

## REAGENTS:

Each **BD ProbeTec ET Urine Processing Kit** contains 4 storage bags - each containing 25 UPPs and 1 desiccant pouch. Each UPP contains >1.4 grams of resin.

**Materials Required But Not Provided:** Sterile, plastic, preservative-free, specimen collection cups; gloves; disinfectant.

**Storage Instructions:** Store kit at 15 – 30°C. Do not use past the expiration date.

**Precautions:** For *in vitro* Diagnostic use.

## Warnings:

1. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>1-4</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
2. The UPP must be in contact with the sample for at least 2 h prior to processing.
3. The UPP is for use with male and female urine specimens to be tested using the **BD ProbeTec ET system**.
4. Optimal performance of the **BD ProbeTec ET** amplified DNA assays requires proper specimen collection, handling and transport.
5. The UPP must be added to the collection cup **AFTER** the specimen has been collected.
6. Storage bags containing unused UPPs **MUST** be carefully resealed after opening. Verify that the desiccant is present prior to resealing the storage bag.
7. Proper labeling should accompany each specimen to the laboratory.

## SAMPLE COLLECTION AND TRANSPORT

The Urine Processing Pouch (UPP) can be added at either the Test Site or at the Collection Site. Instructions are provided for each option.

### Addition of UPP at Test Site

#### Urine Specimen Collection

1. The patient should not have urinated for at least 1 h prior to specimen collection.
2. Collect specimen in a sterile, plastic, preservative-free specimen collection cup.
3. The patient should collect the first 15 – 20 mL of voided urine (the first part of the stream – NOT midstream).  
**NOTE:** During the clinical evaluation, testing urine volumes up to 60 mL was included in the performance estimates.
4. Label with patient identification and date/time collected.

#### Urine Storage and Transport (Addition of UPP at Test Site):

**NOTE:** Specimens must be shipped in an insulated container with ice using either an overnight or 2-day delivery vendor. Storage up to 4 days has been validated with clinical specimens; storage up to 6 days has been demonstrated with seeded specimens.

1. Store and transport urine specimens to the test site at 2 – 8°C within 4 – 6 days of collection.
2. Add the UPP to the urine specimen collection cup. Wear clean gloves when handling the UPP and urine specimen.

**NOTE:** Do not place UPP on any work surface. Remove UPP from pouch with a freshly gloved hand or with clean, sterile forceps.

**NOTE:** Add UPP carefully to avoid splashing.

3. Cap the collection cup and swirl gently to ensure that the UPP is completely submerged in urine.
4. The UPP must be in contact with the urine specimen for at least 2 h prior to processing.
5. Do not freeze the urine specimen.

#### Addition of UPP at Collection Site

#### Urine Specimen Collection

1. The patient should not have urinated for at least 1 h prior to specimen collection.
2. Collect specimen in a sterile, plastic, preservative-free specimen collection cup.
3. The patient should collect the first 15 – 20 mL of voided urine (the first part of the stream – NOT midstream).

**NOTE:** During the clinical evaluation, testing urine volumes up to 60 mL were included in the performance estimates.

4. Immediately add the UPP to the specimen collection cup. Wear clean gloves when handling the UPP and urine specimen.  
**NOTE:** Do not place UPP on any work surface. Remove UPP from pouch with a freshly gloved hand or with clean, sterile forceps.  
**NOTE:** Add UPP carefully to avoid splashing.
5. Cap collection cup and swirl gently to ensure that the UPP is completely submerged in urine.
6. Label with patient identification and date/time collected.

#### Urine Storage and Transport (Addition of UPP at Collection Site):

**NOTE:** If specimens cannot be transported directly to the testing laboratory under ambient temperatures (15 – 27°C) and must be shipped, an insulated container with ice should be used with either an overnight or 2-day delivery vendor. Storage up to 4 days has been validated with clinical specimens; storage up to 6 days has been demonstrated with seeded specimens.

1. Store and transport the urine specimens containing a UPP to the laboratory or test site at 2 – 8°C within 4 – 6 days of collection or at 15 – 27°C within 2 days of collection.
2. Do not freeze the urine specimen.
3. The UPP must be in contact with the urine specimen for at least 2 h prior to processing.

For domestic and international shipments, specimens should be labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

#### REFERENCES

1. National Committee for Clinical Laboratory Standards. 1997. Approved Guideline, M29-A. Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue. NCCLS, Wayne, PA.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
3. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories. HHS Publication (CDC), 4th edition. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 021-0045.

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BENEX Limited