BD Control Set for the BD ProbeTec™

Chlamydia trachomatis/Neisseria gonorrhoeae (CT/GC) Q^X Amplified DNA Assays

| | | | | | | | |
|--------------------------------|--------------|----------|-----------------|----------------|----------|-----------------------------------------------------------------------------------|------------------------------------------------|
| English: | pages | 1 | Deutsch: | Seiten | 3 |  | 8090434(04) |
| Français : | pages | 2 | Español: | páginas | 4 | 0086 | 441125 |

Свържете се с местния представител на BD за инструкции. / Kontaktiraj lokalnog predstavnika BD za upute. / Pokynu vám poskytné místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласңыз. / Naudojite se svom lokalnom predstavniku kompanije BD įgaliojoto atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Instrukcie ziskate u miestneho zástupcu spoločnosti BD. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin.

INTENDED USE

The Control Set for the **BD ProbeTec™** *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Q^X Amplified DNA Assays contains Positive and Negative Controls that are intended for Quality Control use with the **BD ProbeTec™** *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Q^X Amplified DNA Assays when tested with the **BD Viper™** System.

SUMMARY AND EXPLANATION

Quality control must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory’s standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices. The **BD ProbeTec** *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Positive and Negative Controls must be included with each **BD Viper** run. Controls must be positioned according to the **BD Viper** Instrument User’s Manual. The CT/GC Q^X Positive Control monitors for substantial reagent failure only. The CT/GC Q^X Negative Control monitors for reagent and/or environmental contamination. The CT/GC Q^X Positive Control contains *in vitro* cloned plasmids containing target regions for the CT and GC Q^X Amplified DNA Assays. These controls may be used for internal quality control or users may develop their own internal quality control material.¹ Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations. Refer to CLSI C24-A3 for additional guidance on appropriate internal quality control testing practices.¹ The Positive Control contains approximately 2400 copies per mL of pCTB4 and pGCint3 linearized plasmids.

REAGENTS

Materials Provided:

Each Control Set for the **BD ProbeTec** *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Q^X Amplified DNA Assays contains: 24 CT/GC Q^X Positive Control Tubes containing approximately 2400 copies per mL of pCTB4 and pGCint3 linearized plasmids in carrier nucleic acid, and 24 CT/GC Q^X Negative Control Tubes containing carrier nucleic acid alone. The concentrations of the pCTB4 and pGCint3 plasmids are determined by UV spectrophotometry.

Materials Required But Not Provided:

BD ProbeTec *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Q^X Amplified DNA Assays and accessories, Nitrile gloves.

Storage and Handling Requirements: Reagents may be stored at 2 – 33 °C. Do not freeze.

WARNINGS AND PRECAUTIONS

- For *in vitro* Diagnostic Use.
- The Control Set for the **BD ProbeTec** *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Q^X Amplified DNA Assays is used to evaluate the performance of the **BD ProbeTec** *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Q^X Amplified DNA Assays and should not be used with other methods.
- For additional warnings, cautions, and notes specific to the **BD Viper**, consult the **BD Viper** Instrument User’s Manual.
- Because the CT/GC Q^X Positive Control is used for both CT Q^X and GC Q^X testing, correct positioning of the microwell strips is important for final results reporting.
- Do not re-hydrate the controls prior to loading in the **BD Viper** Lysing Rack.
- Do not eat, drink or smoke in areas where specimens or kit reagents are handled.

QUALITY CONTROL PREPARATION

The CT/GC Q^X Positive and Negative Controls (i.e., assay controls) do not require addition of fluid by the user prior to loading in the **BD Viper** Specimen Rack.

Interpretation of Quality Control Results

The CT/GC Q^X Positive and Negative Controls must test as positive and negative, respectively, in order to obtain patient results. If controls do not perform as expected, the run is considered invalid and patient results will not be reported by the instrument. If either of the controls does not provide the expected results, repeat the entire run using a new set of controls, new extraction tubes, new extraction reagent trough, and new microwells. If the repeat QC does not provide the expected results, contact BD Technical Services. Refer to the **BD ProbeTec** *Chlamydia trachomatis* and/or *Neisseria gonorrhoeae* (CT/GC) Q^X Amplified DNA Assays package inserts (Cat. No. 441124, 441126, 442959 and 442842) for further interpretation of the quality control results.

LIMITATIONS OF THE PROCEDURE

The CT/GC Q^X Control Set may not adequately control for specimen processing. To prepare and test specimen processing controls, refer to the **BD ProbeTec** *Chlamydia trachomatis* or *Neisseria gonorrhoeae* (CT/GC) Q^X Amplified DNA Assays package inserts (Cat. No. 441124, 441126, 442959 and 442842).

REFERENCE

- Clinical and Laboratory Standards Institute. 2006. Approved Guideline C24-A3. Statistical quality control for quantitative measurement procedures: principles and definitions, 3rd ed. CLSI, Wayne, PA.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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BD, BD Logo, BD ProbeTec, and BD Viper are trademarks of Becton, Dickinson and Company. © 2014 BD.
EC REP Benex Limited