
Revisions

BALTSO0191 Version 7.0 Template 4

Rev from	Rev to	Job #
01	02	7015-04

NOTES:

1. BD Catalog Number: 241050, 240939
2. Blank (Sheet) Size: Length: 18" Width: 30"
3. Number of Pages: 10 Number of Sheets: 1
4. Page Size: Length: 18" Width: 6" Final Folded Size: 2.25" x 4"
5. Ink Colors: No. of Colors: 3 PMS#: 2755 Blue, PMS 032 Red, PMS Standard Black
6. Printed two sides: Yes  No
7. Style (see illustrations below): # 4

8. See specification control no. 8810131 for material information.
9. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

REVISED BY  
By Bryan Smith at 9:45 am, Nov 03, 2014

Label Design <i>[Signature]</i> Date <i>11/3/14</i>	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	<b>BD</b> Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer <i>[Signature]</i> Date <i>11/03/14</i>		
Checked By <i>[Signature]</i> Date <i>11-04-14</i>	Category and Description	Sheet: 1 of 11
Part Number: 8810131JAA	Package Insert, Francisella Tularensis Antigen & Antiserum	Scale: N/A
		<b>A</b>

# BD Francisella Tularensis Antigen and Antiserum for Febrile Antigen Agglutination Tests

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## INTENDED USE

Febrile Antigens are used in agglutination tests as an aid in the diagnosis of certain febrile diseases such as tularemia. The patient's serum is tested directly for homologous antibodies by either a slide or tube agglutination test. These tests are qualitative and semi-quantitative. The rapid slide test is used primarily as a screening procedure especially useful when large numbers of sera must be examined. The tube test should be used to confirm positive results obtained by the slide test.

Francisella Tularensis Control Antiserum is used in agglutination tests in conjunction with the Febrile Antigen as a positive control.

## SUMMARY AND EXPLANATION

Serological diagnosis of patients suspected of having infectious diseases characterized by persistent fever is dependent upon demonstration of an agglutination reaction between the appropriate antigen and the patient's serum.

The natural response to the invasion by pathogenic organisms is the production of antibodies. This immune response is highly individualized and in addition to the host's physiological status and genetic capabilities, a number of other factors are involved with the production of antibodies to a particular stimulus. These include the antigenicity of the organism, the total amount introduced to the host and the route of introduction, and whether the host has had previous exposure to the organism. These factors will determine the rate of antibody formation, the amount of antibodies produced and their persistence in the circulatory system.

## PRINCIPLES OF THE PROCEDURE

When the human body is invaded by any pathogenic microbiological agent, a variety of antibodies are formed. Among these antibodies are the agglutinins. Serum containing specific agglutinins in combination with homologous antigens under properly controlled conditions is capable of causing a visible agglutination. The degree to which agglutination occurs depends upon the concentration of the antigen, the number of antibodies present, the composition of the salt solution and the temperature.<sup>1</sup>

Sera from normal patients may show positive agglutination with febrile antigens due to previous immunization, past infection or the presence of antibodies to related antigens. In general, the titers found in these instances will be lower and remain at a constant level. Titers detected as a result of active infection or recent immunization with an organism containing homologous antigens will usually be higher and tend to rise over a period of time.

It is, therefore, necessary to evaluate two or more serum samples taken at 3- to 5-day intervals after the onset of the disease. A progressive increase in titer is the prime evidence of recent infection or immunization.<sup>2</sup>

## REAGENTS

Francisella Tularensis Antigen is a suspension of nonviable bacterial cells containing 0.5% formaldehyde.

Francisella Tularensis Control Antiserum is a polyclonal rabbit antiserum containing 30% – 50% glycerin as a preservative.

## Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>3-6</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.

## Cat. No. Description

Febrile Antigen

241050 BD Francisella Tularensis Antigen, one 5 mL vial.

## WARNING

**H317** May cause an allergic skin reaction. **H351** Suspected of causing cancer. **P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P302 + P352** IF ON SKIN: Wash with plenty of water. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

**Storage Instructions:** On receipt, store at 2 – 8 °C. Francisella Tularensis Antigen and Francisella Tularensis Control Antiserum are ready for use as supplied. Do not open until ready to use. Protect from evaporation and contamination. The expiration date on each vial applies to product in intact container stored as directed.

## SPECIMEN COLLECTION

Using aseptic techniques, collect 5 – 10 mL of blood in a blood tube without anticoagulant (i.e., **BD Vacutainer™** Brand Glass Evacuated Tube). This type of tube is recommended because the collection, clotting and centrifugation of the specimen can all be done without transferring the blood. Allow the blood to clot for 20 – 30 min. Loosen the clot from the side of the container with a clean glass rod or applicator stick, stopper the tube and centrifuge. Draw off the serum and refrigerate at 2 – 8 °C until testing, or immediately freeze the serum at -15 °C or below if the test is delayed more than 4 h.

## PROCEDURE

**Materials Provided:** Francisella Tularensis Antigen and Francisella Tularensis Control Antiserum.

### Materials Required But Not Provided:

- 3" x 6" glass slide (circled)
- Pipettes
- Toothpicks, glass rods or other applicator sticks
- 0.85% saline
- Water bath (37 °C)
- Tubes

### Test Procedure

**NOTE: ALL MATERIALS AND EQUIPMENT MUST BE AT ROOM TEMPERATURE AT THE TIME OF TEST PERFORMANCE.**

Directions for use should be followed carefully.

### Slide Test

- Using a 3" x 6" slide with circles (circles can be created using a wax pencil), dispense serum (patient or control) onto slide. Dispense serum samples into separate circles as follows:

Circle	Serum
1	0.08 mL
2	0.04 mL
3	0.02 mL
4	0.01 mL
5	0.005 mL

- Negative Control: Into a separate circle, dispense 1 drop of 0.85% sterile saline solution onto the 3" x 6" slide.
- Invert the bottle several times to mix antigen.
- Add one drop of antigen to each circle using the dropper supplied with the vial.
- Using a new toothpick or applicator for each circle, mix circle contents thoroughly.
- Rock the slide for 1 min near a light source and observe for agglutination.

### Tube Test

- Prepare dilutions of patient and control sera.

#### Patient Sera

- Add 0.5 mL of 0.85% saline to each of 8 tubes.
- Place 0.8 mL of saline in a separate tube and add 0.2 mL of serum. Mix well with a separate pipette and add 0.5 mL from this tube to the first tube of 8 tubes prepared in step "a".
- Using a clean pipette for each patient serum, mix and transfer 0.5 mL from tube 1 to tube 2 and serial dilute through tube 7. Discard 0.5 mL from tube 7. Tube 8 is the antigen control.

**Note:** Reading References may be set up as described in step 3 of RESULTS, "Reading References."

#### Control Serum

Prepare serial dilutions of the control serum in 0.5 mL amounts in tubes in the following manner:

- Place 8 tubes in a rack.
  - Pipette 0.9 mL of 0.85% saline into the first tube and 0.5 mL into each of the remaining tubes.
  - Add 0.1 mL of the serum to the first tube containing 0.9 mL of saline solution.
  - Mix well with a pipette and transfer 0.5 mL of the first tube to the second tube. Mix thoroughly.
  - Continue carrying the 0.5 mL of the serum dilution through tube 7. Discard 0.5 mL from tube 7 after mixing thoroughly. Tube 8 is the antigen control tube.
- Prepare Antigen Dilution
    - Dilute Francisella antigen to one part antigen in forty-nine parts of 0.85% saline (1:50), e.g., 0.25 mL antigen and 12.25 mL saline. Mix well.
  - Mix Antigen and Serum

Add 0.5 mL of diluted antigen to each tube containing patient sera or control serum, as well as to the antigen control tube. Shake tubes for a period of 10 sec. The resultant dilutions are 1:20 through 1:1280, respectively.
  - Incubate all tubes for 20 h in a 37 °C water bath.

## User Quality Control

At time of use, apply both positive and negative controls to check performance of the antigen, techniques and methodology.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

### Slide test

1. Record Agglutination		
Degree of Agglutination		
100 %	background clear to slightly hazy	4+
75 %	background slightly cloudy	3+
50 %	background moderately cloudy	2+
25 %	background very cloudy	1+
None		Negative

### 2. Determination of Antibody Titer

The titer is the reciprocal of the highest dilution producing at least 2+ agglutination.

Serum, mL	Correlation Dilution
0.08	1:20
0.04	1:40
0.02	1:80
0.01	1:160
0.005	1:320

### Tube Test

#### 1. Record Agglutination

Very gently agitate each tube and record agglutination as follows:

100 %	Clear liquid above agglutinins	4+
75 %	Slightly cloudy	3+
50 %	Moderately cloudy	2+
25 %	Very Cloudy	1+
Negative	Very cloudy with no agglutinins	

#### 2. Measure Titer

The titer is the reciprocal of the highest final dilution producing at least 2+ agglutination.

Tube Number	Dilution
1	1:20
2	1:40
3	1:80
4	1:160
5	1:320
6	1:640
7	1:1280

#### 3. Reading References

The following Reading References table will be helpful in determining how the liquid portion of the tube test should appear for differing degrees of agglutination. This is for the liquid portion only.

Antigen, mL	Saline, mL	Liquid Reading
0.5 mL	0.5 mL	Negative
0.25 mL	0.75 mL	2+
0.125 mL	0.875 mL	3+

**NOTE:** If auto-agglutination of the antigen is suspected, prepare the "Reading References" mixtures exactly as indicated above. Agglutinated material in any or all of the tubes indicates that the antigen is unstable and should be discarded.

As in other agglutination tests, no single test or determination can be considered diagnostic. Many variables are involved (See "Limitations of the Procedure"). However, a fourfold rise in antibody titer demonstrable in paired acute and convalescent phase sera collected at approximately five-day intervals is strongly suggestive of a specific presumptive diagnosis. The chart below gives an approximate indication of the significance of serum titers.<sup>2</sup>

Disease	Febrile Antigen	Serum Agglutinins		
		Appear	Maximum	Titer and Significance
Tularemia	Francisella tularensis	7 to 14 days	4 to 8 weeks	1:160 = strongly suggestive

### LIMITATIONS OF THE PROCEDURE

Prozone reactions are possible. Care should be taken to observe agglutination at the highest dilutions if no clumping is seen in lower dilutions.

Patients occasionally fail to develop any serum agglutinins.

Cross reactions may occur between Brucella and Francisella antigens and antisera. Therefore, parallel tests should be run with those antigens. Generally a higher titer is obtained with the homologous antigen.

Cross reactions, previous vaccinations, anamnestic responses,<sup>7</sup> antibiotic therapy, other diseases known or unknown, prozones, and autoagglutinins, as well as other factors, may affect these tests.

Consider the stage of disease when the specimen is collected. It is possible that peak titers occur only during convalescence.

For more complete information, consult appropriate references.<sup>8-10</sup>

These products are not intended as a substitute for culture. An appropriate attempt should be made to recover and identify the etiologic organism.

Only the product catalog numbers listed in this insert are to be used together. Do not use with other product catalog numbers.

### PERFORMANCE CHARACTERISTICS

- The positive control antiserum should produce 2+ or greater agglutination at 1:80 in the slide or tube test when tested with the corresponding homologous febrile antigen.

<b>Febrile Antigen</b>	<b>Homologous Control Antiserum</b>
Francisella tularensis	Francisella tularensis

- The negative control should show no agglutination with the Francisella Antigen.

### AVAILABILITY

#### Cat. No. Description

Febrile Antigen

241050 **BD** Francisella Tularensis Antigen, one 5 mL vial.

Febrile Control Antiserum

240939 **BD** Francisella Tularensis Antiserum, one 5 mL vial.

### REFERENCES

- Davidsohn, I, and J.B. Henry (ed.). 1969. Todd-Sanford clinical diagnosis by laboratory methods, 14th ed. W.B. Saunders Company, Philadelphia.
- Sonnenwirth, A.C. 1970. Gradwohl's clinical laboratory methods and diagnosis, 7th ed., p. 1482. The C.V. Mosby Co., St. Louis.
- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
- U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
- Foomen, J., and H.R. Morgan. 1954. Amer. J. Med. Soc. 228:520.
- Ormsbee, R.A. 1985. Rickettsiae, p. 845-855. In E.H. Lennette, A. Balows, W. J. Hausler, Jr., and H.J. Shadomy (ed.), Manual of clinical microbiology, 4th ed. American Society for Microbiology, Washington, D.C.
- Kolmer, J.A., E.H. Spaulding, and H.W. Robinson. 1951. Approved laboratory technic, 5th ed. Appleton Century-Crofts, New York.
- Sack, R.B. 1986. Serologic tests for the diagnosis of enterobacterial infections, p. 359-362. In N.R. Rose, H. Friedman, and J.L. Fahey (ed.), Manual of clinical laboratory immunology, 3rd ed. American Society for Microbiology, Washington, D.C.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Benex Limited