BD TB Stain Kits and Reagents

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TB Stain Kit K TB Carbofuchsin KF TB Decolorizer TB Brilliant Green K	For staining mycobacteria by the Kinyoun (cold) acid-fast procedure.	1 1 x 250 mL 1 x 250 mL 1 x 250 mL	Cat. No. 212522
TB Stain Kit ZN TB Carbofuchsin ZN TB Decolorizer TB Methylene Blue	For staining mycobacteria by the Ziehl-Neelsen (hot) acid-fast procedure.	1 1 x 250 mL 1 x 250 mL 1 x 250 mL	212520
TB Fluorescent Stain Kit M TB Auramine M TB Decolorizer TM TB Potassium Permanganate	For staining mycobacteria by the Morse, Blair, Weiser and Sproat fluorescent procedure.	1 1 x 250 mL 1 x 250 mL 1 x 250 mL	212519
TB Fluorescent Stain Kit T TB Auramine-Rhodamine T TB Decolorizer TM TB Potassium Permanganate	For staining mycobacteria by the Truant, Brett and Thomas fluorescent procedure.	1 1 x 250 mL 1 x 250 mL 1 x 250 mL	212521
TB Auramine M	For staining mycobacteria by the Morse, Blair, Weiser and Sproat fluorescent procedure.	4 x 250 mL	212514
TB Auramine-Rhodamine T	For staining mycobacteria by the Truant, Brett and Thomas fluorescent procedure.	4 x 250 mL	212515
TB Brilliant Green K	For staining mycobacteria by the Kinyoun (cold) and Ziehl-Neelsen (hot) acid-fast procedures.	4 x 250 mL	212523
TB Carbofuchsin KF	For staining mycobacteria by the Kinyoun (cold) acid-fast procedure.	4 x 250 mL	212518
TB Carbofuchsin ZN	For staining mycobacteria by the Ziehl-Neelsen (hot) acid-fast procedure.	4 x 250 mL	212511
TB Decolorizer	For staining mycobacteria by the Kinyoun (cold) and Ziehl-Neelsen (hot) acid-fast procedures.	4 x 250 mL	212517
TB Decolorizer TM	For staining mycobacteria by the Truant, Brett and Thomas and the Morse, Blair, Weiser and Sproat fluorescent procedures.	4 x 250 mL	212512
TB Methylene Blue	For staining mycobacteria by the Kinyoun (cold) and Ziehl-Neelsen (hot) acid-fast procedures.	4 x 250 mL	212516
TB Potassium Permanganate	For staining mycobacteria by the Truant, Brett and Thomas and the Morse, Blair, Weiser and Sproat fluorescent procedures.	4 x 250 mL	212513

INTENDED USE

TB Stain Kits are used to stain smears prepared from specimens and cultures suspected of containing mycobacteria for early presumptive diagnosis of mycobacterial infection and characterization of bacterial isolates.

SUMMARY AND EXPLANATION

The microscopic staining technique is one of the earliest methods devised for detecting the tubercle bacillus and it remains a standard procedure.¹ The unique acid-fast characteristic of mycobacteria makes the staining technique valuable in early presumptive diagnosis as well as in providing information about the number of acid-fast bacilli present.

Fluorescent microscopy offers numerous advantages over classic methods for detecting mycobacteria because of its speed and simplicity, the ease of examining the slide and the reliability and superiority of the method.²

The TB Stain Kit K uses the Kinyoun (cold) acid-fast procedure.³

The TB Stain Kit ZN uses the Ziehl-Neelsen (hot) acid-fast procedure.⁴

The TB Fluorescent Stain Kit M uses the auramine O acid-fast fluorescent procedure.⁵

The TB Fluorescent Stain Kit T uses the auramine-rhodamine acid-fast fluorescent procedure described by Truant, Brett and Thomas.⁶

PRINCIPLES OF THE PROCEDURE

A major characteristic of mycobacteria is that they are "acid-fast," i.e., once stained with the aniline dye, basic fuchsin, they are difficult to decolorize, and retain the red color even when treated with a mixture of acid and alcohol.

In the Ziehl-Neelsen stain, methylene blue is the counterstain. In the Kinyoun stain, brilliant green provides the necessary contrast.

Fluorescent dyes possess the property of emitting visible radiation when stimulated by shorter wavelength ultra violet light, i.e., a mercury vapor lamp or a strong, blue appropriately filtered light source.⁷

Auramine O and rhodamine B are reported to bind to mycolic acids in the mycobacterial cell wall.⁸ Unbound dye is removed by the acid-alcohol decolorizer solution. Potassium permanganate serves as a counterstain to suppress nonspecific background fluorescence.

REAGENTS

Approximate Formula* Per Liter

TB Stain Kit K

TB Carbofuchsin KF

Basic Fuchsin	15.0 g
Phenol, USP	45.0 g
Isopropanol	200.0 mL
Ethanol	50.0 mL
Distilled Water	750.0 mL

TB Decolorizer

Hydrochloric Acid	30.0 mL
Ethanol/Methanol	970.0 mL

TB Brilliant Green K

Brilliant Green	2.0 g
Sodium Hydroxide	0.02 g
Distilled Water	1000.0 mL

TB Stain Kit ZN

TB Carbofuchsin ZN

Basic Fuchsin	1.7 g
Phenol, USP	50.0 g
Isopropanol	95.0 mL
Distilled Water	905.0 mL

TB Decolorizer

Hydrochloric Acid	30.0 mL
Ethanol/Methanol	970.0 mL

TB Methylene Blue

Methylene Blue, USP	2.4 g
Ethanol/Methanol	300.0 mL
Distilled Water	700.0 mL

TB Fluorescent Stain Kit M

TB Auramine M

Auramine O	2.0 g
Phenol, USP	4.0 g
Glycerine, USP	100.0 mL
Isopropanol	250.0 mL
Distilled Water	650.0 mL

TB Decolorizer TM

Hydrochloric Acid	5.0 mL
Isopropanol	700.0 mL
Distilled Water	300.0 mL

TB Potassium Permanganate

Potassium Permanganate	5.0 g
Distilled Water	1000.0 mL

TB Fluorescent Stain Kit T

TB Auramine-Rhodamine T

Auramine O	12.0 g
Rhodamine B	6.0 g
Phenol, USP	80.0 g
Glycerine, USP	600.0 mL
Isopropanol	140.0 mL
Distilled Water	260.0 mL

TB Decolorizer TM

Hydrochloric Acid	5.0 mL
Isopropanol	700.0 mL
Distilled Water	300.0 mL

TB Potassium Permanganate

Potassium Permanganate	5.0 g
Distilled Water	1000.0 mL

*Adjusted and/or supplemented as required to meet performance criteria.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁹⁻¹² and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Biosafety Level 2 practices and procedures, containment equipment and facilities are required for non-aerosol-producing manipulations of clinical specimens such as preparation of acid-fast smears. All aerosol-generating activities must be conducted in a Class I or II biological safety cabinet. Biosafety Level 3 practices, containment equipment and facilities are required for laboratory activities in the propagation and manipulation of cultures of *M. tuberculosis* and *M. bovis*. Animal studies also require special procedures.¹¹

These reagents are harmful or fatal if swallowed and can cause eye irritation or blindness. In the event of eye contact, flush thoroughly with an eye wash system or tap water for at least 15 min and seek medical advice. The dye solutions and decolorizer solution are highly flammable and vapors may be harmful; use away from open flame. Avoid contact with skin, eye or mucous membranes. Use all reagents in a well ventilated area, avoiding prolonged inhalation, skin contact and ingestion.

TB Carbofuchsin ZN: Toxic in contact with skin and if swallowed. Causes burns. Keep container tightly closed. Avoid contact with skin and eyes. Wear suitable protective clothing.

TB Decolorizer TM: Causes burns. Keep container tightly closed and in a well ventilated place. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately.

TB Potassium Permanganate: Irritating to eyes, respiratory system and skin. Keep container tightly closed. Avoid contact with skin and eyes. Wear suitable protective clothing.

TB Auramine M: Flammable. Irritating to eyes, respiratory system and skin. After contact with skin, wash immediately with plenty of water. Keep container tightly closed. Keep away from sources of ignition - No Smoking. Wear suitable protective clothing.

TB Auramine-Rhodamine T: Flammable. Toxic in contact with skin and if swallowed. Causes burns. Possible risks of irreversible effects. Keep container tightly closed. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately.

TB Decolorizer: Harmful by inhalation and if swallowed. Keep container tightly closed and in a well ventilated place. Keep away from sources of ignition - No smoking. Wear suitable protective clothing.

TB Methylene Blue: Harmful by inhalation and if swallowed. Irritating to eyes, respiratory system and skin. Possible risks of irreversible effects. Keep container tightly closed. Keep away from sources of ignition - No Smoking. Wear suitable protective clothing and gloves.

TB Carbofuchsin KF: Harmful in contact with skin and if swallowed. Irritating to eyes and skin. Possible risks of irreversible effects. Keep container tightly closed. Keep away from sources of ignition - No Smoking. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing and gloves.

TB Carbofuchsin ZN

Phenol

Basic Fuchsin

DANGER

H314 Causes severe skin burns and eye damage. **H341** Suspected of causing genetic defects. **H350** May cause cancer.

P260 Do not breathe dust/fume/gas/mist/vapours/spray. **P303/361/353** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. **P305/351/338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

TB Carbofuchsin KF

Phenol

Basic Fuchsin

DANGER

H226 Flammable liquid and vapour. **H303** May be harmful if swallowed. **H314** Causes severe skin burns and eye damage. **H333** May be harmful if inhaled. **H341** Suspected of causing genetic defects. **H350** May cause cancer.

P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking. **P260** Do not breathe dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303/361/353** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. **P304/312** IF INHALED: Call a POISON CENTER or doctor/physician if you feel unwell. **P403/235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

TB Decolorizer Stain Set K

Methanol

DANGER

H225 Highly Flammable liquid and vapour. **H303** May be harmful if swallowed. **H333** May be harmful if inhaled. **H371** May cause damage to organs.

P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking. **P260** Do not breathe dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303/361/353** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. **P304/312** IF INHALED: Call a POISON CENTER or doctor/physician if you feel unwell. **P403/235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

TB Decolorizer Stain Set ZN

Methanol

ANGER

H225 Highly Flammable liquid and vapour. **H371** May cause damage to organs.

P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking. **P260** Do not breathe dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303/361/353** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. **P403/235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Gram Decolorizer

ANGER

H225 Highly Flammable liquid and vapour. **H319** Causes severe eye irritation. **H336** May cause drowsiness or dizziness.

P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking. **P303/361/353** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. **P403/235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Methylene Blue

ARNING

H226 Flammable liquid and vapour.

P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303/361/353** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. **P403/235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

TB Decolorizer TM

ANGER

H225 Highly Flammable liquid and vapour. **H314** Causes severe skin burns and eye damage. **H336** May cause drowsiness or dizziness.

P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking. **P260** Do not breathe dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303/361/353** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. **P403/235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

TB Auramine M

WARNING

H226 Flammable liquid and vapour. **H319** Causes severe eye irritation. **H333** May be harmful if inhaled. **H336** May cause drowsiness or dizziness.

P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303/361/353** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. **P403/235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

TB Auramine – Rhodamine T

Phenol

Auramine O

DANGER

H226 Flammable liquid and vapour. **H303** May be harmful if swallowed. **H314** Causes severe skin burns and eye damage. **H333** May be harmful if inhaled. **H341** Suspected of causing genetic defects. **H351** Suspected of causing cancer.

P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking. **P260** Do not breathe dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303/361/353** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. **P304/312** IF INHALED: Call a POISON CENTER or doctor/physician if you feel unwell. **P403/235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Brilliant Green

WARNING

H302 Harmful if swallowed.

TB Potassium Permanganate

WARNING

H315 Causes skin irritation. **H319** Causes serious eye irritation.

P280 Wear protective gloves/protective clothing/eye protection/face protection. **P305/351/338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Storage: Store TB Stain Kits and Reagents at 15 - 30 °C. Reagents that have been removed from the packing carton should be stored in the dark. The expiration date applies to the product in its intact container when stored as directed.

SPECIMEN COLLECTION AND PREPARATION

Apply a thin smear of the specimen directly on a clean microscope slide. Sputum (direct specimen or concentrated sputum), spinal or other body fluid (catheterized urine), homogenized tissue and impression slides are satisfactory specimens for smears.¹

Allow smear to dry.

Fix the smear to the slide by passing the slide through a low flame 2 - 3 times, avoiding excessive heat, or use an electric slide warmer at 65 to 75°C for at least 2 h.¹³

PROCEDURE

Materials Provided*: TB Stain Kit K, TB Stain Kit ZN, TB Quick Stain Kit, TB Fluorescent Stain Kit M, TB Fluorescent Stain Kit T.

**Individual reagents available separately.*

Materials Required But Not Provided: Microscope slides - new or cleaned in acid dichromate solution, staining rack, filter paper and microscope with oil immersion lens or fluorescent microscope.

Test Procedure

KINYOUN STAIN

TB Stain Kit K

Instructions:

1. Place slides on a staining rack and flood with TB Carbofuchsin KF for 4 min. Do not heat.
2. Wash gently in running water.
3. Decolorize with TB Decolorizer for 3 - 5 s.
4. Wash gently in running water.
5. Counterstain with either TB Brilliant Green K or TB Methylene Blue (available separately) for 30 s.
6. Wash gently in running water.
7. Air dry. If using TB Methylene Blue, dry over gentle heat.

ZIEHL-NEELSEN STAIN

TB Stain Kit ZN

Instructions:

1. Place slides on a staining rack and place a piece of filter paper, larger than the size of the smear, on each slide. Flood with TB Carbofuchsin ZN. Heat gently to steaming and allow to steam for 5 min.
2. Remove the filter paper. Wash gently in running water.
3. Decolorize with TB Decolorizer with two changes of reagent for 1 - 2 min until no more red color appears in washing.
4. Wash gently in running water.
5. Counterstain with either TB Methylene Blue or TB Brilliant Green K for 30 s.
6. Wash gently in running water.
7. Dry over gentle heat.

MORSE STAIN

TB Fluorescent Stain Kit M

Instructions:

1. Place slides on a staining rack and flood with TB Auramine M for 15 min.
2. Wash gently in running water.
3. Decolorize with TB Decolorizer TM for 30 - 60 s.
4. Wash slides gently in running water.
5. Counterstain with TB Potassium Permanganate for 2 min.
6. Wash gently in running water.
7. Air dry.
8. Examine under a microscope fitted with filter sets as described by Chapin and Murray.¹⁴

TRUANT STAIN

TB Fluorescent Stain Kit T

Instructions:

1. Place slides on a staining rack and flood with TB Auramine-Rhodamine T that has been thoroughly shaken prior to use. Leave undisturbed for 20 - 25 min at room temperature.
2. Wash gently in running water.
3. Decolorize with TB Decolorizer TM for 2 - 3 min.
4. Wash gently in running water.
5. Counterstain with TB Potassium Permanganate for 4 - 5 min.
6. Wash gently in running water.
7. Blot lightly. Dry in air or very gently over a flame.
8. Examine under a microscope fitted as described by Truant et al.,⁶ or Chapin and Murray.¹⁴

User Quality Control: It is recommended that controls be run each day of use using **BBL™** AFB Slide (Cat. No. 231391) or known acid-fast and nonacid-fast microorganisms.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

Organism	Expected Results
<i>Mycobacterium tuberculosis</i> H37Ra ATCC™ 25177	Positive
<i>Bacillus subtilis</i> ATCC 6633	Negative

LIMITATIONS OF THE PROCEDURE

A positive staining reaction provides presumptive evidence of the presence of mycobacteria in the specimen. A negative staining reaction does not necessarily indicate that the specimen will be culturally negative for mycobacteria. For positive identification of mycobacteria, other methods must be employed (cultural, molecular, etc.).

EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS^{1-8,14}

	Acid-fast Bacilli	Nonacid-fast Bacilli
TB Stain Kit K	dark red to pink	green
TB Stain Kit ZN	dark pink to red	blue
	Mycobacteria	
TB Fluorescent Stain Kit M	bright yellow-green fluorescence	
TB Fluorescent Stain Kit T	reddish-orange fluorescence	

REFERENCES

- Nolte, F.S. and B. Metchock. 1995. *Mycobacterium*, p. 400 – 436. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
- Tenover, F.C., J.T. Crawford, R.E. Huebner, L.J. Geiter, C.R. Horsburgh, Jr., and R.C. Good. 1993. The resurgence of tuberculosis: is your laboratory ready? J. Clin. Microbiol. 31: 767 – 770.
- Ebersole, L.L. 1994. Acid-fast stain procedures, p.3.5.1. – 3.5.11. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
- Murray, R.G.E., R.N. Doetsch and C.F. Robinow. 1994. Determinative and cytological light microscopy, p. 21 – 41. In P. Gerhardt, R.G.E. Murray, W.A. Wood and N.R. Krieg (ed.), Methods for general and molecular bacteriology, American Society for Microbiology, Washington, D.C.
- Heifets, L.B. and R.C. Good. 1994. Current laboratory methods for the diagnosis of tuberculosis, p. 85 – 110. In B.R. Bloom (ed.), Tuberculosis: pathogenesis, protection, and control, American Society for Microbiology, Washington, D.C.
- Truant, J.P., W.A. Brett and W. Thomas, Jr. 1962. Fluorescence microscopy of tubercle bacilli stained with auramine and rhodamine. Henry Ford Hosp. Med. Bull. 10: 287 – 296.
- Cernoch, P.L., R.K. Enns, M.A. Saubolle and R.J. Wallace, Jr. 1994. Cumitech 16A, Laboratory diagnosis of the mycobacterioses. Coord. ed. A.S. Weissfeld. American Society for Microbiology, Washington, D.C.
- Baron, E.J. and S.M. Finegold. 1990. Bailey & Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Company, Baltimore, MD.
- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17: 53-80.
- U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
- Kent, P.T. and G.P. Kubica, 1985. Public Health Mycobacteriology. A guide for the level III laboratory, p. 58, U.S. Dept. of Health and Human Services, Centers for Disease Control, Atlanta Ga.
- Chapin, K.C. and P.R. Murray. 1999. Stains, p. 1674 – 1686. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Becton Dickinson Pty Ltd.
Benex Limited