# **BD ProbeTec™ Q<sup>x</sup> Collection Kit for Endocervical or Lesion Specimens**

English: pages 1 – 2    Français : pages 8 – 9  
Dansk: side 3 – 4    Norsk: sider 9 – 11  
Deutsch: Seiten 4 – 6    Svenska: sidan 11 – 12  
Español: páginas 6 – 7

L009847(01)  
2012-12

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcję użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндің жергілікті БД өкіліне жүгініп нұсқа алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USES

The **BD ProbeTec™ Q<sup>x</sup> Collection Kit for Endocervical or Lesion Specimens** is used to collect and transport patient endocervical specimens to the laboratory for testing with the **BD ProbeTec CT/GC Q<sup>x</sup> Amplified DNA Assays** and the **BD ProbeTec™ TV Q<sup>x</sup> Amplified DNA Assay** on the **BD Viper™ System** in Extracted Mode. The Collection Kit also may be used to collect and transport external anogenital lesion specimens to the laboratory for testing with the **BD ProbeTec HSV1&2 Q<sup>x</sup> Amplified DNA Assays** on the **BD Viper System** in Extracted Mode.

## SUMMARY AND EXPLANATION

The **BD ProbeTec Q<sup>x</sup> Collection Kit for Endocervical or Lesion Specimens** consists of a sterile cleaning swab with a white shaft, a sterile collection swab with a pink shaft that is scored for breaking and a tube of Q<sup>x</sup> Swab Diluent for specimen transport. The polyester fiber-tipped cleaning swab with the white shaft (not to be used for lesion collection for the HSV1&2 Q<sup>x</sup> Amplified DNA Assays) is provided to clean the cervical os of blood and mucus prior to endocervical specimen collection.† The polyurethane foam-tipped swab with the pink-colored shaft is used to collect the specimen for testing. After the specimen is collected, the cap of the diluent tube is removed, the swab is placed into the Q<sup>x</sup> Swab Diluent tube and the shaft is broken at the score mark. The tube is recapped and transported to the testing facility.

†For the CT/GC Q<sup>x</sup> and TV Q<sup>x</sup> collection

## REAGENTS

**Materials Provided:** Each **BD ProbeTec Q<sup>x</sup> Collection Kit for Endocervical or Lesion Specimens** contains:

- one sterile polyester fiber-tipped cleaning swab with a white shaft\*
- one sterile polyurethane foam-tipped collection swab with a pink shaft\*
- one tube of Q<sup>x</sup> Swab Diluent

The Swab Diluent contains approximately 2 mL of potassium phosphate/potassium hydroxide buffer with DMSO and preservative.

**Materials Required But Not Provided:** Gloves

**Storage Instructions:** Store collection kit at 2 – 33 °C. Do not use beyond the expiration date. Swab sterility is guaranteed if primary swab container is intact.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

1. **Do not pre-wet the swab with the Q<sup>x</sup> Swab Diluent prior to collecting specimen.**
2. Pathogenic microorganisms including hepatitis viruses and Human Immunodeficiency Virus may be present in clinical specimens. "Standard Precautions"<sup>1-4</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
3. Optimal performance of the **BD ProbeTec Q<sup>x</sup> Amplified DNA Assays** requires proper specimen collection, handling and transport.
4. Use only collection and transport systems labeled for use with the **BD ProbeTec Q<sup>x</sup> Amplified DNA Assays**.
5. The polyester fiber-tipped swab with the white shaft is for cleaning the cervical os and should not be used for collecting the patient specimen.
6. Proper labeling should accompany each specimen to the laboratory.
7. The Q<sup>x</sup> Swab Diluent contains dimethyl sulfoxide (DMSO). DMSO is harmful by inhalation, in contact with skin and if swallowed. Keep container in a well-ventilated place. Avoid contact with eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. Wear suitable protective clothing and gloves.
8. When breaking the shaft of the swab, take care to avoid splashing of contents.
9. Dispose of unused Q<sup>x</sup> Swab Diluent tubes according to the Material Safety Data Sheet.
10. Take care to avoid cross-contamination during the specimen handling steps. Specimens can contain high levels of organisms. Ensure that specimen containers do not contact one another, and discard used materials without passing over open containers. If gloves come in contact with specimen, change gloves to avoid cross-contamination.

11. Specimens must be collected and tested before the expiration date of the Q<sup>x</sup> Swab Diluent tube.
12.  **BD ProbeTec** Q<sup>x</sup> Collection Kit for Endocervical or Lesion Specimens is for single use only; reuse may cause a risk of infection and/or inaccurate results.

#### **SPECIMEN COLLECTION AND TRANSPORT**

##### **Endocervical Specimen Collection**

1. Remove the cleaning swab from packaging.
2. Using the polyester fiber-tipped cleaning swab with the white shaft, remove excess blood and mucus from the cervical os.
3. Discard the used cleaning swab.
4. Remove the pink collection swab from the packaging.
5. Insert the collection swab into the cervical canal and rotate for 15 – 30 s.
6. Withdraw the swab carefully. Avoid contact with the vaginal mucosa.
7. Uncap the Q<sup>x</sup> Swab Diluent tube.
8. Fully insert the collection swab into the Q<sup>x</sup> Swab Diluent tube.
9. Break the shaft of the swab at the score mark. Use care to avoid splashing of contents.
10. **Tightly** recap the tube.
11. Label the tube with patient information and date/time collected.
12. Transport to the laboratory.

##### **Endocervical Specimen Swab Storage and Transport**

The Q<sup>x</sup> Swab Diluent tube with swab must be stored and transported to the laboratory and/or test site at 2 – 30 °C within 30 days of collection. Alternatively, the Q<sup>x</sup> Swab Diluent tube may be stored and transported frozen (-20 °C) to the laboratory and/or test site within 180 days of collection.

##### **Lesion Specimen Collection**

1. Open the inner Q<sup>x</sup> Swab packaging and **dispose of the swab with the white shaft**.
2. Remove the pink collection swab from the packaging.
3. Swab the base of the exposed lesion firmly to absorb exudates and cellular material **with the pink swab**.
4. Uncap the Q<sup>x</sup> Swab Diluent tube.
5. Fully insert the pink collection swab into the Q<sup>x</sup> Swab Diluent tube.
6. Break the shaft of the swab at the score mark. Use care to avoid splashing of contents.
7. **Tightly** recap the tube.
8. Label the tube with patient information and date/time collected.
9. Transport to the laboratory.

##### **Lesion Specimen Swab Storage and Transport**

The Q<sup>x</sup> Swab Diluent tube with swab must be stored and transported to the laboratory and/or test site at 2 – 30 °C within 14 days of collection. Alternatively, the Q<sup>x</sup> Swab Diluent tube may be stored and transported frozen (-20 °C) to the laboratory and/or test site within 180 days of collection.

For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

#### **REFERENCES**

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guide for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
3. U. S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5<sup>th</sup> ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/31/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or [www.bd.com/ds](http://www.bd.com/ds).

# **BD ProbeTec Q<sup>x</sup> Collection Kit for Endocervical or Lesion Specimens**

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Becton Dickinson France S.A.S.