
Rev From	Rev To	ECO #	Date	Appr.
1082	0903	1905-03		

Notes

- BD Cat. No. 224051 & 235161
- Blank (Sheet) Size : Length: 14" Width: 21"
 Number of Pages: 12 Number of Sheets: 1
 Page Size: Length 14" Width 3.5" Final Folded Size: 2.33" X 3.5"
- Style (see illustrations below): #4

- See Specification Control No. S1238 for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755 (blue)
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 250 Schilling Circle Cockeysville, MD. 21030-0243 USA	
Proofer	Date			
Checked By	Date			
Part Number: S1238		Category and Description Package Insert Difco USR Antigen and Test Control Serum Set	Sheet: 1 of 13 Scale: 1:1	A

BD Difco™ USR Antigen Difco™ USR Test Control Serum Set

CE S1238
2003/09

English: pages 1 - 3 Italiano: pagine 7 - 10
Français: pages 3 - 5 Español: páginas 10 - 12
Deutsch: Seiten 5 - 7

See symbol glossary at end of insert.
Voir le glossaire des symboles à la fin de la notice.
Consulte el glosario de símbolos al final del prospecto.
Vedere il glossario dei simboli alla fine del foglio illustrativo.
Siehe Symbol-Erklärungen am Ende der Packungsbeilage.
Se symbolförteckningen vid slutet av bipacksedeln.
Se symbolglossaret i slutningen af indlægssedlen.
Consulte o glossário de símbolos no fim do folheto informativo.
Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου.

Kontakt den lokale BD repræsentant for at få instruktioner.
Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες.
Contacte o seu representante local da BD para obter instruções.
Kontakt lokal Becton Dickinson-representant för anvisningar.

INTENDED USE

Difco™ Unheated Serum Reagin (USR) Antigen is recommended for use in the USR test for the detection of reagin, an antibody-like substance, by the qualitative and quantitative slide flocculation tests, as an aid in the diagnosis of syphilis.

The USR Test Control Serum Set is standardized human sera used to control the USR test.

SUMMARY AND EXPLANATION

Treponema pallidum is the causative agent of syphilis. Syphilis is a chronic infection with clinical manifestations that occur in distinct stages. Specific laboratory tests are recommended for the detection of each stage of the disease.

During the primary stage, treponemes present in the characteristic lesion, a chancre, are detectable by dark-field microscopy¹ or by the Direct Fluorescent Antibody Test for *T. pallidum* (DFA-TP). During the secondary stage, most serological tests for syphilis are reactive and treponemes may be found in the lesions by using dark-field microscopy. The latent period, which is asymptomatic, may last for years. Serological tests are usually reactive in the early latent period, but the reactivity of nontreponemal tests decreases during the late latent period. Symptoms of tertiary or late stage syphilis may occur 10-20 years after initial infection. Approximately 71% of patients in the tertiary stage of syphilis have reactive non-treponemal tests.^{2,3} In the tertiary stage, treponemal tests will usually be reactive and are the only basis for diagnosis. The lesions in tertiary syphilis will have few treponemes. Neurosyphilis and late cardiovascular syphilis are complications of tertiary syphilis.

Since the clinical manifestations of syphilis can be confused with other infectious or noninfectious conditions, proper diagnosis must include microscopic examination of lesion material and serological results.²

Difco USR Antigen is a nontreponemal antigen composed of cardiolipin, cholesterol and lecithin. The antigen is a modification of the Venereal Disease Research Laboratory (VDRL) Antigen emulsion in a USR suspending solution. The antigen suspension contains choline chloride, which enhances the reactivity of reagin in unheated serum.⁴ The USR Test is suitable for qualitative as well as quantitative determinations.

Nontreponemal tests measure reagin, an antibody-like substance that can be detected in syphilitic serum. Reagin is also occasionally found in the serum of persons with other acute or chronic diseases. Reactive nontreponemal tests aid in the diagnosis of latent subclinical syphilis and are effective tools for detecting cases in epidemiological investigations. Nontreponemal tests are superior to treponemal tests for following the response to therapy.²

Nontreponemal antigen tests are not entirely specific for syphilis, nor do they have satisfactory sensitivity in all stages of syphilis. Whenever the results of a nontreponemal antigen test disagree with the clinical impression, a treponemal antigen test such as the Fluorescent Treponemal Antibody-Absorption (FTA-ABS)^{2,3} test should be performed. Nontreponemal tests such as the USR, are used to screen patient serum and provide only presumptive evidence for syphilis. Treponemal tests such as the FTA-ABS are used for confirmation.

The likelihood of obtaining a reactive USR test result in various stages of untreated syphilis has been reported as follows:⁵

Stage of Untreated Syphilis	% Reactive
Primary	80
Secondary	100
Latent	95

PRINCIPLES OF THE PROCEDURE

In the USR test procedure, the patient's unheated serum is mixed with a buffered saline suspension of USR Antigen containing cardiolipin, lecithin and cholesterol. The combination of reagin and USR Antigen forms microscopic clumping called flocculation.

REAGENTS

USR Antigen is 0.03% cardiolipin and 0.9% cholesterol dissolved in absolute alcohol with sufficient lecithin (approximately 0.2%) to produce standard reactivity. The antigen is suspended in a solution containing EDTA, choline chloride and phosphate with 0.2% thimerosal as a preservative.^{5,6}

USR Test Control Serum Set contains 3 mL each, of the following lyophilized human sera with 0.02% thimerosal as a preservative:

- Nontreponemal Antigen Reactive Serum
- USR Weakly Reactive Serum
- Nontreponemal Antigen Nonreactive Serum

These reagents are standardized to provide reactive, weakly reactive and nonreactive readings, respectively, when tested according to the USR test procedure.

Warnings and Precautions

1. For *in vitro* Diagnostic Use.
2. **WARNING!** POTENTIAL BIOHAZARDOUS REAGENTS. Each donor unit used in preparation of USR Test Control Serum Set was tested by an FDA-licensed method for the presence of the antibody to Human Immunodeficiency Virus (HIV) and for hepatitis B surface antigen (HbsAg) and found negative (were not repeatedly reactive).
3. Because no test method can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent, specimens and control solutions should be handled as though capable of transmitting an infectious disease. "Standard Precautions"⁷⁻¹⁰ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
4. **USR Test Control Serum Set**
The Packaging of This Product Contains Dry Natural Rubber.

Storage: Store USR Antigen at 2-8°C. If the original 3 mL quantity exceeds what is needed for one testing period, transfer the remainder from the first day's use to one or more aliquot vials and store at 2-8°C.

Store the lyophilized sera in the USR Test Control Serum Set at 2-8°C. Store the reconstituted control sera at 2-8°C or divide into aliquots sufficient for one day of testing and store at -20°C. Do not thaw and refreeze.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Product Deterioration: Expiration date applies to product in its intact container when stored as directed. Do not use if the product is caked, discolored or shows other signs of deterioration.

SPECIMEN COLLECTION AND PREPARATION

Collect a blood specimen by aseptic venipuncture into a clean, dry tube without anticoagulant. After the specimen has clotted, centrifuge the specimen at 1,500-2,000 rpm for 5 min to obtain test serum. Store serum specimens at room temperature for no longer than 4 h; for prolonged storage, keep at 2-8°C for up to 5 days or maintain below -20°C. Serum specimens must be clear, free of hemolysis and show no visible evidence of bacterial contamination, such as turbidity or particulate matter. Refer to appropriate references for additional information on collection of specimens.^{2,11,12}

PROCEDURE

Materials Provided: Difco USR Antigen, Difco USR Test Control Serum Set.

Materials Required But Not Provided: 0.9% saline, nondisposable glass syringe, 1-2 cc, nondisposable calibrated 18-gauge needles without bevel, micropipettor, 50 µL, pipettes, serological, graduated to tip; 1.0 mL, graduated in 1/100 mL; 5.0 mL, graduated in 1/10 mL; 10.0 mL, graduated in 1/10 mL.

Slides, 2 x 3 inches with paraffin or ceramic rings approximately 14 mm in diameter and high enough to prevent spillage during rotation.

Slide holder for 2 x 3 inch slides

Mechanical rotator, adjustable to 180 ± 2 rpm circumscribing a circle 19 mm in diameter on a horizontal plane.

Light microscope with 10X ocular and 10X objective, sterile purified water, absolute alcohol, acetone, timer, applicator sticks.

Reagent Preparation: Equilibrate all materials to room temperature (23-29°C) before performing the tests. Ensure that all glassware and pipettes are clean and free of detergent residues.

USR Antigen is ready to use.

USR Test Control Serum Set: To reconstitute the control sera, add 3 mL sterile purified water and rotate gently to completely dissolve the contents.

Preparation of Specific Glassware:

Syringes with needles:

1. Prerinse with tap water.
2. Soak and hand wash thoroughly in a glassware detergent solution.
3. Rinse with tap water 6-8 times.
4. Rinse with unused purified water.
5. Rinse with absolute alcohol.
6. Rinse with acetone.
7. Air dry until the acetone odor is completely eliminated.
8. Remove needles from syringes for storage.

Ceramic-ringed slides:

1. Prerinse with tap water.
2. Wash with a glassware detergent solution. Avoid prolonged soaking of ceramic-ringed slides in detergent solution because the ceramic rings will become brittle and flake off.

3. Rinse with tap water 3-4 times.
4. Rinse with unused purified water.
5. Wipe dry with a clean lint-free cloth. If cleaned slides do not allow serum to spread evenly within the inner surface of the circle, proceed as follows.
6. Scrub the slides with a nonscratching cleanser.
7. Rinse, dry and polish with a clean, lint-free cloth.

Testing the Accuracy of the Antigen Suspension Needle

1. The accuracy of the test depends on the amount of antigen suspension used. Check the calibration of the needle periodically to ensure delivery of the correct volume of USR Antigen suspension.
2. For the qualitative and quantitative tests on serum, dispense the antigen suspension from a syringe fitted with an 18-gauge needle without bevel that will deliver 45 ± 1 drops of antigen suspension per mL when held vertically.
3. Place the needle on a 1 mL syringe. Fill the syringe with 0.5 mL of USR Antigen suspension. Holding the syringe in a vertical position, count the number of drops delivered in 0.5 mL. The needle is correctly calibrated if 22-23 drops are delivered in 0.5 mL.
4. Adjust or replace the needle if it does not meet this specification. Repeat calibration on the new or adjusted needle.
5. After each day of use, clean the dispensing needle, bottle and syringe by rinsing with water, alcohol and acetone, in that order. Remove the needle from the syringe after cleaning.

USR Qualitative Slide Test on Serum

For reliable and reproducible test results, the USR Antigen suspension, controls and test specimens must be at 23-29°C when tests are performed.

1. Pipette 50 μ L of unheated serum into one ring of a paraffin- or ceramic-ringed slide using a safety pipetting device. Do not use a glass slide with concavities, wells or glass rings. Spread the serum with a circular motion of the pipette tip so that the serum covers the entire inner surface of the paraffin or ceramic ring.
2. Gently swirl to resuspend the USR Antigen and withdraw the desired quantity with a syringe and needle.
3. Hold the syringe and needle containing the USR Antigen suspension in a vertical position. Dispense several drops to clear the needle of air. Add exactly 1 free-falling drop (22 μ L) of antigen suspension to each circle containing serum. Do not allow the needle to touch the serum.
4. Place the slide on the mechanical rotator. Rotate the slide for 4 min at 180 ± 2 rpm. If the environment is dry, cover the slides with a moist, humidifying cover during rotation to prevent excessive evaporation.
5. Immediately after rotating the slide, remove the slide from the rotator and read the test results microscopically, using a 10X ocular and a 10X objective.

USR Quantitative Test on Serum

1. To quantitate serum samples to an endpoint titer, prepare serum dilutions on the slide at 1:1, 1:2, 1:4 and 1:8, as follows.
2. Dispense 50 μ L of 0.9% saline in circles 2-4. Do not spread the saline.
3. Dispense 50 μ L of serum in circles 1 and 2.
4. Mix the saline and the serum in circle 2 by drawing the mixture up and down in the pipette approximately 8 times. Mix gently to prevent bubbles.
5. Transfer 50 μ L from circle 2 (1:2) to circle 3 and mix.
6. Transfer 50 μ L from circle 3 (1:4) to circle 4 (1:8), mix, and then discard 50 μ L from circle 4.
7. Spread the contents of each circle to fill the inner surface of the ring using a clean applicator stick for each circle.
8. Holding the syringe and needle containing the USR Antigen suspension in a vertical position, dispense several drops to clear the needle of air. Then, add exactly 1 free-falling drop (22 μ L) of antigen suspension to each circle containing serum. Do not allow the needle to touch the serum.
9. Place the slide on the mechanical rotator. Rotate the slide for 4 min at 180 ± 2 rpm. If the environment is dry, cover the slides with a moist, humidifying cover during rotation to prevent excessive evaporation.
10. Immediately after rotating the slide, remove the slide from the rotator and read the test results microscopically using a 10X ocular and a 10X objective.
11. If the highest dilution tested (1:8) is reactive, prepare a 1:8 dilution of the test specimen by adding 0.1 mL of serum to 0.7 mL of 0.9% saline. Mix thoroughly. Retest as in steps 1-9, above.

User Quality Control

Test antigen suspension reactivity with control sera (Reactive, Weakly Reactive, and Nonreactive). Use the antigen suspension only if it produces the expected reactivity with the control sera (Reactive, Weakly Reactive, and Nonreactive).

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Qualitative Test

1. Read and record results as follows:
 - Medium to large clumps - Reactive (R)
 - Small clumps - Weakly Reactive (WR)
 - No clumping or very slight roughness - Nonreactive (N)
2. Verify that the control sera results are as expected. If reactions are not as expected, the test is invalid and results cannot be reported.
3. Perform a quantitative test on all serum specimens that produce Reactive, Weakly Reactive or "rough" Nonreactive results, since prozone reactions are occasionally encountered.

Quantitative Test

Report the titer as the highest dilution that produces a Reactive result.

Table 1. Sample quantitative USR test result.

If Reactive results are obtained through dilution 1:32, prepare further two-fold serial dilutions in 0.9% saline (1:64, 1:128 and 1:256) and retest using the quantitative test procedure.

Undiluted (1:1)	Serum Dilutions					Report
	1:2	1:4	1:8	1:16	1:32	
R	WR	N	N	N	N	Reactive, undiluted
R	R	WR	N	N	N	Reactive, 1:2 dilution
R	R	R	WR	N	N	Reactive, 1:4 dilution
WR	WR	R	R	WR	N	Reactive, 1:8 dilution
N (rough)	WR	R	R	R	N	Reactive, 1:16 dilution
WR	N	N	N	N	N	Weakly Reactive, undiluted

INTERPRETATION OF TEST RESULTS

If control results are not as expected, patient results should not be reported.

1. The results of the serum USR test must be confirmed by a treponemal test.
2. The diagnosis of syphilis depends on the results of the USR test, treponemal confirmatory test, clinical signs and symptoms, and risk factors.
3. A Reactive USR test may indicate past or present infection with a pathogenic treponeme. However, it may be a false-positive reaction. A false positive is determined if the confirmatory treponemal test is negative.
4. A Nonreactive USR test with clinical evidence of syphilis may indicate early, primary syphilis, a prozone reaction in secondary syphilis, or late syphilis.
5. A Nonreactive USR test with no clinical evidence of syphilis indicates no current infection, an effectively treated infection, or latent syphilis.
6. A quantitative USR test detects changes in reagin titer. Therefore, a serum specimen showing a four-fold increase in titer on a repeat specimen may indicate an infection, a reinfection or a treatment failure. Likewise, a four-fold decrease during treatment indicates adequate syphilis therapy.

LIMITATIONS OF THE PROCEDURE

1. Positive results obtained by using USR Antigen should not be considered as conclusive evidence that the patient is syphilitic. Conversely, a nonreactive USR test, by itself, does not rule-out the diagnosis of syphilis.
2. A prozone reaction may occur in which reactivity with true positive undiluted serum is inhibited. The prozone phenomenon often gives Weakly Reactive or "rough" Nonreactive results in the qualitative test. Specimens with such nonreactive results must be quantitatively tested.
3. Biological false-positive reactions can occur with nontreponemal tests in persons who abuse drugs, have diseases such as lupus erythematosus, mononucleosis, malaria, leprosy or viral pneumonia, or who have recently been immunized.¹²
4. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
5. If the temperature of the testing area, specimens or reagents is less than 23°C, test reactivity is decreased. If the temperature is greater than 29°C, test reactivity is increased.¹²
6. Test results are unpredictable when testing hemolyzed, contaminated or extremely turbid serum specimens.
7. Test results may be erroneous if speed and time of rotation are not correct.

PERFORMANCE CHARACTERISTICS^{13,14}

The performance of the USR slide test was compared to a rapid plasma reagin (RPR) card test and the VDRL slide test in a study conducted by Pettit, Larsen, Pope, Perryman and Adams.¹³ Quantitative tests were performed on 664 syphilitic sera from treated and untreated patients in different stages of the disease. Also included were 95 false-positive sera and 36 sera from normal donors. The VDRL slide test and the RPR card test were used as a basis for comparison with the USR test. Results of the USR quantitative tests paralleled the VDRL and RPR tests with approximately 84% agreement.¹³

In a second study by Saxena and Ghosh,¹⁴ the performance of the USR slide test was compared to the standard VDRL slide test and the complement fixation (CF) test. Tests were performed on 476 samples of blood. The results obtained are recorded in the table below and show that the USR test can detect 93% of the syphilitic reactive sera.

Number of sera	USR Test Results	VDRL Test Results	CF Test Results
43	+	+	+
2	-	+	-
1	+	-	-
430	-	-	-

AVAILABILITY

Cat. No.	Description
224051	Difco™ USR Antigen, 6 x 3 mL
235161	Difco™ USR Test Control Serum Set, 1 set Set contains: Nontreponemal Antigen Reactive Serum, 3 mL USR Weakly Reactive Serum, 3 mL Nontreponemal Antigen Nonreactive Serum, 3 mL Aliquot Vials, 3 vials

REFERENCES

- Kennedy, E.J. and E.T. Creighton. 1998. Dark field microscopy for the detection and identification of *Treponema pallidum*, p. 119-134. In S.A. Larsen, V. Pope, R.E. Johnson and S.J. Kraus (ed.), Manual of tests for syphilis, 9th ed. American Public Health Association, Washington, D.C.
- Janda, W.M. (ed.). 1994. Immunology, p. 9.7.1-9.7.20. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.
- Norris, S.J., V. Pope, R.E. Johnson and S.A. Larsen. 2003. *Treponema* and other human host-associated spirochetes, p. 955-971. In P.R. Murray, E.J. Baron, J.H. Jorgensen, M.A. Pfaller and R.H. Tenover (ed.), Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Perine, P.L., A.L. Wallace, J.H. Blount and S.T. Brown. 1981. Syphilis, p. 631-673. In A. Ballows and W.J. Hausler, Diagnostic procedures for bacterial, mycotic and parasitic infections. American Public Health Association, Washington, D.C.
- Portnoy, J. and W. Garson. 1960. New and improved antigen suspension for rapid reagin test for syphilis. Public Health Rep. 75:985-988.
- Portnoy, J., H.W. Bossak, V.H. Falcone and A. Harris. 1961. A rapid reagin test with unheated serum and new improved antigen suspension. Public Health Rep. 76:933-935.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
- Thomas, R.B. and J.M. Miller. 2003. Specimen collection, transport and processing: bacteriology. In P.R. Murray, E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover, Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Larsen, S.A., V. Pope, R. E. Johnson and S.J. Kraus. 1998. A manual of tests for syphilis, 9th ed. American Public Health Association, Washington, D.C.
- Pettit, D.E., S.A. Larsen, V. Pope, M.D. Perryman, and M.R. Adams. 1982. Unheated serum reagin test as a quantitative test for syphilis. J. Clin Microbiol. 15:238-242.
- Saxena, R.S. and B.N. Ghosh. 1969. A comparative study of unheated serum reagin (USR), VDRL and CF tests for syphilis. Ind. Jour. Med. Res. 57:2203-2204.

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Difco is a trademark of Difco Laboratories, a subsidiary of Becton, Dickinson and Company. BD and BD Logo are trademarks of Becton, Dickinson and Company. © 2003 BD.
BENEX Limited