
Rev From	Rev To	ECO #	Date	Appr.
0803	1103	2404-03		

Notes


- BD Cat. No. 224491, 224492, 223441, 232591, 232592, 232661, 224391, 224401
- Blank (Sheet) Size : Length: 18" Width: 36"
 Number of Pages (panels): 12 Number of Sheets: 1
 Page Size: Length 18" Width 6" Final Folded Size: 2.25" x 6"
- Style (see illustrations below): #4


#1

#2


#3


#4


#5

#6

#7

- See Specification Control No. S1333 for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 250 Schilling Circle Cockeysville, MD. 21030-0243 USA	
Proofer	Date			
Checked By	Date			
Part Number:	S1333	Category and Description Package Insert FA Human Globulin Antiglobulin (Rabbit)	Sheet: 1 of 13 Scale: 1:1	A

BD Difco™ FTA Antigen
Difco™ FTA Serum Reactive
Difco™ FTA Serum Non-Reactive
Difco™ FTA Sorbent
Difco™ FTA Sorbent Control
Difco™ FA Human Globulin Antiglobulin (Rabbit)
Difco™ Tween™ 80
Difco™ FA Buffer, Dried
Difco™ FA Mounting Fluid pH 7.2

English: pages 1-2 Italiano: pagine 7-8
 Français : pages 3-4 Espagnol: páginas 9-10
 Deutsch: Seiten 5-6

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INTENDED USE

These reagents are used in the FTA-ABS (Fluorescent Treponemal Antibody-Absorption) Test. The FTA-ABS test is an indirect fluorescent antibody procedure for detecting human antibody against *Treponema pallidum*, the causative agent of syphilis.

Difco FTA Reagents are not FDA cleared (approved) for use in testing (i.e., screening) blood or plasma donors.

SUMMARY AND EXPLANATION

Treponema pallidum is the causative agent of syphilis, a chronic infection with many clinical manifestations. These manifestations occur in distinct stages and detection of each stage requires different laboratory tests.

Since the clinical manifestations of syphilis can be confused with other infectious diseases or with noninfectious conditions that cause skin lesions, proper diagnosis must be based on microscopic examination of lesion material and serological test results.²

The FTA-ABS test is a standard diagnostic test for syphilis as defined by the Centers for Disease Control and Prevention (CDC). Other standard treponemal tests include the Fluorescent Treponemal Antibody-Absorption Double Staining Test (FTA-ABS DS) and the Microhemagglutination Assay for Antibodies to *Treponema pallidum* (MHA-TP).

Treponemal antigen tests, such as the FTA-ABS test, are used as confirmatory tests in diagnostic cases when the initial nontreponemal test is reactive. The FTA-ABS test also has utility with patients for whom the clinical, historical or epidemiological evidence of syphilis disagrees with nontreponemal tests. The persistent reactivity of the FTA-ABS test to a treated case of syphilis, sometimes for life, minimizes its use for following response to therapy. Therefore, the FTA-ABS test is also unreliable in detecting new untreated cases in epidemiological investigations. The test should not be used as a routine screening procedure.^{3,4} The likelihood of obtaining a reactive FTA-ABS test result in various stages of untreated syphilis has been reported as follows:²

Stage of Untreated Syphilis	% Reactive
Primary	84
Secondary	100
Latent	100
Tertiary (Late)	96

PRINCIPLES OF THE PROCEDURE

Patient serum is diluted 1:5 in sorbent and layered on a microscope slide fixed with *T. pallidum*. If the patient's serum contains antibodies, these antibodies will coat the treponemes on the slide. Fluorescein-labeled antihuman immunoglobulin is added. It combines with the patient antibodies already adhering to the *T. pallidum* and produces fluorescein-stained spirochetes that can be observed with a fluorescence microscope.^{5,6}

REAGENTS

Difco FTA Antigen (also known as *T. pallidum* antigen) is a lyophilized, standardized, killed suspension of *T. pallidum* (Nichols strain).

Difco FTA Serum Reactive is lyophilized, standardized syphilitic human sera containing 0.02% thimerosal as a preservative. It is used to prepare:

- Reactive Control Serum (4+) – Unabsorbed
- Reactive Control Serum (4+) – Absorbed
- Minimally Reactive Control Serum (1+)

It is used as a positive control in the FTA-ABS test.

Difco FTA Serum Non-Reactive is lyophilized, standardized, non-syphilitic human sera containing 0.02% thimerosal as a preservative. It is used to make Nonreactive Control Serum (N). It is used as a negative control in the FTA-ABS test.

Difco FTA Sorbent is a lyophilized, standardized extract of the nonpathogenic Reiter treponeme (*Treponema phagedenis*) prepared from broth culture. It is used to remove nonpathogenic treponeme antibodies during preparation of the test specimen.

Difco FTA Sorbent Control is lyophilized, standardized, non-syphilitic human sera containing 0.02% thimerosal as a preservative. It is used to make Non-specific Control Serum – Unabsorbed, which demonstrates at least 2+ nonspecific reactivity at a 1:5 dilution in FA Buffer, and Non-specific Control Serum – Absorbed, which demonstrates essentially no reactivity at a 1:5 dilution in FTA Sorbent.

Difco FA Human Globulin Antiglobulin (Rabbit) is lyophilized, fluorescein-conjugated (FITC) antihuman globulin containing 0.02% thimerosal as a preservative. It is used to show the presence of human syphilitic antibodies on the treponemal antigen.

Difco Tween™ 80 is Polysorbate 80, U.S.P. It is used to prepare 2% Tween 80, which acts as a dispersing agent.

Difco FA Buffer, Dried is phosphate-buffered saline (PBS) which, upon rehydration, yields a 0.85% NaCl solution buffered to pH 7.2. **Difco FA Buffer, Dried** is used in preparing:

- Reactive Control Serum (4+) – Unabsorbed
- Minimally Reactive Control Serum (1+)
- Nonreactive Control Serum (N)
- Non-specific Staining Control – Unabsorbed

Difco FA Mounting Fluid pH 7.2 is standardized, reagent grade glycerin adjusted to pH 7.2 for use in mounting specimens on slides to be viewed under the fluorescence microscope.

Warnings and Precautions:

- For *in vitro* Diagnostic Use.
- Difco FTA Serum Reactive**
Difco FTA Serum Non-Reactive
Difco FTA Sorbent Control
WARNING! POTENTIAL BIOHAZARDOUS REAGENTS. Each donor unit used in the preparation of these reagents was tested by a FDA-licensed method for the presence of the antibody to Human Immunodeficiency Virus (HIV) as well as for hepatitis B surface antigen (HBsAg) and found to be negative (were not repeatedly reactive).

Because no test method can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health Manual *Biosafety in Microbiological and Biomedical Laboratories*, 1999.

- Difco FTA Antigen**
Difco FTA Serum Reactive
Difco FTA Serum Non-Reactive
Difco FTA Sorbent
Difco FTA Sorbent Control
Difco FA Human Globulin Antiglobulin (Rabbit)
 The Packaging of This Product Contains Dry Natural Rubber.
- Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁷⁻¹⁰ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Storage: Store unopened products as specified below:

Difco FTA Antigen	2-8°C
Difco FTA Serum Reactive	2-8°C
Difco FTA Serum Non-Reactive	2-8°C
Difco FTA Sorbent	2-8°C
Difco FTA Sorbent Control	2-8°C
Difco FA Human Globulin Antiglobulin (Rabbit)	2-8°C in the dark
Difco Tween 80	15-30°C
Difco FA Buffer, Dried	Below 30°C
Difco FA Mounting Fluid pH 7.2	15-30°C

Rehydrated **Difco FTA Antigen** stored at 2-8°C is stable for 1 week.

Rehydrated **Difco FA Buffer** showing turbidity or mold growth should be discarded. Discard 2% Tween 80 that exhibits a precipitate or pH change.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Product Deterioration: Expiration date applies to product in its intact container when stored as directed. Do not use if the product is caked, discolored or shows other signs of deterioration.

SPECIMEN COLLECTION AND PREPARATION

Test serum: Collect a blood specimen by aseptic venipuncture into a clean, dry tube without anticoagulant. After the specimen has clotted, centrifuge the specimen at 1,500-2,000 rpm for 5 min to obtain test serum. Store specimens at room temperature for up to 4 h or at 2-8°C for up to 5 days; serum specimens may be frozen at or below -20°C.

Test and control sera: Equilibrate the sera to room temperature, then heat at 56°C for 30 min. Reheat previously heated sera for 10 min on the day of testing. Cool to room temperature before testing. Bacterial contamination or excessive hemolysis may render a specimen unsuitable for testing. Such specimens should not be tested.

PROCEDURE

Materials Provided

- Difco FTA Antigen
- Difco FTA Serum Reactive
- Difco FTA Serum Non-Reactive
- Difco FTA Sorbent
- Difco FTA Sorbent Control
- Difco FA Human Globulin Antiglobulin (Rabbit)
- Difco Tween 80
- Difco FA Buffer, Dried
- Difco FA Mounting Fluid, pH 7.2

Materials Required but Not Provided

- Timer
- Serological pipettes, 0.2 mL, 1 mL, 5 mL
- Micropipettors delivering 10-200 µL
- Test tubes, 12 x 75 mm
- Water bath (56°C)
- Vortex mixer
- Platinum loop, 2 mm, 26 gauge

- Slides, plain or frosted, 1 x 3 inch, 1 mm thick, inscribed with 2 x 1 cm circles
- Staining dish with removable slide carriers
- Slide board or holder
- Moisture chamber
- Acetone
- Bibulous paper
- Purified water
- Incubator, 35-37°C
- Oil, Immersion
- Cover slips, No. 1, 22 mm square
- Fluorescence microscope assembly:
- Lamps: HBO-50, HBO-100, HBO-200 or Xenon XBO-150; 6X 5A Tungsten
- Ocular: 10X
- Objective: 10X, 40X (Fluorite)
- Filters: BG-12 or KP490, K515 or K530
- Condenser: Dark-field D1.20-1.40

Reagent Preparation

FTA Antigen: Rehydrate with 1 mL purified water and rotate to completely dissolve the contents. This solution will yield approximately 3.5 x 10⁷ treponemes per mL. Mix thoroughly with a disposable pipette and rubber bulb, drawing the suspension into and expelling it from the pipette 8-10 times to break treponemal clumps and ensure an even distribution of treponemes. Confirm the even distribution by dark-field examination. Use FTA Antigen in its entirety to prepare antigen smears on the day it is rehydrated. Approximately 200-300 slides may be prepared with 1 mL of antigen.

To prepare FTA Antigen smears:

- Wipe inscribed slides with clean gauze and, if necessary, alcohol to remove dust particles.
- Using a platinum wire loop (2 mm, 26 gauge), smear 1 loopful of reconstituted FTA Antigen within the 2 circles. Air dry at room temperature for at least 15 min.
- Immerse the dry slide into acetone for 10 min to fix the treponemal antigen smear to the slide; air dry. Fix no more than 50 slides per 200 mL of acetone.
- Use slides immediately or store at or below -20°C after acetone fixation. Thaw before use; do not refreeze. Use within 1 year, but only if satisfactory results are obtained with test controls.

FTA Serum Reactive: Rehydrate with 5 mL purified water and rotate gently to completely dissolve the contents. Aliquot in 0.4 mL amounts and store at or below -20°C. Do not refreeze thawed aliquot. Approximately 12 tests may be obtained per 5 mL vial. This serum should be heated at 56°C for 30 min before use.

FTA Serum Non-Reactive: Rehydrate with 5 mL purified water and rotate gently to completely dissolve the contents. Aliquot in 0.4 mL amounts and store at or below -20°C. Approximately 90-100 tests may be obtained per 5 mL vial. This serum should be heated at 56°C for 30 min before use.

FTA Sorbent: Rehydrate with 5 mL purified water and rotate gently to completely dissolve the contents. Store at 2-8°C or aliquot and store at -20°C. The quantity of FTA Sorbent used for each test sample or serum is 0.2 mL. The quantity of FTA Sorbent needed for 3 controls is 0.6 mL. Approximately 20-25 tests may be performed with 5 mL of FTA Sorbent.

FTA Sorbent Control: Rehydrate with 0.5 mL purified water and rotate gently to completely dissolve the contents. Aliquot in 0.25 mL amounts and store at or below -20°C. For each test, 0.1 mL of FTA Sorbent Control is needed. Approximately 2 tests may be performed per 0.5 mL vial because of evaporation from heating. This serum should be heated at 56°C for 30 min before using.

FA Human Globulin Antiglobulin (Rabbit): Rehydrate with 1 mL or 5 mL purified water, depending on label directions. Aliquot in 0.5 mL amounts and store at or below -20°C. Each lot is supplied with a dilution titer. Since conditions and equipment differ from one laboratory to another, it is necessary to titer and test a new lot of conjugate with the fluorescence microscope assembly currently in use.^{3,4,11}

- Prepare serial dilutions in 2% Tween 80, including the titer specified on the vial.
- Test each dilution per the Test Procedure with Reactive Control Serum (4+) and Non-specific Staining Control.
- Test a known lot of conjugate using the Reactive Control Serum (4+), Minimally Reactive Control Serum (1+) and Non-specific Staining Control as controls of the reagents and test conditions.
- The dilution to be used in the test is the dilution that produces 1 doubling dilution lower than the 4+ endpoint. The 4+ endpoint is the highest dilution of conjugate yielding 4+ fluorescence with the Reactive Control Serum (4+).

FA Buffer, Dried: Dissolve 10 g in 1 liter of purified water and rotate gently to completely dissolve the contents. Store at 2-8°C. Use the solution if it is free of mold growth and turbidity.

Tween 80: Heat the bottle of Tween 80 and a flask containing 98 mL FA Buffer to 56°C in a water bath. Add 2 mL of Tween 80 to the buffer and rinse the pipette thoroughly in the buffer. Check the pH and adjust to pH 7.2 with 1N NaOH. Discard if a precipitate develops or the pH changes.

Test Procedure

This procedure conforms to those published by the American Public Health Association.^{6,11}

- FTA Antigen smears:** Obtain previously prepared smears, thaw and dry if appropriate, and identify the frosted end of the slides to correspond with each test and control serum to be tested.
- Prepare the following test and control sera in appropriately identified tubes no more than 30 min before testing and mix thoroughly (at least 8 times):

Test Serum (1:5): Dilute 0.05 mL of heated (or reheated) test serum in 0.2 mL FTA Sorbent.

Reactive Control Serum (4+) – Unabsorbed: Dilute 0.05 mL FTA Serum Reactive in 0.2 mL FA Buffer (PBS).

Reactive Control Serum (4+) – Absorbed: Dilute 0.05 mL FTA Serum Reactive in 0.2 mL FTA Sorbent.

Minimally Reactive Control Serum (1+): Dilute FTA Serum Reactive, as indicated on the label, in FA Buffer (PBS) to yield a 1+ fluorescence. The minimal degree of fluorescence that can be reported as reactive is 1+ fluorescence.

Nonreactive Control Serum (N) (1:40): Prepare a 1:40 dilution of FTA Serum Non-Reactive by adding 0.05 mL (50 µL) of serum to 1.95 mL FA Buffer (PBS).

Non-specific Serum Control – Unabsorbed (≥ 2+ nonspecific reactivity): Dilute 0.05 mL FTA Sorbent Control in 0.2 mL FA Buffer (PBS).

Non-specific Serum Control – Absorbed (nonreactive, – to ±): Dilute 0.05 mL FTA Sorbent Control in 0.2 mL FTA Sorbent.

Non-specific Staining Control – Unabsorbed: Use 0.03 mL FA Buffer (PBS) undiluted.

Non-specific Staining Control – Absorbed: Use 0.03 mL FTA Sorbent undiluted.

- FTA Antigen smears:** Cover the previously identified FTA Antigen smears with 0.03 mL of the corresponding test or control serum prepared above, making certain that the entire smear is covered.
- Place the slides in a moist chamber to prevent evaporation and incubate at 35-37°C for 30 min.
- Place the slides in a slide carrier and rinse as follows:
 - Rinse in running FA Buffer for 5 sec.
 - Soak in FA Buffer for 5 min.
 - Agitate by dipping in and out of the buffer 30 times.
 - Repeat the soaking and agitation in fresh buffer.
 - Rinse in running distilled water for 5 sec.
 - Gently blot dry with bibulous paper.
- FA Human Globulin Antiglobulin (Rabbit):** Dilute the conjugate to its working titer (determined above) using 2% Tween 80 in FA Buffer.
- FTA Antigen smears:** Cover each test and control smear with approximately 0.03 mL of diluted FA Human Globulin Antiglobulin (Rabbit). Spread uniformly to cover the entire smear.
- Repeat steps 4 and 5.
- Mount the slides immediately using a small drop of FA Mounting Fluid pH 7.2 and apply a cover slip, being careful not to trap air bubbles in the mounting fluid.
- Immediately examine the slides microscopically for intensity of fluorescence using the microscope assembly described above. If it is necessary to delay reading, store the slides in the dark and read within 4 h. Results are valid only if the quality control pattern is satisfactory.
- Verify the presence of treponemes on the nonreactive control slides by dark-field microscopy.

User Quality Control

Rehydrate and dilute reagents per directions above. Test as described. Tests failing to exhibit the following control results are unsatisfactory and should not be reported.^{4,11}

Serum Tested	Expected Fluorescence	Interpretation
Reactive Control Serum – Unabsorbed	4+	Reactive
Reactive Control Serum – Absorbed	3+ to 4+	Reactive
Minimally Reactive Control Serum	1+	Reactive
Nonreactive Control Serum	N	Nonreactive
Non-specific Serum Control – Unabsorbed	2+ to 4+	Reactive
Non-specific Serum Control – Absorbed	N to ±	Nonreactive
Non-specific Staining Control – Unabsorbed	N	Nonreactive
Non-specific Staining Control – Absorbed	N	Nonreactive

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Using the 1+ serum control as a reading standard, record the intensity of fluorescence of the treponemes and report as follows. Retest all specimens with an initial test fluorescence of 1+.

Intensity of Fluorescence	Initial Test	Retest	Report
Moderate to strong	2+ to 4+	NA	Reactive
	1+	>1+	Reactive
Equivalent to 1+ control	1+	1+	Reactive minimal*
	1+	<1+	Nonreactive
Visible staining but <1+	± to <1+	NA	Nonreactive
None or vaguely visible, not distinct "Moth eaten" or "beaded"	-	NA	Nonreactive Atypical ¹

*Equivalent result.

¹Atypical or beaded pattern of fluorescence has been described in patients with lupus or other autoimmune diseases.¹¹

Specimens with an initial reading of 1+ and a retest reading of 1+ or greater are reported as reactive. All other results are reported as nonreactive.

Retesting nonreactive specimens is not necessary.

Without historical or clinical evidence of treponemal infection, equivocal test results (see table) suggest the need for testing a second specimen obtained 1-2 weeks after the initial specimen.

LIMITATIONS OF THE PROCEDURE

- When the treponemal test results and the clinical opinion disagree, repeat the treponemal test and obtain additional clinical and historical information. If the disagreement persists, send the specimen to a reference laboratory, such as the local state health department, for additional confirmatory tests. The final diagnosis depends on the clinical judgment of a specialist very experienced in sexually transmitted diseases.^{2,3}
- The test should not be used to follow the response to therapy nor can it be relied on to detect new, untreated cases in epidemiological investigations.
- "Atypical" fluorescence and false-positive results have been associated with patients having active systemic, discoid and drug-induced varieties of lupus erythematosus¹²⁻¹⁵ and other autoimmune diseases.
- Elderly patients may exhibit unexplained FTA-ABS reactions.
- At times, deciding whether a reading is weak or vaguely visible may be difficult. The ability to make this distinction is critical, since a nonreactive (vaguely visible to none) serum is not retested.
- The FTA-ABS test may be reactive in persons from areas where yaws or pinta was, or is, endemic.

PERFORMANCE CHARACTERISTICS¹⁶

The performance of the FTA-ABS test was compared to an ELISA test, the VDRL slide test and the MHA-TP test in a study conducted by Pope, Hunter and Feeley¹⁶ at the Centers for Disease Control in Atlanta, Georgia. Two hundred and ninety-seven (297) sera were tested by an enzyme-linked immunosorbent assay (ELISA), VDRL slide test, FTA-ABS test and the microhemagglutination assay for *T. pallidum* antibodies (MHA-TP). Table 1 below shows a comparison of test reactivity by category of disease when testing 75 sera from syphilitic individuals and 222 sera from nonsyphilitic individuals.¹⁶

Table 1

Category	No. Tested	No. with reactivity to:			
		FTA-ABS	ELISA	MHA-TP	VDRL
Primary syphilis					
Untreated	22	22	17	9	17
Treated	2	2	1	2	1
Secondary syphilis					
Untreated	8	8	8	6	8
Treated	12	12	12	12	12
Latent syphilis					
Untreated	13	13	12	12	10
Treated	10	10	10	9	9
Treatment unknown	1	1	1	1	1
Questionable latent	3	3	3	3	0
Neurosyphilis	3	3	2	2	0
Cardiovascular syphilis	1	1	1	1	0
Nonsyphilitic					
Presumed Normal	178	1	1	1	0
Biological false-positive	15	1	0	0	7
Diseases other than syphilis	29	1	1	0	1

Table 2 below lists the sensitivities and specificities of each method published in this study.¹⁶

Table 2

Test	Sensitivity	Specificity
FTA-ABS	100%	97.8%
ELISA	89.3%	98.5%
MHA-TP	76%	98.2%
VDRL	93.3%	92.7%

AVAILABILITY

Cat. No.	Description
223441	Difco™ FTA Antigen, 1 mL
224401	Difco™ FTA Serum Non-Reactive, 5 mL
224391	Difco™ FTA Serum Reactive, 5 mL
232591	Difco™ FTA Sorbent, 5 mL
232592	Difco™ FTA Sorbent, 6 x 5 mL
232661	Difco™ FTA Sorbent Control, 6 x 0.5 mL
223143	Difco™ FA Buffer, Dried, 6 x 10 g
223142	Difco™ FA Buffer, Dried, 100 g
224491	Difco™ FA Human Globulin Antiglobulin (Rabbit), 1 mL
224492	Difco™ FA Human Globulin Antiglobulin (Rabbit), 5 mL
223291	Difco™ FA Mounting Fluid pH 7.2, 6 x 5 mL
231181	Difco™ Tween™ 80, 100 g

REFERENCES

- Johnson, R.M. Letter. July 1, 1994. Department of Health & Human Services, Public Health Service, Food and Drug Administration, Rockville, Md.
- Janda, W.M. (ed.). 1994. Immunology, p. 9.7.1-9.7.20. *In* H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.
- Norris, S.J., V. Pope, R.E. Johnson and S.A. Larsen. 2003. *Treponema* and other human host-associated spirochetes, p. 955-971. *In* P.R. Murray, E.J. Baron, J.H. Tenover, M.A. Tenover and R.H. Tenover (eds.), Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Turgeon, M.L. 1990. Immunology and serology in laboratory medicine. The C.V. Mosby Company, St. Louis, Mo.
- Hunter, E.F., W.E. Deacon and P.E. Meyer. 1964. An improved FTA test for syphilis; the absorption procedure (FTA-ABS). *Publ. Hlth. Report.* 79:410-412.
- Wentworth, B.B., and F.N. Judson. 1984. Laboratory methods for the diagnosis of sexually transmitted diseases. American Public Health Association, Washington, D.C.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000, p. 0021-0045.*
- Larsen, S.A., V. Pope, R.E. Johnson and E.J. Kennedy, Jr. (ed.). 1998. A manual of tests for syphilis, 9th ed. American Public Health Association, Washington, D.C.
- Kraus, S.J., J.R. Haserick and M.A. Lantz. 1970. Fluorescent treponemal antibody absorption test reactions in lupus erythematosus. *N. Engl. J. Med.* 282:1287-1290.
- Shore, R.N., and J.A. Faricelli. 1977. Borderline and reactive FTA-ABS results in lupus erythematosus. *Arch. Dermatol.* 113:37-41.
- Monson, R.A. 1973. Biological false-positive FTA-ABS test in drug induced lupus erythematosus. *J.A.M.A.* 224:1028-1030.
- Anderson, B., and M.T. Stillman. 1978. False-positive FTA-ABS in hydralazine-induced lupus. *J.A.M.A.* 239:1392-1493.
- Pope, V, E.F. Hunter and J.C. Feeley. 1982. Evaluation of the microenzyme-linked immunosorbent assay with *Treponema pallidum* antigen. *J. Clin. Microbiol.* 15:630-634.

FTA-ABS Test Procedure

Abbreviated Schematic

Step 1.	Step 2.	Step 3.	Step 4.	Step 5.
Prepare sera and reagents.	Dilute sera.	Add test/control serum to appropriate FTA Antigen smears.	Add conjugate to FTA Antigen smears.	Record reactions of test and control sera. Verify that the control sera provided the expected results.

FTA Antigen	FTA Antigen Smear
Rehydrate with 1 mL purified water. Prepare smears. Fix with acetone. Use as "FTA Antigen smear."	Thaw, dry and identify sufficient FTA Antigen smears to correspond with each of the test and control sera to be tested.

Test (Patient) Serum				
Heat at 56°C for 30 min or reheat previously heated serum for 10 min.	Dilute 1:5 by adding 0.05 mL serum to 0.2 mL in FTA Sorbent.	Apply 0.03 mL test serum to an FTA Antigen smear. Incubate smear. Rinse.	Apply 0.03 mL conjugate to the smear. Incubate, rinse and mount slide. Examine microscopically.	Dependent on antibody status or test serum.

FTA Serum Reactive				
Rehydrate with 5 mL purified water. Heat at 56°C for 30 min.	Dilute 1:5 by adding 0.05 mL FTA Serum Reactive to 0.2 mL FA Buffer.	Apply 0.03 mL Reactive Control Serum - Unabsorbed to an FTA Antigen Smear. Incubate smear. Rinse.	Apply 0.03 mL conjugate to the smear. Incubate, rinse and mount slide. Examine microscopically.	4+ Reactive
	Dilute 1:5 by adding 0.05 mL FTA Serum Reactive to 0.2 mL FTA Sorbent.	Apply 0.03 mL Reactive Control Serum - Absorbed to an FTA Antigen smear. Incubate smear. Rinse.	Apply 0.03 mL conjugate to the smear. Incubate, rinse and mount slide. Examine microscopically.	3+ to 4+ Reactive

Minimally Reactive Control Serum (1+)				
Dilute FTA Serum Reactive to 1+ fluorescence (labeled titer) in FA Buffer.	Apply 0.03 mL Minimally Reactive Control Serum to an FTA Antigen Smear. Incubate smear. Rinse.	Apply 0.03 mL conjugate to the smear. Incubate, rinse and mount slide. Examine microscopically.		1+ Reactive

FTA Serum Non-Reactive				
Rehydrate with 5 mL purified water. Heat at 56°C for 30 min.	Dilute 1:40 by adding 0.05 mL FTA Serum Non-Reactive to 1.95 mL FA Buffer.	Apply 0.03 mL Nonreactive Control Serum to an FTA Antigen smear. Incubate smear. Rinse.	Apply 0.03 mL conjugate to the smear. Incubate, rinse and mount slide. Examine microscopically.	N - Nonreactive

FTA Sorbent Control				
Rehydrate with 0.5 mL purified water. Heat at 56°C for 30 min.	Dilute 1:5 by adding 0.05 mL FTA Sorbent Control to 0.2 mL FA Buffer.	Apply 0.03 mL Nonspecific Serum Control - Unabsorbed to an FTA Antigen Smear. Incubate smear. Rinse.	Apply 0.03 mL conjugate to the smear. Incubate, rinse and mount slide. Examine microscopically.	2+ to 4+ Reactive
	Dilute 1:5 by adding 0.05 mL FTA Sorbent Control to 0.2 mL FTA Sorbent.	Apply 0.03 mL Nonspecific Serum Control - Absorbed to an FTA Antigen smear. Incubate smear. Rinse.	Apply 0.03 mL conjugate to the smear. Incubate, rinse and mount slide. Examine microscopically.	N to ± Nonreactive

FA Buffer, Dried				
Dissolve 10 g in 1 L purified water.	Use 0.03 mL FA Buffer undiluted as the diluent (above) and as the Nonspecific Staining Control - Unabsorbed.	Apply 0.03 mL Nonspecific Staining Control - Unabsorbed to an FTA Antigen smear. Incubate smear. Rinse.	Apply 0.03 mL conjugate to the smear. Incubate, rinse and mount slide. Examine microscopically.	N - Nonreactive
	Use 0.03 mL FA Sorbent undiluted as the diluent (above) and as the Nonspecific Staining Control - Absorbed.	Apply 0.03 mL Nonspecific Staining Control - Absorbed to an FTA Antigen smear. Incubate smear. Rinse.	Apply 0.03 mL conjugate to the smear. Incubate, rinse and mount slide. Examine microscopically.	N - Nonreactive

FA Human Globulin Antiglobulin (Rabbit)				
Rehydrate with 1 mL or 5 mL purified water. Determine titer if a new lot.	Dilute to labeled titer with 2% Tween. Use as "Conjugate."			

Tween 80				
Heat Tween 80 and FA Buffer to 56°C.				
Add 2 mL Tween 80 to 98 mL FA Buffer. Adjust to pH 7.2. Use as "2% Tween."				

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Tween is a trademark of ICI Americas, Inc.
Difco is a trademark of Difco Laboratories, a subsidiary of Becton, Dickinson and Company.
BD and BD Logo are trademarks of Becton, Dickinson and Company. ©2003 BD.
**BENEX Limited**