BD BACTEC™ Peds Plus™/F Culture Vials (plastic) **Soybean-Casein Digest Broth with Resins in a Plastic Vial**

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INTENDED USE

BD BACTEC Peds Plus™/F culture vials (enriched Soybean-Casein Digest broth with CO₂) are for aerobic blood cultures. Principal use is with the **BD BACTEC** fluorescent series instruments for the qualitative culture and recovery of aerobic microorganisms (mainly bacteria and yeast) from pediatric and non-pediatric blood specimens which are generally less than 5 mL in volume.

SUMMARY AND EXPLANATION

The sample to be tested is inoculated into one or more vials which are inserted into the **BD BACTEC** fluorescent series instrument for incubation and periodic reading. Each vial contains a chemical sensor which can detect increases in CO₂ produced by the growth of microorganisms. The sensor is monitored by the instrument every ten minutes for an increase in its fluorescence, which is proportional to the amount of CO₂ present. A positive reading indicates the presumptive presence of viable microorganisms in the vial. Detection is limited to microorganisms that will grow in a particular type of medium.

Resins have been described for the treatment of blood specimens both prior to and after their inoculation into culture media. Resins have been incorporated into **BD BACTEC** culture media to enhance recovery of organisms without a need for special processing.^{1-3,8}

PRINCIPLES OF THE PROCEDURE

If microorganisms are present in the test sample inoculated into the **BD BACTEC** vial, CO₂ will be produced when the organisms metabolize the substrates present in the vial. Increases in the fluorescence of the vial sensor caused by the higher amount of CO₂ are monitored by the **BD BACTEC** fluorescent series instrument. Analysis of the rate and amount of CO₂ increase enables the **BD BACTEC** fluorescent series instrument to determine if the vial is positive, i.e., that the test sample contains viable organisms.

REAGENTS

The **BD BACTEC** culture vials contain the following reactive ingredients prior to processing:

List of Ingredients	BD BACTEC Peds Plus/F
Processed Water	40 mL
Soybean-Casein Digest Broth	2.75% w/v
Yeast Extract	0.25% w/v
Animal Tissue Digest	0.10% w/v
Sodium Pyruvate	0.10% w/v
Dextrose	0.06% w/v
Sucrose	0.08% w/v

List of Ingredients	BD BACTEC Peds Plus/F
Hemin	0.0005% w/v
Menadione	0.00005% w/v
Sodium Polyanetholsulfonate (SPS)	0.02% w/v
Pyridoxal HCl (Vitamin B ₆)	0.001% w/v
Nonionic Adsorbing Resin	10.0% w/v
Cationic Exchange Resin	0.6% w/v

All **BD BACTEC** media are dispensed with added CO₂.

Warnings and Precautions:

The prepared culture vials are for *in vitro* diagnostic use. This Product Contains Dry Natural Rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁴⁻⁷ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Prior to use, each vial should be examined for evidence of contamination such as cloudiness, bulging or depressed septum, or leakage. DO NOT USE any vial showing evidence of contamination. A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, gas or contaminated culture media could be refluxed into the patient's vein. Vial contamination may not be readily apparent. If a direct draw procedure is used, monitor the process closely to avoid refluxing materials into the patient.

Prior to use, the user should examine the vials for evidence of damage or deterioration. Vials displaying turbidity, contamination, or discoloration (darkening) should not be used. On rare occasions a vial may not be sealed sufficiently. The contents of the vials may leak or spill, especially if the vial is inverted. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling should be performed in a biological safety cabinet if possible, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or **BD Luer-Lok™** tips.

Storage Instructions

The **BD BACTEC** vials are ready for use as received and require no reconstitution or dilution. Store in a cool, dry place (2–25 °C), out of direct light.

SPECIMEN COLLECTION

The specimen must be collected using sterile techniques to reduce the chance of contamination. The range of blood volume which can be cultured is 0.5 to 5.0 mL. If the volume of blood cultured is less than 0.5 mL, recovery of some fastidious organisms, such as *Haemophilus* species, may require the use of an appropriate supplement, as described later in this Package Insert. It is recommended that the specimen be inoculated into the **BD BACTEC** vials at bedside. Most commonly, a syringe with a **BD Luer-Lok** brand tip is used to draw the sample. If appropriate, a **BD Vacutainer®** brand Needle Holder and a **BD Vacutainer** brand Blood Collection Set, **Vacutainer Safety-Lok™** Blood Collection Set or other tubing "butterfly" set may be used. If using a needle and tubing set (direct draw), carefully observe the direction of blood flow when starting sample collection. The vacuum in the vial will usually exceed 5 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. When the recommended 1–3 mL has been drawn, the flow should be stopped by crimping the tubing and removing the tubing set from the **BD BACTEC** vial. **The inoculated BD BACTEC vial should be transported as quickly as possible to the laboratory.**

PROCEDURE

Remove the flip-off cap from the **BD BACTEC** vial top and inspect the vial for cracks, contamination, excessive cloudiness, and bulging or indented septum. DO NOT USE if any defect is noted. Before inoculating, swab the septum with alcohol (iodine is not recommended). Aseptically inject or draw directly a maximum of 5 mL of specimen per vial (see Limitations of the Procedure). Inoculated vials should be placed in the **BD BACTEC** fluorescent series instrument as soon as possible for incubation and monitoring. If placement of an inoculated vial into the instrument has been delayed and visible growth is apparent, it should not be tested in the **BD BACTEC** fluorescent series instrument, but rather it should be subcultured, Gram-stained and treated as a presumptively positive vial.

Vials entered into the instrument will be automatically tested every ten minutes for the duration of the testing protocol period. Positive vials will be determined by the **BD BACTEC** fluorescent series instrument and identified as such (see the appropriate **BD BACTEC** Fluorescent Series Instrument User's Manual). The sensor inside the vial will not appear visibly different in positive and negative vials, however, the **BD BACTEC** fluorescent series instrument can determine a difference in fluorescence.

If at the end of the testing period a negative vial appears visually positive (i.e., chocolate blood, bulging septum, and/or lysed), it should be subcultured and Gram-stained and treated as a presumptively positive.

Positive vials should be subcultured and Gram-stained. In a great majority of cases, organisms will be seen and a preliminary report can be made to the physician. Subcultures to solid media and a preliminary direct antimicrobial susceptibility test may be prepared from fluid in the **BD BACTEC** vials.

Subculturing: Prior to subculturing, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, use an appropriate venting unit (BD Cat. No. 249560, or equivalent). The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

DO NOT USE culture vials past their expiration date.

DO NOT USE culture vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each carton of media. Quality Control Certificates list test organisms, including ATCC® cultures specified in the CLSI Standard M22, *Quality Control for Commercially Prepared Microbiological Culture Media*. The range of time-to-detection in hours was ≤ 72 hours for each of the organisms listed on the Quality Control Certificate for this medium:

Peds Plus Medium Organisms

<i>Streptococcus pyogenes</i> ATCC 19615	<i>Neisseria meningitidis</i> ATCC 13090
<i>Escherichia coli</i> ATCC 25922	<i>Alcaligenes faecalis</i> ATCC 8750
<i>Streptococcus pneumoniae</i> * ATCC 6305	<i>Haemophilus influenzae</i> ATCC 19418
<i>Pseudomonas aeruginosa</i> ATCC 27853	<i>Staphylococcus aureus</i> ATCC 25923
<i>Candida albicans</i> ATCC 18804	

* CLSI-recommended strain

For information on Quality Control for the **BD BACTEC** fluorescent series instrument, refer to the appropriate **BD BACTEC** Fluorescent Series Instrument User's Manual.

RESULTS

A positive sample is determined by the **BD BACTEC** fluorescent series instrument and indicates the presumptive presence of viable microorganisms in the vial.

LIMITATIONS OF THE PROCEDURE

Contamination

Care must be taken to prevent contamination of the sample during collection and inoculation into the **BD BACTEC** vial. A contaminated sample will give a positive reading, but will not indicate a relevant clinical sample. Such a determination must be made by the user based on such factors as type of organism recovered, occurrence of the same organism in multiple cultures, patient history, etc.

Recovery of SPS Sensitive Organisms From Blood Samples

Because blood can neutralize the toxicity of SPS toward organisms sensitive to SPS (such as some *Neisseria* species), the presence of recommended volumes of blood (1–3 mL) can help to optimize recovery of these organisms.

Some organisms may be dependent on having a minimum amount of blood in the medium for optimal growth. Fastidious organisms, such as certain *Haemophilus* species, require growth factors from the blood specimen, such as NAD, or factor V. Optimal growth of these organisms is dependent on having greater than a minimum of 0.5 mL blood in the specimen. If the blood specimen volume is very small (0.5 mL or less), an appropriate supplement may be required for recovery of these organisms. **BD BACTEC FOS™** Fastidious Organism Supplement (catalog number 442153) may be used as a nutritional supplement.

Nonviable Organisms

A Gram-stained smear from a culture medium may contain small numbers of nonviable organisms derived from media constituents, staining reagents, immersion oil, glass slides, and specimens used for inoculation. In addition, the patient specimen may contain organisms that will not grow in the culture medium or in media used for subculture. Such specimens should be subcultured to special media as appropriate.

Antibiotic Activity

Neutralization of the antibiotic activity by resins varies depending on dosage level and timing of specimen collection.

Studies have demonstrated that the resins present in this medium do not adequately neutralize meropenem preparations.

Based on the neutralization study with *Candida albicans* and Fluconazole, the media showed some degree of neutralization; however the results were inconclusive. The adequacy of antifungal neutralization by resins in the **BD BACTEC Peds Plus/F** vial (plastic) is unknown.

Recovery of *Streptococcus pneumoniae*

In aerobic media, *S. pneumoniae* will typically be visually and instrument positive, but in some cases no organism will be seen on Gram stain or recovered on routine subculture. If an anaerobic vial was also inoculated, the organism can usually be recovered by performing an aerobic subculture of the anaerobic vial, since this organism has been reported to grow well under anaerobic conditions.⁹

General Considerations

Recovery of isolates will be achieved by adding the recommended volume 1–3 mL of blood. Use of lower or higher volumes may adversely affect recovery and/or detection. Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms. False negative readings may result when certain organisms are present which do not produce enough CO₂ to be detected by the system or significant growth has occurred before placing the vial into the system. False positivity may occur when the white blood cell count is high. The default 5-day (120 hours) protocol was utilized for all analytical testing with the **BD BACTEC Peds Plus/F** culture media and protocol lengths of >5 days have not been evaluated.

EXPECTED VALUES AND SPECIFIC PERFORMANCE CHARACTERISTICS

Internal studies have demonstrated that antibiotics are effectively neutralized by the resins used in **BD BACTEC** resin media. In these tests, antibiotics were added in clinically relevant concentrations directly to resin media prior to inoculation with susceptible strains. These tests showed equivalent performance in **BD BACTEC Peds Plus** in plastic vial compared to **BD BACTEC Peds Plus** in glass vial.

A total of 984 paired sets at 10–100 CFU per vial were evaluated across the four instruments comprising the **BD BACTEC** fluorescent-series instrument family: **BD BACTEC 9050**, **BD BACTEC 9240**, **BD BACTEC FX** and the **BD BACTEX FX40**. Of the 984 paired sets 953 sets recovered organisms within the instrument series. There were 18 sets with no detection of organisms in either the plastic or glass bottle that included *Candida albicans* (4 sets) *Haemophilus influenzae* (9 sets) and *Haemophilus parainfluenzae* (5 sets). There were 4 sets with no detection in the plastic bottle that included *Candida albicans* (2 sets), *Enterococcus faecalis* (1 set) and *Haemophilus influenzae* (1 set). There were 9 sets with no detection in the glass bottle that included *Candida albicans* (3 sets), *Haemophilus influenzae* (1 set) *Haemophilus parainfluenzae* (4 sets) and *Pediococcus acidilactici* (1 set) with no detection in the glass bottle. The detection rate of *Candida albicans*, *Enterococcus faecalis* and *Pediococcus acidilactici* was 73%, 98% and 98% respectively under these test conditions. The *Haemophilus* species detection rates were between 8%–40% due to the quality (freshness) and volume of blood used in the test (bagged blood with 0.5 mL added to the bottle).

There were five organisms with false negative results (i.e., end of protocol, instrument negative vials with a positive terminal subculture) observed with the **BD BACTEC Peds Plus**/F medium contained in a plastic vial using 0.5 mL of bagged blood: *H. influenzae* inoculated at 54, 65 CFU, *Haemophilus parainfluenzae* inoculated at 4, 58 CFU, *Candida glabrata*, inoculated at 1 CFU, *Micrococcus luteus* inoculated at 0 CFU, and *Cryptococcus neoformans* inoculated at 0 CFU. Three *Haemophilus influenzae* strains were species that were retested using 0.5 and 1 mL fresh instead of bagged blood and were detected in both glass and plastic vials.

The following organisms were evaluated in the analytical studies: *Abiotrophia defectiva*, *Acinetobacter lwoffii*, *Actinobacillus actinomycetemcomitans*, *Aerococcus viridans*, *Alcaligenes faecalis*, *Bacillus subtilis*, *Candida albicans*, *Candida glabrata*, *Candida tropicalis*, *Candida parapsilosis*, *Cardiobacterium hominis*, *Corynebacterium jeikeium*, *Cryptococcus neoformans*, *Eikenella corrodens*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Granulicatella adiacens*, *Haemophilus influenzae*, *Haemophilus influenzae*, type a, *Haemophilus influenzae*, type b, *Haemophilus parainfluenzae*, *Kingella kingae*, *Klebsiella pneumoniae*, *Leuconostoc mesenteroides*, *Micrococcus luteus*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Pediococcus acidilactici*, *Proteus mirabilis*, *Providencia stuartii*, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Stenotrophomonas maltophilia*, *Stomatococcus mucilaginosus*, *Streptococcus agalactiae*, four strains of *Streptococcus pneumoniae*, *Streptococcus pyogenes*, and *Streptococcus sanguinis* (formerly *S. sanguis*).

In microbial detection limit testing, a total of 360 paired sets at target inoculum levels of 0 to 1 and 1 to 10 CFU per vial were evaluated. This study was designed to assess the capability of the **BD BACTEC** blood culture media tested to detect one CFU, when present. Of the 360 paired sets tested, 196 grew and detected in both devices, 42 detected in glass bottles only, 57 detected in plastic bottles only, and 65 did not detect in either. There were a total of 107 pair sets that were not detected in plastic bottles, in which 36 showed organism growth on the inoculum plate: *Neisseria meningitidis* (5 CFU), *Haemophilus parainfluenzae* (4 CFU), *Staphylococcus epidermidis* (2 CFU), 1 CFU each for *Candida albicans*, *Candida glabrata*, *Escherichia coli*, *Staphylococcus aureus*, and *Streptococcus sanguinis*. The remaining 71 pair sets showed no organism growth (0 CFU) on the inoculum plate: *Cryptococcus neoformans*, *Enterococcus faecalis*, *Haemophilus influenzae*, *Haemophilus parainfluenzae*, *Micrococcus luteus*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Staphylococcus epidermidis*, and *Streptococcus pneumoniae*.

AVAILABILITY

Cat. No. Description

442020 **BD BACTEC™ Peds Plus™/F Medium**, Case of 50 vials

REFERENCES

1. Wallis, C. et al. 1980. Rapid isolation of bacteria from septicemic patients by use of an antimicrobial agent removal device. *J. Clin. Microbiol.* 11:462-464.
2. Applebaum, P.C. et al. 1983. Enhanced detection of bacteremia with a new BACTEC resin blood culture medium. *J. Clin. Microbiol.* 7:48-51.
3. Pohlman, J.K. et al. 1995. Controlled clinical comparison of Isolator and BACTEC 9240 Aerobic/F resin bottle for detection of bloodstream infections. *J. Clin. Microbiol.* 33:2525-2529.
4. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
5. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17: 53-80.
6. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
7. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
8. Flayhart, D. et al. 2007. Comparison of BACTEC Plus blood culture media to BacT/Alert FA blood culture media for detection of bacterial pathogens in samples containing therapeutic levels of antibiotics. *J. Clin. Microbiol.* 45:816-821.
9. Howden, R.J. 1976. Use of anaerobic culture for the improved isolation of *Streptococcus pneumoniae*. *J. Clin. Pathol.* 29:50-53.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

BD BACTEC Peds Plus/F Culture Vials (plastic)

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ATCC is a trademark of the American Type Culture Collection.
© 2016 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.
Becton Limited