

Rev from	Rev to	ECO #
0504	0904	3050-04

Notes:


1. BD Cat. Number 270308
2. Blank (Sheet) Size : Length: 22" Width: 13"
Number of Pages: 6 Number of Sheets: 1
Page Size: Length 22" Width 4.25" Final Folded Size: 2.75" x 4.5"
3. Style (see illustrations below): # 7


#1


#2


#3


#4

#5

#6


#7

4. See Specification Control Number 8012590 for Material Information
5. Ink Colors: Printed two sides Yes No
No. of Colors: 1 PMS# Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

VS Controlled By BD Caribe, LTD.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION		 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	Date			
Checked By	Date	Category and Description		Sheet: 1 of 7
Part Number: 8012590		Package Insert BBL GasPak CO ₂ System Envelopes		Scale: 1:1

A

BD BBL™ GasPak™ CO₂ System Envelopes

English: pages 1 Italiano: pagine 3 - 4 8012590
 Français : pages 1 - 2 Español: páginas 4 2004/09
 Deutsch: Seiten 2 - 3

See symbol glossary at end of insert. / Viz popis symbolů na konci příbalového letáku. / Se symbolglossaret i slutningen af indlægssedlen. / Zie lijst met symbolen aan het einde van de bijsluter. / Vaadake sümbole seletust infolehe lõpus. / Katso pakkausseloste lopussa olevaa kuvamerkkien sanastoa. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / A jelmagyarázat a használati utasítás végén található. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Žr. informaciniu lapelio pabaigoje pateikiamą simbolų glosarijų. / Se i symbolforklaringen på slutten av produktvedleget. / Zobacz objaśnienie symboli na końcu ulotki. / Consulte o glossário de símbolos no fim do folheto informativo. / Pozri slovník symbolov na konci letáka. / Consulte el glossario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD representant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviseletétől. / Naudojimo instrukcijų teiraukitės vietas BD įgaliotojo atstovo. / Aby uzyskać instrukcję użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contate o seu representante local da BD para obter instruções. / Instruções získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar.

INTENDED USE

BBL™ GasPak™ CO₂ System Envelopes are disposable carbon dioxide generating envelopes that produce an atmosphere of approximately 4 to 10% carbon dioxide in BBL GasPak jars suitable to support the growth of aerobic and aerotolerant microorganisms and tissue culture that benefit from a CO₂ enriched aerobic atmosphere.

SUMMARY AND EXPLANATION

The stimulatory effects of carbon dioxide on the growth and reproduction of bacteria have been studied extensively. Wherry and Oliver were among the first to report that the primary isolation of the gonococcus was greatly enhanced by an increased carbon dioxide atmosphere.¹ An environment enriched with up to 10% carbon dioxide is essential for optimal growth of *N. gonorrhoeae* and *N. meningitidis*.²⁻⁴ Optimum growth of *Mycobacterium* and *Brucella* with a 10% carbon dioxide-enriched environment has been documented.^{5,6} The stimulatory effect of 5 to 10% carbon dioxide on the growth of pneumococci has also been demonstrated.⁷

PRINCIPLES OF THE PROCEDURE

The GasPak CO₂ System Envelope is activated by the addition of water, which passes through a series of channels to a filter paper wick. The wick regulates the flow of water to the gas-generating tablet in the tablet chamber, providing a controlled release of gas. The wick delays the introduction of water into the tablet chamber while the lid is placed on the jar.

Approximately 4 to 10% carbon dioxide, measured 1 h after activation at room temperature (21 to 24°C), is generated by the sodium bicarbonate plus citric acid tablet. Catalyst is not required for this reaction.

REAGENTS

BBL GasPak CO₂ System Envelopes

Each GasPak carbon dioxide generator envelope contains: 1 Sodium bicarbonate plus citric acid tablet and 1 piece of filter paper.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

ANY UNSATISFACTORY ENVELOPE MUST BE CUT OPEN AND EXAMINED FOR INTACT OR PARTIALLY INTACT TABLETS. ANY REMAINING TABLETS SHOULD BE PLACED IN A CONVENIENT SINK AND DISSOLVED IN FLOWING WATER. DO NOT DISCARD INTACT ENVELOPES.

Storage Instructions: On receipt, store envelopes in a dry environment at 2-30°C. Avoid storage in close proximity to laboratory chemicals.

GasPak CO₂ System Envelopes are ready to use and should be activated according to instructions below and on the envelope. The expiration date is for unopened and intact envelopes.

Do not open until ready to use.

Product Deterioration: Do not use if envelopes show evidence of having been punctured, crushed, distorted, exposed to moisture, or exhibit other signs of deterioration.

PROCEDURE

Material Provided: GasPak CO₂ System Envelopes.

Materials Required But Not Provided: Culture media in Petri dishes, tubes, bottles or other suitable containers, Incubator (35°C), GasPak 100 system or GasPak 150 system, CO₂ Indicator Reagent Droppers, Pipette, 10 mL, or syringe.

Test Procedure:

1. Place inoculated plates or tubes (in GasPak tube holder) in the GasPak 100 or GasPak 150 rack. Cut off the corner of the envelope along the dotted line and place the envelope behind the shielded clip of the rack with the printed side toward the outside. One envelope is to be used in each GasPak 100 system and three envelopes are to be used in each GasPak 150 system. Do not crease, fold or crush the envelope.
2. Break the ampule of the CO₂ Indicator Reagent Dropper (Cat. No. 271055) per label instructions. Shake well. Place CO₂ Indicator on top of the stack of plates. *Do not place the indicator vertically along side the stack of plates.*
3. Add 10 mL of tap, distilled or deionized water through the open/cut corner of the envelope. Use a pipette or syringe inserting only the tip into the open corner. Never push any object down into the envelope.
4. Close the GasPak Jar promptly after envelope is activated and finger tighten the lid clamp in the manner described in the insert directions for the system being used. Incubate at conditions appropriate for the organisms being cultured, but not above 42°C.
5. After use, open the jar and allow it to aerate for approximately 15 s prior to removing the rack with its contents.

User Quality Control:

1. Each system should be tested periodically for its ability to provide adequate conditions for growth of carbon dioxide-requiring aerobic bacteria, e.g., *Neisseria gonorrhoeae*.

2. The following quality control check list will help assure proper performance:

- a. System should be upright at all times while in operation.

- b. **Envelope:** Keep upright and cut open along the dotted line indicated at the upper right-hand corner. The gas-generating tablet must be present and not displaced or crushed, and the foil channels and reaction chambers must not be crushed or distorted and must be able to accept 10 mL of water.

Improper positioning of envelope or use of a damaged envelope may lead to imbalances in reactivity. Damage may result in spattering and "boiling over" of contents of envelope and hence inadequate gas output.

- c. **Plastic Lid and Clamp:** Inspect for cracks, chips, split O-ring grooves and other irregularities. Most cracks (craze or stress lines) develop on prolonged use or excessive tightening of clamps on lids. Early cracks may appear as fine black lines near the center of the lid or along the upper edges of the GasPak 100 clamp near the screw. A cracked lid and clamp exhibiting these fine black lines should not be used.

- d. **O-Ring Gasket:** Examine for cracks, bumps, dirt and improper seating in the O-ring groove. Examine O-ring groove wall for splitting or cracking. The O-ring must be resilient, clean and snugly retained within the O-ring groove for satisfactory operation of the system. Replace as necessary.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

An atmosphere of approximately 4 to 10% carbon dioxide is generally obtained within 1 h.

other irregularities. Most cracks (craze or stress lines) develop on prolonged use or excessive tightening of clamps on lids. Early cracks may appear as fine black lines near the center of the lid or along the upper edges of the GasPak 100 clamp near the screw. A cracked lid and clamp exhibiting these fine black lines should not be used.

d. **O-Ring Gasket:** Examine for cracks, bumps, dirt and improper seating in the O-ring groove. Examine O-ring groove wall for splitting or cracking. The O-ring must be resilient, clean and snugly retained within the O-ring groove for satisfactory operation of the system. Replace as necessary.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

An atmosphere of approximately 4 to 10% carbon dioxide is generally obtained within 1 h.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of GasPak Carbon Dioxide Generator Envelopes are tested for specific product characteristics. Samples are analyzed for the amount of carbon dioxide emitted from the envelope using an Illinois CO₂ and O₂ Analyzer. Carbon dioxide levels are between 4 and 10%.

AVAILABILITY

Cat. No.	Description
270308	BBL™ GasPak™ CO ₂ System Envelopes, Disposable CO ₂ Generator Envelopes, Pkg. of 10 envelopes.
271055	BBL™ CO ₂ Indicator Reagent Droppers, Ctn. of 50 indicators.

REFERENCES

- Wherry, W.B., and W.W. Oliver. 1916. Adaptation to certain tensions of oxygen as shown by gonococcus and other parasitic and saprophytic bacteria. *J. Infect. Dis.* 19:288-298.
- Kellogg, D.S., K.K. Holmes, and G.A. Hill. 1976. Cumitech 4, Laboratory diagnosis of gonorrhea. Coordinating ed., S. Marcus and J.C. Sherris. American Society for Microbiology, Washington, D.C.
- Bovre, K. 1984. *Neisseriaceae*, p. 288-310. In N.R. Krieg, and J.G. Holt (ed.), Bergey's manual of systematic bacteriology, volume 1. Williams & Wilkins, Baltimore.
- Devaux, D.L., G.L. Evans, C.W. Arndt, and W.M. Janda. 1987. Comparison of the Gono-Pak system with the candle extinction jar for recovery of *Neisseria gonorrhoeae*. *J. Clin. Microbiol.* 25:571-572.
- Beam, E.R. and G.P. Kubica. 1968. Stimulatory effects of carbon dioxide on the primary isolation of tubercle bacilli on agar-containing medium. *Am. J. Clin. Pathol.* 50:395-397.
- Alton, G.G., L.M. Jones, and D.E. Pietz. 1975. W.H.O. laboratory techniques in brucellosis. W.H.O. Monogr. Ser. 55.
- Kempner, W., and C. Schlayer. 1942. Effect of CO₂ on the growth rate of the pneumococcus. *J. Bacteriol.* 43:387-396.

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BD, BD Logo, BBL and GasPak are trademarks of Becton, Dickinson and Company. © 2004 BD.
BENEX Limited