BD Difco™ QC Antigens Salmonella O Difco™ QC Antigen Salmonella Vi

CE 8085881(02)
2015-04

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INTENDED USE

BD Difco™ QC Antigens Salmonella O and **BD Difco QC Antigen Salmonella Vi** are used in the quality control testing of **BD Difco Salmonella O Antisera** and **BD Difco Salmonella Vi Antiserum** by slide agglutination tests.

SUMMARY AND EXPLANATION

The use of **BD Difco Salmonella O Antisera** and **BD Difco Salmonella Vi Antiserum** in the serologic identification of *Salmonella* requires the use of quality control test suspensions to verify that the antisera are performing as expected. Most laboratories are required to test antisera with positive and negative controls prior to use.^{1,2} **BD Difco QC Antigens Salmonella O** and **BD Difco QC Antigen Salmonella Vi** are chemically stabilized and inactivated suspensions of known strains of the genus *Salmonella* designed for use as homologous controls for testing the efficacy of the *Salmonella* typing antisera employed in routine laboratory procedures.

The chart below includes the **BD Difco QC Antigens Salmonella O** or **BD Difco QC Antigen Salmonella Vi** recommended as the homologous (positive) control antigen. The homologous control antigen has in common certain identifying antigen(s) with the antiserum.

The **BD Difco QC Antigens Salmonella O** and **BD Difco QC Antigen Salmonella Vi** may also be used to demonstrate a heterologous (negative) control reaction. For use as a negative control, use a QC Antigen Salmonella that contains antigens unrelated to those in the homologous control.

BD Difco Salmonella Antisera	BD Difco QC Antigen Homologous Control
Poly A-I & Vi	A, B, C, D, E ₁ , E ₂ , E ₄ , F, G ₁ , H, I, Vi
Poly A	A, B, C, D, E ₁ , E ₂ , E ₄
Poly B	C ₁ , C ₂ , F, G ₁ , H
Poly C	I
Group A Factors 1, 2, 12	A
Group B Factors 1, 4, 5, 12, 27	B
Group C ₁ Factors 6, 7, 14	C ₁
Group C ₂ Factors 6, 8, 20	C ₂
Group D ₁ Factors 1, 9, 12	D
Group E Factors 1, 3, 10, 15, 19, 34	E ₁ , E ₂ , E ₄
Group F Factor 11	F
Group G Factors 13, 22, 23	G ₁
Group H Factors 1, 6, 14, 24, 25	H
Group I Factor 16	I
Vi	Vi
Factor 2	A
Factor 4	B
Factors 4, 5	B
Factor 5	B
Factor 7	C ₁
Factor 8	C ₂
Factor 9	D

BD Difco Salmonella Antisera	BD Difco QC Antigen Homologous Control
Factor 10	E ₁
Factor 15	E ₂
Factor 19	E ₄
Factor 22	G ₁
Factor 14	H

Note: For a complete and current explanation of the classification of *Salmonella*, consult appropriate references.³⁻⁵

PRINCIPLES OF THE PROCEDURE

Serologic procedures that confirm the identification of an organism are usually agglutination reactions. Agglutination reactions may be either homologous or heterologous. Homologous reactions occur between a microorganism (antigen) and the corresponding antibody. These reactions occur rapidly and are strong. Heterologous reactions occur when a microorganism (antigen) reacts with an antibody produced in response to some other species or serotype. These reactions occur slowly and are weak.

Heterologous reactions may be unexpected and unpredictable and may lead to confusion in serologic identification. Therefore, only strongly positive homologous agglutination reactions should be regarded as significant.

REAGENTS

BD Difco QC Antigen	Homologous Identifying Antigen
Group A	2
Group B	4,5
Group C ₁	7
Group C ₂	8
Group D	9
Group E ₁	10
Group E ₂	15
Group E ₄	19
Group F	11
Group G ₁	22
Group H	14
Group I	16
Vi	Vi

When used as described (see PROCEDURE) each vial of antigen is sufficient for 20 slide tests.

BD Difco QC Antigens Salmonella O are preserved with 0.5% phenol USP (v/w) except for **BD Difco QC Antigen Salmonella Vi** which contains approximately 0.1% sodium azide.

Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Observe aseptic technique and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

Difco QC Antigens Salmonella O or **Difco QC Antigen Salmonella Vi** are not to be used for immunization of man or animals.

WARNING: BD Difco QC Antigen Salmonella Vi contains sodium azide. Sodium azide is toxic by inhalation, by skin contact, and if swallowed. Contact with acid liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Warning

H302 Harmful if swallowed. **P264** Wash thoroughly after handling. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage

Store **BD Difco QC Antigens Salmonella O** and **BD Difco QC Antigen Salmonella Vi** at 2 – 8 °C.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Expiration date applies to product in its intact container when stored as directed. Do not use if the product is discolored or shows other signs of deterioration.

BD Difco QC Antigens Salmonella O and **BD Difco** QC Antigen Salmonella Vi are stable through the expiration date on the label when stored as described. Exposure to temperatures below 2 °C can result in autoagglutination. Antigens must be smooth uniform suspensions; examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded.

PROCEDURE

Materials Provided: **BD Difco** QC Antigens Salmonella O, **BD Difco** QC Antigen Salmonella Vi.

Materials Required But Not Provided: Agglutination slides, Applicator sticks and Sterile 0.85% saline.

Reagent Preparation

BD Difco QC Antigens Salmonella O and **BD Difco** QC Antigen Salmonella Vi are ready for use. Shake well before use to suspend the organisms.

Before using the **BD Difco** QC Antigens Salmonella O or **BD Difco** QC Antigen Salmonella Vi, examine the **BD Difco** Salmonella Antisera Poly, Group or Factor chosen for use. The antisera must meet all specifications for the product (see product information for **BD Difco** Salmonella Antisera).

Equilibrate all materials to room temperature prior to performing the tests.

Ensure that all glassware and pipettes are clean and free of residues such as detergents.

Test Procedure

1. Positive control – Dispense 1 drop (approximately 35 µL) of the **BD Difco** Salmonella Antisera to be tested on an agglutination slide. After shaking well, add 1 drop of the appropriate **BD Difco** QC Antigens Salmonella O or **BD Difco** QC Antigen Salmonella Vi and mix thoroughly.
2. Negative control – Dispense 1 drop of sterile 0.85% saline on the agglutination slide. After shaking well, add 1 drop of the appropriate **BD Difco** QC Antigens Salmonella O or **BD Difco** QC Antigen Salmonella Vi and mix thoroughly.
3. Rotate the slides for 1 min and read for agglutination.

User Quality Control

At the time of use, apply both homologous and heterologous controls to check performance of the antigen, techniques and methodology.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Interpretation of Results

Read and record results as follows:

4+	100% agglutination (background clear to slightly hazy)
3+	75% agglutination (background slightly cloudy)
2+	50% agglutination (background moderately cloudy)
1+	25% agglutination (background cloudy)
–	No agglutination

LIMITATIONS OF THE PROCEDURE

1. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False positive reactions can occur.
2. Exposure to temperatures below 2 °C can cause autoagglutination. Antigens must be smooth, uniform suspensions. Examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded.
3. Allow the QC Antigen, the antisera and all equipment used to be at room temperature at the time of testing. The test reagents, if cold may cause false-negative reactions. Shake the antigen well before use to suspend the organisms.
4. **BD Difco** QC Antigens Salmonella O will react with their corresponding homologous **BD Difco** Salmonella Antisera Polyvalent, Group or Factor Antisera. Some reactions with single factor antisera may give reactions somewhat weaker than the polyvalent or grouping antisera due to specificity of the single factor antisera for the identifying antigen(s).
5. **BD Difco** Salmonella Antisera Polyvalent, Group or Factor Antisera may have a particular avidity for a given factor and agglutinate with another **BD Difco** QC Antigen Salmonella having the common factor.
6. Rough reactions can occur when using the **BD Difco** QC Antigens Salmonella O or **BD Difco** QC Antigen Salmonella Vi and 0.85% saline. If this occurs, product is not useable and should be discarded.

EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS

1. The positive control should show 3+ or greater agglutination.
2. The negative control should show no agglutination.

AVAILABILITY

Cat. No.	Description
221301	BD Difco™ QC Antigen Salmonella O Group A, 1 mL
221311	BD Difco™ QC Antigen Salmonella O Group B, 1 mL
221321	BD Difco™ QC Antigen Salmonella O Group C ₁ , 1 mL
221331	BD Difco™ QC Antigen Salmonella O Group C ₂ , 1 mL
221341	BD Difco™ QC Antigen Salmonella O Group D, 1 mL
211750	BD Difco™ QC Antigen Salmonella O Group E ₁ , 1 mL
211751	BD Difco™ QC Antigen Salmonella O Group E ₂ , 1 mL
221371	BD Difco™ QC Antigen Salmonella O Group E ₄ , 1 mL
211753	BD Difco™ QC Antigen Salmonella O Group F, 1 mL
211754	BD Difco™ QC Antigen Salmonella O Group G ₁ , 1 mL
221401	BD Difco™ QC Antigen Salmonella O Group H, 1 mL
211757	BD Difco™ QC Antigen Salmonella O Group I, 1 mL
221421	BD Difco™ QC Antigen Salmonella Vi, 1 mL

REFERENCES

1. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller and R.H. Tenover (ed.). 2003. Manual of clinical microbiology, 8th ed., vol. 1. American Society for Microbiology, Washington, D.C.
2. Isenberg, H.D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
3. Ewing, W.H. (ed.). 1986. Edwards and Ewing's identification of *Enterobacteriaceae*, 4th ed. Elsevier Science Publishing Co., Inc., New York, NY.
4. McWhorter-Murlin, A.C., and F.W. Hickman Brenner. 1994. Identification and serotyping of *Salmonella* and an update of the Kauffmann-White scheme. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, Atlanta, GA.
5. Popoff, M.Y., and L. LeMinor. 1997. Antigenic formulas of the *Salmonella* serovars. WHO Collaborating Centre for Reference and Research on *Salmonella*. Institut Pasteur, Paris, France.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Difco is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company. BD and BD Logo are trademarks of Becton, Dickinson and Company. © 2015 BD.
Benex Limited