BD Difco™ Haemophilus Influenzae Antisera

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INTENDED USE

Difco™ Haemophilus Influenzae Antisera are recommended for use in slide agglutination tests for serotyping *Haemophilus influenzae*.

SUMMARY AND EXPLANATION

H. influenzae constitutes part of the normal respiratory flora of humans and many animal species and can cause a variety of diseases from chronic respiratory infections to meningitis. Most of the *H. influenzae* isolates associated with meningitis belong to serotype b.¹

H. influenzae is a nonmotile, facultative anaerobe requiring both factor X (hemin) and factor V (nicotinamide adenine dinucleotide, NAD) for *in vitro* growth. In microscopic morphology, the organism is a pleomorphic gram-negative coccobacillus that sometimes forms threads or filaments.

The presence of a polysaccharide capsule is a major virulence factor for strains of *H. influenzae* that cause systemic infections. *H. influenzae* is divided into serologic groups a,b,c,d,e and f based on capsular polysaccharides. Most encapsulated strains that cause infection belong to serotype b.² The encapsulated strains are referred to as "typeable strains." Unencapsulated or "nontypeable" strains may also cause infection, but infections caused by unencapsulated strains are usually limited to the upper respiratory tract.²

Antigenic similarities exist between *H. influenzae* and a number of unrelated bacteria. For example, *H. influenzae* serotype b shares cross-reacting antigens with *Streptococcus pneumoniae* serotypes 6, 15a, 29 and 35a (see "LIMITATIONS OF THE PROCEDURE").

PRINCIPLES OF THE PROCEDURE

Identification of *H. influenzae* includes the isolation of the microorganism, biochemical identification and serologic confirmation.

Serologic confirmation involves the reaction in which the microorganism (antigen) reacts with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction should be rapid, tightly binding (high affinity) and not dissociate (high avidity).

Because a microorganism (antigen) may agglutinate with an antibody produced in response to some other species, heterologous reactions are possible. Such unexpected and perhaps unpredictable reactions may lead to some confusion in serologic identification. A positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism.

REAGENTS

Difco Haemophilus Influenzae Antisera are stable, lyophilized, polyclonal rabbit antisera containing approximately 0.1% sodium azide as a preservative.

When properly rehydrated and used as described (see PROCEDURE) each vial of **Difco** Haemophilus Influenzae Antiserum is sufficient reagent for 20 slide tests.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

This product contains Dry Natural Rubber.

Observe aseptic technique and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

WARNING: This product contains sodium azide. Sodium azide is toxic by inhalation, by skin contact, and if swallowed. Contact with acid liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Warning

H302 Harmful if swallowed. **P264** Wash thoroughly after handling. **P301 + P312** IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage: Store lyophilized and rehydrated **Difco™** Haemophilus Influenzae Antisera at 2–8 °C.

Expiration date applies to product in its intact container when stored as directed.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Discard any antiserum that is cloudy or has a precipitate after rehydration or storage unless it can be clarified by centrifugation or filtration and demonstrates proper reactivity with validated positive and negative controls. Gross turbidity indicates contamination and such antisera should be discarded.

SPECIMEN COLLECTION AND PREPARATION

H. influenzae can be recovered from clinical specimens on chocolate agar. For specific recommendations consult appropriate references.^{2,3} Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *H. influenzae*. After these criteria are met, serologic identification can be performed.

PROCEDURE

Materials Provided: **Difco** Haemophilus Influenzae Antiserum Poly, **Difco** Haemophilus Influenzae Antiserum Type a, **Difco** Haemophilus Influenzae Antiserum Type b, **Difco** Haemophilus Influenzae Antiserum Type c, **Difco** Haemophilus Influenzae Antiserum Type d, **Difco** Haemophilus Influenzae Antiserum Type e and **Difco** Haemophilus Influenzae Antiserum Type f.

Materials Required But Not Provided: Agglutination slides, applicator sticks, sterile purified water, sterile 0.85% saline and bacteriological loops.

Reagent Preparation

Equilibrate all materials to room temperature prior to performing the tests. Ensure that all glassware and pipettes are clean and free of residues, such as detergent.

To rehydrate **Difco** Haemophilus Influenzae Antisera, add 1 mL sterile purified water and rotate to completely dissolve the contents.

User Quality Control

At the time of use, test both positive and negative control cultures to check performance of the antisera, techniques and methodology.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Test isolate for autoagglutination

1. From the test culture on chocolate agar, transfer a loopful of growth to a drop (approximately 35 µL) of sterile 0.85% saline on a clean slide and emulsify the organism.
2. Rotate the slide for 1 min, then observe for agglutination.
3. If agglutination (autoagglutination) occurs, the culture is rough and cannot be tested. Subculture to chocolate agar, incubate and test the organism again as described in steps 1 and 2.

If no agglutination occurs, proceed with testing the organism as described (see Procedural Steps).

Procedural Steps

Culture isolates are tested with **Difco** Haemophilus Influenzae Antiserum Poly for presumptive identification followed by testing with monospecific antisera.

1. Put 1 drop of the **Difco** Haemophilus Influenzae Antiserum to be tested on an agglutination slide.
2. Transfer a loopful of growth of the test organism to the drop of antiserum and mix thoroughly.
3. Rotate the slide for 1 min and read for agglutination.

RESULTS

1. Read and record results as follows:
 - 4+ 100% agglutination (background clear to slightly hazy)
 - 3+ 75% agglutination (background slightly cloudy)
 - 2+ 50% agglutination (background moderately cloudy)
 - 1+ 25% agglutination (background cloudy)
 - No agglutination
2. The positive control should produce a 3+ or greater agglutination.
3. The negative control should show no agglutination.
4. For testing isolates, a 3+ or greater agglutination within 1 min is a positive result.

Positive agglutination using **Difco** Haemophilus Influenza Antiserum Poly and one of the individual typing sera combined with adequate biochemical results confirms the serotype of *H. influenzae* present.

LIMITATIONS OF THE PROCEDURE

1. Correct interpretation of serologic reactions depends upon culture purity and morphological characteristics and biochemical reactions that are consistent with identification of the microorganism as *H. influenzae*.
2. Serologic methods alone cannot identify the isolate as *H. influenzae*.
3. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False-positive reactions may occur.
4. Rough culture isolates do occur and will agglutinate spontaneously causing agglutination of negative controls (autoagglutination). Smooth colonies must be selected and tested in serologic procedures.
5. *H. influenzae* has antigenic similarities with several unrelated bacteria. Cross-reactions can occur between *H. influenzae* and strains of *S. pneumoniae*, *Escherichia coli* and several species of *Staphylococcus* and *Bacillus*.
6. **Difco** Haemophilus Influenzae Antisera have been tested using undiluted cultures taken from agar media. These antisera have not been tested using antigen suspensions in saline or other diluents. If the user applies variations in the recommended steps (see PROCEDURE), each lot of antiserum must be tested with known control cultures to verify expected reactions under the modified procedure.

PERFORMANCE CHARACTERISTICS⁴

The performance of **Difco**™ Haemophilus Influenzae Antiserum Type b, manufactured by Difco Laboratories, was compared in a slide agglutination test to four other tests for serotyping isolates of *Haemophilus influenzae* in a study by Himmelreich, Barenkamp and Storch.⁴ Slide agglutination using another commercial antiserum, coagglutination, latex agglutination and counterimmunoelectrophoresis with multiple antisera were also used to serotype isolates. Eighty (80) reference strains were used in this study. Of these isolates, 39 were *Haemophilus influenzae* type b.⁴

The table below lists the number of correctly identified isolates obtained for each test system.

Strain	Difco™ Haemophilus Influenzae Antiserum Type b Slide Agglutination Test	Commercial Slide Agglutination Test	Commercial Latex Agglutination Test	Commercial Coagglutination Test	Counter Immuno- electrophoresis
<i>H. influenza</i> type b	39/39	39/39	39/39	39/39	9/39
<i>H. influenza</i> not type b	37/37*	37/37*	37/37**	41/41	41/41
Organisms not tested	4	4	4	–	–

* 4 samples out of the 41 *H. influenzae* not type b autoagglutinated and could not be tested by the slide agglutination method.
 **4 samples out of the 41 *H. influenzae* not type b gave nonspecific reactions, agglutination of the isolate with both the *H. influenzae* type b and control latex preparations.

AVAILABILITY

Cat. No.	Description
222371	BD Difco™ Haemophilus Influenzae Antiserum Poly, 1 mL
222501	BD Difco™ Haemophilus Influenzae Antiserum Type a, 1 mL
222361	BD Difco™ Haemophilus Influenzae Antiserum Type b, 1 mL
227891	BD Difco™ Haemophilus Influenzae Antiserum Type c, 1 mL
227901	BD Difco™ Haemophilus Influenzae Antiserum Type d, 1 mL
227911	BD Difco™ Haemophilus Influenzae Antiserum Type e, 1 mL
227921	BD Difco™ Haemophilus Influenzae Antiserum Type f, 1 mL

REFERENCES

1. Insel, R. and P. Anderson. 1986. *Haemophilus influenzae* type b: assays for the capsular polysaccharide and for antipolysaccharide antibody, p. 379-384. In N.R. Rose, H. Friedman, and J.L. Fahey (ed.), Manual of clinical laboratory immunology, 3rd ed. American Society for Microbiology, Washington, D.C.
2. Campos, J.M. 1999. *Haemophilus*, p. 604-611. In P. R. Murray, E.J. Baron, M.A. Tenover, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
3. Pezzlo, M. 1992. Aerobic bacteriology, p.1.0.1-1.20.47. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol.1. American Society for Microbiology, Washington, D.C.
4. Himmelreich, C.A., S.J. Barenkamp and G.A. Storch. 1985. Comparison of methods for serotyping of *Haemophilus influenzae*. J. Clin. Microbiol. 21: 158-160.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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© 2016 BD
Difco is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company. BD and BD Logo are trademarks of Becton, Dickinson and Company.
Benex Limited