BD Difco™ Shigella Antisera Poly

English: pages 1–4 Italiano: pagine 11–14

Français : pages 4–7 Español: páginas 14–17

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 8085884(05)
2017-03

R_x Only

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD Difco™ Shigella Antisera Poly are used for the identification of *Shigella* species by the slide agglutination test. BD Difco Alkaescens-Dispar Antiserum Poly is used for the identification of Alkaescens-Dispar Group of microorganisms by the slide agglutination test.

SUMMARY AND EXPLANATION

Shigella species cause the human diarrheal disease shigellosis (classic bacillary dysentery). The range of illness is from mild diarrhea to severe dysentery characterized by abdominal cramps and frequent passage of bloody, mucoid stools. While the disease is usually self-limiting, it can be life threatening to the young, the elderly and malnourished persons. *Shigella* species are carried primarily in humans and are not generally distributed in nature. While transmission is usually direct person-to-person and through contaminated water supplies, foodborne outbreaks do occur.

The genus *Shigella* belongs to the family *Enterobacteriaceae*. *Shigella* species are facultatively anaerobic, gram-negative bacilli that typically are oxidase negative, lactose negative, H₂S negative and non-gas producing. *Shigella* and *Escherichia* are genetically related. Certain strains of *E. coli* may resemble *Shigella* biochemically because both can be lactose-negative, nonmotile or non-gas producing. These anaerogenic, nonmotile types have historically been called the Alkaescens-Dispar group and are presently classified as biotypes of *E. coli*. Consult appropriate references for biochemical tests specific for differentiating these strains from *Shigella*.¹⁻⁶

Serological testing with polyvalent and group specific antisera should be used to confirm the identification of isolates that are morphologically and biochemically identified as *Shigella* species. *Shigella* species are nonmotile, so serological identification is based on somatic (O) antigens. However, some strains have envelope antigens that prevent agglutination in somatic antisera. Heating the suspension at 100 °C for 30–60 minutes destroys these interfering antigens. The four named species or serotypes of *Shigella* are *S. dysenteriae* (10 serovars), *S. flexneri* (six serovars), *S. boydii* (15 serovars) and *S. sonnei*. For a complete and current explanation of the classification of *Shigella*, consult appropriate references.^{1,2}

PRINCIPLES OF THE PROCEDURE

Identification of *Shigella* species includes the isolation of the microorganism, biochemical identification and serological confirmation. Serological confirmation involves the reaction in which the microorganism (antigen) reacts with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, does not dissociate (high avidity) and binds strongly (high affinity).

Because a microorganism (antigen) may agglutinate with antibodies produced in response to other species, heterologous reactions are possible. Such unexpected and, perhaps, unpredictable reactions may lead to some confusion in serological identification. Therefore, a positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism. Homologous reactions occur rapidly and are strong. Heterologous reactions form slowly and are weak.

REAGENTS

BD Difco Shigella Antisera Poly and BD Difco Alkaescens-Dispar Antiserum Poly are lyophilized, polyclonal rabbit antisera containing approximately 0.2% sodium azide as a preservative.

BD Difco Shigella Antisera Poly are absorbed when necessary to render each lot of serum as specific as practical. Antisera are absorbed to a certain point without reducing homologous reactions to an unsatisfactory level. They have been absorbed inter- and intra-specifically except that *Shigella* antisera are not prepared from or tested for:

- S. dysenteriae* provisional serotypes,
- S. flexneri* X and Y variants, or
- Alkaescens-Dispar Groups

Antiserum**Reacts with**

Shigella Antiserum Poly Group A	<i>S. dysenteriae</i> types 1–7
Shigella Antiserum Poly Group A1	<i>S. dysenteriae</i> type 8ab, 8ac, 9, 10
Shigella Antiserum Poly Group B	<i>S. flexneri</i> types 1–6
Shigella Antiserum Poly Group C	<i>S. boydii</i> types 1–7
Shigella Antiserum Poly Group C1	<i>S. boydii</i> types 8–13
Shigella Antiserum Poly Group C2	<i>S. boydii</i> types 14–18
Shigella Antiserum Poly Group D	<i>S. sonnei</i> I and II
Alkalescens-Dispar Antiserum Poly	Alkalescens-Dispar Groups 1,2,3 and 4

When rehydrated and used as described, each vial of **BD Difco** Shigella Antisera Poly and **BD Difco** Alkalescens-Dispar Antiserum Poly contains sufficient reagent for 60 slide tests.

Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

WARNING: This product contains sodium azide. Sodium azide is toxic by inhalation, by skin contact, and if swallowed. Contact with acid liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Warning

H302 Harmful if swallowed.

P264 Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage

Store lyophilized and rehydrated **BD Difco** Shigella Antisera Poly and **BD Difco** Alkalescens-Dispar Antiserum Poly at 2–8 °C.

The expiration date applies to the product in its intact container when stored as directed.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Discard any antiserum that is cloudy or has a precipitate after rehydration or storage unless it can be clarified by centrifugation or filtration and demonstrates proper reactivity with validated positive and negative controls. Gross turbidity indicates contamination and such antisera should be discarded.

SPECIMEN COLLECTION AND PREPARATION

Clinical specimens: *Shigella* can be recovered on selective differential media such as Hektoen Enteric Agar or XLD Agar. For specific recommendations, consult appropriate references.²⁻⁴ Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as a *Shigella* species. After these criteria are met, serological identification can be performed.

Food samples: *Shigella* can be recovered from various types of foods when samples are processed to recover injured microorganisms and to prevent overgrowth of competing microorganisms. Consult appropriate references for recommended procedures when testing food samples.^{5,6} After following an established protocol, determine that a pure culture of the microorganism has been obtained. Biochemical test reactions should be consistent with the identification of the organism as a *Shigella* species. After these criteria are met, serological identification can be performed.

PROCEDURE

Materials Provided: **BD Difco** Shigella Antisera Poly, **BD Difco** Alkalescens-Dispar Antiserum Poly

Materials Required But Not Provided: Agglutination slides with 1-inch squares, Applicator sticks, Water Bath, boiling, Sterile 0.85% NaCl solution.

Reagent Preparation

Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents. To rehydrate, add 3 mL of sterile 0.85% NaCl solution and rotate gently to completely dissolve the contents.

Test isolate for autoagglutination

1. From the test culture on nonselective media, transfer a loopful of growth to a drop of sterile 0.85% saline on a clean slide and emulsify the organism.
2. Rotate the slide for one minute, then observe for agglutination.
3. If agglutination (autoagglutination) occurs, the culture is rough and cannot be tested. Subculture to nonselective agar, incubate and test the organism again as described in steps 1 and 2.
4. If no agglutination occurs, proceed with testing the organism.

Test Procedure

Use this procedure to test the isolate with each selected antiserum.

1. Dispense 1 drop (35 µL) of the antiserum to be tested on an agglutination slide.

- Negative control:** Dispense 1 drop of 0.85% sterile NaCl solution on an agglutination slide.
- From a solid agar medium, transfer a portion of a loopful of isolated colonies to each reaction area above and mix thoroughly.
- Quality control testing:** Dispense 1 drop of each **BD Difco** Shigella Antiserum to be tested on an agglutination slide. Add a portion of a loopful of an appropriate stock culture of known serological identification.
- Mix each reaction area with a separate applicator stick and rock for 1 minute. Read for agglutination.

User Quality Control

At the time of use, apply both homologous and heterologous controls to check performance of the antiserum, techniques and methodology. Stock cultures of known serological identification may be used as homologous and heterologous controls.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

- Read and record results as follows:
 - 4+ 100% agglutination (background is clear to slightly hazy)
 - 3+ 75% agglutination (background is slightly cloudy)
 - 2+ 50% agglutination (background is moderately cloudy)
 - 1+ 25% agglutination (background is cloudy)
 - No agglutination
- The positive control should show 3+ or greater agglutination.
- The negative control should show no agglutination.
- For test isolates, a 3+ or greater agglutination after 1 minute is a positive result.
- If no agglutination occurs or agglutination is weak, follow this procedure to remove blocking envelope antigens:
 - Prepare a dense suspension of the isolate from an agar medium in 3–5 mL of sterile 0.85% NaCl solution.
 - Heat in a boiling water bath for 30–60 minutes and cool. The suspension should not show precipitation after heating. If this occurs, select another colony for testing.
 - Centrifuge at 1,000 rpm for 10–15 minutes.
 - Aspirate and discard the supernatant.
 - Resuspend the sediment in 0.5 mL sterile 0.85% NaCl solution.
 - Use a drop of the suspension and perform the slide agglutination test as outlined above.
- A partial (less than 3+) or delayed agglutination reaction should be considered negative.
- If test results for either the quality controls or negative control are not as described, the test is invalid and results cannot be reported.

LIMITATIONS OF THE PROCEDURE

- Correct interpretation of serological reactions depends on culture purity, morphological characteristics and biochemical reactions that are consistent with identification of the microorganism as a *Shigella* species.
- Serological methods alone cannot identify the isolate as a *Shigella* species.
- Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent making a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False-positive reactions may occur.
- Rough culture isolates do occur and will agglutinate spontaneously, causing agglutination of the negative control (autoagglutination). Smooth colonies must be selected and tested in serological procedures.
- BD Difco** Shigella Antiserum Poly and **BD Difco** Alkalescens-Dispar Antiserum Poly have been tested using cultures taken directly from agar media. These antisera have not been tested using antigen suspensions in NaCl solution or other diluents. If the user applies variations in the recommended steps, each lot of antiserum must be tested with known control cultures to verify expected reactions under the modified procedure.

PERFORMANCE CHARACTERISTICS

The performance of **BD Difco** Shigella Antisera, manufactured by Difco Laboratories, was compared in slide agglutination tests to five other commercially available antisera and one commercially available latex agglutination test for serological grouping of *Shigella* in a study by Lefebvre, Gosselin, Ismail, Lorange, Lior and Woodward.⁷ Forty-six (46) *Shigella* reference strains, 50 clinical *Shigella* strains and 9 non-*Shigella* strains were used in this study.⁷

The table below lists the results obtained.

Product	Sensitivity ^a	Specificity ^b			
		Group A	Group B	Group C	Group D
BD Difco Antisera	94 (85/90)	99 (90/91)	96 (67/70)	100 (86/86)	98 (81/83)
Commercial Antisera #1	85 (66/78)	100 (94/94)	97 (68/70)	100 (95/95)	99 (82/83)
Commercial Antisera #2	93 (88/95)	97 (86/89)	97 (68/70)	95 (79/83)	95 (79/83)
Commercial Antisera #3	85 (71/84)	99 (96/97)	97 (68/70)	88 (76/86)	98 (81/83)
Commercial Antisera #4	47 (42/90)	100 (91/91)	100 (70/70)	100 (86/86)	99 (82/83)
Commercial Antisera #5	91 (77/85)	98 (91/93)	99 (69/70)	98 (87/89)	96 (80/83)
Latex Agglutination Test	91 (84/92)	100 (89/89)	100 (70/70)	100 (86/86)	99 (82/83)

a Denominators varied according to the number of homologous strains recognized by the commercial antisera.

b Total number of non-*Shigella* and heterologous *Shigella* strains tested.

In this study, for the strains evaluated, the overall sensitivity of the **BD Difco** Shigella Antisera was 94%. The specificity for **BD Difco** Shigella Group A was 99%, **BD Difco** Shigella Group B was 96%, **BD Difco** Shigella Group C was 100% and **BD Difco** Shigella Group D was 98%.⁷

In this study, cross-reactivity was observed between:

- *P. shigelloides* and group D antisera
- *H. alvei* and *S. flexneri* and group B antisera
- *E. coli* ATCC® 13706 and group D antisera
- *E. coli* 0157:H19 and group B antisera

AVAILABILITY

Cat. No.	Description
228341	BD Difco ™ Shigella Antiserum Poly Group A, 3 mL
227761	BD Difco ™ Shigella Antiserum Poly Group A1, 3 mL
228351	BD Difco ™ Shigella Antiserum Poly Group B, 3 mL
228361	BD Difco ™ Shigella Antiserum Poly Group C, 3 mL
227771	BD Difco ™ Shigella Antiserum Poly Group C1, 3 mL
227781	BD Difco ™ Shigella Antiserum Poly Group C2, 3 mL
228371	BD Difco ™ Shigella Antiserum Poly Group D, 3 mL
228381	BD Difco ™ Alkaescens-Dispar Antiserum Poly, 3 mL

REFERENCES

1. Ewing, W.H. (ed.). 1986. Edwards and Ewing's identification of *Enterobacteriaceae*, 4th ed. Elsevier Science Publishing Co., Inc., New York, N.Y.
2. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller and R.H. Tenover (ed.). 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
3. Forbes, B.A., D.F. Sahm, A.S. Weissfeld and E.A. Trevino. 2002. Bailey & Scott's diagnostic microbiology, 11th ed. Mosby-Year Book, Inc., St. Louis, Mo.
4. Pezzlo, M. (ed.). 1994. Aerobic bacteriology, p. 1.0.1-1.20.47. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
5. Association of Official Analytical Chemists. 2001. FDA Bacteriological analytical manual online. <<http://www.cfsan.fda.gov/~ebam/bam-mm.html>>.
6. Downes, F.P. and K. Ito (eds.). 2001. Compendium of methods for the microbiological examination of foods, 4th edition. American Public Health Association, Washington, D.C.
7. Lefebvre, J., f. Gosselin, J. Ismail, M. Lorange, H. Lior and D. Woodward. 1995. Evaluation of commercial antisera for *Shigella* serogrouping. J. Clin. Microbiol. 33:1997-2001.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najnižja dovoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrense / Dolna granica temperatury / Limite minimo de temperatura / Limitá minimá de temperaturá / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sicaklık alt sınırı / Мінімальна температура / 温度下限

Control / Контролно / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Controllo / Контроль / kontrol / Контроль / 对照

Positive control / Положителен контрол / Pozitivní kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positiive kontroll / Contrôle positif / Pozitivna kontrola / Pozitiv kontroll / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiamą kontrolę / Pozitívá kontrolé / Positiive controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контрол / Pozitif kontrol / Позитивный контрол / 阳性对照试剂

Negative control / Отрицателен контрол / Negativní kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negatiive kontroll / Contrôle négatif / Negativna kontrola / Negativ kontroll / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiamą kontrolę / Negatívá kontrolé / Negatiive controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контрол / Negatif kontrol / Негативный контрол / 阴性对照试剂

Method of sterilization: ethylene oxide / Метод на стерилизация: этиленов оксид / Způsob sterilizace: etylenoxid / Steriliseringmethode: ethylenoxid / Sterilisationsmethode: Ethylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξειδίο / Método de esterilización: óxido de etileno / Steriliseerimismetod: etüleenoksiid / Méthode de stérilisation: oxyde d'éthylène / Metoda sterilizacije: etilen oksid / Sterilizálás módszere: etilén-oxid / Metodo di sterilizzazione: ossido di etilene / Стерилизация әдісі – этилен тотығы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizēšanas metode: etilēnoksīds / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringsmetode: etylenoksid / Metoda sterylizacji: tienek etylu / Método de esterilizaçõ: óxido de etileno / Metodá de sterilizare: oxid de etilená / Метод стерилизації: этиленоксид / Metodá sterilizácie: etylénoxid / Metoda sterilizacije: etilen oksid / Steriliseringsmetod: etenoxid / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: этиленоксидом / 灭菌方法: 环氧乙烷

Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringsmetode: bestrålning / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseerimismetod: kiirgus / Méthode de stérilisation: irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Стерилизация әдісі – сәуле түсіру / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringsmetode: bestrålning / Metoda sterylizacji: napromienianie / Método de esterilizaçõ: irradiação / Metodá de sterilizare: iradiere / Метод стерилизації: облучение / Metodá sterilizácie: ožiarenie / Metoda sterilizacije: ozračavanje / Steriliseringsmetod: strålning / Sterilizasyon yöntemi: ıradasyon / Метод стерилизації: опроміненням / 灭菌方法: 辐射

Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Bioloogilised riskid / Risques biologiques / Biološki rizik / Biológiallag veszélyes / Rischio biologico / Биологиялық төуекелдер / 생물학적 위험 / Biologinis pavojus / Biologiskie riski / Biologisk risiko / Biologisk risiko / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Биологическая опасность / Biologické riziko / Biološki rizici / Biologisk risk / Biyolojik Riskler / Биологична небезпека / 生物学风险

Caution, consult accompanying documents / Внимание, направте справка в придружаващите документи / Pozor! Prostudujte si příloženou dokumentaci! / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα / Precaución, consultar la documentación adjunta / Ettevaatust! Lugeda kaasnevat dokumentatsiooni / Attention, consulter les documents joints / Upozorenje, koristiti prateću dokumentaciju / Figyelem! Olvassa el a mellékelt tájékoztatót / Attensjon, se consultare la documentazione allegata / Абайланыз, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dmesio, žiūrėkite priedamus dokumentus / Piesardzība, skatīt pavaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsiktig, se vedlagt dokumentasjon / Należy zwracać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Vůstraha, pozri sprievodné dokumenty / Pažnja! Pogledajte priložena dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlikle verilen belgelerle bagvurun / Увага: див. супутню документацію / 小心, 请参阅附带文档。

Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Øvre temperaturgrænse / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite superior de temperatura / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dovoljena temperatura / Felső hőmérsékleti határ / Limite superiore di temperatura / Температураның рұқсат етілген жоғарғы шегі / 상한 온도 / Aukščiausia laikymo temperatūra / Augšējā temperatūras robeža / Hoogste temperatuurlimiet / Øvre temperaturgrense / Górná granica temperatury / Limite máximo de temperatura / Limitá maximá de temperaturá / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Øvre temperaturgräns / Sicaklık üst sınırı / Мінімальна температура / 温度上限

Keep dry / Παзете сухо / Skladujte v suchém prostředí / Opbevares tørt / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivais / Conserver au sec / Držati na suhom / Száraz helyen tartandó / Tenere all'asciutto / Κυρρακ κүйінде ұста / 건조 상태 유지 / Laikykite sausiai / Uzglabāt sausu / Droog houden / Holdes tørt / Przechowywać w stanie suchym / Manter seco / A se feri de umezeală / Не допускать попадания влаги / Uchovávať v suchu / Držite na suvom mestu / Förvaras tørt / Kuru bir şekilde muhafaza edin / Беретти від вологи / 请保持干燥

Collection time / Време на събиране / Čas odběru / Opsamlingsstidspunkt / Entnahmeuhrzeit / Ωρα συλλογής / Hora de recogida / Kogumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жинау уақыты / 수집 시간 / Paémimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de colheita / Ora colectării / Время сбора / Doba odberu / Vreme prikupljanja / Uppsamlingsstid / Toplama zamanı / Час забору / 采集时间

Peel / Обелете / Otevířete zde / Åbn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Ўстиғи қабатын алып таста / 벗기 / Plésti čia / Atlīmēt / Schillen / Trekk av / Oderwaç / Destacar / Se dezlipeste / Отклеить / Odtrhnite / Oljuštiti / Dra isär / Ayırma / Відклеїти / 撕下

Perforation / Перфорация / Perforace / Perforering / Διείρση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Тесик тесу / 결취선 / Perforacja / Perforăcija / Perforatje / Perforacja / Perforação / Perforare / Перфорация / Perforácia / Perforasyon / Перфорация / 穿孔

Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nerouživejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packungnicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιό. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használnia, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Ереп пакер бұзылған болса, пайдаланба / 패키지가 손상된 경우 사용 금지 / Jei pakuotė pažeista, nenaudoti / Nelletot, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не использовать при повреждении упаковки / Nepoužívať, ak je obal poškodený / Не користіть ако je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar gölmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包装破损, 请勿使用

Keep away from heat / Παзете от топлина / Nevystavujte přílišnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od izvora topline / Övja a melegtől / Tenere lontano dal calore / Салқын жерде сақта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargăt no karstuma / Beschermen tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródła ciepła / Manter ao abrigo do calor / A se feri de căldură / Не нагревать / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Беретти від дії тепла / 请远离热源

Cut / Срежете / Odstrhňte / Klip / Schneiden / Κόψτε / Cortar / Lögata / Découper / Rezi / Vágja ki / Tagliare / Kesçiz / 잘라내기 / Kirpti / Nogriet / Knippen / Kutt / Odciąć / Cortar / Decupați / Отрезать / Odstrhните / Iseći / Klipp / Kesme / Pozřizati / 剪下

Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskuupeäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тзбекүні / 수집 날짜 / Paémimo data / Savākšanas datums / Verzameldatum / Dato prøvetaking / Data pobrania / Data de colheita / Data colectării / Дата сбора / Datum odberu / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ATCC is a trademark of the American Type Culture Collection.

Difco is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company.

© 2017 BD. BD, and the BD Logo are trademarks of Becton, Dickinson and Company.