BD Difco™ Listeria Antisera and Antigens

English: pages 1 – 4 Italiano: pagine 11 – 14

Français : pages 4 – 7 Español: páginas 14 – 17

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 8085890(04)
2016-03

R_x Only

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Inštrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD Difco™ Listeria O Antisera Types 1, 4 and Poly are used for identifying *Listeria monocytogenes* in the slide agglutination test.

BD Difco Listeria O Antigens Types 1 and 4 (Slide) are used as positive controls in the slide agglutination test.

SUMMARY AND EXPLANATION

First described in 1926 by Murray, Webb and Swann,¹ *Listeria monocytogenes* is a widespread problem in public health and the food industries. This organism can cause human illness and death, particularly in immunocompromised individuals and pregnant women.² The first reported foodborne outbreak of listeriosis was in 1985³ and, since then, microbiological and epidemiological evidence from both sporadic and epidemic cases of listeriosis has shown that the principal route of transmission is via the consumption of foodstuffs contaminated with *L. monocytogenes*.⁴

The organism has been isolated from commercial dairy and other food processing plants and is ubiquitous in nature, being present in a wide-range of unprocessed foods and in soil, sewage, silage and river water.⁵

Listeria spp. grow over a pH range of 4.4–9.6 and survive in food products with pH levels outside these parameters.⁶ *Listeria* spp. are microaerophilic, gram-positive, asporogenous, non-encapsulated, non-branching, regular, short, motile rods. Motility is most pronounced at 20 °C.

Strains of *Listeria* spp. are divided into serotypes based on cellular (O) and flagellar (H) antigens.⁷ Thirteen serotypes of *L. monocytogenes* are known. Most human disease is caused by serotypes 1/2a, 1/2b and 4b.⁸

PRINCIPLES OF THE PROCEDURE

Identification of *L. monocytogenes* includes both biochemical and serological confirmation.

Serological confirmation requires that the microorganism (antigen) react with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, does not dissociate (high avidity) and binds strongly (high affinity).

Because a microorganism (antigen) may agglutinate with an antibody produced in response to another species, heterologous reactions are possible. Such unexpected and perhaps unpredictable reactions may lead to some confusion in serological identification. A positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism.

Agglutination of the somatic antigen in the slide test appears as a firm granular clumping. Homologous reactions occur rapidly and are strong (3+). Heterologous reactions form slowly and are weak.

REAGENTS

BD Difco Listeria O Antisera Types 1, 4 and Poly are lyophilized, polyclonal rabbit antisera containing approximately 0.1% sodium azide as a preservative. Listeria O Antisera Types 1 and 4 are specific for the respective serotypes of *L. monocytogenes* while Listeria O Antiserum Poly contains agglutinins for *L. monocytogenes* serotypes 1 and 4.

BD Difco Listeria O Antigens Types 1 and 4 (Slide) are suspensions of appropriate *L. monocytogenes* serotypes containing 0.3% formaldehyde as a preservative. When used according to the suggested procedure, the reagents will yield the following:

Reagent	Vial	Number of tests
Listeria O Antiserum	1 mL	400 slide tests
Listeria O Antigen (Slide)	5 mL	100 slide tests

Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

WARNING: This product contains sodium azide. Sodium azide is toxic by inhalation, by skin contact, and if swallowed. Contact with acid liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Warning

Cat. No.	Description
223001	BD Difco™ Listeria O Antiserum Type 1, 1 mL
223011	BD Difco™ Listeria O Antiserum Type 4, 1 mL
223021	BD Difco™ Listeria O Antiserum Poly, 1 mL

H302 Harmful if swallowed.

P264 Wash thoroughly after handling. **P301+P312** IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Danger

Cat. No.	Description
223031	BD Difco™ Listeria O Antigen Type 1 (Slide), 5 mL
223041	BD Difco™ Listeria O Antigen Type 4 (Slide), 5 mL

H317 May cause an allergic skin reaction. **H350** May cause cancer.

P261 Avoid breathing dust/fume/gas/mist/vapors/spray. **P280** Wear protective gloves.

P302+P352 IF ON SKIN: Wash with plenty of water. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage: Store lyophilized and rehydrated **BD Difco** Listeria O Antisera at 2–8 °C.

Store **BD Difco** Listeria O Antigen (Slide) at 2–8 °C.

The expiration date applies to the product in its intact container when stored as directed.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Discard any antiserum that is cloudy or has a precipitate after rehydration or storage unless it can be clarified by centrifugation or filtration and demonstrates proper reactivity with validated positive and negative controls. Gross turbidity indicates contamination and such antisera should be discarded.

Examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded.

SPECIMEN COLLECTION AND PREPARATION

Clinical specimens: *Listeria* can be recovered on selective differential media such as Oxford Agar, Modified Oxford Agar, LPM Agar or PALCAM Medium. For specific recommendations on isolation of *Listeria* from clinical specimens, consult appropriate references.⁸⁻¹⁰ Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *L. monocytogenes*. After these criteria are met, serological identification can be performed.

Food or dairy samples: *Listeria* can be recovered when samples are processed to recover injured microorganisms and prevent overgrowth of competing microorganisms. Consult appropriate references for recommended procedures for the isolation of *Listeria* from foods.^{6,11,12} Having followed an established protocol, isolate a pure culture of the microorganism and confirm that biochemical test reactions are consistent with the identification of the organism as *L. monocytogenes*. After these criteria are met, serological identification can be performed.

PROCEDURE

Materials Provided: **BD Difco** Listeria O Antiserum Type 1, **BD Difco** Listeria O Antiserum Type 4, **BD Difco** Listeria O Antiserum Poly, **BD Difco** Listeria O Antigen Type 1 (Slide), **BD Difco** Listeria O Antigen Type 4 (Slide).

Materials Required But Not Provided:

Slide Test: **BD Difco** FA Buffer (Dried), agglutination slides, applicator sticks, water bath (80–100 °C), droppers.

Reagent Preparation: Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents.

Listeria O Antisera: To rehydrate, add 1 mL sterile purified water to each vial. Rotate gently to dissolve contents completely.

Listeria O Antigens (Slide) are ready to use.

Test Procedure

Slide Test

1. Rehydrate **BD Difco** FA Buffer as per label directions.
2. Suspend a loopful of the test isolate growth from a solid noninhibitory agar medium in approximately 5 mL of **BD Difco** FA Buffer.
3. Heat the organism suspension at 80–100 °C (in a water bath) for 1 h.
4. Centrifuge the suspension at 2,000–5,000 RCF for 15–20 min and remove the bulk of the supernatant fluid.
5. Resuspend the organism in the remaining portion of supernatant fluid.
6. Prepare a 1:20 dilution of the desired **BD Difco** Listeria Antiserum in 0.85% NaCl solution.
7. On an agglutination slide, dispense 1 drop of the desired diluted **BD Difco** Listeria Antiserum onto each of two separate areas. Dispense 1 drop of **BD Difco** FA Buffer onto a third area of the same slide. The first drop of antiserum will be used for the test isolate and the second for the positive control. The drop of **BD Difco** FA Buffer will be used for the negative control.
8. Dispense 1 drop of organism suspension from step 5 to the first drop of antiserum.
9. **Positive control:** Dispense 1 drop of homologous **BD Difco** Listeria O Antigen (Slide) to the second drop of antiserum.
10. **Negative control:** Dispense 1 drop of organism suspension from step 5 to the drop of **BD Difco** FA Buffer on the agglutination slide.
11. Rotate the slide for 1–2 min and read for agglutination.

User Quality Control: At the time of use, test both positive and negative antigen controls to check performance of the antiserum, techniques and methodology. If results of the positive control or negative control are not as described, the test is invalid and results cannot be reported.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Slide Test

1. Read and record results as follows:

4+	100% agglutination; background is clear to slightly hazy.
3+	75% agglutination; background is slightly cloudy.
2+	50% agglutination; background is moderately cloudy.
1+	25% agglutination; background is cloudy.
–	No agglutination
2. For the test isolate, a 3+ or greater agglutination within 1–2 min is a positive result.
3. A partial (less than 3+) or a delayed agglutination reaction should be considered negative.

LIMITATIONS OF THE PROCEDURE

1. Serological techniques employing **BD Difco** Listeria O Antisera serve as corroborative evidence for the identification of *L. monocytogenes*. Final identification cannot be made without consideration of morphological, serological and biochemical characterization.
2. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent making a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False positive reactions may occur.
3. Rough culture isolates occur and will agglutinate spontaneously, causing agglutination of the negative control (autoagglutination). If autoagglutination occurs, the culture is rough and cannot be tested. Subculture to a noninhibitory medium, incubate and test the organism again.
4. Agglutination reactions of 3+ or greater in the slide test are interpreted as positive reactions. Cross-reactions resulting in a 1+ or 2+ agglutination are likely since there are somatic antigens shared among different organisms such as staphylococci, enterococci and *Bacillus* species.¹⁰
5. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
6. Exposure of Listeria O Antigens to temperatures below 2 °C can result in autoagglutination. Antigens must be smooth uniform suspensions. Examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded.
7. Discard any **BD Difco** Listeria O Antiserum that is cloudy or has a precipitate after rehydration or storage.

PERFORMANCE CHARACTERISTICS

1. The positive control, using the appropriate Listeria O Antigen, should show 3+ or greater agglutination in the slide test.
2. The negative control should show no agglutination.

Sensitivity of **BD Difco** Listeria O Antisera Types 1, 4 and Poly is determined by demonstrating appropriate reactivity in the slide test, as defined in the "Results" section using the respective

homologous **BD Difco** Listeria O Antigen. Specificity is determined by demonstrating non-reactivity against non-related (heterologous) **BD Difco** Listeria O Antigens.

AVAILABILITY

Cat. No.	Description
223001	BD Difco ™ Listeria O Antiserum Type 1, 1 mL
223011	BD Difco ™ Listeria O Antiserum Type 4, 1 mL
223021	BD Difco ™ Listeria O Antiserum Poly, 1 mL
223031	BD Difco ™ Listeria O Antigen Type 1 (Slide), 5 mL
223041	BD Difco ™ Listeria O Antigen Type 4 (Slide), 5 mL
223143	BD Difco ™ FA Buffer, Dried, 6 x 10 g
223142	BD Difco ™ FA Buffer, Dried, 100 g

REFERENCES

1. Murray, E.G.D., R.A. Webb and M.B.R. Swann. 1926. A disease of rabbits characterized by large mononuclear leucocytosis caused by a hitherto undescribed bacillus *Bacterium monocytogenes* (n. sp.). *J. Path. Bact.* 29:407 439.
2. Monk, J.D., R.S. Clavero, L.R. Beuchat, M.P. Doyle and R.E. Brackett. 1994. Irradiation inactivation of *Listeria monocytogenes* and *Staphylococcus aureus* in low and high fat, frozen and refrigerated ground beef. *J. Food Prot.* 57:969 974.
3. Wehr, H.M. 1987. *Listeria monocytogenes* – a current dilemma special report. *J. Assoc. Off. Anal. Chem.* 70:769 772.
4. Bremer, P.J., and C.M. Osborne. 1995. Thermal death times of *Listeria monocytogenes* in green shell mussels (*Perna canaliculus*) prepared for hot smoking. *J. Food Prot.* 58:604 608.
5. Patel, J.R., C.A. Hwang, L.R. Beuchat, M.P. Doyle and R.E. Brackett. 1995. Comparison of oxygen scavengers for their ability to enhance resuscitation of heat injured *Listeria monocytogenes*. *J. Food Prot.* 58:244 250.
6. Ryser, E.T., and C.W. Donnelly. 1992. *Listeria*, p. 343–356. In F.P. Downes and K. Ito (ed.), *Compendium of methods for the microbiological examination of foods*, 4th ed. American Public Health Association, Washington, D.C.
7. Seeliger, H.P.R., and K. Hohn. 1979. Serotyping of *Listeria monocytogenes* and related species, p. 31–49. In T. Bergen and J.R. Norris (ed.), *Methods in microbiology*, vol. 13. Academic Press, London, England.
8. Bille, J., J. Rocourt and B. Swaminathan. 2003. *Listeria* and *Erysipelothrix*, p. 461–471. In P.R. Murray, E.J. Baron, J.H. Jorgensen, M.A. Pfaller and R.H. Tenover (ed.), *Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.
9. Pezzlo, M. 1994. Aerobic bacteriology, p. 1.0.1.–1.20.47. In H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol. 1. American Society for Microbiology, Washington, D.C.
10. Forbes, B.A., D.F. Sahm, A.S. Weissfeld and E.A. Trevino. 2002. *Bailey & Scott's diagnostic microbiology*, 11th ed. Mosby-Year Book, Inc., St. Louis, Mo.
11. Association of Official Analytical Chemists. 2001. *FDA Bacteriological analytical manual online*. <<http://www.cfsan.fda.gov/~ebam/bam-mm.html>>.
12. Flowers, R. S., W. Andrews, C. W. Donnelly, and E. Koenig. 1993. Pathogens in milk and milk products. In Marshall, R. T., (ed.), *Standard methods for the examination of dairy products*, 16th ed. American Public Health Association, Washington, D.C.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Difco is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company.
BD and BD Logo are trademarks of Becton, Dickinson and Company. © 2016 BD
Benex Limited