# **BD BBL™ Sensi-Disc™ Antimicrobial Susceptibility Test Discs**

English: pages 1 – 2      Italiano: pagine 5 – 6  
Français: pages 2 – 3      Español: páginas 6 – 8  
Deutsch: Seiten 4 – 5

 8085891(02)  
2015-04

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкциите. / Pokynu vám poskytné místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítás kérje a BD helyi képviselőtől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o reprezentante locale da BD para instruções. / Pentru instrucțiunile, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

**INTENDED USE** – These discs are used for semi-quantitative *in vitro* susceptibility testing by the agar disc diffusion test procedure of fungal pathogens.

**SUMMARY AND EXPLANATION** – Agar diffusion methods employing dried filter paper discs impregnated with specific concentrations of antimicrobial agents were developed in the 1940s. In order to eliminate or minimize variability in this testing, Bauer et al. developed a standardized procedure in which Mueller Hinton Agar was selected as the test medium.<sup>1,2</sup>

Various regulatory agencies and standards-writing organizations subsequently published standardized reference procedures based on the Bauer-Kirby method. Among the earliest and most widely accepted of these standardized procedures were those published by the U.S. Food and Drug Administration (FDA)<sup>3</sup> and the World Health Organization (WHO).<sup>4,5</sup> The procedure was adopted as a consensus standard by the Clinical and Laboratory Standards Institute (CLSI) and is periodically updated.<sup>6,7</sup>

A disk diffusion method for testing *Candida* species was developed and in 2004, CLSI Approved Guideline M44-A, *Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts*, was released.<sup>8</sup> The disk diffusion method often provides qualitative results 24 hours sooner than broth dilution, making antifungal susceptibility testing more readily available to some clinical laboratories and providing a reduced cost alternative. **The latest CLSI documents should be consulted for current recommendations.**

**PRINCIPLES OF THE PROCEDURE** – Discs containing antifungal agents are applied to the surface of Mueller Hinton Agar supplemented with 2% Glucose and 0.5 µg/mL Methylene Blue Dye (GMB) plates that have been inoculated with pure cultures of clinical isolates. Following incubation, the plates are examined and the zones of inhibition surrounding the discs are measured and compared with established zone size ranges for individual antifungal agents in order to determine the agent(s) most suitable for use in antimicrobial therapy.<sup>8</sup>

**REAGENTS – BD BBL™ Sensi-Disc™** brand discs are 6-mm discs prepared by impregnating high quality absorbent paper with accurately determined amounts of antibiotic or other chemotherapeutic agents. Discs are clearly marked on both sides with letters and numbers designating the agent and the drug content. (See chart containing concentrations of reactive ingredients). The drug content of discs is assayed by the methods established by the FDA or by methods similar or comparable to those published in the United States *Federal Register*.

**BD BBL Sensi-Disc** agents are furnished in cartridges containing 50 discs each. The last disc in each cartridge is marked "X" and contains the drug as coded. Cartridges are for use in **BD BBL Sensi-Disc** Dispensers; these include a Single Disc Dispenser, an 8-Place Dispenser for 100 mm-style Petri dishes, 6- and 8-Place Self-Tamping Dispensers for 100 mm-style dishes and a Self-Tamping 12-Place Dispenser for 150 mm-style plates.

**Warnings and Precautions:** For *in vitro* Diagnostic Use.

Follow directions for use; disc performance depends not only on disc potency, but on use of proper inoculum and control cultures, functional pretested plates, proper storage temperature and other factors.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. Sterilize cultures, containers and other contaminated materials after use.

#### **Storage Instructions:**

1. On receipt, store discs at -20 – +8°C.
2. Allow containers to come to room temperature before opening. Return unused discs to the refrigerator when application of discs has been completed. Once a cartridge of discs has been removed from its sealed packaging, it should be placed in a tightly sealed, desiccated container.
3. Use the oldest discs first.
4. Discard expired discs. Also, cartridges from which discs have been frequently removed during a week and discs left out overnight in the laboratory should be discarded, or else the discs should be tested for acceptable performance prior to continued use.
5. If the discs form incorrect zones with the recommended control organisms, the entire procedure should be checked; faulty zone size may be due to the disc, the inoculation, the preparation or depth (about 4 mm) of medium, or other factors.

The expiration date applies only to discs in intact containers, stored as directed.

**SPECIMENS** – Direct specimens should not be employed in this test. Test must be used with isolated cultures. See Directions, which include preparation of inoculum. If possible, cultures should be derived from specimens obtained from patients prior to the initiation of antifungal therapy.

#### **PROCEDURE**

**Material Provided:** Sensi-Disc™ susceptibility test discs as labeled.

**Materials Required But Not Provided:** Ancillary culture media, reagents, quality control organisms and laboratory equipment required to perform disc diffusion susceptibility testing by the standardized procedure. Prepare a 0.5 McFarland turbidity standard by adding 0.5 mL of 0.048 M BaCl<sub>2</sub> [1.175% (wt/vol) BaCl<sub>2</sub>·2H<sub>2</sub>O] to 99.5 mL of 0.18 M [0.36N] H<sub>2</sub>SO<sub>4</sub> [1% (vol/vol)]. Verify by using a spectrophotometer with a 1-cm light path and matched cuvette; absorbance at 625 nm should be 0.08 – 0.13.

#### **Directions, Including User Controls:**<sup>8</sup>

1. Preparation of inoculum with test and control cultures.
  - a. Select five distinct colonies of approximately 1 mm in diameter from a 24-hour-old culture of *Candida* species grown on blood agar or Sabouraud dextrose agar. Transfer with inoculation needle or loop into 5 mL of sterile 0.145 M NaCl (8.5g/L NaCl; 0.85% saline). Vortex for 15 s.
  - b. Dilute or add more colonies, if required, to obtain turbidity equivalent to the 0.5 McFarland turbidity standard. Alternatively, standardize the inoculum photometrically.<sup>8</sup>
2. Inoculation.
  - a. Within 15 min, dip a sterile cotton swab into the properly adjusted inoculum and rotate it firmly several times against the upper inside wall of the tube to express excess fluid.
  - b. Streak the entire agar surface of a Mueller Hinton + GMB agar plate three times, turning the plate 60° between streakings to obtain even inoculation.
  - c. The lid may be left ajar for 3 – 5 min, but no more than 15 min, to allow for any surface moisture to be absorbed before applying the drug-impregnated discs.
3. Select appropriate discs (such as recommended in CLSI M44, Tables 1 and 2).
4. Apply the discs by means of a **BD BBL** dispenser, using aseptic precautions. Deposit discs so that the centers are at least 24 mm apart. If discs have been placed on the agar with other than the Self-Tamping Dispensers, press them down with a sterile needle or forceps to make contact with the surface.
5. Within 15 min, place the plates agar side up in a 35 ± 2°C incubator.
6. Examine the plates after 20 – 24 h of incubation. The diameters of the zones of complete inhibition are measured, as determined by gross visual inspection. Zones are measured to the nearest whole millimeter. For further details in measuring zones of inhibition, consult CLSI M44.<sup>8</sup> If only isolated colonies grow, the inoculum is too light and the test should be repeated. Zones around discs containing different drugs are not comparable for the purpose of comparing activity of drugs. Read at 48 h only when insufficient growth is observed after 24 h of incubation.
7. Control tests using recommended cultures should be included each day susceptibility testing is performed or weekly if satisfactory performance can be documented according to the CLSI standard.<sup>8</sup>

**RESULTS**<sup>8</sup>– NOTE: Recommended interpretive criteria are based on usual dosage regimens and routes of administration in the U.S.

The susceptible (S) category implies that isolates are inhibited by the usually achievable concentrations of antifungal agent tested when the recommended dosage is used for the site of infection. The intermediate (I) category implies that an infection due to the isolate may be appropriately treated in body sites where the drugs are physiologically concentrated or when a high dosage of drug is used. The resistant (R) category implies that isolates are not inhibited by the usually achievable concentrations of the agent with normal dosage schedules and clinical efficacy of the agent against the isolate has not been reliably shown in treatment studies.

#### **LIMITATIONS OF THE PROCEDURE**

1. The classifications of resistant (R), intermediate (I) and susceptible (S) vary only by one millimeter, which is within normal laboratory error. Some cultures may give a borderline zone that varies from day to day or from laboratory to laboratory; such cultures are relatively uncommon.
2. Antifungal agents other than those listed in the Chart may be in current use. Susceptibility tests employing these agents should be interpreted on the basis of presence or absence of a definite zone of inhibition and should be considered as only qualitative until such time as interpretive zones have been established. All zone diameters should be recorded.

#### **AVAILABILITY**

Cat. No.	Description
232201	<b>BD BBL™ Sensi-Disc™</b> Fluconazole, 25 µg
232202	<b>BD BBL™ Sensi-Disc™</b> Voriconazole, 1 µg

#### **REFERENCES**

1. Bauer, A.W., W.M.M. Kirby, J.C. Sherris, and M. Turck. 1966. Antibiotic susceptibility testing by a standardized single disk method. *Am. J. Clin. Pathol.* 45:493-496.
2. Ryan, K.J., F.D. Schoenkecht, and W.M.M. Kirby. 1970. Disc sensitivity testing. *Hospital Practice* 5:91-100.
3. Federal Register. 1972. Rules and regulations. Antibiotic susceptibility discs. *Fed. Regist.* 37:20525-20529. Erratum, 38:2756, 1973.
4. Ericsson, H.M., and J.C. Sherris. 1971. Antibiotic sensitivity testing. Report of an international collaborative study. *Acta Pathol. Microbiol. Scand. Sec. B. Suppl.* 217:1-90.
5. World Health Organization Expert Committee on Biological Standardization. 1977. Technical report series 610. W.H.O., Geneva.
6. Clinical and Laboratory Standards Institute. 2006. Approved standard M2-A9. Performance standards for antimicrobial disk susceptibility tests, 9th ed. CLSI, Wayne, Pa.
7. Clinical and Laboratory Standards Institute. 2008. M100-S18 (M2). Disk Diffusion Supplemental Tables, CLSI, Wayne Pa.
8. Clinical and Laboratory Standards Institute. 2004. Approved guideline M44-A. Method for antifungal disk diffusion susceptibility testing of yeasts, CLSI, Wayne Pa.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or [www.bd.com/ds](http://www.bd.com/ds).

Zone Diameter Interpretive Chart†									
			Zone Diameter Interpretive Standards (mm)			Control Zone Diameter Limits (mm)			
Antifungal Agent	Disc Code	Disc Potency	R	I	S	<i>C. albicans</i> ATCC® 90028	<i>C. parapsilosis</i> ATCC 22019	<i>C. tropicalis</i> ATCC 750	<i>C. krusei</i> ATCC 6258
Fluconazole <sup>a</sup>	FCA-25	25 µg	≤14	15-18	≥19	28-39	22-33	26-37	– <sup>b</sup>
Voriconazole	VOR-1	1 µg	≤13	14-16	≥17	31-42	28-37	– <sup>b</sup>	16-25

† Interpretive criteria and quality control recommendations from the FDA-approved drug labeling. Interpretive criteria and quality control adapted in part from CLSI Document M44-A: Disk Diffusion Tables, Performance Standards for Susceptibility Testing of Yeasts, with permission. The complete standard may be obtained from the Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA

<sup>a</sup> Isolates of *C. krusei* are assumed to be intrinsically resistant to fluconazole, and their MICs should not be interpreted using this scale.

<sup>b</sup> A quality control range has not been established for this strain/antimicrobial agent combination due to its extensive interlaboratory variation during initial quality control studies.

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Benex Limited