# **BD BACTEC™ Lytic/10 Anaerobic/F Culture Vials** **Soybean-Casein Digest Broth in a Plastic Vial**

English: pages 1 – 4      Italiano: pagine 10 – 13  
Français : pages 4 – 7      Español: páginas 13 – 16  
Deutsch: Seiten 7 – 10      Português: páginas 16 – 19

8089974(06)  
2016-07

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytné místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukciju teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinize temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

## INTENDED USE

**BD BACTEC™ Lytic/10 Anaerobic/F** culture vials (prereduced enriched Soybean-Casein Digest broth with CO<sub>2</sub>) are for anaerobic blood cultures. Principal use is with the **BD BACTEC** fluorescent series instruments for the qualitative culture and recovery of anaerobic microorganisms from blood.

## SUMMARY AND EXPLANATION

The sample to be tested is inoculated into one or more vials which are inserted into the **BD BACTEC** fluorescent series instrument for incubation and periodic reading. Each vial contains a chemical sensor which can detect increases in CO<sub>2</sub> produced by the growth of microorganisms. The sensor is monitored by the instrument every ten minutes for an increase in its fluorescence, which is proportional to the amount of CO<sub>2</sub> present. A positive reading indicates the presumptive presence of viable microorganisms in the vial. Detection is limited to microorganisms that will grow in a particular type of medium.

## PRINCIPLES OF THE PROCEDURE

If microorganisms are present in the test sample inoculated into the **BD BACTEC** vial, CO<sub>2</sub> will be produced when the organisms metabolize the substrates present in the vial. Increases in the fluorescence of the vial sensor caused by the higher amount of CO<sub>2</sub> are monitored by the **BD BACTEC** fluorescent series instrument. Analysis of the rate and amount of CO<sub>2</sub> increase enables the **BD BACTEC** fluorescent series instrument to determine if the vial is positive; i.e., that the test sample contains viable organisms.

## REAGENTS

The **BD BACTEC Lytic/10 Anaerobic/F** culture vials contain the following active ingredients prior to processing:

### List of Ingredients

Processed Water .....	40 mL
Soybean-Casein Digest Broth .....	2.75% w/v
Yeast Extract .....	0.2% w/v
Animal Tissue Digest .....	0.05% w/v
Dextrose .....	0.2% w/v
Hemin .....	0.0005% w/v
Menadione .....	0.00005% w/v
Sodium Citrate .....	0.02% w/v
Thiols .....	0.1% w/v
Sodium Pyruvate .....	0.1% w/v
Saponin .....	0.26% w/v
Antifoaming Agent .....	0.01% w/v
Sodium Polyanetholsulfonate (SPS) .....	0.035% w/v

All **BD BACTEC** media are dispensed with added CO<sub>2</sub>. Anaerobic media are prereduced and dispensed with added CO<sub>2</sub> and N<sub>2</sub>. Composition may have been adjusted to meet specific performance requirements.

## Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

**Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>1-4</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.**

Prior to use, each vial should be examined for evidence of damage, contamination or deterioration. Vials displaying evidence of damage or contamination such as leakage, cloudiness, discoloration (darkening), bulging or depressed septum should not be used.

A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, contaminated culture media could be refluxed into the patient's vein. Vial contamination may not be readily apparent. When using direct draw procedures, monitor the process closely to avoid refluxing materials into the patient.

On rare occasions, a vial may not be sealed sufficiently, which may result in the contents of the vial leaking or spilling. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: Before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling should be performed in a biological safety cabinet if possible, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or **BD Luer-Lok™** brand tips.

### Storage Instructions

The **BD BACTEC** vials are ready for use as received and require no reconstitution or dilution. Store at 2–25 °C in a dry location, **out of direct light**.

### SPECIMEN COLLECTION

The specimen must be collected using sterile techniques to reduce the chance of contamination. Published studies have shown that the recommended specimen volume is 8–10 mL.<sup>5,6</sup> It is recommended that the specimen be inoculated into the **BD BACTEC** vials at bedside. Most commonly, a 10 cc or 20 cc syringe with a **BD Luer-Lok** brand tip is used to draw the sample. If appropriate, a **BD Vacutainer®** brand Needle Holder and a **BD Vacutainer** brand Blood Collection Set, **BD Vacutainer Safety-Lok™** Blood Collection Set or other tubing "butterfly" set may be used. If using a needle and tubing set (direct draw), carefully observe the direction of blood flow when starting sample collection. The vacuum in the vial will usually exceed 10 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. When the desired 8–10 mL has been drawn, the flow should be stopped by crimping the tubing and removing the tubing set from the **BD BACTEC** vial. Sample volumes as low as 3 mL can be used, however, recovery will not be as great as with larger volumes. **The inoculated BD BACTEC vial should be transported as quickly as possible to the laboratory.**

### PROCEDURE

Remove the flip-off cap from **BD BACTEC** vial top and inspect the vial for cracks, contamination, excessive cloudiness, and bulging or indented stoppers. **DO NOT USE** if any defect is noted. Before inoculating, swab the septum with alcohol (iodine is **not** recommended). Aseptically inject or draw directly 8–10 mL of specimen per vial. If sample volumes of 3–4 mL are used, recovery will not be as great as with larger volumes (see Limitations of the Procedure). **Inoculated anaerobic vials should be placed in the BD BACTEC fluorescent series instrument as soon as possible** for incubation and monitoring. If placement of an inoculated vial into the instrument has been delayed and visible growth is apparent, it should not be tested in the **BD BACTEC** fluorescent series instrument, but rather it should be subcultured, Gram-stained and treated as a presumptively positive bottle.

Vials entered into the instrument will be automatically tested every ten minutes for the duration of the testing protocol period. Positive vials will be determined by the **BD BACTEC** fluorescent series instrument and identified as such (see the appropriate **BD BACTEC** fluorescent series instrument User's Manual). The sensor inside the bottle will not appear visibly different in positive and negative vials, however the **BD BACTEC** fluorescent series instrument can determine a difference in fluorescence.

Blood will lyse immediately upon addition to the **BD BACTEC** Lytic/10 Anaerobic/F Medium. The blood will appear chocolized or very dark initially. If at the end of the testing period a **BD BACTEC** Lytic/10 Anaerobic/F vial is observed to have a bulging septum, it should be subcultured, Gram-stained or treated as presumptive positive.

Positive vials should be subcultured and a Gram-stained slide prepared. In a great majority of cases, organisms will be seen and a preliminary report can be made to the physician. Subcultures to selective media and a preliminary direct antimicrobial susceptibility test may be prepared from fluid in the **BD BACTEC** vials.

**Subculturing:** Prior to subculturing, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, insert a sterile needle with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

For maximum yield of isolates, negative cultures may be checked by stain and/or subcultured at some point prior to discarding as negative.

### QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

**DO NOT USE** culture vials past their expiration date.

**DO NOT USE** culture vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each carton of media. Quality Control Certificates list test organisms, including ATCC® cultures specified in the CLSI Standard, *Quality Control for Commercially Prepared Microbiological Culture Media*.<sup>7</sup>

The range of time-to-detection in hours was ≤ 72 hours for each of the organisms listed on the Quality Control Certificate for this medium:

*Clostridium perfringens* ATCC 13124  
*Streptococcus pneumoniae* ATCC 6305  
*Clostridium histolyticum* ATCC 19401

*Bacteroides fragilis*\* ATCC 25285  
*Escherichia coli* ATCC 25922

*Bacteroides vulgatus* ATCC 8482  
*Staphylococcus aureus* ATCC 25923

\*CLSI Strain

For information on Quality Control for the **BD BACTEC** fluorescent series instrument, refer to the appropriate **BD BACTEC** fluorescent series instrument User's Manual.

## LIMITATIONS OF THE PROCEDURE

### Contamination

Care must be taken to prevent contamination of the sample during collection and inoculation into the **BD BACTEC** vial. A contaminated sample will give a positive reading, but will not indicate a relevant clinical result. Such a determination must be made by the user based on such factors as type of organisms recovered, occurrence of the same organism in multiple cultures, patient history, etc.

### Recovery of SPS Sensitive Organisms from Blood Samples

Because blood can neutralize the toxicity of SPS toward organisms sensitive to SPS (such as *P. anaerobius*), the presence of maximum volumes of blood (i.e., up to 10 mL) can help to optimize recovery of these organisms. To enhance the growth of SPS sensitive organisms when less than 8 mL of blood is inoculated, additional whole human blood may be added.

Some fastidious organisms, such as certain *Haemophilus* species, require growth factors, such as NAD, or factor V, which are provided by the blood specimen. If the blood specimen volume is 3.0 mL or less, an appropriate supplement may be required for recovery of these organisms.

**BD BACTEC FOS™** Fastidious Organism Supplement may be used as a nutritional supplement.

### Nonviable Organisms

A Gram-stained smear from culture medium may contain small numbers of nonviable organisms derived from media constituents, staining reagents, immersion oil, glass slides, and specimens used for inoculation. In addition, the patient specimen may contain organisms that will not grow in the culture medium or in media used for subculture. Such specimens should be subcultured to special media as appropriate.<sup>9</sup>

### General Considerations

Optimum recovery of isolates will be achieved by adding 8–10 mL of blood.<sup>5,6</sup> Use of lower or higher volumes may adversely affect recovery and/or detection times. Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms. False negative readings may result when certain organisms are present which do not produce enough CO<sub>2</sub> to be detected by the system or significant growth has occurred before placing the vial into the system. False positivity may occur when the white blood cell count is high. The default 5 day protocol was utilized for all analytical testing with this device and protocols longer than 5 days have not been evaluated.

## EXPECTED RESULTS

Performance of the **BD BACTEC** Lytic/10 Anaerobic/F medium in glass vials has been established by a number of published external clinical studies.<sup>9,10</sup> Seeded laboratory studies performed by BD have shown equivalent performance of the **BD BACTEC** Lytic/10 Anaerobic/F medium in plastic vials compared to **BD BACTEC** Lytic/10 Anaerobic/F medium in glass vials.<sup>11</sup>

A total of 342 paired sets at 10 to 100 CFU per vial were evaluated for recovery, with 100% recovering in both the **BD BACTEC** Lytic/10 Anaerobic/F medium contained in a plastic vial and the **BD BACTEC** Lytic/10 Anaerobic/F medium contained in a glass vial. This study included a diverse set of anaerobic and aerobic microorganisms frequently isolated in blood. The median time to detection (TTD) difference between the paired sets was 10 minutes, in favor of the **BD BACTEC** Lytic/10 Anaerobic/F medium contained in a plastic vial. Ninety-five percent of the TTD differences between the paired sets were between -1.68 hours faster for the glass vial and 3 hours faster for the plastic vial.

The following anaerobes were evaluated in the analytical studies: *Bacteroides fragilis*, *B. ovatus*, *B. thetaiotaomicron*, *B. vulgatus*, *Clostridium histolyticum*, *C. novyi*, *C. perfringens*, *Fusobacterium nucleatum*, *Porphyromonas asaccharolytica* (formerly *Bacteroides melaninogenicus* subsp. *asaccharolyticus*) and *Veilonella parvula*. The facultative anaerobe *S. pneumoniae* was also tested.

A subset of organisms, including *Finegoldia magna* (formerly *Peptostreptococcus magnus*) and *Peptoniphilus asaccharolyticus* (formerly *Peptostreptococcus asaccharolyticus*) were evaluated on the **BD BACTEC** FX instrument at 10 to 100 CFU per vial and demonstrated 100% recovery in both the **BD BACTEC** Lytic/10 Anaerobic/F medium contained in a plastic vial and the **BD BACTEC** Lytic/10 Anaerobic/F medium contained in a glass vial.

In microbial detection limit testing, a total of 312 paired sets at inoculum levels of 0 to 1 and 1 to 10 CFU per vial were evaluated. This study was designed to assess the capability of the **BD BACTEC** blood culture media tested to detect one CFU, when present. Of the 312 paired sets tested, 191 grew and detected in both devices and 44 did not detect in either. Twenty-nine (29) cultures grew and detected only in the **BD BACTEC** Lytic/10 Anaerobic/F medium contained in a glass vial. Forty-eight (48) cultures grew and detected only in **BD BACTEC** Lytic/10 Anaerobic/F medium contained in a plastic vial. One of 12 replicates of *Porphyromonas asaccharolytica* (ATCC 25260, 4 CFU per bottle) failed to detect in the **BD BACTEC** Lytic/10 Anaerobic/F medium contained in a plastic vial. Signal analysis demonstrated no evidence of growth in the replicate and a terminal subculture yielded no growth; indicating that there were likely no viable organisms inoculated into the vial.

## AVAILABILITY

Cat. No.	Description
442021	<b>BD BACTEC™</b> Lytic/10 Anaerobic/F Culture Vials

## REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.

3. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
5. Lin, H-H, et al. 2012. Evaluation of the blood volume effect on the diagnosis of bacteremia in automated blood culture systems. Doi:10.1016/j.jmii. 2012. 03.2012.
6. Reimer, L.G. et al. 1997. Update on detection of bacteremia and fungemia. Clin. Micro. Rev. 10:444-465.
7. Clinical and Laboratory Standards Institute. 2004. Approved Standard M22-A3. Quality control of commercially prepared microbiological culture media, 3rd ed., CLSI, Wayne, Pa.
8. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller (ed.). 2007. Manual of clinical microbiology, 9th ed. American Society for Microbiology, Washington, D.C.
9. Hollick, GE, et al. 1996. Diagnostic Microbiology and Infectious Disease Journal. 24:191-196.
10. Rohner, P. et al. 1997. Advantage of combining resin with Lytic BACTEC blood culture media. J. Clin. Micro. 35:2634-2638.
11. Data available from BD Life Sciences.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA

Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

#### Australian Sponsor:

Becton Dickinson Pty Ltd.  
4 Research Park Drive  
Macquarie University Research Park  
North Ryde, NSW 2113  
Australia

ATCC is a trademark of the American Type Culture Collection.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.