88-0890-1
Revised: July 1999

BBL™ Antimicrobial Supplement for Selective Isolation of *Listeria monocytogenes*

<i>Listeria</i> Selective Supplement	Pkg. of 10 vials, lyophilized; each reconstitutes to 2 mL.	Cat. No. 212402
--------------------------------------	--	------------------------

English:	pages	1 – 2	Italiano:	pagine	4 – 5
Français:	pages	2 – 3	Español:	páginas	5 – 6
Deutsch:	Seiten	3 – 4			

INTENDED USE – *Listeria* Selective Supplement is an antimicrobial supplement for addition to basal media for the selective isolation of *Listeria monocytogenes*, e.g., LPM Agar Base and Oxford Agar Base, Modified.

SUMMARY AND EXPLANATION – *L. monocytogenes* is the etiologic agent for a variety of animal and human infections.¹ *L. monocytogenes* has also been implicated in outbreaks of human listeriosis in which the vehicle of infection was milk.² Various foods have served as sources of *Listeria* in epidemics of listeriosis.³

LPM Agar was developed by Lee and McClain⁴ as a modification of McBride Agar for the isolation of *L. monocytogenes*. Oxford Agar Base, Modified⁵ is a modification of the Oxford *Listeria* selective medium developed by Curtis et al.⁶ The selective supplement contains the antimicrobial agent, moxalactam. When the selective supplement is added to these basal media, the complete media are suitable for use in the isolation of *L. monocytogenes*, particularly from beef⁴ and processed meat and poultry products.⁵

PRINCIPLES OF THE PROCEDURE – *Listeria* Selective Supplement (moxalactam) is added to sterile basal medium to inhibit both gram-positive and gram-negative species, including *Staphylococcus*, *Proteus* and *Pseudomonas*.

REAGENT

Approximate Formula* per 1 mL of Reconstituted Solution: Moxalactam..... 0.01 g

*Adjusted and/or supplemented as required to meet performance criteria.

Precautions: For Laboratory Use.

This product contains dry natural rubber.

Listeria Selective Supplement is for use in culture media and not for use in human or animal therapy. Observe aseptic techniques in the reconstitution and use of this medium supplement.

Storage Instructions and Reconstitution: On receipt, store at -20 to +8°C. Protect from light. After reconstitution, use immediately.

Reconstitute each lyophilized vial by aseptically adding 2.0 mL of sterile purified water with a sterile syringe and needle. Invert the vial several times to assure complete solution.

The expiration date applies to the lyophilized product, stored as directed.

Product Deterioration: Examine reconstituted reagents at time of use for evidence of contamination, evaporation or other signs of deterioration.

PROCEDURE

Material Provided: *Listeria* Selective Supplement.

Materials Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Instructions: Prepare the basal agar medium according to the label directions. Cool the sterile basal medium to 45 to 50°C. Aseptically add 2.0 mL of the reconstituted selective supplement to the 1 L of sterile cooled basal agar medium. Mix gently but thoroughly and distribute into sterile Petri dishes. Inoculate the medium for isolation by streaking with test samples, preferably previously enriched by incubation in UVM Modified *Listeria* Enrichment Broth and then in Fraser Broth, Modified. Incubate the plates at 35 ± 2°C for 24 to 48 h in an aerobic atmosphere. Refer to the reference for additional information.⁵

User Quality Control: Examine the reconstituted supplement as described under "Product Deterioration." Check performance of the finished medium by inoculation with pure cultures of stable control organisms. The following cultures are recommended:

TEST STRAIN	INCUBATION 24 to 48 h at 35 ± 2°C	
	LPM Agar	Oxford Agar, Modified
<i>Listeria monocytogenes</i> ATCC® 19112	Growth; colonies will show blue to blue-green iridescence when held obliquely under intense light.	Growth; colonies black with blackening of the medium.
<i>Staphylococcus aureus</i> ATCC 25923	Inhibition (partial to complete).	Inhibition (partial to complete); no black zones.
<i>Escherichia coli</i> ATCC 25922	Inhibition (partial to complete).	Inhibition (partial to complete); no black zones.
<i>Pseudomonas aeruginosa</i> ATCC 10145	Inhibition (partial to complete).	Inhibition (partial to complete); no black zones.

RESULTS – Observe plates for the presence of growth. After 24 and 48 h of incubation, colonies of *L. monocytogenes* on LPM Agar are small, circular, smooth, entire, raised, translucent and exhibit blue-green iridescence when examined with oblique transmitted light. After 24 h of incubation, *Listeria* strains on Oxford Agar, Modified form black colonies with blackening of the medium. After 48 h, colonies are black with sunken centers and surrounded by black zones.

LIMITATIONS OF THE PROCEDURE – A single medium is rarely adequate for detecting all organisms of potential significance in a specimen. The agents in selective media may inhibit some strains of the desired species or permit growth of a species they were designed to inhibit, especially if the species is present in large numbers in the specimen. Specimens cultured on selective media should, therefore, also be cultured on nonselective media to obtain additional information and help ensure recovery of potential pathogens.

REFERENCES

1. Bille, J., and M.P. Doyle. 1991. *Listeria* and *Erysipelothrix*, p. 287-295. In A. Balows, W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.), Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
2. Hayes, P.S., J.C. Feeley, L.M. Graves, G.W. Ajello, and D.W. Fleming. 1986. Isolation of *Listeria monocytogenes* from raw milk. *Appl. Environ. Microbiol.* 51:438-440.
3. Schlech, W.F. III, P.M. Lavigne, R.A. Bortolussi, A.C. Allen, E.V. Haldane, A.J. Wort, A.W. Hightower, S.E. Johnson, S.H. King, E.S. Nicholls, and C.V. Broome. 1983. Epidemic listeriosis — evidence for transmission by food. *N. Engl. J. Med.* 308:203-206.
4. Lee, W.H., and D. McClain. 1986. Improved *Listeria monocytogenes* selective agar. *Appl. Environ. Microbiol.* 52:1215-1217.
5. McClain, D., and W.H. Lee. May 24, 1989. FSIS method for the isolation and identification of *Listeria monocytogenes* from processed meat and poultry products. Laboratory communication no. 57. Microbiology Division, Food Safety and Inspection Service, U.S. Department of Agriculture, Beltsville, Md.
6. Curtis, G.D.W., R.G. Mitchell, A.F. King, and E.J. Griffin. 1989. A selective differential medium for the isolation of *Listeria monocytogenes*. *Appl. Microbiol.* 8:95-98.

TECHNICAL INFORMATION: In the United States telephone Technical Services, toll free (800) 638-8663.

© 1999 Becton Dickinson and Company

BBL is a trademark of Becton Dickinson and Company.
ATCC is a trademark of the American Type Culture Collection.

**BECTON
DICKINSON**

Becton Dickinson Microbiology Systems
Becton Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA

Becton Dickinson France S.A.
11 rue Aristide Bergès
38800 Le Pont de Claix, France

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U.S.A.