88-0933-1  
Revised: July 1999**BBL™ SEPTI-CHEK™ AFB**  
**Mycobacteria Culture System**English: page 1 Italiano: pagin 4  
Français : page 2 Español: página 5  
Deutsch: Seite 3 Japanese: page 6**INTENDED USE**

The **BBL™ SEPTI-CHEK™** AFB Culture Bottle, **BBL™ SEPTI-CHEK™** AFB Supplement and **BBL™ SEPTI-CHEK™** AFB Slide are intended for use as an integrated, contained, in vitro diagnostic system for the detection and isolation of mycobacteria from sputum, bronchial washings or aspirate, body fluids (pleural, CSF, ascites, synovial), urine, stool, biopsy tissues or wounds and skin when mycobacteriosis is suspected.

**SUMMARY AND EXPLANATION**

Mycobacteria are infectious agents responsible for persistent and chronic clinical conditions most commonly of a respiratory nature, and some are spread by coughing and expectoration.<sup>1,2</sup> As a group, these organisms are capable of infecting nearly every organ of the body. This has become increasingly evident and frequent in cases of underlying immune dysfunction associated with conditions such as HIV infection and drug abuse.<sup>1,2</sup>

Historically, the most frequently encountered human pathogenic *Mycobacterium* is *M. tuberculosis*.<sup>1,2</sup> In some populations, this has been exceeded by the *Mycobacterium avium* complex.<sup>1,2</sup>

The mycobacteria have several characteristics which distinguish them from other bacteria. The most notable of these are resistance to decolorization after staining, relatively slow growth, and resistance to digesting agents.<sup>1</sup> Conventional culture media used for mycobacteria isolation are most often Lowenstein-Jensen (LJ), or Middlebrook 7H9, 7H10 or 7H11.<sup>1</sup> Growth in these media is observed in a range of 3 – 56 days, depending on the species isolated and concentration of viable bacteria. The **SEPTI-CHEK** AFB system is designed to promote rapid growth in a broth medium, while providing simultaneous subculture on three different media, a presumptive identification, and the ability to perform identification and susceptibility testing.

**PRINCIPLES OF THE PROCEDURE**

An appropriate patient sample is obtained, digested, concentrated and decontaminated as described in the "Specimen Collection" section of this insert. A measured portion of the treated specimen is then transferred aseptically to the **SEPTI-CHEK** AFB Culture Bottle to which **SEPTI-CHEK** AFB Culture Supplement has been added. The **SEPTI-CHEK** AFB Slide is immediately attached to the top of the bottle in place of the bottle cap in order to maintain the CO<sub>2</sub>-enriched atmosphere. The system is inverted and rotated initially, then periodically during incubation until colonies are observed on the slide or until the specimen is judged to be negative for mycobacteria (8 weeks). Both the broth and slides should be observed for growth. Acid-fast and Gram stains, additional subcultures, or probe assays may be performed from slide colonies and early broth growth can be used for additional testing, if desired.

**REAGENTS**

**SEPTI-CHEK AFB Mycobacteria Culture Bottle:** Modified Middlebrook 7H9 Broth.

Each bottle contains approximately 20 mL of medium in a CO<sub>2</sub>-enriched atmosphere.

Approximate Formula\* Per L of Processed Water

Casein Hydrolysate .....	1.053 g	Magnesium Sulfate Heptahydrate .....	0.053 g
Dipotassium Hydrogen Phosphate .....	1.579	L-glutamic Acid .....	0.526
Sodium Chloride .....	0.895	Ferrous Ammonium (II) Citrate .....	0.042
Ammonium Sulfate .....	0.526	Biotin .....	0.0005
Potassium Dihydrogen Phosphate .....	1.579	Malachite Green .....	0.001
Tri-Sodium Citrate Dihydrate .....	0.421		

\*Adjusted and/or supplemented as required to meet performance criteria.

**SEPTI-CHEK AFB Slide:** A slide with three numbered sides containing media –

Side 1 – Non-selective Middlebrook 7H11 Agar, modified, suitable for growth of most mycobacteria.

Side 2 – Egg based medium, suitable for growth of most mycobacteria.

Side 3 – Chocolate Agar for initial isolation of bacteria other than mycobacteria, which will grow on this agar only after extended incubation.

**BBL™ SEPTI-CHEK™ AFB Slide:** Approximate Formula\* Per L of Processed Water

Side 1: 7H11	Side 2: Egg	Side 3: Chocolate
Middlebrook 7H11 Agar Base .....	Middlebrook 7H11 Agar Base .....	GC Agar Base .....
21.0 g	21.0 g	36.0 g
Glycerol .....	Glycerol .....	Hemoglobin Powder .....
6.2	6.2	15.0
OADC .....	OADC .....	IsoVitalEX™ Enrichment ...
40.0 ml	40.0 ml	10.0
	Malachite Green .....	Granulated Agar .....
	0.05 g	4.0
	Egg Powder .....	
	50.0 (from whole eggs)	

\*Adjusted and/or supplemented as required to meet performance criteria.

**SEPTI-CHEK AFB Mycobacteria Culture Supplement:** Antibiotic and enrichment supplement for culture of mycobacteria. Approximately 2.2 g per bottle containing:

Approximate Formula\* Per L of Processed Water

HCl .....	0.16 g	Sodium Desoxycholate .....	0.01 g
Oleic Acid .....	0.45	Polymixin B Sulfate .....	0.13
Bovine Albumin Fract. V .....	100.00	Amphotericin B .....	0.13
Pyridoxal .....	0.02	Nalidixic Acid .....	0.40
Catalase .....	0.40	Trimethoprim Lactate.....	0.13
Glucose Monohydrate .....	40.00	Adiocillin .....	0.21
Glycerin .....	63.05		

\*Adjusted and/or supplemented as required to meet performance criteria.

**Precautions:** *in vitro* Diagnostic

As with any diagnostic procedure, clinical considerations and professional judgement should be used in interpreting any test result, particularly a single test result.

The **SEPTI-CHEK** AFB System is intended for use in laboratories experienced in handling mycobacteria or mycobacteria-containing specimens.<sup>1</sup> Since *M. tuberculosis* is highly infectious, specimens and **SEPTI-CHEK** AFB components should be handled by qualified personnel in compliance with "Universal Precautions"<sup>1</sup> and institutional guidelines. Bottles and slides should be autoclaved prior to being discarded.

**Instructions for Reagent Handling:** **SEPTI-CHEK** AFB Culture Bottles must be inoculated with 1 mL of reconstituted **SEPTI-CHEK** AFB Culture Supplement prior to use.

To reconstitute **SEPTI-CHEK** AFB Culture Supplement, add 9 mL of sterile distilled water. Gently swirl vial to dissolve contents. Dispense 1 mL to each **SEPTI-CHEK** AFB Culture Bottle to be used. One bottle of reconstituted **SEPTI-CHEK** AFB Culture Supplement is enough for 10 bottles of culture medium.

**Storage and Stability:**

**SEPTI-CHEK AFB Mycobacteria Culture Bottle** – Store at 15 – 25°C. Protect from light. Do not use if turbid. Broth should appear pale green to colorless. Broth to which 1 mL of supplement has been added is stable at 2 – 8°C for 14 days. DO NOT FREEZE. Do not inoculate after expiration date shown on bottle label.

**SEPTI-CHEK AFB Slide** – Store at 2 – 8°C. DO NOT FREEZE. Do not use dehydrated or contaminated slides. Do not attach to bottle after expiration date shown on package.

**SEPTI-CHEK AFB Mycobacteria Culture Supplement** – Store at 2 – 8°C. Protect from light. Supplement reconstituted with sterile distilled water can be stored at 2 – 8°C or at -20°C for up to 14 days prior to use. Do not use to supplement culture bottle after the expiration date shown on the bottle label.

**SPECIMEN COLLECTION**

All specimens should be collected and transported as recommended by the CDC, the *Manual of Clinical Microbiology* or your laboratory procedure manual.<sup>1,2</sup>

Specimen Type	Specific Specimen Handling Recommendations	Pretreatment before Culture Inoculation <sup>1,2</sup>
Sputum Bronchial Washing Bronchial Aspirate	3 early morning samples on successive days, not to exceed 10 mL each.	Decontaminate with N-acetyl-L-cysteine/NaOH (NALC-NaOH) or equivalent method. Resuspend pellet to final volume of 1 – 2 mL. Inoculate.
Body Fluid (Pleural, CSF, Ascites, Synovial)		Specimens which are collected aseptically and are expected to have no other bacteria can be inoculated without decontamination. If specimen is suspected of containing other bacteria, decontaminate with NALC-NaOH or equivalent method and proceed as for sputum.
Stool		Take small sample from specimen. Proceed as for sputum.
Tissues and Biopsy	Specimens may be frozen when a delay in processing is necessary.	Grind in a sterile tissue homogenizer in sterile saline. Bring volume to 2.0 mL. Inoculate. If specimen is suspected of containing other bacteria, decontaminate with NALC-NaOH or equivalent method prior to inoculation.
Urine	3 clean catch or catheterized morning urine specimens on successive days. 24 h pooled specimen is not recommended due to potential contamination.	Centrifuge (≥ 3,000 x g, 30 min). Decant and resuspend in 2 – 5 mL sterile distilled water. Decontaminate with NALC-NaOH or equivalent method. Resuspend pellet to final volume of 1 – 2 mL. Inoculate.
Wounds and Skin		For swabs, suspend in sterile saline, vortex 30 sec, and decontaminate with NALC-NaOH or equivalent method. Resuspend to 2.0 mL. Inoculate.

**PROCEDURE**

**Materials Provided:** Ten 20-mL **SEPTI-CHEK** AFB Mycobacteria Culture Bottles, Ten **SEPTI-CHEK** AFB Slides, One vial **SEPTI-CHEK** AFB Mycobacteria Culture Supplement sufficient for ten slides/bottles (See "Availability").

**Materials Not Provided:** Sterile specimen collection container and centrifuge tubes (non-breakable, sealed), sterile pipettes or needle and syringe, sterile distilled water, N-acetyl-L-cysteine powder, 4% NaOH, tissue homogenizer, sterile swab (tissue, wound, biopsy), centrifuge (2000 – 3500 x g), incubator (30 – 37°C), autoclave, vertical laminar flow hood, and bottle racks (See "Availability").

**Instructions for Use When Using a SEPTI-CHEK AFB Slide:**

- All manipulations should be carried out aseptically under a vertical laminar flow hood using standard precautions for mycobacteria.
- Reconstitute the **SEPTI-CHEK** AFB Mycobacteria Culture Supplement with 9 mL sterile distilled water. Add 1 mL supplement to each **SEPTI-CHEK** AFB Mycobacteria Culture Bottle. (Remove screw cap or inject through stopper.)
- Prepare specimen according to the chart found in "Specimen Collection" section of this insert.
- Add 0.5 – 1.0 mL of treated specimen to culture bottle by removing the screw cap.  
NOTE: The larger the volume of specimen used, the greater the sensitivity of the system. Particulate or turbid specimens may decrease the ability to recognize mycobacterial growth.
- Hold the **SEPTI-CHEK** AFB Slide upright while unscrewing the white bottom cap of the slide container. Do not touch the inner thread. Immediately attach the slide to the culture bottle. Inspect for tightness of fit. Label slide so as not to obscure any bacterial growth.
- Slowly rotate the system 360° around the vertical axis twice. (The **BBL™ SEPTI-CHEK™** Small Bottle Inversion Rack can be used to manipulate multiple systems. See "Availability.")  
NOTE: Delaying the first inversion until 48 h after incubation may decrease contamination.
- Incubate in a vertical position at 37°C (except for skin specimens, which require one system at 37°C and another at 30°C). A CO<sub>2</sub> incubator is **not** required, but may be used.
- Invert and rotate daily during the first week and then weekly for eight weeks or until growth is observed. In the event of condensation on the inside of the slide housing, carefully rotate the bottle so the broth washes the condensation clear but does not contact the slide. Inversion may wash off colonies and is not recommended once any colonies are seen.

**Instructions for Use When NOT Using a SEPTI-CHEK AFB Slide:**

- All manipulations should be carried out aseptically under a vertical laminar flow hood using standard precautions for mycobacteria.
- Reconstitute the **SEPTI-CHEK** AFB Mycobacteria Culture Supplement with 9 mL sterile distilled water. Add 1 mL supplement to each **SEPTI-CHEK** AFB Mycobacteria Culture Bottle. (Remove screw cap or inject through stopper.)
- Prepare specimen according to the chart found in "Specimen Collection" section of this insert.
- Add 0.5 – 1.0 mL of treated specimen to culture bottle by either removing the screw cap or by injecting the specimen directly into the bottle through the rubber stopper using a sterile syringe with a permanently attached needle or **Luer-Lok®** brand tip. Prior to injection of specimen, disinfect rubber stoppers using a 70% isopropyl alcohol solution.  
NOTE: The larger the volume of specimen used, the greater the sensitivity of the system. Particulate or turbid specimens may decrease the ability to recognize mycobacterial growth. If the specimen is injected, briefly loosen screw cap to increase the oxygen content of the bottle. Close cap tightly.

TRANSLINGUA  
ENG – 5 pt on 5.75 pt leading  
MLL – 4.75 pt on 5.25 pt leadingEdited: BDMS Package Design  
Date: 05/27/97  
MAP Format & revisionsReformatted: Translingua  
Date: 07/15/97  
MAP Format & revisionsEdited: BDMS Package Design  
Date: 09/10/97  
Re-sizeEdited: Translingua  
Date: 10/22/97  
AA-L Jap &  
Lasers/Negs/Disk-MLL & JapTRIM: 4.5" W x 18.5" L  
COLOR: Standard BLACK  
Lasers: MLL 85%  
Jap. 100%

5. If specimen is injected into culture bottle, disinfect rubber stopper immediately using a mycobactericidal disinfectant followed by cleaning with a 70% isopropyl alcohol solution.
6. Observe the contents of the bottle daily during the first week to detect rapidly growing mycobacteria or potential contaminants. After one week, examine negative cultures weekly for up to eight weeks. Observed changes may include turbidity, flakes or granules in broth.
7. Make subcultures to appropriate media and stains of broth when growth is suspected.

**Quality Control:** Test procedures for the SEPTI-CHEK AFB Mycobacteria Culture Bottle with supplement added, and the SEPTI-CHEK AFB Slide have been prepared in compliance with National Committee for Clinical Laboratory Standards (NCCLS) Vol. 5, No. 20, *Quality Assurance for Commercially Prepared Microbiological Culture Media*.

Upon receipt of a new shipment or lot number, one of the following organisms is recommended for inoculation of SEPTI-CHEK AFB Mycobacteria Culture Bottle or SEPTI-CHEK AFB Slide at 10<sup>6</sup> CFU/bottle:

a)	<i>M. tuberculosis</i> ..... ATCC <sup>®</sup> 25177	<i>M. intracellulare</i> ..... ATCC 13950
	<i>M. kansasii</i> ..... ATCC 12478	<i>M. fortuitum</i> ..... ATCC 6841
	<i>M. scrofulaceum</i> ..... ATCC 19811	

and one of the following organisms is recommended for inoculation at 10<sup>6</sup> CFU/bottle:

b)	<i>Escherichia coli</i> ..... ATCC 25922	<i>Staphylococcus aureus</i> ..... ATCC 25923
	<i>Streptococcus pneumoniae</i> ..... ATCC 6305	

The growth in the bottle of the (a) group is checked after 21 and 28 days at 37°C with the following anticipated results:

(+) = Weak Growth.  
+ = Moderate-Strong Growth.

Mycobacterium	Growth After	
	21 Days	28 Days
<i>M. tuberculosis</i>	+	+
<i>M. kansasii</i>	+	+
<i>M. scrofulaceum</i>	(+)	(+)
<i>M. intracellulare</i>	(+)	(+)
<i>M. fortuitum</i>	+	+

Growth of the (a) group on slides should be observed after 21 days with the following expected results:

+ = Growth.  
- = No Growth.

The (b) group should show no or weak growth in the bottle after 7 days at 37°C.

Mycobacterium	Side 1: 7H11	Side 2: Egg	Side 3: Choc.
<i>M. tuberculosis</i>	+	+	+/-
<i>M. kansasii</i>	+	+	+/-
<i>M. scrofulaceum</i>	+	+	+/-
<i>M. intracellulare</i>	+	+	-
<i>M. fortuitum</i>	+	+	+/-

If the SEPTI-CHEK AFB Mycobacteria Culture Bottle with supplement added or the SEPTI-CHEK AFB Slide do not give the expected results, do not use the media until Technical Services (U.S.A.) has been contacted, 1-800-638-8663.

**RESULTS**

Growth on Slide			Interpretation of Growth on the Slide: Conclusion
1 7H11	2 Egg	3 Choc.	
+	-	-	High probability of mycobacteria.
-	+	-	High probability of mycobacteria.
+	+	-	High probability of mycobacteria.
+	+	+	High probability of mycobacteria. Moderate probability of mycobacteria and contaminant.
-	-	+	High probability of non-mycobacterial contaminant.

Confirmation or differentiation of non-mycobacterial contaminants should always be done by microscopic examination for acid-fast bacilli.

**Interpretation of Growth in Broth:**

Observation	Conclusion
Broth only	High probability of mycobacteria.
Broth & Side 1 & 2 or 1 or 2, or 1, 2 & 3	High probability of mycobacteria.
Not in Broth or on Slide	Uninoculated system or no viable organisms in sample.
Broth and Side 3 only	Probable mycobacteria and non-mycobacterial contaminant. Low probability of non-mycobacterial contaminant alone. Low probability of mycobacteria alone.

NOTE: Mycobacteria may sometimes generate a non-homogeneous turbidity in the growth culture medium, but often give rise to merely grains or flakes of varying diameter and quantity. Microscopic inspection of positive cultures is recommended for detection of acid-fast bacilli.

**LIMITATIONS OF PROCEDURE**

Mycobacteria are often fastidious, and a single negative culture should not be used to rule out their presence.

Cultures with large amounts of non-mycobacterial contaminants should be repeated and/or subcultured by transfer of broth from an inoculated system to a new system. Mycobacterial growth can often be inhibited or obscured by contaminants but enhanced by subculture. If mixture of contaminant and mycobacteria is present, digestion and decontamination may allow recovery of mycobacteria. The BBL™ SEPTI-CHEK™ AFB system is intended for use in the detection, isolation and subculture of mycobacteria. An acid-fast stain from the culture broth and/or slide is required for confirmation. Species identification from growth on the slide may be performed using either conventional methods or probe technologies.

**PERFORMANCE CHARACTERISTICS**

The SEPTI-CHEK AFB system permits the recovery of mycobacteria from clinical specimens. Three independent published studies have evaluated the performance of the SEPTI-CHEK AFB system as shown in the following tables:<sup>9,10,11</sup>

**STUDY A<sup>9</sup>**

Organism	# of Isolates	% Recovery			SEPTI-CHEK TTD <sup>‡</sup> (days)
		SEPTI-CHEK AFB	BACTEC™ 460TB	Conventional Media*	
<i>M. tuberculosis</i> complex	75	100	N/A	90	23
<i>M. avium</i> complex	40	100	N/A	60	18 – 21
<i>Mycobacterium</i> species <sup>†</sup>	17	100	N/A	71	24 – 26

\* Conventional media reflects combination of Lowenstein-Jensen and 7H11.  
† *M. fortuitum* (4), *M. chelonae* (3), *M. goodii* (5), *M. kansasii* (4) and *M. xenopi* (1).  
‡ TTD (Time-To-Detection) based on SEPTI-CHEK AFB broth-slide combination.

**STUDY B<sup>10</sup>**

Organism	# of Isolates	% Recovery			SEPTI-CHEK TTD <sup>‡</sup> (days)
		SEPTI-CHEK AFB	BACTEC™ 460TB	Conventional Media*	
<i>M. tuberculosis</i> complex	34	100	100	82	19 – 37
<i>M. avium</i> complex	25	92	92	68	18 – 20
<i>Mycobacterium</i> species <sup>†</sup>	2	50	50	50	36

\* Lowenstein-Jensen. † *M. goodii* (2).

**STUDY C<sup>11</sup>**

Organism	# of Isolates	% Recovery			SEPTI-CHEK TTD <sup>‡</sup> (days)
		SEPTI-CHEK AFB	BACTEC™ 460TB	Conventional Media*	
<i>M. tuberculosis</i> complex	9	100	89	67	20
<i>M. avium</i> complex	47	94	74	62	14
<i>Mycobacterium</i> species <sup>†</sup>	10	90	40	40	17

\* Lowenstein-Jensen. † *M. chelonae* (4), *M. fortuitum* (2), *M. scrofulaceum* (1) and *M. goodii* (3).

**AVAILABILITY**

Cat. No.	Description	Description
243558	BBL™ SEPTI-CHEK™ AFB Mycobacteria Culture Bottles, 10 x 20 mL bottles.	BBL™ SEPTI-CHEK™ Small Bottle Inversion Rack.
	BBL™ SEPTI-CHEK™ AFB Mycobacteria Culture Supplement, 5 x 2.2 g bottles.	BBL™ SEPTI-CHEK™ Slide, cpu 10.
	BBL™ SEPTI-CHEK™ AFB Slide, 10 per carton.	BBL™ SEPTI-CHEK™ Slide, cpu 50.

For specific catalog number information, visit our website <http://www.bd.com/microbiology>, or contact the nearest Becton Dickinson Microbiology Systems office.

**REFERENCES**

1. Baron, E.J., L.R. Peterson, and S.M. Finegold. *Bailey and Scott's diagnostic microbiology*, 9th ed. St Louis: Mosby-Year Book, Inc., 1994: 590-633.
2. Wentworth, B.B., coordinating ed. *Diagnostic procedures for bacterial infections*, 7th ed. Washington, DC: American Public Health Association, Inc., 1987: Chapter 30.
3. Joklik, W., H.P. Willett, D.B. Amos, C. Wilfert eds., *Zinsser microbiology*, 19th ed. New York: Appleton and Lange, 1988: Chapter 33.
4. Kiehn, T.E., et al., *J. Clin. Microbiol.*, 21:168-73, 1985.
5. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover, eds. *Manual of clinical microbiology*, 6th ed. Washington, DC: American Society for Microbiology, 400-437, 1995.
6. Mac Donald, N., Paper C-24 Abstracts American Soc. Microb. Annual Meeting 1988:336.
7. Kent, P.T., and G.P. Kubica. *Public health mycobacteriology: A guide for the level III laboratory*. U.S. Department of Health and Human Services, PHS CDC, 1985.
8. *Bloodborne Pathogens. Code of Federal Regulations, Title 29, Part 1910.1030, Fed. Regist.*, 1991, 56:64175-64182.
9. D'Amato, R.F., Isenberg, H.D., Hochstein, L., Mastellone, A.J., and Alperstein, P.J. *Clin. Microbiol.*, 29:2906-2908, 1991.
10. Sewell, D.L., Rashad, A.L., Rourke, W.J., Poor, S.L., McCarthy, J.C. and Pfaller, M.A. *J. Clin. Microbiol.*, 31:2689-2691, 1993.
11. Whittier, P.S., Westfall, K., Setterquist, S. and Hopfer, R.L. *Eur. J. Clin. Microbiol. and Infect. Dis.*, 11:915-917, 1992.

TECHNICAL INFORMATION: In the United States, telephone Technical Services, toll free (800) 638-8663.

© 1999 Becton Dickinson and Company  
BACTEC, BBL, IsoVistaX, Luer-Lok and SEPTI-CHEK are trademarks of Becton Dickinson and Company.  
ATCC is a trademark of the American Type Culture Collection.

**BECTON DICKINSON**

Becton Dickinson Microbiology Systems  
Becton Dickinson and Company  
7 Loveton Circle  
Sparks, Maryland 21152 USA

Becton Dickinson France S.A.  
11 rue Aristide Bergès  
38800 Le Pont de Claix, France

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