
Revisions

SO 0191-5

Rev from	Rev to	ECO #
2009/08	2010/06	5361-10

Notes:


- BD Cat. Number 240862, 240863
- Blank (Sheet) Size: Length: 14" Width: 25"
 Number of Pages: 10 Number of Sheets: 1
 Page Size: Length 14" Width 5" Final Folded Size: 7" x 5"
- Style (see illustrations below): # 4


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No
 No. of Colors: 1 PMS# 2755 Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: <div>8809541</div>		Category and Description <div>Package Insert, MycoPrep</div>	Sheet: 1 of 11	A
			Scale: N/A	

BBL™ MycoPrep™

Specimen Digestion/Decontamination Kit For Processing of Mycobacterial Specimens

English:	pages	1 – 2	Italiano:	pagine	5 – 6
Français :	pages	2 – 3	Español:	páginas	6 – 7
Deutsch:	Seiten	3 – 4	Português:	páginas	7 – 8

 8809541
2010/06

Pokyny vám poskytnú miestni zástupcov spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instruccie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіңіздің жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

The **BBL™ MycoPrep™** Specimen Digestion/Decontamination Kits are used for the digestion and decontamination of clinical specimens suspected to contain mycobacteria, especially *Mycobacterium tuberculosis*.

SUMMARY AND EXPLANATION

The majority of clinical specimens sent to the mycobacteriology laboratory for cultural confirmation of suspected mycobacterial infection (e.g., sputum, bronchial or tracheal lavage) are contaminated by rapidly growing normal flora. To maximize the mycobacterial yield, contaminated specimens require treatment with a digestion and decontamination procedure. N-acetyl-L-cysteine-sodium hydroxide (NALC-NaOH) solution is recommended as a gentle but effective digesting and decontaminating agent.¹

PRINCIPLES OF THE PROCEDURE

Sodium hydroxide (NaOH) can be used as both a digestant and decontaminant. As a mucolytic agent, it is most effective at a final specimen concentration of 2%. However, as a decontaminating agent, this concentration is toxic to both contaminants and to some mycobacteria.¹

N-acetyl-L-cysteine (NALC) is also a mucolytic agent. In the **BBL™ MycoPrep™** reagent bottle, the NALC is combined with 2% NaOH. When the reagent is diluted with an equal volume of specimen, it provides effective digestion and decontamination with a final concentration of 1% NaOH, which is less toxic to mycobacteria.¹

Sodium citrate is included in the reagent to bind heavy metal ions that may be present in the specimen and which can inactivate NALC.¹

Because NALC loses mucolytic activity on standing, the **BBL™ MycoPrep™** reagent contains the NALC component in a sealed glass ampule within the NaOH-citrate solution. The ampule is broken and the reagents gently mixed before use.

Packets of pre-weighed, powdered phosphate buffer, pH 6.8, are included in the **BBL™ MycoPrep™** kit for use in washing the digested-decontaminated specimen. The phosphate buffer decreases the activity of the NALC-NaOH solution and lowers the specific gravity of the specimen before the mycobacteria are recovered by centrifugation.

REAGENTS

BBL™ MycoPrep™ Reagent

Approximate Formula* Per L Purified Water
NaOH..... 20.0 g
Trisodium Citrate (Na₃C₆H₅O₇•2H₂O) 14.5 g
Each sealed glass ampule within the bottle
contains a minimum of 0.370 g NALC (C₅H₉NO₃S).

*Adjusted and/or supplemented as required to meet performance criteria.

BBL™ MycoPrep™ Phosphate Buffer

Approximate Formula* Per 500 mL Purified Water
Disodium Phosphate (Na₂HPO₄) 2.37 g
Monopotassium Phosphate (KH₂PO₄) 2.27 g
Final pH 6.8

Warnings and Precautions: For *in vitro* Diagnostic Use.

The NALC-NaOH reagent contains strong alkali and causes severe burns. Take off all contaminated clothing immediately. Gloves and eye/face protection must be worn. NaOH is irritating to the eyes and skin. In the event of eye or skin contact, rinse immediately with an eye wash system or tap water for at least 15 min and seek medical advice. If ingested, give milk, egg white or large amounts of water and seek medical advice. Keep out of reach of children.

Caution: Break ampule close to its center *one time only*. Do not manipulate ampule further as the plastic bottle may puncture and injury may occur.

Storage Instructions: On receipt, store at 15 to 25°C. Do not freeze. Do not open until ready to use.

Product Deterioration: Do not use reagents if ampules are broken, missing or do not appear to contain powder. Do not use if there is visible evidence of deterioration (e.g., reagent develops a yellow color). Do not use phosphate buffer if packages are torn or unsealed.

SPECIMEN COLLECTION AND HANDLING

Biosafety Level 2 practices and procedures, containment equipment and facilities are required for non-aerosol-producing manipulations of clinical specimens such as preparation of acid-fast smears. All aerosol-generating activities must be conducted in a Class I or II biological safety cabinet. Biosafety Level 3 practices, containment equipment and facilities are required for laboratory activities in the propagation and manipulation of cultures of *M. tuberculosis* and *M. bovis*. Animal studies also require special procedures.²

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"²⁻⁵ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Refer to appropriate texts for details of specimen collection and handling procedures.^{1,6-8}

PROCEDURE

Materials Provided: See "Availability."

Materials Required But Not Provided: Ancillary laboratory equipment required for this procedure.

Test Procedure:

1. Prepare the **BBL™ MycoPrep™** Phosphate Buffer as needed, by pouring contents of one packet into a 500 mL volumetric flask and fill to line with purified water. Transfer the buffer solution to a screw-capped container and, with cap loosened, autoclave at 121°C for 15 min. Cool to room temperature and tighten cap.
2. Using caution not to spill, loosen screw-cap on the **MycoPrep** Reagent bottle. Locate ampule in bottle, squeeze excess air from the bottle and tighten cap. With bottle in the upright position, squeeze the bottle until the ampule breaks. (*Note: the 150 mL bottle contains two ampules that must be broken.*) Shake gently to dissolve the NALC. Avoid excessive agitation. **ONCE AMPULE IS BROKEN, USE REAGENT WITHIN 24 H.**¹
3. In a biological safety cabinet, using a sterile, aerosol-free 50 mL centrifuge tube with screw cap, add equal amounts of specimen and activated NALC-NaOH solution (approximately 10 mL of each).
4. Cap the centrifuge tube and mix on a Vortex-type mixer until specimen is liquified. If specimen is especially viscous, add more NALC-NaOH solution and repeat mixing.
5. Allow mixture to stand at room temperature for 15 min with occasional gentle shaking. Avoid over-treating the specimen.
6. Add the prepared phosphate buffer to the 50 mL mark on the centrifuge tube and mix. Centrifuge for 15 – 20 min at 3000 x g.
7. Carefully decant all of the supernatant fluid.
8. Add a small quantity of phosphate buffer of pH 6.8 (e.g., 0.5 to 2.0 mL) and resuspend the sediment. Use the suspension for the preparation of smears and the performance of mycobacteriological procedures.^{1,6-8}

User Quality Control: For each lot or shipment, examine kit components as described under "Product Deterioration." Process a culture containing acid-fast organisms as per established laboratory quality control procedures for mycobacterial specimens.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

The **MycoPrep** Specimen Digestion/Decontamination Kit is used for the digestion and decontamination of clinical specimens (e.g., sputum, bronchial or tracheal lavage) suspected to contain mycobacteria. If procedures are followed correctly, viscous specimens will be liquefied and contamination by normal flora will be reduced or eliminated.

LIMITATIONS OF THE PROCEDURE

No one method of digestion-decontamination is suitable for all clinical specimens in all situations. When selecting a procedure, choose the mildest procedure that will reduce contamination.

PERFORMANCE CHARACTERISTICS

The **MycoPrep** kit was used in support of clinical trials for the **BACTEC™ MGIT™ 960** System. One clinical trial site used **MycoPrep** for the digestion/decontamination of clinical specimens prior to testing. Of the 63 specimens positive for mycobacteria, none were contaminated and required reprocessing.⁹

AVAILABILITY

Cat. No.	Description
240862	BBL™ MycoPrep™ Specimen Digestion/Decontamination Kit, consisting of ten 75 mL bottles of Reagent (NALC-NaOH Solution) and 5 packages of Phosphate Buffer (pH 6.8).
240863	BBL™ MycoPrep™ Specimen Digestion/Decontamination Kit, consisting of ten 150 mL bottles of Reagent (NALC-NaOH Solution) and 10 packages of Phosphate Buffer (pH 6.8).

REFERENCES

1. Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. USDHHS. Centers for Disease Control, Atlanta.
2. U.S. Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. 2007. Biosafety in microbiological and biomedical laboratories, 5th ed. HHS Publication No. (CDC) 93-8395. U.S. Government Printing Office, Washington, D.C.
3. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI, Wayne, Pa.
4. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
5. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
6. Cernoch, P.L., R.K. Enns, M.A. Saubolle, and R.J. Wallace, Jr. 1994. Cumitech 16A, Laboratory diagnosis of the mycobacterioses. Coord. ed., A.S. Weissfeld. American Society for Microbiology, Washington, D.C.
7. Pfyffer, G.E., B.A. Brown-Elliott, and R.J. Wallace, Jr. 2003. *Mycobacterium*: general characteristics, isolation, and staining procedures, p. 532-559. In P.R. Murray, E.J. Baron, J.H. Jorgensen, M.A. Pfaller and R.H. Tenover. (eds.), Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
8. Isenberg, H.D. (ed). 2004. Clinical microbiology procedures handbook, vol. 1, 2 and 3, 2nd ed. American Society for Microbiology, Washington, D.C.
9. Data on file, BD Diagnostics.

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© 2010 BD.
BD, BD Logo, BBL, BACTEC, MGIT and MycoPrep are trademarks of Becton, Dickinson and Company.
Benex Limited