BD BBL™ MGIT™ AST SIRE System

For the Antimycobacterial Susceptibility Testing of *Mycobacterium tuberculosis*

English: pages 1 – 5 Italiano: pagine 14 – 18
Français: pages 5 – 9 Español: páginas 19 – 23
Deutsch: Seiten 10 – 14

8809591JAA(03)
2016-09

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нусқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukciju teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

The **BBL™ MGIT™ AST SIRE System** is a rapid manual qualitative procedure for susceptibility testing of *Mycobacterium tuberculosis*, from culture, to streptomycin, isoniazid, rifampin and ethambutol.

SUMMARY AND EXPLANATION

Antimycobacterial susceptibility testing is necessary for the proper treatment of patients with tuberculosis. The treatment of tuberculosis is commonly through a multiple drug regimen using the primary antimycobacterial drugs, streptomycin, isoniazid, rifampin and/or ethambutol. It is important that the antimycobacterial drugs utilized show appropriate activity against *Mycobacterium tuberculosis*, i.e., susceptibility.

Multidrug resistant *Mycobacterium tuberculosis* (MDR-TB) has recently become a serious public health problem.¹ Resistance to any of the four primary drugs, streptomycin (STR), isoniazid (INH), rifampin (RIF), and ethambutol (EMB), makes the disease more difficult and expensive to treat. The rapid detection of these strains is critical to the effective treatment of the patient.

Two methods have been widely used for antimycobacterial susceptibility testing. The first method, known as the Method of Proportion,² uses Middlebrook and Cohn 7H10 or 7H11 Agar. It compares colony counts on drug-containing and drug-free media. Resistance to a drug is detected when 1% or more of the population is resistant to the drug concentration under test. Results are generally available after 21 days of incubation. The second method is based on growth in liquid culture and generally takes from 3 to 14 days.

The **BBL MGIT AST System** provides the susceptibility result within 14 days and allows appropriate antibiotic therapy to be implemented sooner than with the Method of Proportion.

PRINCIPLES OF THE PROCEDURE

The **BBL MGIT Mycobacteria Growth Indicator Tube** is a tube containing a modified Middlebrook 7H9 Broth which, when supplemented with **BBL MGIT OADC enrichment**, supports the growth and detection of mycobacteria (see **BBL MGIT Products** package insert). The **MGIT** tube contains a fluorescent compound embedded in silicone on the bottom of a 16 x 100 mm round-bottom tube. The fluorescent compound is sensitive to the presence of oxygen dissolved in the broth. The initial concentration of dissolved oxygen quenches the emissions from the compound, and little fluorescence can be detected. Later, actively respiring microorganisms consume the oxygen and allow more fluorescence to be observed using a 365 nm UV transilluminator or longwave UV light.

The **MGIT AST System** is a three to fourteen day qualitative test. The test is based on growth of the *Mycobacterium tuberculosis* strain in a drug-containing tube compared to a drug-free tube. **MGIT** tubes are observed daily from the third day after inoculation. The absence of fluorescence in the drug-containing tube two days beyond the appearance of fluorescence in the Growth Control tube is indicative of susceptibility of the organism to that drug. Fluorescence of a drug-containing tube on or within two days of fluorescence of the Growth Control tube is indicative of resistance of the organism to that drug.

REAGENTS

MGIT AST SIRE Kit contains two each lyophilized vials of streptomycin, isoniazid, rifampin and ethambutol.

Approximate Formula* Per Vial Lyophilized Streptomycin: Streptomycin 160 µg.
Approximate Formula* Per Vial Lyophilized Isoniazid: Isoniazid..... 20 µg.
Approximate Formula* Per Vial Lyophilized Rifampin: Rifampin..... 200 µg.
Approximate Formula* Per Vial Lyophilized Ethambutol: Ethambutol..... 700 µg.

*Adjusted and/or supplemented as required to meet performance criteria.

PRODUCT DETERIORATION

Some variation in the appearance of the lyophilized SIRE drugs may occur. This results from the lyophilization process and does not affect performance of the products.

Directions For Use: Reconstitute each **BBL MGIT** Streptomycin lyophilized vial with 4 mL of sterile distilled/deionized water to make a stock solution of 40 µg/mL.

Reconstitute each **BBL MGIT** Isoniazid lyophilized vial with 4 mL of sterile distilled/deionized water to make a stock solution of 5 µg/mL.

Reconstitute each **BBL MGIT** Rifampin lyophilized vial with 4 mL of sterile distilled/deionized water to make a stock solution of 50 µg/mL.

Reconstitute each **BBL MGIT** Ethambutol lyophilized vial with 4 mL of sterile distilled/deionized water to make a stock solution of 175 µg/mL.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Laboratory procedures involving mycobacteria require special equipment and technique to minimize biohazard.³ Biosafety Level 2 practices and procedures, containment equipment and facilities are required for non-aerosol-producing manipulations of clinical specimens such as preparation of acid-fast smears. All aerosol-generating activities must be conducted in a Class I or II biological safety cabinet. Biosafety Level 3 practices, containment equipment and facilities are required for laboratory activities in the propagation and manipulation of cultures of *M. tuberculosis* and *M. bovis*. Animal studies also require special procedures.⁴

BBL MGIT AST SIRE

BBL MGIT AST SIRE-Rifampin, Lyophilized

Warning

H302+H332 Harmful if swallowed or if inhaled

P261 Avoid breathing dust/fume/gas/mist/vapours/spray. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P271** Use only outdoors or in a well-ventilated area. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Read and follow directions contained in all appropriate package inserts including the **BBL MGIT** Mycobacteria Growth Indicator Tube (see "Availability").

Wear UV protective glasses when observing fluorescence and use only longwave illumination (365 nm). DO NOT USE SHORTWAVE UV LIGHT FOR VIEWING.

Prior to use, the user should examine the tubes for evidence of contamination or damage. Discard any tubes if they appear unsuitable or exhibit fluorescence prior to use. Dropped tubes should be examined carefully. If damage is seen, the tube should be discarded.

A nephelometer must be used for the preparation of isolate suspensions from solid media (e.g., Lowenstein-Jensen Medium).

Autoclave all inoculated **MGIT** tubes prior to disposal.

Storage Instructions: On receipt, store the lyophilized vials at 2–8 °C. Once reconstituted, the antibiotic solutions may be frozen and stored at -20 °C up to 6 months, not to exceed the original expiration date. Once thawed, use immediately. Discard unused portions.

SPECIMEN PREPARATION

All preparations detailed below are from cultures of *Mycobacterium tuberculosis*. The laboratory should confirm, by appropriate identification techniques, that the isolate to be tested is a pure culture.

Preparation of the Isolate from Solid Media:

1. Add 4 mL of **BBL** Middlebrook 7H9 Broth to a 16.5 x 128 mm sterile tube with cap containing 8 – 10 glass beads.
2. Scrape with a sterile loop as many colonies as possible from growth up to 14 days old, trying not to remove any solid medium. Suspend the colonies in the Middlebrook 7H9 Broth. The suspension should exceed in turbidity a 1.0 McFarland standard.
3. Vortex the suspension for 2–3 min to break up the larger clumps.
4. Let the suspension sit for 20 min without disturbing.
5. Transfer the supernatant fluid to another 16.5 x 128 mm sterile tube with cap (avoid transferring any of the sediment) and let sit for another 15 min.
6. Transfer the supernatant fluid (it should be smooth, free of any clumps) to a third 16.5 x 128 mm sterile tube.
7. Adjust the suspension to a 0.5 McFarland standard using a nephelometer.
8. Dilute 1.0 mL of the 0.5 McFarland-adjusted suspension in 4 mL of sterile saline (1:5 dilution). The inoculum is now ready. Proceed to "Inoculation Procedure for Susceptibility Test."

Preparation from a Positive MGIT Tube:

1. For the preparation of the test inoculum, a positive **MGIT** tube should be used the day after it first becomes positive, up to and including three days following initial positivity. A tube which has been positive longer than four days should be subcultured to a fresh **MGIT** tube and used from one to three days following positivity in that **MGIT** tube. Vortex the **MGIT** tube for 10 s.
2. Pipet 1.0 mL of the suspension from the **MGIT** tube into 4 mL of sterile saline (1:5 dilution). The inoculum is now ready. Proceed to "Inoculation Procedure for Susceptibility Test."

PROCEDURE

Materials Provided: **BBL MGIT AST SIRE** Kit containing two each lyophilized vials of streptomycin, isoniazid, rifampin and ethambutol.

Materials Required But Not Provided: **BBL MGIT** Mycobacteria Growth Indicator Tubes, **BBL MGIT** OADC, ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Negative Control: An unopened, uninoculated **MGIT** tube is used for the Negative Control.

Preparation of the Positive Control:

1. Empty broth from an uninoculated **MGIT** tube.
2. Label the tube as a Positive Control and record the date.
3. Prepare a 0.4% sodium sulfite solution (0.4 g in 100 mL sterile distilled or deionized water).
4. Add 5 mL of sodium sulfite solution to the empty **MGIT** tube. Replace the cap, tighten and allow to stand at room temperature for 1 h. Do not incubate.
5. Positive Control tubes can be used many times. Each Positive Control tube can be used for up to 4 weeks when stored at room temperature.

Inoculation Procedure for MGIT Susceptibility Test:

1. Label five **MGIT** tubes for each test isolate. Label one as **MGIT** GC (Growth Control), one as **MGIT** STR, one as **MGIT** INH, one as **MGIT** RIF, and the last as **MGIT** EMB.
2. Aseptically add 0.5 mL of **MGIT** OADC to each tube.
3. Aseptically pipet, using a micropipet, 100 μ L of 40 μ g/mL **MGIT** STR solution to the appropriately labelled **MGIT** tube. Aseptically pipet 100 μ L of 5 μ g/mL **MGIT** INH solution to the appropriately labelled **MGIT** tube. Aseptically pipet 100 μ L of 50 μ g/mL **MGIT** RIF solution to the appropriately labelled **MGIT** tube. Aseptically pipet 100 μ L of 175 μ g/mL **MGIT** EMB solution to the appropriately labelled **MGIT** tube. No antibiotics should be added to the **MGIT** GC tube.

Drug	Concentration of Drug after Reconstitution	Volume Added to MGIT Tubes for Test	Final Concentration In MGIT Tubes
MGIT STR	40 μ g/mL	100 μ L	0.8* μ g/mL
MGIT INH	5 μ g/mL	100 μ L	0.1* μ g/mL
MGIT RIF	50 μ g/mL	100 μ L	1.0* μ g/mL
MGIT EMB	175 μ g/mL	100 μ L	3.5* μ g/mL

*Equivalent to CDC³ recommended critical drug concentrations.

4. Inoculate using a pipet, 0.5 mL of the 1:5 organism suspension (see "Specimen Preparation") into each of the five **MGIT** tubes. Wipe the tubes with a tuberculocidal disinfectant. Tightly recap the tubes and mix well.
5. Incubate the labelled **MGIT** tubes at 37 °C.
6. Streak 0.1 mL of the 1:5 organism suspension to a **Trypticase™** Soy Agar with 5% Sheep Blood (TSA II) plate. Enclose in a plastic bag. Incubate at 35–37 °C.
7. Check the blood agar plate at 48 h for bacterial contamination.
8. If the blood agar plate shows no growth, then proceed to "Reading **MGIT** Tubes."
9. If the blood agar plate shows growth, discard **MGIT** tubes and repeat testing with pure culture.

Reading MGIT Tubes:

1. Remove the **MGIT** tubes from the incubator on the third day after inoculation and read with a 365 nm UV transilluminator or longwave UV light.
NOTE: It is important to read the AST tube every day beginning on Day 3, until results can be interpreted.
2. Compare the **MGIT** GC tube to the Positive and Negative Control tubes. The Positive Control should strongly fluoresce (very bright orange color on the bottom of the tube and also an orange reflection on the meniscus). The Negative Control tube should have very little fluorescence.
3. If fluorescence of the **MGIT** GC tube looks more like the Positive Control than the Negative Control, it is positive for growth. Once the **MGIT** GC tube is positive, it is used to interpret the drug-containing tubes. The drug-containing tubes are interpreted on the same day the **MGIT** GC is positive and for up to two additional days, according to the "Interpretation of Test Results" section, not to exceed fourteen days.
4. If the GC tube has no fluorescence and looks more like the Negative Control, reincubate the tubes and continue to read daily until twelve days after inoculation of all tubes. If the GC result is equivocal (difficult to determine whether an orange fluorescence is present), then the tube should be considered negative and reincubated.
5. If the GC tube is not positive by the twelfth day of the test, the test is invalid.

Interpretation of Test Results: Interpret the **MGIT** result as Susceptible if the drug-containing tube does NOT fluoresce within two days of onset of fluorescence in the GC tube. Interpret the **MGIT** result as Resistant if the drug-containing tube fluoresces on or within two days of the day of onset of fluorescence in the GC tube. When interpreting resistance, finalize the result as soon as the **MGIT** GC and the drug-containing tubes fluoresce.

User Quality Control: Upon receipt of a new shipment or lot number of **BBL MGIT AST SIRE** Kit vials, it is suggested that the control organism shown below be inoculated into tubes with drugs (see "Inoculation Procedure for Susceptibility Test"). Upon observing the proper results, as shown below, the **MGIT AST SIRE** drugs are ready for use in testing patient isolates. If the proper results are not observed, repeat the test. If, after repeating the test, the proper results are still not observed, do not use the media until you have contacted Technical Services at (800) 638-8663 (United States only).

Strain	GC	MGIT STR	MGIT INH	MGIT RIF	MGIT EMB
<i>M. tuberculosis</i> ATCC® 27294	Fluorescence within 3–7 days	No fluorescence within 2 days of GC	No fluorescence within 2 days of GC	No fluorescence within 2 days of GC	No fluorescence within 2 days of GC

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

LIMITATIONS OF THE PROCEDURE

Suspensions made from solid media must be allowed to settle for the prescribed times prior to standardization. Inoculum preparations made from solid media without the use of a nephelometer may give inaccurate results due to incorrect biomass.

The test is uninterpretable if the Growth Control does not fluoresce within twelve days of inoculation.

Use only pure cultures of *M. tuberculosis*. Cultures which are contaminated or which contain multiple strains of mycobacteria may give erroneous results.

PERFORMANCE CHARACTERISTICS

The performance of the **MGIT AST SIRE** System was established in two clinical evaluations conducted at regional reference centers for mycobacterial susceptibility testing and university hospital-based laboratories in areas of high prevalence of *M. tuberculosis* resistance to isoniazid and/or rifampin. The **MGIT AST** System was compared to the Method of Proportion. The clinical evaluation at four sites comparing INH and RIF results between the **MGIT AST** System and Method of Proportion using 7H10 Media (INH at 0.2 µg/mL, RIF at 1.0 µg/mL) included 259 clinical isolates. The clinical evaluation comparing STR and EMB results between the **MGIT AST** System and Method of Proportion included 138 clinical isolates: 103 isolates from two sites using 7H10 Media (STR at 2.0 µg/mL, EMB at 5.0 µg/mL), and 35 isolates from one site using LJ Media (STR at 4.0 µg/mL, EMB at 1.0 µg/mL).

Data was analyzed and interpreted qualitatively for category agreement (S/S or R/R), and overall (site combined) percent agreements were as follows: STR = 94.9%, INH = 93.1%, RIF = 98.5%, and EMB = 93.5%.⁵

Tables 1 and 2 show comparative performance between the **MGIT AST SIRE** System and the Method of Proportion.

Table 1

Number of Isolates with Indicated Susceptibility Results					
Drug	MGIT and MOP R	MOP S MGIT R	MOP R MGIT S	MGIT and MOP S	Total Isolates Tested
STR	24	5	2	107	138
INH	70	13	5	171	259
RIF	61	1	3	194	259
EMB	12	5	4	117	138

S = Susceptible R = Resistant

Reproducibility of **MGIT AST SIRE** System test results were compared to expected results for a panel of 5 ATCC strains and 16 qualified strains which included several strains resistant to each of the drugs. The reproducibility results were: 97% for STR, 94% for INH, 98% for RIF, and 94% for EMB. Individual site reproducibility results ranged from 92% to 100% for combined drug results.

Table 2

Performance Characteristics (%)					
Drug	Sensitivity	Specificity	Predictive Value of Susceptibility	Predictive Value of Resistance	Category Agreement
STR	92.3	95.5	98.2	82.8	94.9
INH	93.3	92.9	97.2	84.3	93.1
RIF	95.3	99.5	98.5	98.4	98.5
EMB	75.0	95.9	96.7	70.0	93.5

AVAILABILITY

Cat. No.	Description	Cat. No.	Description
245119	BD BBL™ MGIT™ AST SIRE Kit, carton of 8 lyophilized vials.	221818	BD BBL™ MGIT™ Normal Saline, 5 mL, carton of 10.
245111	BD BBL™ MGIT™ Mycobacteria Growth Indicator Tubes, 4 mL, carton of 25 tubes.	221819	BD BBL™ MGIT™ Normal Saline, 5 mL, carton of 100.
245113	BD BBL™ MGIT™ Mycobacteria Growth Indicator Tubes, 4 mL, carton of 100 tubes.	295939	BD BBL™ MGIT™ Middlebrook 7H9 Broth, 8 mL, carton of 10.
245116	BD BBL™ MGIT™ OADC, 15 mL, carton of 6 vials.	297345	BD BBL™ MGIT™ Water, 5 mL, carton of 100.

REFERENCES

1. Barenfanger, J. 1993. Making your lab safe against multi-drug resistant *Mycobacterium tuberculosis*. Clin. Microbiol. Newsl. 15: 76-80.
2. Clinical and Laboratory Standards Institute (CLSI). Susceptibility Testing of Mycobacteria, Nocardiae, and Other Aerobic Actinomycetes; Approved Standard-Second Edition. CLSI document M24-A2. Clinical and Laboratory Standards Institute, Wayne, PA, USA
3. Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. USDHHS. Centers for Disease Control, Atlanta.
4. U.S. Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. 1999. Biosafety in microbiological and biomedical laboratories, 4th ed. HHS Publication No. (CDC) 93-8395. U.S. Government Printing Office, Washington, D.C.
5. Data on file, BD Diagnostic Systems.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113 Australia

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.