BD BBL™ Crystal™ Identification Systems **Rapid Gram-Positive ID Kit**

8809711JAA(02)
2015-01

For Export Use Only

English: pages 1 – 5 Italiano: pagine 15 – 20
Français: pages 6 – 10 Español: páginas 20 – 24
Deutsch: Seiten 11 – 15

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie ziskate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

The **BD BBL™ Crystal™** Rapid Gram-Positive (RGP) Identification (ID) System is a miniaturized identification method employing modified conventional, fluorogenic and chromogenic substrates. It is intended for the identification of frequently isolated aerobic gram-positive bacteria.^{1,2,13,16}

SUMMARY AND EXPLANATION

Micromethods for the biochemical identification of microorganisms were reported as early as 1918.³ Several publications reported on the use of the reagent-impregnated paper discs and micro-tube methods for differentiating enteric bacteria.^{3,4,7,17,19} The interest in miniaturized identification systems led to the introduction of several commercial systems in the late 1960s, and they provided advantages in requiring little storage space, extended shelf life, standardized quality control and ease of use.

In general, many of the tests used in the **BD BBL Crystal** ID Systems are modifications of classical methods. These include tests for fermentation, oxidation, degradation and hydrolysis of various substrates. In addition, there are chromogen and fluorogen linked substrates, as in the **BD BBL Crystal** RGP ID panel, to detect enzymes that microbes use to metabolize various substrates.^{5,7,8,9,11,12,14,15}

The **BD BBL Crystal** RGP ID kit is comprised of (i) **BD BBL Crystal** RGP ID panel lids, (ii) **BD BBL Crystal** bases and (iii) **BD BBL Crystal** ANR, GP, RGP, N/H ID Inoculum Fluid (IF) tubes. The lid contains 29 dehydrated substrates and a fluorescence control on tips of plastic prongs. The base has 30 reaction wells. Test inoculum is prepared with the inoculum fluid and is used to fill all 30 wells in the base. When the lid is aligned with the base and snapped in place, the test inoculum rehydrates the dried substrates and initiates test reactions.

Following an incubation period, the wells are examined for color changes or presence of fluorescence that result from metabolic activities of the microorganisms. The resulting pattern of the 29 reactions is converted into a ten-digit profile number that is used as the basis for identification.¹⁸ Biochemical and enzymatic reaction patterns for the 29 **BD BBL Crystal** RGP ID substrates for a wide variety of microorganisms are stored in the **BD BBL Crystal** RGP ID data base. Identification is derived from a comparative analysis of the reaction pattern of the test isolate to those held in the database. A complete list of taxa that comprises the current database is provided in Table 1.

PRINCIPLES OF THE PROCEDURE

The **BD BBL Crystal** RGP ID panels contain 29 dried biochemical and enzymatic substrates. A bacterial suspension in the inoculum fluid is used for rehydration of the substrates. The tests used in the system are based on microbial utilization and degradation of specific substrates detected by various indicator systems. Enzymatic hydrolysis of fluorogenic substrates containing coumarin derivatives of 4-methylumbelliferone (4MU) or 7-amino-4-methylcoumarin (7-AMC), results in increased fluorescence that is easily detected visually with a UV light source.^{11,12,14,15} Chromogenic substrates upon hydrolysis produce color changes that can be detected visually. In addition, there are tests that detect the ability of an organism to hydrolyze, degrade, reduce or otherwise utilize a substrate in the **BD BBL Crystal** ID Systems.

Reactions employed by various substrates and a brief explanation of the principles employed in the system are described in Table 2. Panel location in referred tables indicates the row and column where the well is located (example: 1J refers to Row 1 in column J).

REAGENTS

The **BD BBL Crystal** RGP ID panel contains 29 enzymatic and biochemical substrates. Refer to Table 3 for a list of active ingredients.

Precautions: *in vitro* Diagnostic

After use, all infectious materials including plates, cotton swabs, inoculum fluid tubes, and panels must be autoclaved prior to disposal or incineration.

STORAGE AND HANDLING/SHELF LIFE

Lids: Lids are individually packaged and must be stored unopened in a refrigerator at 2 – 8 °C. **DO NOT FREEZE.** Visually inspect the package for holes or cracks in the foil package. Do not use if the packaging appears to be damaged. Lids in the original packaging, if stored as recommended, will retain expected reactivity until the date of expiration.

Bases: Bases are packaged in two sets of ten, in **BD BBL Crystal** incubation trays. The bases are stacked facing down to minimize air contamination. Store in a dust-free environment at 2 – 30 °C, until ready to use. Store unused bases in the tray, in plastic bag. Empty trays should be used to incubate inoculated panels.

Inoculum Fluid: **BD BBL Crystal** ANR, GP, RGP, N/H ID Inoculum Fluid (IF) is packaged in two sets of ten tubes. Visually inspect the tubes for cracks, leaks, etc. Do not use if there appears to be a leak, tube or cap damage or visual evidence of contamination (i.e., haziness, turbidity). Store tubes at 2 – 25 °C. Expiration dating is shown on the tube label. Only **BD BBL Crystal** ANR, GP, RGP, N/H Inoculum Fluid should be used with **BD BBL Crystal** RGP ID panel. On receipt, store the **BD BBL Crystal** RGP ID kit at 2 – 8 °C. Once opened, only the lids need to be stored at 2 – 8 °C. The remaining components of the kit may be stored at 2 – 25 °C. If the kit or any of the components are stored refrigerated, each should be brought to room temperature prior to use.

SPECIMEN COLLECTION AND PROCESSING

BD BBL Crystal ID Systems are not for use directly with clinical specimens. Use isolates from media such as **Trypticase™** Soy Agar with 5% Sheep Blood (TSA) or Columbia Agar with 5% Sheep Blood (Columbia). Use of selective media such as Phenylethyl Alcohol Agar with 5% Sheep Blood (PEA) or Columbia CNA Agar with 5% Sheep Blood (CNA) is also acceptable. Media containing esculin should not be used. The test isolate must be a pure culture, no more than 18 – 24 h old for most genera; for some slow growing organisms up to 48 h may be acceptable. When swabs are utilized, only cotton-tipped applicators should be used to prepare the inoculum suspensions. Some polyester swabs may cause problems with inoculation of the panels. (See "Limitations of the Procedure".) Once lids are removed from the sealed pouches, they must be used within 1 h to ensure adequate performance. The plastic cover should remain on the lid until used.

The incubator used should be humidified to prevent evaporation of fluid from the wells during incubation. The recommended humidity level is 40 – 60%. The usefulness of **BD BBL Crystal** ID Systems or any other diagnostic procedure performed on clinical specimens is directly influenced by the quality of the specimens themselves. It is strongly recommended that laboratories employ methods discussed in the *Manual of Clinical Microbiology* for specimen collection, transport and inoculation onto primary isolation media.^{1,16}

TEST PROCEDURE

Materials Provided: BD BBL Crystal RGP ID Kit –


- 20 **BD BBL Crystal** RGP ID Panel Lids,
- 20 **BD BBL Crystal** Bases,
- 20 **BD BBL Crystal** ANR, GP, RGP, N/H ID IF Tubes. Each tube has approximately 2.3 ± 0.15 mL of Inoculum Fluid containing: KCl 7.5 g, CaCl_2 0.5 g, Tricine N-[2-Hydroxy-1, 1-bis (hydroxymethyl)methyl] glycine 0.895 g, purified water to 1000 mL.
- 2 incubation trays,
- 1 **BD BBL Crystal** RGP ID Report Pad.

Materials Required But Not Provided: Sterile cotton swabs (*do not use polyester swabs*), incubator (35 – 37 °C) non- CO_2 (40 – 60% humidity), McFarland No. 2 standard, **BD BBL Crystal** Panel Viewer, **BD BBL Crystal** ID System Electronic Codebook or **BD BBL Crystal** RGP Manual Codebook, and appropriate culture media.


Also required are the necessary equipment and labware used for preparation, storage and handling of clinical specimens.

Test Procedure: **BD BBL Crystal** RGP ID System requires a Gram stain.


1. Remove lids from pouch. Discard desiccant. Once removed from the pouch, covered lids should be used within 1 h. Do not use the panel if there is no desiccant in the pouch.
2. Take an inoculum fluid tube and label with patient's specimen number. Using aseptic technique, with the tip of a sterile cotton swab (*do not use a polyester swab*) or a wooden applicator stick or disposable plastic loop, pick colonies of the same morphology from one of the recommended media (see section "Specimen Collection and Processing").
3. Suspend colonies in a tube of **BD BBL Crystal** ANR, GP, RGP, N/H ID Inoculum Fluid.
4. Recap tube and vortex for approximately 10 – 15 sec. The turbidity should be equivalent to a McFarland No. 2 standard. If the inoculum suspension concentration is in excess of the recommended McFarland standard, one of the following steps is recommended:
 - a. Use a fresh tube of inoculum fluid to prepare a new inoculum suspension equivalent to a McFarland No. 2 standard.
 - b. If additional colonies are unavailable for preparation of a new inoculum suspension, using aseptic techniques, dilute the inoculum by adding the minimum required volume (not to exceed 1.0 mL) of 0.85% sterile saline or inoculum fluid to bring down the turbidity equivalent to a McFarland No. 2 standard. Remove the excess amount added to the tube with a sterile pipet so that the final volume of inoculum fluid is approximately equivalent to that of the original volume in the tube (2.3 ± 0.15 mL). Failure to make this adjustment in volume will result in spilling of the inoculum suspension over the black portion of the base rendering the panel unusable.


5. Take a base, and mark the patient's specimen number on the side wall.
6. Pour entire contents of the inoculum fluid tube into target area of the base.

7. Hold base in both hands and roll inoculum gently along the tracks until all of the wells are filled. Roll *back* any excess fluid to the target area and place the base on a bench top. Due to the high cell concentrations used in **BD BBL Crystal** RGP ID panels, the inoculum should be carefully rolled across the tracks to ensure a proper fill of all wells. Make sure there is no excess fluid between the wells before the lid is aligned.


8. Align the lid so that the labeled end of the lid is on top of the target area of the base.


9. Push down until a slight resistance is felt. Place thumb on edge of lid towards middle of panel on each side and push downwards simultaneously until the lid snaps into place (listen for two "clicks").

Purity Plate: Using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the base and inoculate an agar slant or plate (any appropriate medium) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the slant or plate for 24 – 48 h at 35 – 37 °C under appropriate conditions. The purity plate or slant may also be used for any supplementary tests or serology, if required.


Incubation: Place inoculated panels in incubation trays. Ten panels can fit in one tray (5 rows of 2 panels). All panels should be incubated **face down** (larger windows facing up; label facing down) in a non-CO₂ incubator with 40 – 60% **humidity**. Trays should not be stacked more than two high during incubation. The incubation time for panels is **4 h** at 35 – 37 °C. **NOTE:** The incubator door should not be opened repeatedly during the incubation period (preferably less than 3 times). Panels should be read within 30 min after removing from incubator.

Reading: After the recommended period of incubation, remove the panels from the incubator. All panels should be read **face down** (larger windows up; label facing down) using the **BD BBL Crystal** Panel Viewer. Refer to the color reaction chart and/or Table 3 for an interpretation of the reactions. Use the **BD BBL Crystal** RGP Report Pad to record reactions.

- a. Read columns E thru J first, using the regular (white) light source.
- b. Read columns A thru D (fluorescent substrates) using the UV light source in the panel viewer. A fluorescent substrate well is considered positive *only if* the intensity of the fluorescence observed in the well is *greater* than the Negative Control well (4A).

Calculation of BD BBL Crystal Profile Number: Each test result (except 4A, which is used as a fluorescence negative control) scored positive is given a value of 4, 2, or 1, corresponding to the row where the test is located. A value of 0 (zero) is given to any negative result. The values resulting from each positive reaction in each column are then added together. A 10-digit number is generated; this is the profile number.

Example:	A	B	C	D	E	F	G	H	I	J
4	+	+	–	–	+	+	+	–	+	–
2	–	+	+	+	–	+	–	+	+	–
1	+	–	+	–	+	–	–	+	+	–
Profile	1	6	3	2	5	6	4	3	7	0

*(4A) = fluorescent negative control

The resulting profile number and cell morphology, if known should be entered on a PC in which the **BD BBL Crystal** ID System Electronic Codebook has been installed to obtain the identification. A manual codebook is also available. If a PC is not available contact, BD Diagnostics Technical Services for assistance with the identification.

User Quality Control: Quality control testing is recommended for each lot of panels as follows –

1. Inoculate a panel with *Streptococcus pyogenes* ATCC™ 19615 per recommended procedure (refer to "Test Procedure").
2. Prior to incubation, let panel remain at room temperature for 1 min (not more than 2 min).
3. Read and record reactions with the aid of the panel viewer and color reaction chart.
4. If any of the wells (except 1J) are positive per color reaction chart (after 1 – 2 min), DO NOT USE PANELS from this lot. Contact BD Diagnostics Technical Services.
5. If all wells are negative, then incubate panel for 4 h at 35 – 37 °C.
6. Read panel with the panel viewer and color reaction chart; record reactions using the report pad.
7. Compare recorded reactions with those listed in Table 4. If discrepant results are obtained, confirm purity of quality control strain before contacting BD Diagnostics Technical Services.
8. The incubator door should not be opened repeatedly during the incubation period (preferably less than 3 times).

Expected test results for additional quality control test strains are listed in Table 5.

LIMITATIONS OF THE PROCEDURE

The **BD BBL Crystal** RGP ID System is designed for the taxa provided. Taxa other than those listed in Table 1 are not intended for use in this system.

During clinical evaluations, interpretation errors were seen at one site with *Enterococcus faecium* and *Streptococcus mitis*. It is therefore recommended that identification of these two species be confirmed whenever the user finds it appropriate.

The **BD BBL Crystal** RGP ID database was developed with **BBL™** brand media. Reactivity of some substrates in miniaturized identification systems may be dependent upon the source media used in inoculum preparations. We recommend the use of the following media for use with the **BD BBL Crystal** RGP ID System: **TSA II** or Columbia Blood Agar. Use of selective media, such as PEA or CNA is also acceptable. Media containing esculin should not be used.

BD BBL Crystal Identification Systems use a modified microenvironment; therefore, expected values for its individual tests may differ from information previously established with conventional test reactions. The accuracy of the **BD BBL Crystal** RGP ID System is based on statistical use of specially designed tests and an exclusive database.

While **BD BBL Crystal** RGP ID System aids in microbial differentiation, it should be recognized that minor variations may exist in strains within species. Use of panels and interpretation of results require a competent microbiologist. The final identification of the isolate should take into consideration the source of the specimen, aerotolerance, cell morphology, colonial characteristics on various media as well as metabolic end products as determined by gas-liquid chromatography, when warranted.

Only cotton-tipped applicator swabs should be used to prepare the inoculum suspension as some polyester swabs may cause the inoculum fluid to become viscous. This may result in insufficient inoculum fluid to fill the wells. Once lids are removed from the sealed pouches they must be used within 1 h to ensure adequate performance. The plastic cover should remain on the lid until used.

The incubator where panels are placed should be humidified to prevent evaporation of inoculum fluid from the wells during incubation. The recommended humidity level is 40 – 60%.

The panels, after inoculation, should only be incubated **face down** (larger windows facing up; label facing down) to maximize the effectiveness of substrates.

If the **BD BBL Crystal** test profile yields a "No identification" result and culture purity has been confirmed, then it is likely that (i) the test isolate is producing *atypical BD BBL Crystal reactions* (which may also be caused by procedural errors), (ii) the test species is not part of the intended taxa or (iii) the system is unable to identify the test isolate with the required level of confidence. Conventional test methods are recommended when user error has been ruled out.

PERFORMANCE CHARACTERISTICS

Reproducibility: In an external study involving three clinical laboratories, (total of three evaluations), the reproducibility of **BD BBL Crystal** RGP ID substrate (29) reactions was studied by replicate testing. The reproducibility of the individual substrate reactions ranged from 94.7% – 100%. The overall reproducibility of the **BD BBL Crystal** RGP ID panel was determined to be 99.4%.²⁰

Accuracy of Identification: The performance of **BD BBL Crystal** RGP ID System was compared to currently available commercial systems using **clinical isolates and stock cultures**. A total of three studies were conducted in three independent laboratories. Fresh, routine isolates arriving in the clinical laboratory, as well as previously identified isolates of the clinical trial sites' choice were utilized to establish performance characteristics.

Out of 604 total isolates tested from the studies, 550 (91.1%) were correctly identified (including isolates that required supplemental testing) by the **BD BBL Crystal** RGP Identification System. A total of 53 (8.8%) isolates were incorrectly identified, and a message of "No Identification" was obtained for one (0.2%) isolate.²⁰

AVAILABILITY

Cat. No.	Description	Cat. No.	Description
245150	BD BBL™ Crystal™ RGP Rapid Gram-Positive ID System, 1 Kit.	221165	BD BBL™ Columbia Agar with 5% Sheep Blood, pkg. of 20.
245038	BD BBL™ Crystal™ ANR, GP, RGP, N/H ID Inoculum Fluid.	221263	BD BBL™ Columbia Agar with 5% Sheep Blood, ctn. of 100.
245031	BD BBL™ Crystal™ Panel Viewer, Domestic model, 110 V, 60 Hz.	221352	BD BBL™ Columbia CNA Agar with 5% Sheep Blood, pkg. of 20.
245032	BD BBL™ Crystal™ Panel Viewer, European model, 220 V, 50 Hz.	221353	BD BBL™ Columbia CNA Agar with 5% Sheep Blood, ctn. of 100.
245033	BD BBL™ Crystal™ Panel Viewer, Japanese model, 100 V, 50/60 Hz.	221179	BD BBL™ Phenylethyl Alcohol Agar with 5% Sheep Blood, pkg. of 20.
245034	BD BBL™ Crystal™ Panel Viewer, Longwave UV Tube.	221277	BD BBL™ Phenylethyl Alcohol Agar with 5% Sheep Blood, ctn. of 100.
245036	BD BBL™ Crystal™ Panel Viewer, White Light Tube.	221239	BD BBL™ Trypticase™ Soy Agar with 5% Sheep Blood (TSA II), pkg. of 20.
245001	BD BBL™ Crystal™ ID System Electronic Codebook.	221261	BD BBL™ Trypticase™ Soy Agar with 5% Sheep Blood (TSA II), ctn. of 100.
245041	BD BBL™ Crystal™ Identification Systems Rapid Gram-Positive Manual Codebook.	212539	BD BBL™ Gram Stain Kit, pkg. of 4 x 250 mL bottles.

REFERENCES

1. Balows, A., W.J. Hausler, Jr., K.L. Hermann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
2. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey and Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
3. Bronfenbrenner, J., and M.J. Schlesinger. 1918. A rapid method for the identification of bacteria fermenting carbohydrates. Am. J. Public Health. 8:922-923.
4. Cowan, S.T., and K.J. Steel. 1974. Manual for the identification of medical bacteria. 2nd ed. Cambridge University Press, Cambridge.
5. Edberg, S.C., and C.M. Konnick. 1986. Comparison of b-glucuronidase-based substrate systems for identification of *Escherichia coli*. J. Clin. Microbiol. 24:368-371.
6. Ferguson, W.W., and A.E. Hook. 1943. Urease activity of *Proteus* and *Salmonella* organisms. J. Lab. Clin. Med. 28:1715-1720.
7. Hartman, P.A. 1968. Miniaturized microbiological methods. Academic Press, New York.
8. Kamper, P., O. Rauhoff, and W. Dott. 1991. Glycosidase profiles of members of the family *Enterobacteriaceae*. J. Clin. Microbiol. 29:2877-2879.
9. Killian, M., and P. Bulow. 1976. Rapid diagnosis of *Enterobacteriaceae* 1: detection of bacterial glycosidases. Acta Pathol. Microbiol. Scand. Sect. B. 84:245-251.
10. MacFaddin, J.F. 1980. Biochemical tests for identification of medical bacteria, 2nd ed. Williams & Wilkins, Baltimore.
11. Maddocks, J.L., and M. Greenan. 1975. Rapid method for identifying bacterial enzymes. J. Clin. Pathol. 28:686-687.
12. Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial diagnostics. Microbiol. Rev. 55:335-348.
13. Mandell, G.L., R.G. Douglas, Jr. and J.E. Bennett. 1990. Principles and practice of infectious diseases, 3rd ed. Churchill Livingstone Inc., New York.
14. Mangels, J., I. Edvalson, and M. Cox. 1993. Rapid Identification of *Bacteroides fragilis* group organisms with the use of 4-methylumbelliferone derivative substrates. Clin. Infect. Dis. 16(54):5319-5321.
15. Moncla, B.J., P. Braham, L.K. Rabe, and S.L. Hiller. 1991. Rapid presumptive identification of black-pigmented gram-negative anaerobic bacteria by using 4-methylumbelliferone derivatives. J. Clin. Microbiol. 29:1955-1958.
16. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
17. Sanders, A.C., J.E. Faber, and T.M. Cook. 1957. A rapid method for the characterization of enteric pathogen using paper discs. Appl. Microbiol. 5:36-40.
18. Sneath, P.H.A. 1957. The application of computers to taxonomy. J. Gen. Microbiol. 17:201-221.
19. Soto, O.B. 1949. Fermentation reactions with dried paper discs containing carbohydrate and indicator. Puerto Rican J. Public Health. Trop. Med. 25:96-100.
20. Data on file at Becton Dickinson Microbiology Systems.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Panel Location Pos. sur panel / Panel-Position / Posizione nel pannello / Posición en el panel	Test Feature Caractéristiques du test Testsubstrat Caratteristica del test Característica de la prueba	Code Codice Código	Principle (Reference) Principe (référence) Prinzip (Literaturnachweis) Principio (riferimento) Principio (Referencia)
4A	Fluorescent negative control Contrôle négatif fluorescent Fluoreszierende negative Kontrolle Controllo negativo fluorescente Control fluorescente negativo	FCT	Control to standardize fluorescent substrate results. Contrôle pour standardiser les résultats des substrats fluorescents. / Kontrolle zur Standardisierung der fluoreszierenden Substratergebnisse. / Controllo per standardizzare i risultati del substrato fluorescente. / Control para estandarizar los resultados del substrato fluorescente.
2A	4MU-β-D-glucoside 4MU-β-D-Glukosid 4MU-β-D-glucósido	FGC	Enzymatic hydrolysis of the amide or glycosidic bond results in the release of a fluorescent coumarin derivative. ^{5,8,11,12,14,15}
1A	L-proline / Prolin / prolina-AMC	FPR	L'hydrolyse enzymatique des liens amides ou glycosidiques libère un dérivé coumarinique fluorescent.
4B	L-arginine / Arginin / arginina-AMC	FAR	Enzymatische Hydrolyse der Amid- oder
2B	L-methionine / Méthionine / Methionin / metionina-AMC	FME	Glykosidbindungen setzt fluoreszierende Cumarinderivate frei.
1B	4MU-β-D-cellobioside 4MU-β-D-Cellobiosid 4MU-β-D-celobiósido	FCE	L'idrolisi enzimatica del legame ammidico o glicosidico causa il rilascio di un derivato coumarinico fluorescente. La hidrólisis enzimática del enlace amídico o glicosídico resulta en la producción de un derivado coumarínico fluorescente.
4C	4MU-phosphate / Phosphat / fosfato	FHO	
2C	L-pyroglutamic acid-AMC Acide L-pyroglutamique-AMC L-Pyroglutaminsäure-AMC L-acido piroglutamico-AMC L-ácido piroglutámico-AMC	FPY	
1C	L-tryptophan / tryptophane / Tryptophan / triptofano / triptófano-AMC	FTR	
4D	L-valine-AMC / L-Valin-AMC / L-valina-AMC	FVA	
2D	L-phenylalanine-AMC / L-phénylalanine-AMC / L-Phenylalanin-AMC / L-fenilalanina-AMC	FPH	
1D	4MU-α-D-glucoside / Glukosid / glucósido	FGS	
4E	Arabinose / Arabinosio / Arabinosa	ARA	Utilization of carbohydrate results in lower pH and change in indicator (Phenol red). ^{1,2,3,4,7,16}
2E	Maltose / Maltosio / Maltosa	MAL	La métabolisation des carbohydrates entraîne une baisse du pH, ce qui fait virer l'indicateur coloré (rouge de phénol).
1E	Dextrin / Dextrine / Dextrina / Dextrina	DXT	
4F	Mannitol / Mannitolo / Manitol	MNT	
2F	Galactose / Galaktose / Galatlósido / Galactosa	GAL	Die Verwendung von Kohlenhydrat resultiert in niedrigeren pH-Werten und einem Indikatorumschlag (Phenolrot)
1F	N-acetyl-D-glucosamine N-Azetyl-D-Glukosamin N-acetil-D-glucosamina	AGN	L'utilizzo del carboidrato fa abbassare il pH e cambiare colore all'indicatore (rosso fenolo).
4G	Trehalose / Tréhalose / Trealosio / Trehalosa	TRE	La utilización de carbohidratos resulta en una caída del pH y un cambio en el indicador (rojo de fenol).
2G	Mannose / Mannosio / Manosa	MNS	
1G	Maltotriose / Maltotriosio / Maltotriosa	MTT	
4H	o-nitrophenyl-β-D-galactoside (ONPG) & p-nitrophenyl-β-D-glucoside / o-nitrophényl-β-D-galactoside (ONPG) & p-nitrophényl-β-D-glucoside / o-Nitrophenyl-β-D-Galaktosid (ONPG) & p-Nitrophenyl-β-D-Glukosid / o-nitrofenil-β-D-galatlósido (ONPG) & p-nitrofenil-β-D-glucoside / o-nitrofenil-β-D-galactósido (ONPG) y p-nitrofenil-β-D-glucósido	POG	
2H	p-nitrophenyl-α-D-glucoside p-nitrophényl-α-D-glucoside p-Nitrophenyl-α-D-Glukosid p-nitrofenil-α-D-glucoside p-nitrofenil-α-D-glucósido	AGL	

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road
Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia