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INTENDED USE

BD BBL™ Coagulase Plasma, Rabbit and **BD BBL Coagulase Plasma, Rabbit with EDTA** are used to qualitatively determine the pathogenicity of staphylococci using the direct tube method.

SUMMARY AND EXPLANATION

Identification of staphylococci is based on microscopic examination, colonial morphology and cultural and biochemical characteristics. Staphylococci associated with acute infection (*Staphylococcus aureus* in humans and *S. intermedius* and *S. hyicus* in animals) can clot plasma. The most widely used and generally accepted criterion for identification of these pathogenic organisms is based on the presence of the enzyme coagulase.¹ The ability of *Staphylococcus* to produce coagulase was first reported by Loeb² in 1903.

Coagulase binds plasma fibrinogen, causing the organisms to agglutinate or plasma to clot. Two different forms of coagulase can be produced, free and bound. Free coagulase is an extracellular enzyme produced when the organism is cultured in broth. Bound coagulase, also known as clumping factor, remains attached to the cell wall of the organism. The tube test can detect the presence of both bound and free coagulase. Isolates that do not produce clumping factor must be tested for the ability to produce extracellular coagulase (free coagulase).

BD BBL Coagulase Plasma, Rabbit and **BD BBL Coagulase Plasma Rabbit with EDTA** are recommended for performing the direct tube test. The inoculum used for testing must be pure because a contaminant may produce false results after prolonged incubation. For the coagulase test, **BD BBL Coagulase Plasma, Rabbit with EDTA** is superior to citrated plasma because citrate-utilizing organisms such as *Pseudomonas* species, *Serratia marcescens*, *Enterococcus faecalis* and strains of *Streptococcus* will clot citrated plasma in 18 h.³

PRINCIPLES OF THE PROCEDURE

S. aureus produces two types of coagulase, free and bound. Free coagulase is an extracellular enzyme produced when the organism is cultured in broth. Bound coagulase, also known as the clumping factor, remains attached to the cell wall of the organism.

In the direct tube test, free coagulase liberated from the cell acts on prothrombin in the coagulase plasma to give a thrombin-like product. This product then acts on fibrinogen to form a fibrin clot.⁴

The tube test is performed by mixing an overnight broth culture or colonies from a non-inhibitory agar plate into a tube of rehydrated coagulase plasma. The tube is incubated at 37 °C. The formation of a clot in the plasma indicates coagulase production.

REAGENTS

Coagulase Plasma, Rabbit is lyophilized rabbit plasma with 0.85% sodium citrate and 0.85% sodium chloride, approximately.

Coagulase Plasma, Rabbit with EDTA is lyophilized rabbit plasma with 0.15% EDTA (ethylenediaminetetraacetic acid) and 0.85% sodium chloride, approximately.

Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Observe aseptic technique and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving.

Directions for use should be followed carefully.

Storage

Store unopened lyophilized **BD BBL Coagulase Plasma, Rabbit** and **BD BBL Coagulase Plasma, Rabbit with EDTA** at 2–8 °C.

Store reconstituted plasma at 2–8 °C for up to 14 days, or aliquot and freeze promptly at -20 °C for up to 30 days. Do not thaw and refreeze.

Expiration date applies to product in its intact container when stored as directed. Do not use if the product is caked, discolored or shows other signs of deterioration. Examine reconstituted reagents for evidence of contamination, evaporation or other signs of deterioration, such as cloudiness or partial clotting.

SPECIMEN COLLECTION AND PREPARATION

Collect specimens or samples in sterile containers or with sterile swabs and transport immediately to the laboratory according to recommended guidelines.^{1,4-9}

Process each specimen using procedures appropriate for that sample.^{1,4-9}

Select well-isolated colonies. The test described below requires the use of a pure test culture.

Suspicious growth, such as black colonies on Vogel and Johnson Agar or Tellurite Glycine Agar, or golden, hemolytic colonies from **BD Trypticase™ Soy Blood** agar plates should be selected for testing.

Using a bacteriological loop, transfer a well-isolated colony from a pure culture into a tube of sterile Brain Heart Infusion Broth or **BD Trypticase Soy Broth**. Incubate for 18–24 h or until a dense growth is observed. Alternatively, several colonies (1 loopful) taken directly from a non-inhibitory agar plate such as **BD Trypticase Soy Agar** may be used as an inoculum instead of a broth culture.

Test	MSSA* (222 isolates)		MRSA**(78 isolates)		Total (300 isolates)	
	No. of false-negative results	Test Sensitivity (%)	No. of false-negative results	Test Sensitivity (%)	No. of false-negative results	Test Sensitivity (%)
Free Coagulase	0	100	6	92.3	6	98.0
Bound Coagulase	0	100	3	96.1	3	99.0
Latex Agglutination Test #1	0	100	0	100	0	100
Latex Agglutination Test #2	2	99.1	12	84.6	14	95.3
Latex Agglutination Test #3	0	100	0	100	0	100

*Methicillin-susceptible *S. aureus*.

**Methicillin-resistant *S. aureus*.

In a second study by McDonald and Chapin,¹¹ the performance of **BD BBL** Coagulase Plasma was compared in a 2-h tube coagulase test (TCT) to two commercial latex agglutination tests for identifying *S. aureus* directly from blood culture broths and pellets obtained from supernatants of **BD BACTEC™** bottles. One hundred twelve (112) clinical blood culture isolates and 68 negative blood culture bottles seeded with a variety of gram-positive organisms were evaluated.

The table below gives the results of the coagulase test and the latex agglutination tests for both the seeded and clinical specimens.

Cultures and organism(s)	No. aerobic/ No. anaerobic (total)	No. of positive results					
		Direct			Pellet		
		Latex Test #1	Latex Test #2	TCT	Latex Test #1	Latex Test #2	TCT
Seeded							
<i>Staphylococcus aureus</i>	11/8 (19)	0	0	19	2	2	19
Coagulase-negative staphylococci	9/9 (18)	0	0	0	0	0	0
<i>Streptococcus pneumoniae</i>	5/5 (10)	0	0	0	0	0	0
<i>Enterococcus</i> spp.	6/5 (11)	5	5	0	5	5	0
<i>Streptococcus agalactiae</i>	3/5 (8)	0	0	0	0	0	0
<i>Streptococcus pyogenes</i>	1/1 (2)	0	0	0	0	0	0
Clinical							
Coagulase-negative staphylococci	70/3 (73)	0	0	0	0	0	0
<i>Staphylococcus aureus</i>	35/4 (39)	5	4	31	5	4	30

Among the 68 seeded blood culture bottles, the 2-h tube coagulase test using **BD BBL** Coagulase Plasma correctly identified 19 of 19 blood cultures seeded with *S. aureus*. There were no false positives with the 2-h tube coagulase test. Among the 112 clinical specimens tested, the 2-h tube coagulase test correctly identified 31 out of 39 *S. aureus* isolates directly in blood culture broth and 30 out of 39 *S. aureus* isolates in pelleted supernatants, for sensitivities of 79.5 and 76.9%, respectively. The specificity with both seeded and clinical isolates was 100% for the tube coagulase test.

AVAILABILITY

Cat. No. Description

240658	BD BBL™ Coagulase Plasma, Rabbit, 10 X 3.0 mL
240661	BD BBL™ Coagulase Plasma, Rabbit, 10 X 15.0 mL
240827	BD BBL™ Coagulase Plasma, Rabbit with EDTA, 10 X 3.0 mL
240826	BD BBL™ Coagulase Plasma, Rabbit with EDTA, 10 X 15.0 mL

REFERENCES

- Kloos, W. E., and T. L. Bannerman. 1999. *Staphylococcus* and *Micrococcus*, p. 264-282. In P.R. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover, Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
- Loeb, L. 1903. The influence of certain bacteria on the coagulation of the blood. *J. Med. Res.* 10:407-419.
- Bayliss, B.G. and E.R. Hall. 1965. Plasma coagulation by organisms other than *Staphylococcus aureus*. *J. Bacteriol.* 89:101-104.
- Pezzlo, M. (ed.). 1994. Aerobic bacteriology, p. 1.0.0.-1.20.47. In H. D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
- Baron, E.J., L.R. Peterson and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis, MO.
- Association of Official Analytical Chemists. 2000. Official methods of analysis of AOAC International, 17th ed. AOAC International, Arlington, VA.
- Association of Official Analytical Chemists. 2001. FDA Bacteriological analytical manual online. <<http://www.cfsan.fda.gov/~ebam/bam-mm.html>>.
- Downes, F.P. and K. Ito (ed.). 2001. Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.
- Flowers, R.S., W. Andrews, C.W. Donnelly and E. Koenig. 1993. Pathogens in milk and milk products, p. 103-212. In R.T. Marshall (ed.), Standard methods for the examination of dairy products, 16th ed. American Public Health Association, Washington, D.C.
- Luijendijk, A., A. van Belkum, H. Verbrugh and J. Kluytmans. 1996. Comparison of five tests for identification of *Staphylococcus aureus* from clinical samples. *J. Clin. Microbiol.* 34:2267-2269.
- McDonald, C.L. and K. Chapin. 1995. Rapid identification of *Staphylococcus aureus* from blood culture bottles by a classic 2-hour tube coagulase test. *J. Clin. Microbiol.* 33:50-52.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najnižja dovoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrense / Dolna granica temperatury / Limite mínimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sicaklık alt sınırı / Минимальная температура / 温度下限

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ATCC is a trademark of The American Type Culture Collection.

Bacto is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company.

© 2017 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.