# BBL™ Prepared Tubed Media for Differentiation of *Mycobacterium paratuberculosis*

Herrold's Egg Yolk Agar Slants with Mycobactin J and ANV  
Herrold's Egg Yolk Agar Slants (w/o Mycobactin J) with ANV

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## INTENDED USE

Herrold's Egg Yolk Agars are used for the selective isolation and differentiation of *Mycobacterium paratuberculosis*. Mycobacteria other than *M. paratuberculosis* will grow on Herrold's Egg Yolk Agar (w/o Mycobactin), while *M. paratuberculosis* will not grow unless the growth factor, Mycobactin, is provided.

## SUMMARY AND EXPLANATION

Paratuberculosis (Johne's disease) is a chronic granulomatous enteritis of ruminants first described in 1826.<sup>1</sup> Johne and Frothingham further described the disease in 1895 and demonstrated the presence of acid-fast bacilli in affected intestine.<sup>1</sup> The first successful isolation of this organism occurred in 1910.<sup>2</sup> Early nutritional studies demonstrated the stimulatory effect of an ethanol soluble substance obtained from other mycobacteria (e.g., *M. phlei*) called "mycobactin."<sup>3</sup> The organism fails to grow on the medium lacking the growth factor.

Cultivation of *M. paratuberculosis* from fecal specimens is the most frequently used method to diagnose paratuberculosis. Specific procedures are described elsewhere for the cultivation of *M. paratuberculosis* from bovine fecal specimens.<sup>4-7</sup>

## PRINCIPLES OF THE PROCEDURE

Enzymatic digest of casein provides amino acids and other nitrogenous substances. Beef extract provides additional nitrogenous nutrients, vitamins and minerals required for microbial growth. Sodium chloride maintains osmotic equilibrium. Sodium pyruvate is a source of energy for bacterial metabolism. Egg yolk and glycerol provide fatty acids and other nutrients required for the metabolism of mycobacteria.

Most mycobacteria produce endogenous siderophores (iron binding compounds) called mycobactin and will grow on Herrold's Egg Yolk Agar whether or not mycobactin has been added. *M. paratuberculosis* is unable to grow on media lacking mycobactin. Mycobactin J permits earlier and more abundant growth of *M. paratuberculosis* than other mycobactins.<sup>8,9</sup>

Amphotericin B enhances the selectivity of the medium by inhibiting contaminating fungi. Nalidixic acid inhibits contaminating gram-negative organisms and vancomycin inhibits contaminating gram-positive organisms. Malachite green is included to help control contaminants and enhance the visibility of colonies.

## REAGENTS

### Formulae:

Herrold's Egg Yolk Agar with Mycobactin J and ANV

Approximate Formula\*

Pancreatic Digest of Casein	9.0	g
Beef Extract	2.7	g
Sodium Chloride	4.0	g
Sodium Pyruvate	4.1	g
Malachite Green	0.1	g
Amphotericin B	0.05	g
Mycobactin J	0.002	g
Agar	15.3	g
Nalidixic Acid	0.05	g
Vancomycin	0.05	g
Egg Yolks	100.0	mL
Glycerol	27.0	mL
Purified Water	873.0	mL

\*Adjusted and/or supplemented as required to meet performance criteria.

Herrold's Egg Yolk Agar (without Mycobactin J) with ANV contains the same ingredients except that Mycobactin J is omitted.

### Warnings and Precautions:

For Veterinary Use.

Tubes with tight caps should be opened carefully to avoid injury due to breakage of glass.

**Storage Instructions:** On receipt, store media in the dark at 2 – 8°C. Avoid freezing and overheating. Do not open until ready to use. Minimize exposure to light. Media stored as labelled until just prior to use may be inoculated up to the expiration date and incubated for up to 16 weeks.

After storage in the upright position, tubes of egg media may appear to have a large amount of water at the bottom of the slants. The fluid will be reabsorbed into the medium if the tubes are placed on their sides so that the fluid covers the surface of the slant.

Allow the medium to warm to room temperature before inoculation.

**Product Deterioration:** Do not use medium if it shows evidence of microbial contamination, discoloration, drying or other signs of deterioration.

## SPECIMEN COLLECTION AND HANDLING

Refer to appropriate texts for details of specimen collection and handling procedures.<sup>10</sup>

Observe established precautions against microbiological hazards throughout all procedures. All specimens should be handled according to CDC-NIH recommendations, NCCLS guidelines or local institution guidelines for any potentially infectious materials. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.

## PROCEDURE

**Material Provided:** Depending upon which product is ordered, one of the prepared media listed below (see "AVAILABILITY") is provided.

**Materials Required But Not Provided:** Ancillary culture media, reagents, quality control organisms and laboratory equipment are required for this procedure.

**Test Procedure:** Observe aseptic techniques.

1. Weigh out 1 – 2 g of fecal sample.
2. Process specimen using a recommended procedure<sup>4-7</sup> to remove or kill contaminating organisms other than mycobacteria.
3. Inoculate tubes of Herrold's Egg Yolk Agar with and without Mycobactin J with no more than 0.25 mL (5 – 7 drops from a sterile transfer pipette) of the final processed and decontaminated material.
4. Incubate tubes at 35 – 37°C in a slanted position with the caps loose.
5. Tighten caps when medium surface is dry (1 – 2 weeks) and place in an upright position in the incubator.
6. Read and evaluate tubes for growth and contamination every week for up to 16 weeks.

### User Quality Control:

1. Examine medium for signs of deterioration as described under "Product Deterioration."
2. Check performance by inoculating a representative sample of medium with 0.1 mL of a 1:100 dilution of a No. 1 McFarland suspension in water of pure cultures of stable control organisms that produce known, desired reactions. The following test strains are recommended.

Medium	Test Strain	Expected Results
Herrold's Egg Yolk Agar with Mycobactin J and ANV	<i>Mycobacterium paratuberculosis</i> ATCC™ 19698	Growth within 6 weeks
	<i>Mycobacterium intracellulare</i> ATCC 13950	Growth within 3 weeks
Herrold's Egg Yolk Agar (without Mycobactin J) with ANV	<i>Mycobacterium paratuberculosis</i> ATCC 19698	No growth in 6 weeks
	<i>Mycobacterium intracellulare</i> ATCC 13950	Growth within 3 weeks

## RESULTS

Colonies appearing should be evaluated for typical acid-fastness and morphological appearance of *M. paratuberculosis*. *M. paratuberculosis*-like colonies should not appear on the medium lacking Mycobactin J.

To confirm Mycobactin J dependency, suspend suspect colonies in 0.5 mL autoclaved distilled or deionized water so that the turbidity is equivalent to a No. 1.0 McFarland nephelometer standard. Dilute the suspension 100-fold in autoclaved distilled or deionized water. Inoculate 0.1 mL of the suspension onto a single tube of Herrold's Egg Yolk Agar with Mycobactin J and 0.1 mL onto a single tube of the medium without Mycobactin J. Incubate tubes for 1 week with loosened caps at 35 – 37°C. Tighten caps and reincubate. Examine tubes weekly for the presence of growth. Acid-fast cultures that grow only in the presence of Mycobactin J are identified as *M. paratuberculosis*.

## LIMITATIONS OF THE PROCEDURE

Results for individual laboratories may differ depending on specimen storage conditions, the time between collection and receipt in the laboratory, the efficiency of decontamination procedures, the types and amount of spore-forming organisms in the sample, frequency of reading cultures, organism distribution within the inoculum and other conditions which may vary from laboratory to laboratory.

Subcultures of *M. paratuberculosis* may occasionally lose their dependency on Mycobactin. Reduced growth of the organism may be observed after multiple transfers.

A single medium is rarely adequate for detecting all organisms of potential significance in a specimen. The agents in selective media may inhibit some strains of the desired species or permit growth of a species they were designed to inhibit, especially if the species is present in large numbers in the specimen. Specimens cultured on selective media should, therefore, also be cultured on nonselective media to obtain additional information and help ensure recovery of potential pathogens.

## PERFORMANCE CHARACTERISTICS<sup>11</sup>

Cultures were performed in two different laboratories: a state-supported veterinary diagnostic lab (Lab A) and a veterinary school-based research lab (Lab B). Each laboratory processed fecal samples from 25 different animals that were known to be culture positive for *M. paratuberculosis*. Specimens represented animals characterized as high (>50 CFU/slant), moderate (10 – 50), and low (<10) shedders. Samples were removed from storage at –70°C, thawed rapidly, and processed using the double decontamination method as described by Whitlock.<sup>5</sup> Approximately six drops of processed sediment were inoculated onto each Herrold's Egg Yolk Agar (HEYA) slant (with and without added mycobactin). Each lab prepared their own Herrold's Egg Yolk Agar according to the same formulation<sup>1</sup> using different sources of raw materials.

In Lab A, a single tube of BBL HEYA and one tube each from four different lots of lab-prepared HEYA were inoculated using a portion of the same decontaminated specimen. Tubes were observed biweekly, beginning on week 2, for a total of 12 weeks. *M. paratuberculosis* was recovered on all five media for 16 of the specimens. Of the nine remaining cultures, growth was recovered on all nine BBL HEYA slants while one

or more tubes of lab-prepared media did not grow the organism. Growth was evident by week 4 on seven BBL HEYA slants. After 6 weeks of incubation, *M. paratuberculosis* was recovered on the remaining BBL HEYA media and on 15 of the lab-prepared lots. The remaining cultures on lab-prepared media were positive in 8 – 12 weeks.

In Lab B, duplicate tubes of a single lot of BBL HEYA and four different lots of lab-prepared HEYA were inoculated. Cultures were observed weekly for 14 weeks. No growth was observed on any medium for four samples. Two cultures from low shedders were recovered on lab-prepared media and not on BBL HEYA: one on all four lab-prepared lots of media, and one on two of four lab-prepared lots. Similarly, two cultures recovered on BBL HEYA were only recovered on two of four lots of lab-prepared media. This variation in recovery was attributable to the low numbers of organisms in the inoculum (no more than two colonies on any slant). One specimen that was recovered on all lots of lab-prepared media was overgrown with contaminating organisms before growth of *M. paratuberculosis* was detected. Of the remaining 15 cultures, 8 were positive one or more weeks earlier on BBL HEYA than on any lab-prepared lot. Four were positive on BBL HEYA and one lot of lab-prepared medium, while the other lab-prepared media were delayed a week or more. Three cultures were recovered a week earlier on lab-prepared media than on BBL HEYA.

The results of these studies indicate that recovery of *M. paratuberculosis* on BBL HEYA is substantially equivalent to freshly prepared lab media of similar formulation. These results suggest that some cultures may be detectable a week or more earlier on BBL HEYA.

#### AVAILABILITY

Cat. No.	Description
222232	BBL™ Herrold's Egg Yolk Agar Slants with Mycobactin J and ANV, Pkg. of 10 size C tubes.
222233	BBL™ Herrold's Egg Yolk Agar Slants with Mycobactin J and ANV, Ctn. of 100 size C tubes.
222240	BBL™ Herrold's Egg Yolk Agar Slants (without Mycobactin J) with ANV, Pkg. of 10 size C tubes.
222241	BBL™ Herrold's Egg Yolk Agar Slants (without Mycobactin J) with ANV, Ctn. of 100 size C tubes.
245158	BBL™ Herrold's Egg Yolk Agar Slants with Mycobactin J and ANV, Pkg. of 10 Mycoflask™.

#### REFERENCES

- Chiodini, R.J., H.J. van Kruijning, and R.S. Merkal. 1984. Ruminant paratuberculosis (Johne's disease): the current status and future prospects. *Cornell Vet.* 74: 218-262.
- Twort, F.W. 1910. A method for isolating and growing the lepra bacillus of man (preliminary note). *Proc. R. Soc. Lond. Ser. B.* 83: 156-158.
- Twort, F.W., and G.L.Y. Ingram. 1912. A method for isolating and cultivating *Mycobacterium enteritidis chronicae pseudotuberculosis bovis*, Johne, and some experiments on the preparation of a diagnostic vaccine for pseudotuberculosis enteritis of bovines. *Proc. R. Soc. Lond. Ser. B.* 84: 517-542.
- Whipple, D.L., D.R. Callihan, and J.L. Jarnagin. 1991. Cultivation of *Mycobacterium paratuberculosis* from bovine fecal specimens and a suggested standardized procedure. *J. Vet. Diagn. Invest.* 3: 368-373.
- Whitlock, R.H., A.E. Rosenberger, and P.A. Spencer. 1989. Laboratory culture techniques for Johne's disease: a critical evaluation of contamination and incubation times. *Proc. Annu. Meet. U.S. Animal Health Assn.* 93: 382-386.
- Whitlock, R.H. and R.W. Sweeney. 1990. Johne's disease: current aspects of transmission and the art of fecal culture. *Proc. Annu. Meet. Livestock Conservation Inst.*, p. 24-30.
- Stabel, J.R. 1997. An improved method for cultivation of *Mycobacterium paratuberculosis* from bovine fecal samples and comparison to three other methods. *J. Vet. Diagn. Invest.* 9: 375-380.
- Merkal, R.S., and W.G. McCullough. 1982. A new mycobactin, mycobactin J from *Mycobacterium paratuberculosis*. *Curr. Microbiol.* 7: 333-335.
- Thoen, C.O., and K.H. Baum. 1988. Current knowledge of paratuberculosis. *J. Am. Vet. Med. Assoc.* 192: 1609-1611.
- U.S. Department of Agriculture. 1985. Laboratory methods in veterinary mycobacteriology. National Veterinary Services Laboratories, Animal and Plant Health Inspection Service, Ames, Iowa.
- Data on file, BD Diagnostics.

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©2004 BD.
ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL and Mycoflask are trademarks of Becton, Dickinson and Company.
Becton, Dickinson France S.A.