BD BBL™ Taxo™ Differentiation Discs Urea

English: pages 1 – 2 Italiano: pagine 4 – 5
Français: pages 2 – 3 Español: páginas 5 – 6
Deutsch: Seiten 3 – 4

 8820311JAA(02)
2015-04

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BBL™ Taxo™ Differentiation Discs Urea are paper discs containing urea for differentiating microorganisms based on urease production.

SUMMARY AND EXPLANATION

The urease test is a qualitative procedure for determining the ability of microorganisms to hydrolyze urea.¹ This analysis is used for the identification and classification of *Enterobacteriaceae*, *Pseudomonas* spp., *Haemophilus* spp., mycobacteria and other gram-negative bacteria.²⁻⁴ Gray and Roberts found the urease test to be useful in distinguishing between clinically significant yeast isolates from humans.⁵

In 1987, Mills, Grimes and Gherna reported that Differentiation Disc Urea was a rapid test method that detected urea hydrolysis in a larger number of anaerobes than the conventional microbiological test tube assays. Conventional tube assays were more costly and time consuming. They also found the rapid disc method to generate more consistent results that were easier to interpret than the spot test method.⁶

PRINCIPLES OF THE PROCEDURE

The urea in the discs is broken down in the presence of urease. The urease catalyzes the hydrolysis of the carbon-nitrogen bonds producing carbon dioxide, water, ammonia and ammonium carbonate. This results in an alkaline end-product, causing the phenol red indicator in the test solution to turn cerise/red in color.

REAGENTS

BBL Taxo Differentiation Discs Urea are 1/4" yellowish to slightly orange paper discs. The discs contain a sufficient amount of Urea R Broth to yield a positive reaction with organisms producing sufficient urease to hydrolyze the substrate.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow established laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store **BBL Taxo Differentiation Discs Urea** at 2 – 8 °C.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

This product is recommended for use only with pure cultures.^{5,6}

PROCEDURE

Material Provided: **BBL Taxo Differentiation Discs Urea**.

Materials Required But Not Provided: 35 °C incubator or water bath, inoculating loop, Bunsen burner or incinerator, quality control organisms, 12 x 75 mm or 13 x 100 capped test tubes and distilled water.

Test Procedure:

Rapid Urease Test

For microorganisms other than mycobacteria:

1. Dispense 0.5 mL sterile distilled or deionized water into a 12 x 75 mm test tube.
2. Completely emulsify a large loopful of test organism in the water and shake or vortex to form a smooth suspension.
3. Aseptically add one **BBL Taxo Urea Disc** to the tube and cap.
4. Incubate the tube at 35 ± 2 °C for up to 4 h.
5. Read the tube after 10 min, 30 min, 1, 2, 3 and 4 h.
6. Observe and record results.

For mycobacteria:³

1. Dispense 0.5 mL sterile distilled or deionized water into a 13 x 100 mm screw capped test tube.
2. Emulsify a heavy loopful of organism from a 3 week old, egg medium culture.
3. Aseptically add one **BBL Taxo Urea Disc** to the tube and cap.
4. Incubate tubes at 35 °C.
5. Read results after 1 h and daily for up to 3 days.

User Quality Control:

Identity Specifications – Yellowish to slightly orange paper discs. The discs are imprinted with U on both sides.

Cultural Response –

For organisms other than mycobacteria:

Organism	ATCC®	Reaction within 4 h at 35 °C
<i>Escherichia coli</i>	25922	–
<i>Proteus mirabilis</i>	12453	+

For mycobacteria:

Organism	ATCC®	Reaction within 3 days at 35 °C
<i>Mycobacterium scrofulaceum</i>	19981	+
<i>Mycobacterium fortuitum</i>	6841	+
<i>Mycobacterium goodii</i>	14470	–

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

A positive urease reaction is indicated by formation of a cerise/red color in the solution. A negative urease reaction is indicated by no color change or slight orange color in the solution.

LIMITATIONS OF THE PROCEDURE

The rapid urease test may be used as a presumptive diagnostic aid for identifying those organisms that hydrolyze urea. However, additional biochemical and serological tests should be performed for complete identification.

Broth cultures in nutrient media are not usually satisfactory for urease testing, as residual alkali or acid from the culture may give false positive or false negative reactions, respectively.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of **BBL Taxo Differentiation Discs Urea** are tested for performance characteristics.

For microorganisms other than mycobacteria:

Tubes containing 0.5 mL of autoclaved, purified water are inoculated with a large loopful of *Proteus mirabilis* (ATCC 25933), *P. vulgaris* (ATCC 13315), *Salmonella typhimurium* (ATCC 14028) and *Escherichia coli* (ATCC 25922) and mixed gently. Using sterile forceps, one **Taxo Differentiation Discs Urea** is inserted into each inoculated tube. Tubes are incubated at 35 ± 2 °C and examined after 10 min and 1, 2, 3 and 4 h. *P. mirabilis* and *P. vulgaris* give a positive reaction for urease (a cerise/red color in the solution) within 4 h incubation. *S. typhimurium* and *E. coli* give a negative reaction for urease (no color change or a slightly orange color in the solution) after 4 h incubation.

For mycobacteria:

Tubes containing 0.5 mL of autoclaved, purified water are inoculated with a large loopful of 3-week-old egg medium cultures of *Mycobacterium scrofulaceum* (ATCC 19981), *M. fortuitum* (ATCC 6841) and *M. goodii* and emulsified. Using sterile forceps, one **Taxo Differentiation Discs Urea** is inserted into each inoculated tube. Tubes are incubated at 35 ± 2 °C and examined after 1 h and up to 3 days. *M. scrofulaceum* and *M. fortuitum* give a positive reaction for urease (a cerise/red color in the solution) within 3 days of incubation. *M. goodii* gives a negative reaction for urease (no color change or a slightly orange color in the solution) after 3 days incubation.

AVAILABILITY

Cat. No.	Description
231737	BD BBL™ Taxo™ Differentiation Discs Urea

REFERENCES

1. MacFaddin, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed., p. 298-308. Williams and Wilkins, Baltimore.
2. Qadri, S.M.H., S. Zubairi, H.P. Hawley, and E.G. Ramirez. 1984. Simple spot test for rapid detection of urease activity. J. Clin. Microbiol. 20:1198-1199.
3. Kent, P.T., and C.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory, p. 117-119. Centers for Disease Control, Atlanta.
4. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
5. Gray, L.D., and G.D. Roberts. 1988. Identification of medically important yeasts. Clin. Microbiol. Newsl. 10:73-80.
6. Mills, C.K., B.Y. Grimes, and R.L. Gherna. 1987. Three rapid methods compared with a conventional method for detection of urease production in anaerobic bacteria. J. Clin. Microbiol. 11:2209-2210.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

BD BBL Taxo Differentiation Discs Urea

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Benex Limited