# **BD BBL™ Port-A-Cul™ Specimen Collection and Transport Products**

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## INTENDED USE

**BD BBL™ Port-A-Cul™** tubes, jars and vials contain a reduced transport medium and are intended to maintain the viability of anaerobic, facultative and aerobic microorganisms during transport from the patient to the laboratory. Sterile packages are for collection of specimens in clean areas; e.g., surgical suites.

## SUMMARY AND EXPLANATION

Anaerobic microorganisms require special transport methods to prevent even brief exposure to oxygen.<sup>1,2</sup> Maintaining aerobic, anaerobic and facultative organisms in an anaerobic environment provides an ideal method of transport.<sup>3-6</sup>

**BD BBL Port-A-Cul Tubes** are for swab specimens. Swab specimens are inserted into a reduced solidified holding medium, in contrast to other methods in which swabs are placed into fluid media where organisms may be diluted or into an empty tube where organisms may dry out. **BD BBL Port-A-Cul Transport Jars** are for tissue and biopsy specimens. The wide-mouth jar allows for easier insertion of the specimen into the reduced solidified holding medium. **BD BBL Port-A-Cul Vials** are for fluid specimens. Fluid specimens are injected directly onto the solid agar surface.

## PRINCIPLES OF THE PROCEDURE

Swab and tissue specimens are protected from air and remain moist in a column of solid medium under an anaerobic atmosphere. Agar inhibits oxygen diffusion. Reducing agents combine with free oxygen. Salts and buffers provide a protective environment. Resazurin indicates the presence of oxygen (pink to blue).

Fluid specimens are protected from air in a **BD BBL Port-A-Cul** vial because of the anaerobic atmosphere and exchange of reducing agents from the solid medium to the specimen.

## REAGENTS

**BD BBL Port-A-Cul** medium - balanced formula of reducing agents and resazurin in a buffered isotonic agar base.

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

 **BD BBL Port-A-Cul** is for single use only; reuse may cause a risk of infection and/or inaccurate results.

Exercise care in opening container with tight cap.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>7-10</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Autoclave contaminated materials after use and discard in appropriate biohazard waste container.

## BD BBL Port-A-Cul Transport Jar Sterile Pack

Warning

**H302** – Harmful if swallowed.

**P101** – If medical advice is needed, have product container or label at hand. **P264** – Wash thoroughly after handling. **P301 + P312** – IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** – Rinse mouth. **P501** – Dispose of contents/container in accordance with local/regional/national/ international regulations.

## BD BBL™ Port-A-Cul™ Tube

Warning

**H335** – May cause respiratory irritation.

**P101** – If medical advice is needed, have product container or label at hand. **P261** – Avoid breathing dust/fume/gas/mist/vapours/spray. **P304 + P340** – IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P312** – Call a POISON CENTER or doctor/physician if you feel unwell. **P403 + P233** – Store in a well-ventilated place. Keep container tightly closed. **P501** – Dispose of contents/container in accordance with local/regional/national/ international regulations.

## BD BBL Port-A-Cul Tube and Swabs Sterile Pack

Warning

**H302** – Harmful if swallowed.

**P264** – Wash thoroughly after handling. **P301 + P312** – IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** – Rinse mouth. **P501** – Dispose of contents/container in accordance with local/regional/national/international regulations.

## BD BBL Port-A-Cul Vial

Warning

**H335** – May cause respiratory irritation.

**P101** – If medical advice is needed, have product container or label at hand. **P261** – Avoid breathing dust/fume/gas/mist/vapours/spray. **P304 + P340** – IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P312** – Call a POISON CENTER or doctor/physician if you feel unwell. **P403 + P233** – Store in a well-ventilated place. Keep container tightly closed. **P501** – Dispose of contents/container in accordance with local/regional/national/ international regulations.

## BD BBL Port-A-Cul Vial Sterile Pack

Warning

**H335** – May cause respiratory irritation.

**P261** – Avoid breathing dust/fume/gas/mist/vapours/spray. **P304 + P340** – IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P312** – Call a POISON CENTER or doctor/physician if you feel unwell. **P403 + P233** – Store in a well-ventilated place. Keep container tightly closed. **P501** – Dispose of contents/container in accordance with local/regional/national/ international regulations.

**Storage Instructions:** On receipt, store at 20 – 25 °C. Media stored as labeled until just prior to use may be inoculated up to the expiration date and transported for the recommended transportation times.

**Product Deterioration:** Do not use if medium shows evidence of contamination, drying, cracking or other signs of deterioration (flocculation and precipitation are not indications of deteriorated medium) or if integrity of Sterile Pack packaging has been compromised. Do not use if a pink or blue band at the top of the medium is deeper than 3 mm, as this is indicative of excessive oxidation.

## SPECIMEN COLLECTION AND TRANSPORT

Follow accepted specimen collection practices.<sup>11</sup> Fluid specimens are preferred. Do not culture body sites normally contaminated with indigenous flora. Transport to the laboratory at ambient temperatures (20 – 25 °C) within 72 h.

## PROCEDURE

**Material Provided:** BBL Port-A-Cul products (see "Availability").

**Materials Required But Not Provided:** Syringes and needles for Cat. Nos. 221608/221609; swabs for Cat. No. 221606; sterile forceps and materials required for organism isolation, identification and other microbiological procedures.

## Test Procedure:

**SWAB SPECIMENS** – Obtain specimen, loosen screw cap, insert swab into medium within approximately 5 mm from bottom, break shaft evenly with tube lip and quickly replace and tighten cap. (H<sub>2</sub>S may be smelled upon opening tube; this does not adversely affect medium.) Transport to laboratory promptly (within 72 h).

Remove swab by grasping applicator shaft with sterile forceps; streak onto an appropriate growth medium.

**FLUID SPECIMENS** – Remove green flip cap exposing rubber stopper. Swab rubber stopper with disinfectant. Expel air from syringe and needle. Obtain specimen, push needle through stopper and **slowly** inject on surface of agar. Transport to laboratory promptly (within 72 h).

Prior to aspirating specimen for culture, swab stopper with disinfectant.

**TISSUE SPECIMENS** – Obtain specimen, loosen screw cap, insert specimen into medium within approximately 5 mm from bottom and quickly replace and tighten cap. Transport to laboratory promptly (within 72 h).

To remove specimen, loosen cap and remove specimen with sterile forceps.

**STERILE PACK** – Peel envelope apart outside clean area without touching contents of envelope. Wearing sterile gloves, remove sterile tube and swabs, jar or vial. Take into clean area. See directions for use on envelope for swab, fluid, or tissue transport.

#### **Quality Control:**

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

#### **RESULTS**

Little or no growth of microorganisms should occur within the specimen while it is being transported in the **BD BBL™ Port-A-Cul™** medium. However, organisms in logarithmic growth phase or in the nutritious environment of the specimen itself may continue to grow until endogenous or exogenous nutrients are depleted.

#### **LIMITATIONS OF THE PROCEDURE**

For best results, specimens should be transported to the laboratory as quickly as possible, but no longer than 72 h. Overgrowth may occur with polymicrobial infections. Avoid temperature extremes during transport. Organisms in small numbers (< 100 CFU/mL) may not survive longer than 24 h.

Viability of *Neisseria gonorrhoeae* decreases rapidly at temperatures below 35 °C and may not survive more than 8 h even in high concentration.

Caution should be exercised in reporting direct Gram stain and/or other direct microbiological stain results on specimens processed with this medium due to the possible presence of nonviable organisms in the transport medium.

Transport media, staining reagents, immersion oil, glass slides and specimens themselves sometimes contain dead organisms visible upon Gram staining.

#### **PERFORMANCE CHARACTERISTICS**

The survival of *Streptococcus pyogenes* and *Haemophilus influenzae* was found to be better in **BD BBL Port-A-Cul** than in an aerobic transport device.<sup>6</sup> **BD BBL Port-A-Cul** was evaluated with stock cultures of a variety of clinically significant anaerobes and found to maintain viability for at least 72 h at 25 °C.<sup>4</sup> In another study, using simulated specimens stored up to 48 h at ambient temperature, recovery of ten specimens of anaerobic bacteria was greater in **BD BBL Port-A-Cul** than in other controlled conditions.<sup>5</sup>

#### **AVAILABILITY**

<b>Cat. No.</b>	<b>Description</b>
221606	<b>BD BBL™ Port-A-Cul™</b> Tube, 10.
221607	<b>BD BBL™ Port-A-Cul™</b> Tube and Swabs Sterile Pack, 10.
221608	<b>BD BBL™ Port-A-Cul™</b> Vial, 10.
221609	<b>BD BBL™ Port-A-Cul™</b> Vial Sterile Pack, 10.
221602	<b>BD BBL™ Port-A-Cul™</b> Transport Jar Sterile Pack, 10 X 30 mL.

#### **REFERENCES**

1. Loesche, W.J. 1969. Oxygen sensitivity of various anaerobic bacteria. *Appl. Microbiol.* 18:723-727.
2. Chow, A.W., R.D. Leake, T. Yamauchi, B.F. Anthony, and L.B. Guze. 1974. The significance of anaerobes in neonatal bacteremia: analysis of 23 cases and review of the literature. *Pediatrics* 54:736-745.
3. Finegold, S.M., V.L. Sutter, H.R. Attebery, and J.E. Rosenblatt. 1974. Isolation of anaerobic bacteria, p. 365-375. *In* E.H. Lennette, E.H. Spaulding, and J.R. Truant (ed.), *Manual of clinical microbiology*, 2nd ed. American Society for Microbiology, Washington, D.C.
4. Abramson, I.J., G.L. Evans, and P. Bathurst. 1977. Laboratory evaluation of Port-A-Cul™ with stock cultures of anaerobes, aerobes, and facultative organisms, abstr. C 197, p. 68. *Abstr. 77th Annu. Meet. Am. Soc. Microbiol.* 1977.
5. Mena, E., F.S. Thompson, A.Y. Armfield, V.R. Dowell, Jr., and D.J. Reinhardt. 1978. Evaluation of Port-A-Cul transport system for protection of anaerobic bacteria. *J. Clin. Microbiol.* 8:28-35.
6. Hamilton, R.L., G.L. Evans, and T.E. O'Neill. 1979. Survival of facultatively anaerobic bacteria in Port-A-Cul™ tubes and an aerobic transport device, abstr. C 173, p. 338. *Abstr. 79th Annu. Meet. Am. Soc. Microbiol.* 1979.

- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
- U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000, p. 0021-0045.*
- Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover (ed.). 2003. *Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or [www.bd.com/ds](http://www.bd.com/ds).

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Benex Limited
Becton Dickinson Pty Ltd.