# **BD BBL™ Sensi-Disc™ Antimicrobial Discs** For Use in Culture Media

Antimicrobial Agent	Code	Conc.
Streptomycin	*S-50	50 µg
Isoniazid (Isonicotinyl hydrazine)	INH-1 INH-5	1 µg 5 µg
P-Aminosalicylic Acid	PAS-10 PAS-50	10 µg 50 µg
Ethambutol**	EM-25 EM-50	25 µg 50 µg
Ethionamide***	EA-25	25 µg
Rifampin	*RA-25	25 µg

English: pages 1 – 5      Italiano: pagine 13 – 17      8840231JAA(03)  
 Français: pages 5 – 9      Español: páginas 17 – 21      2015-07  
 Deutsch: Seiten 9 – 13

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyttä lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Inštrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

## INTENDED USE

These discs are used in qualitative susceptibility testing procedures in culture media. They serve as a convenient method for addition of antimicrobial agents to culture media, especially for qualitative studies of mycobacteria and related organisms.

Certain susceptibility test discs, i.e., streptomycin, 10 µg, kanamycin, 30 µg and rifampin, 5 µg, are also recommended for use in this procedure.<sup>1,2</sup>

The discs are added to measured amounts of liquified agar medium, e.g., Middlebrook and Cohn 7H10 Agar or Seven H11 Agar, to prepare recommended concentrations of the antimicrobial agent in solidified agar.<sup>2</sup>

Falcon™ four compartment X Plate™ Petri dishes may be used with these discs. Three sections may be employed for various agents or differing concentrations of the same agent and the fourth section as a drug-free control.

\* Not for use in conventional (Bauer-Kirby or standardized) disc susceptibility tests.

\*\* Myambutol™: Wyeth Laboratories

\*\*\* Trecator™: Wyeth Laboratories

## SUMMARY AND EXPLANATION

Since 1882 when Koch showed that the tubercle bacillus was the cause of tuberculosis, a variety of media and methods have been investigated in an effort to reduce the time required for isolation and drug susceptibility testing.

Acid-fast bacilli from human specimens demonstrated by microscopic examination or culture are usually *Mycobacterium tuberculosis*; however, other pathogenic and possibly nonpathogenic acid-fast organisms may be found. These "atypical" or "unclassified" mycobacteria must be recognized and characterized in the laboratory because the prognosis and response to drug therapy differ greatly for disease caused by these organisms.

This disc susceptibility method employing discs in culture media has been shown to be a valuable and reliable tool for many laboratories.<sup>3-5</sup> The simplicity of preparation of disc medium, the convenience and economy of preparing any number of test plates, the stability of drugs in dry discs, the elimination of possible errors in diluting and measuring drug solutions and the immediate recognition of the drug and concentration present in each medium by means of coded discs make this a desirable technique. Only one medium is prepared; the different drugs and their concentrations are provided by the discs.

## PRINCIPLES OF THE PROCEDURE

These discs are employed to impregnate the entire agar bed with a specific concentration of antimicrobial-agent. This is accomplished by first affixing the disc to the center of the plate quadrant and then adding a measured amount of agar. The plate is placed in the refrigerator and held overnight to permit diffusion throughout the medium. The plate is then inoculated with the test organism. Growth or no growth will occur, depending on whether the organism is resistant or susceptible to the agent and concentration employed in the quadrant. Since zones are not formed, the numbers of colonies or lack of growth on each quadrant is used to indicate drug susceptibility.

## REAGENTS

Sensi-Disc brand discs are prepared by impregnating high quality absorbent paper with accurately determined amounts of antimicrobial agents. Discs are clearly marked on both sides with letters and numbers designating the agent and the drug content. Sensi-Disc agents are furnished in cartridges containing 50 discs each. The last disc in each cartridge is marked "X" and contains the drug as coded. The cartridges are for use in Single Disc Dispensers.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures.<sup>1,6</sup> Perform all manipulations in a biological safety cabinet. Sterilize cultures, containers and other contaminated materials after use.

Biosafety Level 2 practices and procedures, containment equipment and facilities are required for non-aerosol-producing manipulations of clinical specimens such as preparation of acid-fast smears. All aerosol-generating activities must be conducted in a Class I or II biological safety cabinet. Biosafety Level 3 practices, containment equipment and facilities are required for laboratory activities in the propagation and manipulation of cultures of *M. tuberculosis* and *M. bovis*. Animal studies also require special procedures.<sup>6</sup>

## Storage Instructions:

1. On receipt, store containers of discs at -20 to +8 °C. If the laboratory refrigerator is frequently opened and closed, and a suitable temperature is not maintained, place there a supply sufficient only for use within a week.
2. Allow containers to come to room temperature before opening. Return unused discs to the refrigerator when application of discs has been completed.
3. Use the oldest discs first.

**Product Deterioration:** Discard expired discs. Also discard cartridges from which discs have been frequently removed during a week. Discard containers left out overnight in the laboratory, or else test the discs for performance. The expiration date applies to discs in intact container stored as directed. Do not open until ready to use.

## SPECIMEN COLLECTION AND HANDLING<sup>1,2</sup>

Samples and digested, concentrated specimens, e.g., sputum, gastric lavage, urine, pleural, spinal, joint and other fluids, exudates and tissues, in which acid-fast bacilli are seen in smears, may be used in the direct susceptibility test method. Cultures derived from samples and specimens may be tested by the indirect method.

## PROCEDURE

**Materials Provided:** In this package are **Sensi-Disc** Antimicrobial discs for Use in Culture Media as labeled.

**Materials Required But Not Provided:**

1. Middlebrook and Cohn 7H10 Agar Base, or Seven H11 Agar Base and Middlebrook OADC Enrichment.
2. Water bath, approximately 45 °C.
3. 100 x 15 mm style four compartment plates.
4. Single Place Dispensers or forceps.
5. Bunsen burner, connected to gas outlet.
6. Inoculating loop made from 24 to 26 gauge platinum or nichrome wire fitted into a standard bacteriological wire loop holder. The wire is formed into a "closed" loop two or three mm in diameter with a shaft about 2 1/2 inches long.
7. Incubator, 35 °C or other appropriate temperature.
8. **Cornwall™** pipette, sterile 10 mL, or other type of pipette for delivering approximately 5 mL of medium.
9. Acid-fast stain reagents and equipment.
10. Microscopes, low power (or dissecting) and oil immersion.
11. Sterile purified water or 0.85% saline, 4.5 mL volumes in screw-capped tubes.
12. Pipettes, 1 mL and Pasteur capillary.
13. 16 x 125 mm screw-cap tubes.
14. Glass or plastic beads, 1 to 2 mm in diameter or high heat polystyrene plastic beads, such as \*Styrone™ 700, 27 clear, No.7
15. Middlebrook 7H9 Broth.
16. Test tube mixer.
17. MacFarland No. 1 standard (prepared by adding 0.1 mL of 1% BaCl<sub>2</sub> to 9.9 mL of 1% H<sub>2</sub>SO<sub>4</sub>).
18. For incubation in a CO<sub>2</sub>-enriched aerobic atmosphere, a **GasPak™** System operated with a **GasPak** Disposable Carbon Dioxide Generator Envelope, **GasPak** vented System, carbon dioxide incubator or other suitable system.

\* Trademark of the Dow Chemical Co.

## Test Procedure:

- A. Preparation of Plates
  1. Prepare the required amount of complete sterile medium. It should be fluid and at approximately 50 °C.
  2. With Single Place Dispensers or sterile forceps, place discs centrally into three of the plate quadrants; the fourth quadrant serves as the drug-free control.
  3. Affix the discs to the centers of the quadrants by adding a drop or two of the medium to each. If the medium is added via a needle, as with a **Cornwall** pipette, the discs may be rimmed with medium to hold them in place. The small amount of agar should gel at once.
  4. Add the remainder of 5 mL of the medium to each quadrant. The resulting concentration of antimicrobial agent per mL of medium will then be one-fifth that of the disc content.
  5. Place the plates in the refrigerator and hold overnight to permit diffusion throughout.
  6. Inoculate the plates as described below.
- B. Direct Drug Susceptibility Test<sup>1</sup>

The direct drug susceptibility test is performed on a specimen in which acid-fast bacilli (AFB) can be demonstrated on a smear of the digested, concentrated specimen.

  1. Stain and read smears of concentrated clinical specimens, recording the average number of AFB found per oil immersion field (observe approximately 20 positive fields).

2. Select dilutions of the specimen digest on the basis of positive smear results as follows:

Microscopy	Inoculate
Less than 1 AFB/field	Undilute and 10 <sup>-2</sup> (1:100)
1 – 10 AFB/field	10 <sup>-1</sup> and 10 <sup>-3</sup> (1:10 and 1:1000)
More than 10 AFB/field	10 <sup>-2</sup> and 10 <sup>-4</sup> (1:100 and 1:10,000)

3. Make the indicated dilutions with 0.5 mL of the specimen digest and 4.5 mL of sterile water or saline. Mix and transfer successively 0.5 mL to 4.5 mL for the required dilutions.
4. With a capillary pipette, inoculate 3 drops onto each quadrant of the drug and control media. Two dilutions of inocula, usually one-hundredfold apart, are planted onto duplicate plates of medium.

C. Indirect Test<sup>1</sup>

Whenever possible, perform the Direct Test. The Indirect Test employs as the inoculum mycobacterial growth from a primary culture on a drug free liquid or solid medium and is indicated when:

1. The smear was negative and the primary culture has grown.
2. The growth on the control of the direct test is inadequate for a reliable report.
3. A reference culture is sent in by another laboratory.

For the test to be valid, the inoculum must be large enough (50 – 100 colonies on the control quadrant) for good growth, but not large enough to stimulate drug resistance due to the growth of spontaneously occurring drug resistant mutants.

1. Remove approximately 2 to 5 mg of growth from the drug-free medium. Pick a small portion of each colony if possible.
2. Transfer to a sterile 16 x 125 mm screw cap tube containing 6 to 8 glass or plastic beads and 3 mL of Polysorbate 80-albumin liquid medium (Middlebrook 7H9 Broth).
3. Homogenize the cells in a test tube mixer for 5 to 10 min.
4. Allow the large particles to settle, withdraw the supernate and adjust the density to approximately MacFarland No. 1 standard with sterile water or 0.85% sterile saline.
5. Dilute to 10<sup>-2</sup> and 10<sup>-4</sup> by making serial tenfold dilutions in sterile distilled water or 0.85% saline.
6. Proceed as for "Direct Test" for inoculation.

D. Incubation

Following inoculation, keep plates shielded from light and place plates, medium side down, in a **GasPak™** CO<sub>2</sub> jar or other suitable system providing an aerobic atmosphere enriched with approximately 3 to 10% CO<sub>2</sub> and incubate at 35 to 37 °C.

Where moisture loss may occur during incubation, place plates in CO<sub>2</sub> permeable plastic (polyethylene) bags. After squeezing air out of bags, seal bags with tape, a heat sealing apparatus, or fold over two or three times and staple.

Stack plates in bags not more than 6 plates in a pile; to prevent the accumulation of drops of moisture in the plates during incubation, place plastic plates on an insulating surface, e.g., styrofoam or several thicknesses of cardboard, etc., and not directly on a metal shelf, tray or basket.

Incubate bags of plates in an aerobic atmosphere enriched with approximately 3 to 10% CO<sub>2</sub> at 35 to 37 °C.

Drug susceptibility tests can be examined after 5 to 7 days with a low-power dissecting microscope and a preliminary of positive culture findings given at the time of initial observation.<sup>7</sup> Final reporting should not be made until after 3 weeks of incubation because growth of resistant mutants may be slower on drug containing media than on the control.<sup>1,2</sup>

NOTE: Cultures from skin lesions suspected to be *M. marinum*, *M. ulcerans* or *M. haemophilum* must be incubated at a temperature not exceeding 33°C for primary isolation; cultures suspected to contain *M. avium* or *M. xenopi* exhibit optimum growth at 40 to 42 °C.<sup>1</sup> Incubate a duplicate culture at 35 to 37 °C.

**User Quality Control:** Check drug-free medium for performance by inoculation with pure cultures of stable control strains of mycobacteria producing known, desired reactions. A drug-susceptible strain of *M. tuberculosis* is satisfactory as a control to demonstrate the activity of drugs in the media used.<sup>2</sup>

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

**RESULTS<sup>1</sup> – Direct and Indirect Tests**

A test is considered valid when there are 50 or more colonies on at least one of the controls.

1. Record bacterial growth on the control and drug medium quadrants as indicated below:

Confluent (500 or more colonies)	4+
Almost Confluent (200 – 500 colonies)	3+
100 – 200 colonies	2+
50 - 100 colonies	1+
Less than 50 colonies	Actual colony count

2. Determine the percent of resistant bacilli by the following formula:

$$\frac{\text{Number of colonies on the drug quadrant}}{\text{Number of colonies on the control quadrant}} \times 100 = \% \text{ Resistance}$$

Sample Calculation	Growth		% Resistant
	Undilute	10 <sup>-2</sup>	
Control	4+	100	
INH 0.2 µg/ mL	1+	10	10%
SM 2.0 µg/ mL	0	0	0%
PAS 2.0 µg/mL	0	0	0%

$$\frac{10 \text{ colonies}}{100 \text{ colonies}} \times 100\% = 10\% \text{ resistance}$$

NOTE: The colony size of organisms growing on drug containing media may be markedly smaller than on the control, especially with PAS; therefore, readings are based on numbers of colonies seen and not their size.

#### LIMITATIONS OF THE PROCEDURE

1. Since the observed drug susceptibility is a comparison of the amount of growth on control and drug-containing media, the inoculum for each culture must be of demonstrated uniformity. Uniform suspension of the inoculum and avoidance of large clumps are essential.
2. Since modification of the proportion of resistant bacteria may occur with *in vitro* subculture, direct susceptibility testing using the processed specimen as inoculum is preferred whenever smear examination reveals acid-fast bacilli.
3. If indirect susceptibility tests are made from cultures having discrete colonies, the inoculum must be prepared from a fully representative selection; all colony types seen should be sampled.

#### PERFORMANCE CHARACTERISTICS

##### Ethambutol Discs

Prior to release, all lots of Ethambutol discs (EM-25 or EM-50) are tested for specific product characteristics. Representative samples of the lot are analyzed for the concentration of Ethambutol. Ethambutol discs (EM-25 or EM-50) are placed in centrifuge tubes along with 5 mL of processed water, 2 mL of bromthymol blue and 10 mL of benzene. The tubes are centrifuged for 5 min. The benzene layer is then transferred to a 1 cm cell and the absorbance is read on an ultraviolet spectrophotometer. The absorbency of the test sample is compared to a control solution consisting of 5 mL of processed water, 2 mL of bromthymol blue and 10 mL of benzene. The potency of the disc is then calculated.

##### Isoniazid Discs

Prior to release, all lots of Isoniazid discs (INH-1 or INH-5) are tested for specific product characteristics. Representative samples of the lot are analyzed for the concentration of Isoniazid. INH discs are placed in a micro centrifuge tube with 500 µL of DI water. After 30 min. discs are vortexed and centrifuged. The extract is Qs to approx. 1 mL with 500 µL of DI water. Each sample is tested for potency by HPLC method. The amount of INH present in the disc is calculated based on a standard curve.

##### Ethionamide Discs

Prior to release, all lots of Ethionamide discs (EA-25) are tested for specific product characteristics. Representative samples of the lot are analyzed for the concentration of Ethionamide. EA discs are placed in a volumetric flask and extracted at room temperature with 10 mL of methanol for 18 to 24 h. The absorbency of each extract solvent is determined by ultraviolet spectrophotometry. The reference solution used is methanol. The amount of Ethionamide present in the discs is calculated based on the absorbance of a reference standard solution (within 10%) of the desired concentration.

##### Streptomycin Discs

Prior to release, all lots of Streptomycin discs (S-50) are tested for specific product characteristics. Representative samples of the lot are analyzed for the concentration of Streptomycin by an assay procedure developed by BD Diagnostics based on the protocol for similar antimicrobial agents in the *Code of Federal Regulations*, Part 460- Antibiotic Drugs Intended for Use in Laboratory Diagnosis, Subpart A – Susceptibility Discs.

##### Rifampin Discs

Prior to release, all lots of Rifampin discs (RA-25) are tested for specific product characteristics. Representative samples of the lot are analyzed for the concentration of Rifampin by an assay procedure developed by BD Diagnostics based on the protocol for similar antimicrobial agents in the *Code of Federal Regulations*, Part 460- Antibiotic Drugs Intended for Use in Laboratory Diagnosis, Subpart A – Susceptibility Discs.

##### P-Aminosalicylic Acid Discs

Prior to release, all lots of P-Aminosalicylic Acid discs (PAS-10 or PAS-50) are tested for specific product characteristics. Representative samples of the lot are analyzed for the concentration of P-Aminosalicylic Acid. PAS discs are placed in a volumetric flask and extracted at room temperature with 10 mL of 0.1 N hydrochloric acid (HCl). After 18 to 24 h, the absorbance of each extract is determined by ultraviolet spectrophotometry. The reference solution used is 0.1 N HCl. The amount of P-Aminosalicylic Acid present in the discs is calculated based on the absorbance of a reference standard solution (within 10%) of the desired concentration.

#### AVAILABILITY

Cat. No.	Description
231570	<b>BD BBL™ Sensi-Disc™</b> Streptomycin, 50 µg, Single Cartridge of 50 discs €€
231571	<b>BD BBL™ Sensi-Disc™</b> Isoniazid (Isonicotinyl hydrazine), 1 µg, Single Cartridge of 50 discs €€
231572	<b>BD BBL™ Sensi-Disc™</b> Isoniazid (Isonicotinyl hydrazine), 5 µg, Single Cartridge of 50 discs €€
231573	<b>BD BBL™ Sensi-Disc™</b> P-Aminosalicylic Acid, 10 µg, Single Cartridge of 50 discs €€
231574	<b>BD BBL™ Sensi-Disc™</b> P-Aminosalicylic Acid, 50 µg, Single Cartridge of 50 discs €€
231575	<b>BD BBL™ Sensi-Disc™</b> Ethambutol, 25 µg, Single Cartridge of 50 discs €€
231576	<b>BD BBL™ Sensi-Disc™</b> Ethambutol, 50 µg, Single Cartridge of 50 discs €€
231577	<b>BD BBL™ Sensi-Disc™</b> Ethionamide, 25 µg, Single Cartridge of 50 discs €€
231578	<b>BD BBL™ Sensi-Disc™</b> Rifampin, 25 µg, Single Cartridge of 50 discs €€

## REFERENCES

1. Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. Centers for Disease Control, Atlanta.
2. Hawkins, J.E., R.J. Wallace, Jr., and B.A. Brown. 1991. In A. Balows, W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.), Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
3. Wayne, L.G., and I. Krasnow. Am. J. Clin. Pathol. 45:769.
4. Griffith, M., M.L. Barrett, H.L. Bodily, and R.M. Wood. 1967. Am. J. Clin. Pathol. 47:812.
5. Griffith, M.E., M.L. Matajack, M.L. Bissett, and R.M. Wood. 1971. Am. Rev. Respir. Dis. 103:423.
6. U.S. Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. 1999. Biosafety in microbiology and biomedical laboratories, 4th ed. HHS Publication (CDC) 93-8395. U.S. Government Printing Office, Washington, D.C.
7. Runyon, E.H. 1973. Am. J. Clin. Pathol. 59:817.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Benex Limited