
Revisions

SO 0191-5

Rev from	Rev to	ECO #
1005	2010/07	5384-10

Notes:

1. BD Cat. Number 231248, 231249
2. Blank (Sheet) Size: Length: N/A Width: N/A
 Number of Pages: N/A Number of Sheets: N/A
 Page Size: Length N/A Width N/A Final Folded Size: N/A
3. Style (see illustrations below): # N/A

4. See Specification Control Number VS-88-4066-1 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: N/A PMS# N/A
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Vendor Spec Held By BD Caribe, LTD.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8840661		Category and Description Package Insert, BBL Taxo Discs for Detecting Lactose Fermenters	Sheet: 1 of 13 <hr/> Scale: N/A	A

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INTENDED USE

Taxo™ ONPG Discs are used for the detection of lactose fermenters, especially those that do not promptly ferment lactose in some routine identification media such as Triple Sugar Iron Agar (TSI Agar) or Kligler Iron Agar.

SUMMARY AND EXPLANATION

The fermentation of lactose is a widely used biochemical characteristic in the differentiation of microorganisms, including *Enterobacteriaceae*.

Rapid fermentation is dependent on the presence of a lactose permease to actively transport lactose across the cell membrane and on the presence in the cell of the enzyme, β -galactosidase, that hydrolyzes the lactose. True lactose nonfermenters are deficient in the latter enzyme.

However, some lactose fermenters that lack permease appear to be nonfermenters because diffusion of lactose into the cell is too slow to be detected in carbohydrate fermentation media.¹

The use of a chromogenic substrate hydrolyzed by β -galactosidase but not dependent on the presence of permease permits the rapid determination of lactose fermentation by a color reaction.

PRINCIPLES OF THE PROCEDURE

β -galactosidase is an inducible enzyme possessed by lactose-fermenting organisms that hydrolyzes lactose to form galactose and dextrose.² The enzyme also hydrolyzes the substrate, ortho-nitrophenyl- β -D-galactosidase (ONPG), which is not dependent on lactose permease to enter the cell. Hydrolysis of ONPG produces a yellow end product, orthonitrophenol.

The test may be performed directly on a dense bacterial suspension in saline or on a culture incubated in a peptone broth, e.g., **Trypticase™** Soy Broth or Tryptophan (**Trypticase**) Broth, 1%. Following addition of a Taxo ONPG Disc to the suspension or culture, the tube is incubated at 35°C for 4 to 6 h and observed for the presence or absence of a yellow color.

REAGENTS

Taxo ONPG Discs are 6-mm discs made from high-quality absorbent paper impregnated with approximately 100 μ g of o-nitrophenyl- β -D-galactopyranoside per disc.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, tubes and other contaminated materials must be sterilized by autoclaving before discarding.

Storage Instructions: On receipt, store at -20 to +8°C. After use, store vial at 2 to 8°C to protect product integrity. The expiration date applies to unopened containers stored as directed. Do not open until ready to use.

Use oldest discs first and discard expired discs. Allow vial to come to room temperature before opening. Return the vial of unused discs to the refrigerator. Discard the vial from which discs have been frequently removed during one week. Discard any vial which has been left at room temperature overnight.

SPECIMENS

These discs are not for use directly with clinical specimens or other sources containing mixed flora. The organism to be presumptively identified must first be isolated as separate colonies by streaking the specimen onto appropriate culture media and then handling as described under "Test Procedure."

PROCEDURE

Material Provided: Taxo ONPG Discs.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure: The test is generally applied to organisms after the use of TSI Agar or Kligler Iron Agar slants and is particularly useful for testing cultures showing no fermentation, or only dextrose fermentation, in these two commonly used media.

1. Prepare tubes containing 0.5 mL of **Trypticase Soy Broth** or **Tryptophan (Trypticase) Broth**, 1% solution. Alternatively, prepare tubes containing 0.2 mL of Normal Saline (0.85% sodium chloride).
2. Inoculate tube of broth or saline with a heavy inoculum of the organism to be tested taken from a TSI Agar or Kligler Iron Agar slant or from the surface of a fresh subculture plate.
3. Aseptically add an ONPG disc to each inoculated tube and to one uninoculated control tube.
4. Incubate tubes at 35°C.
5. Examine and record reactions (presence or absence of a yellow color) after 4 to 6 h of incubation. Continue to incubate tubes exhibiting a negative reaction for a total of 24 h before discarding as negative.

User Quality Control: At the time of use, check performance with pure cultures of stable control organisms producing known desired reactions. The following test strains are recommended:

TEST STRAIN	EXPECTED RESULTS
<i>Escherichia coli</i> ATCC™ 25922	+ (yellow color)
<i>Morganella morganii</i> ATCC 8019	- (no color change)

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

The appearance of a yellow color indicates a positive reaction for lactose fermentation.

LIMITATIONS OF THE PROCEDURE^{1,2}

A heavy inoculum is required in order to obtain a high concentration of enzyme and speed the reaction.

The test cannot be used with yellow-pigmented bacteria unless the suspension is first centrifuged before interpreting results.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of Taxo ONPG Discs are tested for performance characteristics. Representative samples of the lot are placed into tubes of Tryptophan Broth, 1% and Normal Saline for each test organism. The tubes are then inoculated with cultures of *Escherichia coli* (ATCC 25922), *Morganella morganii* (ATCC 8019) and *Shigella sonnei* (ATCC 9290) using a heavy inoculum. After inoculation, the tubes are incubated at 35 ± 2°C. Tubes are read for the formation of a yellow color after 4 to 6 h incubation. Tubes showing a negative reaction are reincubated for an additional 24 to 48 h before discarding as negative. *E. coli* and *S. sonnei* give a positive reaction with the ONPG discs as indicated by the presence of a yellow color in the tubes. *M. morganii* gives a negative reaction as indicated by the absence of a yellow color in the tubes.

AVAILABILITY

Cat. No.	Description
231248	Taxo™ ONPG Discs, Single vial of 50 discs.
231249	Taxo™ ONPG Discs, Pkg. of 6 vials.

REFERENCES

1. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1995. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby – Year Book, Inc., St. Louis.
2. MacFaddin, J.F. 2000. Biochemical tests for identification of medical bacteria, 3rd ed. Lippincott Williams & Wilkins, Baltimore.

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EC/REP Benex Limited