# BD BBL™ Taxo™ Nystatin Discs for Inhibition of Yeasts

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## INTENDED USE

BD BBL Taxo Nystatin Discs are for use in a qualitative procedure for the differentiation and isolation of bacteria from old specimens and those likely to contain a mixed or large population of bacteria and yeasts.<sup>1</sup> They are not for susceptibility testing.

## SUMMARY AND EXPLANATION

Nystatin (fungicidin) was isolated from *Streptomyces noursei* by Hazen and Brown and reported in 1951. They found that a crude preparation protected mice from lethal doses of *Candida albicans*. The drug was later refined and used for clinical fungal infections.<sup>2</sup>

Nystatin is an antimicrobial drug which affects cell membranes and is useful for the suppression of yeasts, especially *Candida*, and some bacteria.

The drug is impregnated into paper discs and applied to agar surfaces. The drug diffuses into the medium inhibiting susceptible organisms.

## PRINCIPLES OF THE PROCEDURE

The discs have been found to be particularly valuable in detection of *Neisseria* in specimens. For *Neisseria*, the BD BBL Taxo Nystatin Discs may be placed on the area of initial (heaviest) inoculation of BBL Thayer-Martin media or Transgrow Medium to supplement the nystatin already present in these media, which is insufficient to inhibit some strains of yeasts.<sup>3</sup> Nystatin Discs are also recommended for use in or on other selective and nonselective media for examination of miscellaneous specimens.

Plates are observed for bacterial growth within zones where nystatin discs have inhibited the growth of yeasts or other undesirable microorganisms.

## REAGENTS

BD BBL Taxo Nystatin Discs are made from high quality absorbent paper impregnated with about 100 units of nystatin\* per disc.

\*Mycostatin™: Bristol-Myers-Squibb.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

Directions for use should be read and followed carefully.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures.

**Storage Instructions:** On receipt, store at -20 to +8 °C. After use, store cartridge to protect product integrity at 2 to 8 °C. Apply control tests at time of use.

Use the oldest discs first and discard expired discs. Also discard cartridges from which discs have been frequently removed during a week or so. Discard containers left out overnight in the laboratory, or else test the discs for performance.

The expiration date applies to produce in intact container stored as directed. Do not open until ready to use.

## PROCEDURE

**Material Provided:** BD BBL Taxo Nystatin Discs

### Materials Required But Not Provided:

1. Suitable agar medium, such as BD Trypticase™ Soy Agar with 5% Sheep Blood, Modified Thayer-Martin (MTM) Agar or Transgrow Medium.
2. Sterile forceps or Single Disc Dispenser.

### Test Procedure:

The recommended method is to streak or swab the specimen on an agar medium suitable for isolation of the desired bacteria. Place one, or preferably two, BD BBL Taxo Nystatin Discs on the area of initial (heaviest) inoculation. When two discs are used, place them approximately 15 mm apart.

**User Quality Control:** Control organisms producing known positive and negative results should be included at time of test. Cultures of *Candida albicans* and *Neisseria gonorrhoeae* are recommended.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

Observe for bacterial growth within zones where the nystatin discs have inhibited growth of yeasts or other undesirable microorganisms.

## PERFORMANCE CHARACTERISTICS

Prior to release, all lots of BD BBL Taxo Nystatin Discs are tested to specific product specifications. Samples of the lot are assayed for nystatin content using *Candida (pseudotropicalis) kefyr* (ATCC® 8553) as the test organism. The assay method is based on the protocol for antimicrobial agents in the United States Code of Federal Regulations, Title 21, Part 460, Subpart A (1998). The acceptable average potency of the BD BBL Taxo Nystatin Discs is 80% to 140% of the theoretical disc potency.

## AVAILABILITY

Cat. No. Description

231550 BD BBL Taxo Nystatin Discs, Cartridge of 50 discs

## REFERENCES

1. Vera, H.D. 1971. Quality control in diagnostic microbiology. Health Lab. Sci. 8:175-189.
2. Hazen, E.L., and R. Brown. 1951. Fungicidin, an antibiotic produced by a soil actinomycete. Proc. Soc. Exp. Biol. Med. 76:93.
3. Thayer, J.D., and J.E. Martin, Jr. 1966. Improved medium selective for cultivation of *N. gonorrhoeae* and *N. meningitidis*. Public Health Rep. 81:559-562.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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ATCC is a trademark of the American Type Culture Collection.
© 2017 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.
Benex Limited