
## Revisions

SO 0191-5

Rev from	Rev to	ECO #
2010-01	2011/10	5954-11

**Notes:**

1. BD Cat. Number: 220093, 220097, 220099, 220105, 200109, 220115, 220116, 220129, 220130, 220132, 220133, 220134, 220135, 220146, 220147, 220148, 220149
2. Blank (Sheet) Size: Length: 23 5/8"    Width: 3 15/32"  
 Number of Pages: 2    Number of Sheets: 1  
 Page Size: Length: N/A    Width: N/A    Final Folded Size: N/A
3. Style (see illustrations below): N/A

4. Material specification controlled by Copan.
5. Ink Colors: Printed two sides     Yes     No  
 No. of Colors: 1    PMS # 2755 Blue
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

\*Peel-off insert / Label. Insert perforated 3/8" from bottom with 3/16" notch for easy peel-off. Top insert size 4 3/16" x 3 15/32" includes 3/8" bottom portion of glued area. Insert is folded under with on horizontal fold and followed by two parallel folds to form the size of 3 1/4" x 3 15/32". Insert to be glued on top 5/16" and bottom 3/8" to roll label, size 3 15/32" (web) x 4 3/8".

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	<b>BD</b> Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L000021/H083		Category and Description Peel-Off Insert / Label, CultureSwab	Sheet: 1 of 3 Scale: N/A	A

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste saadetakse kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojmo informacij za pomoč pri uporabi BD izdelkov. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrucció ziskate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свяжитесь с местным представителем на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompaniji BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

**BBL™ CultureSwab™** devices are sterile ready-to-use systems intended for the collection, transport and preservation of clinical specimens for bacteriological examination.

## SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of bacterial infections involves the collection and safe transportation of a clinical specimen from the patient to the laboratory. This can be accomplished using the **BBL CultureSwab** collection and transport device. Each **BBL CultureSwab** unit is comprised of a sterile peel pouch containing a rayon- or polyester-tipped swab applicator used to collect the sample and a tube containing transport medium into which the swab applicator is placed after sampling.

The **BBL CultureSwab** transport media (Amies Liquid Medium, Liquid Stuart Medium and Cary-Blair Transport Medium) are non-nutritious, buffered with phosphate and provide a reduced environment due to their formulation with sodium thioglycollate.<sup>1</sup> Organisms in the sample material are protected from drying by moisture in the transport medium. The medium is designed to maintain the viability of organisms during transit to the laboratory. **BBL CultureSwab** pouches are made of a plastic film which retards the penetration of atmospheric air into the product.

**BBL CultureSwab** media are contained in a constricted (venturi) tube. Nitrogen gas is flushed into the transport tube during the media filling and capping process. During final packaging of the swab and tube, air is removed from the pouch by vacuum and nitrogen gas is flushed inside.

## REAGENTS

Approximate Formula Per Liter Distilled Water

Amies Liquid Medium	Stuart Liquid Medium	Cary-Blair Agar Medium
Calcium Chloride 0.10 g	Calcium Chloride 0.10 g	Bacteriological Agar 5.60 g
Disodium Phosphate 1.15 g	Mercaptoacetic Acid 1.0 mL	Calcium Chloride 0.09 g
Magnesium Chloride 0.10 g	Sodium Glycerophosphate 10.00 g	Disodium Hydrogen Phosphate 1.10 g
Monopotassium Phosphate 0.20 g		Sodium Chloride 5.00 g
Potassium Chloride 0.20 g		Sodium Thioglycollate 1.15 g
Sodium Chloride 3.00 g		
Sodium Thioglycollate 1.00 g		

**Precautions:** For *in vitro* Diagnostic Use.

It must be assumed that all specimens contain infectious microorganisms; therefore, all specimens should be handled with appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste.

Ⓢ **BBL CultureSwab** is for single use only; reuse may cause a risk of infection and/or inaccurate results.

**Storage:** Store **BBL CultureSwab** devices at 5–25°C.

**Product Deterioration:** Contents sterile if unopened or not damaged. Do not use if they show evidence of damage, dehydration or contamination. Do not use if past expiration date.

## SPECIMEN COLLECTION AND HANDLING

The **BBL CultureSwab** devices are available with different applicator shafts which facilitate the collection of specimens from various sites of the patient. For specific recommendations about collection of specimens for microbiological analysis and primary isolation techniques, consult appropriate references.<sup>2-4</sup>

Once a swab sample is collected, it should be placed in the tube of medium, transported to the laboratory as soon as possible and cultured onto appropriate primary isolation media.

## PROCEDURE

**Materials Provided:** Fifty (50) units of sterile **BBL CultureSwab** devices, which contain media, are contained in each Vi-Pak Pouch.

One hundred (100) units of sterile **BBL CultureSwab** swab and transport tubes are contained in each Vi-Pak Pouch.

**Materials Required But Not Provided:** Appropriate materials for isolating, differentiating and culturing aerobic and anaerobic bacteria. These materials include culture media plates or tubes and incubation systems, gas jars or anaerobic workstations.

## Directions For Use:

The directions for use are printed on each **BBL CultureSwab** unit, along with descriptive diagrams, and are summarized as follows:

1. Peel open the **BBL CultureSwab** pouch.
2. Remove cap from transport tube.
3. Remove applicator swab and collect specimen.

During specimen collection, the applicator tip should only touch the area where the infection is suspected to minimize potential contamination.

4. Place applicator swab in transport tube.
5. Record patient's name and information on tube label.
6. Send specimen to the laboratory for immediate analysis.

## EXPECTED RESULTS

The survival of bacteria in a transport medium depends on many factors. These include the types of bacteria, duration of transport, storage temperature, concentration of bacteria in the sample and formulation of the transport medium. **BBL CultureSwab** devices maintain viability of many microorganisms for 24–48 h. For fastidious bacteria, such as *Neisseria gonorrhoeae* and *Streptococcus pneumoniae*, swab specimens should be plated directly onto culture medium or transported immediately to the laboratory and cultured within 24 h.

## LIMITATIONS OF THE PROCEDURE

**BBL CultureSwab** Amies Liquid Medium and Stuart Liquid Medium are intended for the collection and transport of bacteriological samples only. Preferred samples for anaerobic investigations are: tissue samples obtained during surgical procedures, biopsies from tissue or bone, fluid, pus, or aspirates collected using a syringe. For detailed information and recommendations for transporting fluid and tissue specimens for anaerobic culture, refer to specific publications.<sup>3-7</sup> Samples containing viruses or chlamydia should be collected and transported using alternative specific transport systems.

## PERFORMANCE CHARACTERISTICS

Recovery studies were performed using **BBL CultureSwab** products with a variety of aerobic organisms. Swabs were dosed with inoculum and inserted into the transport tube containing transport media. The tubes were stored at room temperature prior to subculturing onto appropriate culture media. Listed are the organisms tested with the various transport media.

Organism	Amies Liquid form	Transport Media Stuart Liquid form	Cary-Blair Agar
<i>Neisseria gonorrhoeae</i> ATCC™ 43069	*	*	*
<i>Streptococcus pyogenes</i> ATCC 19615	*	*	*
<i>Haemophilus influenzae</i> ATCC 19418	*	*	*
<i>Campylobacter jejuni</i> ATCC 33291			*
<i>Shigella flexneri</i> ATCC 12022			*
<i>Yersinia enterocolitica</i> ATCC 9610			*

## REFERENCES

1. Amies CR. 1967. A modified formula for the preparation of Stuart's transport medium. Can. J. Public Health. 58: 296-300.
2. Isenberg H.D., F.D. Schoenkench, and A. von Graevenitz. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Rubin. American Society for Microbiology, Washington, D.C.
3. Balows, A., W.J. Hausler, Jr, K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
4. Isenberg, H.D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
5. Zavala, M.K., D.M. Citron, E.J.C. Goldstein. 1998. Evaluation of a novel specimen transport system for anaerobic bacteria. Clin. Infect. Dis. 25 (supplement 2): S132-133.
6. Perry, J.L. 1997. Assessment of swab transport systems for aerobic and anaerobic organism recovery. J. Clin. Microbiol. 35: 1269-1271.
7. Summanen, P., E.J. Baron, D.M. Citron, C.A. Strong, H.M. Wexler, and S.M. Finegold. 1993. Dismouth anaerobic bacteriology manual, 5th ed. Star Publishing Co., Belmont, Calif.

**TECHNICAL INFORMATION:** In the United States, telephone Technical Services, toll free (800) 638-8663.

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo, BBL and CultureSwab are trademarks of Becton, Dickinson and Company. © 2011 BD.

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