
Revisions


Rev from	Rev to	ECO #
0905	2010/06	5368-10

Notes:


- BD Cat. Number 231723
- Blank (Sheet) Size : Length: 13.75" Width: 16.125"
 Number of Pages: 6 Number of Sheets: 1
 Page Size: Length 13.75" Width 5.375" Final Folded Size: 1.25" x 5.375"
- Style (see illustrations below): # 4


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number VSL-000149-OPR for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Standard Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VS Controlledby BD Caribe, LTD

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L000149		Category and Description Package Insert, Taxo Differentiation Discs Hippurate	Sheet: 1 of 7 <hr/> Scale: N/A	A

BD BBL™ Taxo™ Differentiation Discs Hippurate

English: pages 1 Italiano: pagine 4
Français: pages 2 Español: páginas 5
Deutsch: Seiten 3

CE L000149
2010/06

Pokyny vám poskytnú miestni zástupcov spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukites vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerininize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкции свяжитесь с местным представителем компании BD. / Өзіндік жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

BBL Taxo Differentiation Discs Hippurate are used for detecting the hydrolysis of sodium hippurate by beta-hemolytic group B streptococci, as well as by other organisms.

SUMMARY AND EXPLANATION

The hippurate test is a qualitative procedure for determining the ability of bacteria to enzymatically hydrolyze sodium hippurate.¹ Hwang and Ederer reported on a rapid method using ninhydrin reagent in which the end product, glycine, was detected rather than benzoic acid.² Since the introduction of this test, many authors have reported on its usefulness. Greenwood and Pickett reported that 90% of *Haemophilus (Gardnerella) vaginalis* tested by this procedure yield positive results.³ Harvey found *Campylobacter fetus* subsp. *jejuni* to be hippurate positive and able to be differentiated from *C. fetus* subsp. *intestinalis*, which does not hydrolyze this substrate.⁴

In 1981, Hebert reported that the rapid hippurate test may be useful for differentiating *Legionella pneumophila* from other *Legionella* species.⁵

PRINCIPLES OF THE PROCEDURE

Sodium hippurate, when hydrolyzed by the bacterial enzyme hippuricase, yields benzoic acid and glycine. The rapid hippurate test employing ninhydrin detects the presence of glycine. Ninhydrin and glycine react, removing the amino group and producing water and ammonia. The ammonia passes through sequential steps to produce a purple color.¹

REAGENTS

BBL Taxo Differentiation Discs Hippurate are paper discs, containing sufficient sodium hippurate to yield a positive reaction with organisms producing sufficient hippuricase to hydrolyze the substrate.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow established laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store BBL Taxo Differentiation Discs Hippurate at 2 – 8°C.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Collect specimens in sterile containers or with sterile swabs and transport immediately to the laboratory in accordance with recommended guidelines.⁶⁻⁸

Process each specimen using procedures appropriate for that sample.⁶⁻⁸

This product is recommended for use only with pure cultures.

PROCEDURE

Material Provided: BBL Taxo Differentiation Discs Hippurate.

Materials Required But Not Provided: 35 – 37°C incubator or water bath, inoculating loop, Bunsen burner or incinerator, quality control organisms, test tubes 13 x 100 mm and ninhydrin reagent.

Test Procedure:

Rapid Hippurate Test

1. Dispense 0.4 mL sterile distilled or deionized water into a small tube (13 x 100 mm) and add one BBL Taxo Hippurate Disc.
2. Completely emulsify a large loopful of test organism in the water and shake or vortex to form a smooth suspension.
3. Incubate the tube in a water bath or incubator at 35 – 37°C for 2 h.
4. Following incubation, add 0.2 mL (5 drops) ninhydrin reagent to each tube and shake gently.
5. Reincubate the tube for 10 – 15 min and read reaction.
6. Observe and record results.

User Quality Control:

Identity Specification – Round, off white (ivory), 1/4" paper disks with HIP printed on both sides, no ragged edges.

Cultural Response – Hydrolysis of Sodium Hippurate is checked with the following cultures.

Organism	ATCC™	Reaction
<i>Streptococcus pyogenes</i>	19615	– (no color change)
<i>Streptococcus agalactiae</i>	13813	+ (purple - blue color)

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

A positive hippurate reaction is indicated by formation of a purple-blue color in the solution.

A negative hippurate reaction is indicated by the absence of a color change in the solution.

LIMITATIONS OF THE PROCEDURE

The rapid hippurate test may be used as a presumptive diagnostic aid for identifying those organisms that hydrolyze hippurate; i.e., *Streptococcus agalactiae*. However, additional biochemical and serological tests should be performed for complete identification.

Some strains of group D streptococci may give a weak positive result. Differentiation of group B from group D streptococci may be made by using BBL™ Bile Esculin Agar.⁸ Only group D streptococci hydrolyze esculin yielding a black discoloration of the medium.

Incubation of the test culture for more than 30 min after the ninhydrin reagent is added can yield false-positive results.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of BBL Taxo Differentiation Discs Hippurate are tested for performance characteristics. Representative samples of the lot are placed in 0.4 mL of sterile, purified water. A large loopful of *Streptococcus pyogenes* group A (ATCC 19615), *S. agalactiae* group B (ATCC 13813) or *Campylobacter fetus* subsp. *jejuni* (ATCC 29428) is added to each tube and shaken gently. The tubes are then incubated at 35 – 37°C. After 2 h incubation, 0.2 mL of ninhydrin reagent is added to each tube and the tubes are then re-incubated at 35 – 37°C for an additional 10 to 15 min. *S. agalactiae* and *C. jejuni* give a positive reaction (purple-blue color change in the tubes). *S. pyogenes* gives a negative reaction (no color change in the tube).

AVAILABILITY

Cat. No. Description

231723 BBL Taxo™ Differentiation Discs Hippurate, 50 discs.

231724 BBL Taxo™ Differentiation Discs Hippurate, 10 x 50 discs.

261201 BBL™ Ninhydrin Reagent Droppers, Ctn. of 50.

REFERENCES

1. MacFaddin, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams and Wilkins, Baltimore.
2. Hwang, M., and G.H. Ederer. 1975. Rapid hippurate hydrolysis method for presumptive identification of group B streptococci. J. Clin. Microbiol. 1: 114-115.
3. Greenwood, J.R., and M.J. Pickett. 1979. Salient features of *Haemophilus vaginalis*. J. Clin. Microbiol. 9: 200-204.
4. Harvey, S.M. 1980. Hippurate hydrolysis by *Campylobacter fetus*. J. Clin. Microbiol. 11: 435-437.
5. Hebert, B.A. 1981. Hippurate hydrolysis by *Legionella pneumophila*. J. Clin. Microbiol. 13: 240-242.
6. Ruoff, K.L. 1995. *Streptococcus*, p. 299-307. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D. C.
7. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc. St. Louis.
8. Pezzlo, M. (ed.). 1994. Aerobic bacteriology, p. 1.0.1.-1.20.47. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D. C.

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Benex Limited