
Revisions

Rev from	Rev to	ECO #
0505	2010/06	5368-10

Notes:

1. BD Cat. Number Varies
2. Blank (Sheet) Size : Length: 13.75" Width: 26.875"
 Number of Pages: 10 Number of Sheets: 1
 Page Size: Length 13.75" Width 5.375" Final Folded Size: 1.25" x 5.375"
3. Style (see illustrations below): # 4

4. See Specification Control Number N/A for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VS Controlled by Caribe, LTD

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L000156		Category and Description Package Insert, BBL Taxo Haemophilus	Sheet: 1 of 11 Scale: N/A	A

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INTENDED USE

BBL™ Taxo™ Differentiation Discs are used for presumptively differentiating *Haemophilus* spp. based on their requirements for X or V factors or both.

SUMMARY AND EXPLANATION

Haemophilus spp. are part of the normal respiratory flora of humans and many animal species. *Haemophilus influenzae* may be an opportunistic secondary invader, usually following a viral infection. This organism can cause a variety of diseases from chronic respiratory infections to meningitis. *Haemophilus parainfluenzae* is rarely implicated in pathological processes. In culture, these two organisms are morphologically and colonially similar, making differentiation difficult. Their differing requirements for growth factors are used for differentiating these species.

Davis¹ and Thjotta and Avery² recognized that the type species *H. influenzae* required two growth factors, one from blood that they called X factor and another that they called V factor. The X factor has been identified as iron protoporphyrin. Hemin is used as a source of this factor. Lwoff and Lwoff³ prepared V factor from baker's yeast. V factor is nicotinic adenine dinucleotide (NAD). These exacting growth requirements aid in differentiation and species identification.

Chocolate agar, horse blood agar and rabbit blood agar all contain both growth factors and support the growth of most *Haemophilus* spp. BBL Hemoglobin provides the X factor, while Supplement B or C or VX will provide the V factor. BBL Fildes Enrichment, employed in culture media in a 5% concentration, is also a source of both X and V factors. The *Staphylococcus* satellite method may also be used for isolating *Haemophilus*. *S. aureus* provides the V factor, and when streaked on sheep blood agar, will enhance the growth of *Haemophilus* around the *Staphylococcus* streak. This method is also used for species identification and differentiation.

Parker and Heoprich⁴ prepared discs containing the X and V factors and used them to differentiate *H. influenzae*, which requires X and V factors, from *H. parainfluenzae*, which requires only V factor for growth. Cultures were inoculated on Heart Infusion Agar, Tryptic Soy Agar or Nutrient Agar. Russell⁵ used discs impregnated with X and V factors to identify organisms belonging to the genus *Haemophilus* based on their requirements for these factors and reported the suitability of this procedure.

PRINCIPLES OF THE PROCEDURE

Haemophilus spp. vary in their requirements for the growth factors V and X. This variation is useful in isolating and differentiating these species. The presence or absence of growth around and/or between discs impregnated with factors V, X and VX is evaluated to determine the species of the organism.

REAGENTS

BBL Taxo Differentiation Discs V, X and VX are paper discs impregnated with the substances listed and imprinted with the appropriate code.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow established laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store BBL Taxo Differentiation Discs V, X and VX at 2 – 8°C.

Disc Code	Disc Content	Disc Size
V	NAD	1/4"
X	Hemin	1/4"
VX	NAD and hemin	1/4"

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

- Specimens for the isolation of *Haemophilus* spp. may be obtained from the throat, nasopharynx, conjunctiva, blood, or spinal fluid.
- Primary specimens should be prepared for differentiation by inoculation onto chocolate agar. Alternatively, blood agar may be used with a staphylococcal cross streak.
- Cultures obtained from an enriched medium should be diluted 1:100 in Tryptic Soy Broth to prevent carryover of growth factors from the medium.
- Uniformly distribute organisms over the entire surface of Heart Infusion Agar or Tryptic Soy Agar. Only pure cultures should be used for differentiation of *Haemophilus* spp.

PROCEDURE

Material Provided: BBL Taxo Differentiation Discs V, X or VX.

Materials Required But Not Provided: Forceps, CO₂ incubator (35 ± 2°C), swabs, inoculating loop, Bunsen burner or incinerator, quality control organisms, Tryptic Soy Agar or Heart Infusion Agar and Tryptic Soy Broth.

Test Procedure:

- Using aseptic technique, place the desired discs on the inoculated agar according to one of the following methods. Press gently.

Two disc method: Place BBL Taxo Differentiation Discs V and X approximately 3 – 5 mm apart (Fig. 1).

Three disc method: Place BBL Taxo Differentiation Discs V, X and VX in the form of an equilateral triangle with at least 30 – 35 mm between discs (Fig. 2).

- Incubate at 35 ± 2°C in a 5 – 10% CO₂ atmosphere for 18 – 24 h.
- Examine the pattern of growth around and/or between the discs.

User Quality Control:


Identity Specifications –

BBL Taxo Differentiation Discs V – Round, white 1/4" paper discs with "V" printed on both sides.

BBL Taxo Differentiation Discs X – Round, brown 1/4" paper discs with "X" printed on both sides.

BBL Taxo Differentiation Discs VX – Round, light tan 1/4" paper discs with "VX" printed on both sides.

Cultural Response – V, X, VX: Inoculate BBL Heart Infusion Agar plates with a dilution of test organisms approximating a McFarland standard 0.5. Place BBL Taxo Differentiation Discs V, BBL Taxo Differentiation Discs X, and BBL Taxo Differentiation Discs VX on the surface of the plates as described under "Test Procedure." Incubate plates under CO₂ at 35 ± 2°C for 18 – 24 h.

Organism	ATCC™	Growth Stimulation		
		V	X	VX
<i>Haemophilus influenzae</i> biogroup <i>aegyptius</i>	11116	–	–	+
<i>Haemophilus parahaemolyticus</i>	10014	+	–	+
<i>Haemophilus parainfluenzae</i>	7901	+	–	+

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Two Disc Method: Growth around only one disc indicates a requirement for only that factor. Growth between the two discs indicates that both V and X factors are required.

Three Disc Method: Growth around the V and VX discs or the X and VX discs indicates a requirement for the single growth factor V or X, respectively. Growth around only the VX disc indicates a requirement for both factors.

Interpret the species of *Haemophilus* according to the growth factor requirements listed below.

Results other than those listed in the table may suggest one or both of the following:

- Improper dilution of cultures taken from enriched media. An inoculum that is too dilute may yield no growth. One that has been under-diluted may contain sufficient growth factors to affect growth over the entire plate or in areas around one or more inappropriate discs on the plate.
- Mixed cultures may result in growth around the discs.

Only presumptive identification of the *Haemophilus* species is possible using **BBL Taxo** Differentiation Discs. Determination of hemolysis and characteristics other than hemolysis may be required for confirmation of a presumptive identification. Biochemical and serological testing may also be used for species identification. Consult appropriate references for further identification of *Haemophilus* spp.^{6,7}

<i>Haemophilus</i> Species	Growth Requirements	
	X	V
<i>Haemophilus aphrophilus</i>	-	-
<i>Haemophilus ducreyi</i>	+	-
<i>Haemophilus haemolyticus</i>	+	+
<i>Haemophilus influenzae</i>	+	+
<i>Haemophilus influenzae</i> biogroup <i>aegyptius</i>	+	+
<i>Haemophilus parahaemolyticus</i>	-	+
<i>Haemophilus parainfluenzae</i>	-	+
<i>Haemophilus paraphrophilus</i>	-	+
<i>Haemophilus segnis</i>	-	+

LIMITATIONS OF THE PROCEDURE

Differentiation of *Haemophilus* species using the above procedures is presumptive in that growth factor requirements are identified. However, similarities in these requirements exist between species. Additional biochemical and serological tests should be performed for complete identification.

PERFORMANCE CHARACTERISTICS

In a study by Parker and Hoepflich, 255 nasopharyngeal and 1,317 throat swabs were examined for the presences of *Haemophilus*-like species using X and V factor disks. From the nasopharyngeal swabs, 59 *H. influenzae* were identified (growth only between the X and V disks) and 12 *H. parainfluenzae* were identified (growth only around the V disk). From the 1,317 throat swabs, 129 *H. influenzae* were identified and 302 were identified as *H. parainfluenzae*.⁴

AVAILABILITY

Cat. No.	Description	Cat. No.	Description
231727	BBL™ Taxo™ Differentiation Discs V, 50 discs/cartridge.	231730	BBL™ Taxo™ Differentiation Discs X, 10 x 50 discs.
231728	BBL™ Taxo™ Differentiation Discs V, 10 x 50 discs.	231731	BBL™ Taxo™ Differentiation Discs VX, 50 discs/cartridge.
231729	BBL™ Taxo™ Differentiation Discs X, 50 discs/cartridge.	231732	BBL™ Taxo™ Differentiation Discs VX, 10 x 50 discs.

REFERENCES

- Davis, J.D. 1917. Food accessory factors (vitamins) in bacterial culture with especial reference to hemophilic bacilli. *J. Infect. Dis.* 21: 392-403.
- Thjotta, T., and O.T. Avery. 1921. Studies on bacterial nutrition. II. Growth accessory substances in the cultivation of hemophilic bacilli. *J. Exper. Med.* 34: 97-114.
- Lwoff, A., and M. Lwoff. 1937. Studies on codehydrogenases. I. Nature of growth factor "V." *Proc. Roy. Soc. London. Ser. B.* 722: 352-359.
- Parker, R.H., and P.D. Hoepflich. 1962. Disk method for rapid identification of *Haemophilus* species. *Am. J. Clin. Pathol.* 37: 319-327.
- Russell, J.P. 1965. The identification of gram-negative pleomorphic bacilli. *Am. J. Clin. Pathol.* 44: 88-93.
- Campos, J.M. 1995. *Haemophilus*, p. 556-566. In P.R. Murray, E.J. Baron, M. A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. *Bailey & Scott's diagnostic microbiology*, 10th ed. Mosby, Inc., St. Louis.

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Benex Limited