L-000158
Revised: May 1999

BBL™ Taxo™ Differentiation Discs Spectinomycin

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INTENDED USE

BBL™ Taxo™ Differentiation Discs Spectinomycin are used to detect spectinomycin-resistant *Neisseria gonorrhoeae* (PPNG).

SUMMARY AND EXPLANATION

Spectinomycin is an aminocyclitol antibiotic isolated from *Streptomyces spectabilis*.¹ It acts similarly to the aminoglycosides by binding to the 30S ribosomal subunits and inhibiting protein synthesis. The drug is recommended for treatment of infections due to *N. gonorrhoeae*, including β -lactamase producing strains.²⁻⁵ Resistant gonococci are rare.^{6,7} The drug has also been used in patients with allergies to penicillin.

Biddle⁸ recommended a 100 μ g spectinomycin disc for testing suspect gonococci. It has been recommended that all PPNG isolates of *N. gonorrhoeae* be screened for spectinomycin resistance.^{8,9}

PRINCIPLES OF THE PROCEDURE

Spectinomycin discs are placed onto inoculated chocolate agar plates and incubated under increased carbon dioxide. Spectinomycin resistance is based on the diameter of the zone of inhibition around the disc.

REAGENTS

BBL Taxo Differentiation Discs Spectinomycin are paper discs that contain 100 μ g of spectinomycin HCl.

Precautions: *in vitro* Diagnostic

Follow established laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store **BBL Taxo** Differentiation Discs Spectinomycin at 2 – 8°C.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Collect specimens or samples in sterile containers or with sterile swabs and transport immediately to the laboratory in accordance with recommended guidelines.^{1,4,10}

Process each specimen, using procedures appropriate for that sample.^{1,4,10}

PROCEDURE

Material Provided: **BBL Taxo** Differentiation Discs Spectinomycin.

Materials Required But Not Provided: Sterile saline, Chocolate Agar, Bunsen burner or incinerator, inoculating loop, swab, McFarland 0.5 turbidity standard, sterile forceps, CO₂ gas generating system, incubator (35 \pm 2°C), ruler (mm) or calipers and quality control organisms.

Test Procedure:

1. Obtain a pure culture of the organism to be tested. Select well-isolated colonies.
2. Using a sterile swab or inoculating loop, transfer enough growth of *Neisseria gonorrhoeae* to be tested to a tube of sterile saline to form a suspension comparable to a McFarland 0.5 turbidity standard (approx. 10⁸ CFU/ml). Agitate the suspension thoroughly.
3. Dip a swab into the suspension and express excess fluid by rotating the swab against the inside wall of the test tube.
4. Inoculate the entire surface of a fresh chocolate agar plate, streaking in three directions by rotating the plate 60° after each streaking, creating a confluent lawn of growth.
5. Allow the inoculum to dry 5 – 15 min with the lid in place.
6. Aseptically place one **BBL Taxo** Spectinomycin Disc on the inoculated surface. Using sterile forceps, gently press the disc down to contact the inoculum.
7. Invert and incubate the plates in a CO₂ incubator or other system providing an aerobic atmosphere supplemented with CO₂ at 35 \pm 2°C for 18 – 24 h.
8. Examine and measure any zone of inhibition of growth around the spectinomycin disc. Record results.

User Quality Control:

Identity Specifications – Round, white 1/4" paper discs. The discs are imprinted with PP over NG.

Cultural Response – To test for presumptive resistance to spectinomycin, inoculate enriched Chocolate Agar using only well characterized strains for control organisms and **BBL Taxo** Differentiation Discs Spectinomycin. Incubate under increased CO₂ at 35 \pm 2°C for 18 – 24 h. A zone of inhibition of <18 mm is presumptive evidence that the test culture is resistant to spectinomycin.

Zone Diameter	Interpretation
Less than 18 mm	Resistant
18 mm or greater	Sensitive

RESULTS

A zone of inhibition of <18 mm indicates presumptive resistance to spectinomycin of the *N. gonorrhoeae* culture tested.

LIMITATIONS OF THE PROCEDURE

BBL Taxo Differentiation Discs Spectinomycin are intended for screening of spectinomycin-resistance only. They are not intended for use in determining the susceptibility of *N. gonorrhoeae* to spectinomycin.

The procedure described here applies only when used with pure cultures.

Some strains of *N. gonorrhoeae* are extremely slow growers and are likely to produce large zones. Therefore, it is important to note these strains when performing this test to correctly assess results.¹¹

AVAILABILITY

Cat. No. Description

231722 **BBL™ Taxo™** Differentiation Discs Spectinomycin, 50 discs.

REFERENCES

1. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
2. Centers for Disease Control. 1983. Penicillin resistant gonorrhoea - North Carolina. Morbid. Mortal. Weekly Rep. 32: 273-275.
3. Centers for Disease Control. 1983. Penicillinase-producing *Neisseria gonorrhoeae* - Los Angeles. Morbid. Mortal. Weekly Rep. 32: 181-183.
4. Isenberg, H.D. (ed.). 1994. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
5. National Committee for Clinical Laboratory Standards. 1997. Approved Standard: M2-A6. Performance standards for antimicrobial disk susceptibility tests, 6th ed. National Committee for Clinical Laboratory Standards, Wayne, Pa.
6. Fekete, T. 1993. Antimicrobial susceptibility testing of *Neisseria gonorrhoeae* and implications for epidemiology and therapy. Clin. Microbiol. Rev. 6: 22-33.
7. Zenilman, J.M., L.J. Nims, M.A. Menegus, F. Nolte, and J.S. Knapp. 1987. Spectinomycin-resistant gonococcal infections in the United States, 1985-1986. J. Infect. Dis. 156:1002-1004.
8. Centers for Disease Control. 1981. Spectinomycin-resistant penicillinase-producing *Neisseria gonorrhoeae* - California. Morbid. Mortal. Weekly Rep. 30: 221-222.

9. World Health Organization. 1981. Surveillance of *Neisseria gonorrhoeae* producing beta-lactamase (PPNG); spectinomycin resistance. *Weekly Epidemiol. Rec.* 56: 158.
10. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. *Bailey & Scott's diagnostic microbiology*, 9th ed. Mosby-Year Book, Inc., St. Louis.
11. Biddle, J.W., J.M. Swenson, and C. Thornsberry. 1978. Disc agar diffusion antimicrobial susceptibility tests with beta-lactamase producing *Neisseria gonorrhoeae*. *J. Antibiot.* 31: 352-358.

TECHNICAL INFORMATION: In the United States, telephone Technical Services, toll free (800) 638-8663.

©1999 Becton Dickinson and Company.
BBL and Taxo are trademarks of Becton Dickinson and Company.
ATCC is a trademark of the American Type Culture Collection.

**BECTON
DICKINSON**

Becton Dickinson Microbiology Systems
Becton Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA

Becton Dickinson France S.A.
11 rue Aristide Bergès
38800 Le Pont de Claix, France

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