BD BBL™ Taxo™ Differentiation Discs Novobiocin

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INTENDED USE

BBL™ Taxo™ Differentiation Discs Novobiocin are recommended for the differentiation of coagulase-negative staphylococci (e.g., *Staphylococcus saprophyticus*) based on novobiocin resistance.

SUMMARY AND EXPLANATION

Staphylococcus saprophyticus has recently been recognized as a significant cause of urinary tract infections, especially in young women.¹⁻³ In their scheme for identification of human staphylococci, Kloos and Schleifer⁴ presented resistance to novobiocin as presumptive identification for this organism. They recommended use of a 5 µg novobiocin disc on specially formulated P agar. Later investigators revised the method, however, and Mueller Hinton medium is now the accepted medium for this test.^{5,6}

PRINCIPLES OF THE PROCEDURE

Resistance to novobiocin is the parameter which differentiates the clinically significant coagulase-negative staphylococci from the less significant species. The method is based on the antibiotic disc diffusion test of Bauer and Kirby.⁷

REAGENTS

BBL Taxo Differentiation Discs Novobiocin are paper discs containing 5 µg novobiocin.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow established laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store BBL Taxo Differentiation Discs Novobiocin at 2 – 8°C. Allow container with discs to reach room temperature before use.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use the product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

This product is to be used with pure cultures.

PROCEDURE

Material Provided: BBL Taxo Differentiation Discs Novobiocin.

Materials Required But Not Provided: Sterile cotton swabs, sterile saline or broth, Mueller Hinton medium plates, sterile forceps, McFarland 0.5 turbidity standard, inoculating loop, ruler or caliper, incubator (35 ± 2°C) and quality control organisms.

Test Procedure:

- Using a sterile swab or inoculating loop, transfer growth from a pure, coagulase-negative *Staphylococcus* culture to a tube of sterile saline or broth. Adjust the turbidity to form a suspension comparable to a McFarland 0.5 standard (approximately 10⁸ CFU/mL). Agitate this suspension thoroughly.
- Dip a swab into the suspension and express excess fluid by rotating the swab against the inside wall of the test tube.
- Inoculate the entire surface of a Mueller Hinton medium plate, streaking in three directions by rotating the plate 60° after each streaking. If the inoculum is satisfactory, there will be a confluent lawn of growth.
- Allow the inoculum to dry approximately 5 min with the lid in place.
- Place one BBL Taxo Differentiation Discs Novobiocin on the inoculated surface.
- Using flamed and cooled forceps, gently press the disc down to ensure complete contact with the agar. Do not move a disc once it has touched the agar because the novobiocin diffuses almost immediately.
- Invert and incubate at 35 ± 2°C for 18 – 24 h.
- Measure the zone of inhibition, if present, with a ruler or caliper.

User Quality Control

Identity Specifications - Round 6 mm paper discs. The discs are imprinted with C on both sides.

Cultural Response - To verify product performance, follow the above procedure using the quality control strains specified.

Organism	ATCC™	Zone of Inhibition
<i>Staphylococcus saprophyticus</i>	15305	≤ 16 mm (Resistant)
<i>Staphylococcus epidermidis</i>	12228	> 16 mm (Susceptible)

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

A zone of inhibition of ≤ 16 mm is presumptive evidence that the test culture is *S. saprophyticus*.

LIMITATIONS OF THE PROCEDURE

BBL Taxo Differentiation Discs Novobiocin (5 µg) are intended for screening of novobiocin resistance only. They are not intended for determining the susceptibility of coagulase-negative staphylococci to novobiocin.

The procedure described here applies only when used with pure cultures.

This test is used for the presumptive identification of *Staphylococcus saprophyticus*. Additional biochemical testing may be needed for complete identification. Other less significant coagulase-negative *Staphylococcus* strains are also novobiocin resistant, such as *S. cohnii*, *S. xylosum* and *S. sciuri*.⁷

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of BBL Taxo Differentiation Discs Novobiocin are tested for specific product specifications. Plates of Mueller Hinton Agar are swab-inoculated with cultures of *Staphylococcus aureus* (ATCC 25923), *S. epidermidis* (ATCC 12228) and *S. saprophyticus* (ATCC 15305). Samples of the lot are placed on the inoculated plates and incubated at 35 ± 2°C. After 18 to 24 h incubation, the plates are examined for zones of inhibition around the discs. *S. aureus* and *S. epidermidis* produce zones of inhibition greater than 16 mm. *S. saprophyticus* is resistant to novobiocin and produces a zone of inhibition less than or equal to 16 mm.

AVAILABILITY

Cat. No. Description
231750 BBL Taxo Differentiation Discs Novobiocin, 50 discs.

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EC REP Benex Limited