
Revisions

Rev from	Rev to	ECO #
2010/06	2011/06	5829-11

Notes:

1. BD Cat. Number 231741, 231735
2. Blank (Sheet) Size : Length: 9" Width: 40"
 Number of Pages: 20 Number of Sheets: 1
 Page Size: Length 9" Width 4" Final Folded Size: 2.25" x 4"
3. Style (see illustrations below): # 4

4. See Specification Control Number L000171 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p> <p style="font-weight: bold; font-size: large;">BD</p> <p>Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA</p>	
Proofer	Date		
Checked By	Date		
Part Number: L000171		Category and Description Package Insert, Taxo TB Niacin	Sheet: 1 of 21 Scale: 1:1

A

BD BBL™ Taxo™ TB Niacin Test Reagents

English: pages 1 – 3 Italiano: pagine 9 – 12

Français : pages 4 – 6 Português: páginas 12 – 14

Deutsch: Seiten 6 – 9 Español: páginas 15 – 17

L000171
2011/06

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзініздің жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

BBL™ Taxo™ TB Niacin Test Strips are paper strips containing reagents for the detection of niacin production by mycobacteria.

BBL Taxo TB Niacin Test Control is used as a positive control in the niacin test.

SUMMARY AND EXPLANATION

Konno¹ first devised the standard niacin test, which was modified by Runyon et al.² Niacin produced by the organism reacts with cyanogen bromide and a primary or secondary amine (usually aniline).

Kilburn and Kubica³ further modified the niacin test, using a single strip impregnated with the reagents. Acidified potassium thiocyanate mixed with chloramine T releases cyanogen chloride, which reacts in the presence of niacin to produce a yellow color. If niacin is not present, no color appears.

BBL Taxo Niacin Test Strips are patterned with some modifications after the Kilburn and Kubica Niacin Test Strip. Kilburn and Kubica³ and Morse et al.⁴ found that the strip test was as sensitive and specific as the classical test.

PRINCIPLES OF THE PROCEDURE

All mycobacteria produce nicotinic acid, or niacin, which plays a vital role in oxidation-reduction reactions. Because of a blocked metabolic pathway, *M. tuberculosis* and some isolates of *M. simiae* and *M. chelonae* produce the largest amounts. The niacin accumulates in the culture medium on which the organisms are growing and can be extracted from the medium.

REAGENTS

BBL Taxo TB Niacin Test Strips are absorbent paper strips impregnated with potassium thiocyanate, chloramine T, citric acid and sodium aminosalicylate.

BBL Taxo TB Niacin Test Control paper discs are impregnated with nicotinamide; the discs, when used according to directions, will yield a yellow solution equivalent to approximately 5 µg niacin.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Only a 3 – 4 weeks old culture containing at least 50 – 100 colonies should be used in the test for niacin production. Fewer colonies may yield questionable or negative results.

The usual precautions for handling *M. tuberculosis* must be followed in performing this test.

The stoppered tubes must be autoclaved after completion of the test and before discarding.

Storage Instructions: Store **BBL Taxo TB Niacin Test Reagents** at 2 – 8°C. The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Collect specimens in sterile containers or with sterile swabs and transport immediately to the laboratory in accordance with recommended guidelines.⁵⁻⁷

Process each specimen as appropriate for that specimen.⁵⁻⁷

An extract of the *Mycobacterium* culture is employed in the test and is prepared using the following methods. NOTE: While cultures grown on egg-based media are reported to yield the most consistent results, satisfactory results are achieved with Middlebrook and Cohn Agar supplemented with 0.1% aspartic acid.⁸

Culture on Egg Media:

1. Use a 3 – 4 weeks old culture having at least 50 – 100 colonies on Lowenstein-Jensen, Petragrani or other coagulated egg medium.
2. To prepare an extract of the *Mycobacterium* culture, add 1.5 mL of sterile distilled or deionized water or saline to the culture. Using a pipette, gently scrape off surface growth. Stab a 1 mL sterile pipette through the growth into the medium to permit extraction of the niacin. Do not use contaminated cultures that show “blued or bleached” slants.
3. Slant the tube so that the surface of the medium is in the horizontal position and is covered with the liquid. Allow to remain in this position for 20 – 30 min.
4. Carefully remove approximately 0.6 mL of the extract with a sterile capillary pipette fitted with a bulb and transfer to the bottom of a 13 x 75 mm test tube. Stopper the tube.

Cultures on BBL Middlebrook and Cohn 7H10 Agar Slants with Aspartic Acid:

1. Use a 3 – 4 weeks old culture having at least 50 – 100 colonies.
2. To prepare an extract of the *Mycobacterium* culture, add 1.5 mL of sterile distilled or deionized water or saline to the culture. Using a pipette, gently scrape off surface growth. Stab a 1 mL sterile pipette through the growth into the medium to permit extraction of the niacin.
3. Allow the fluid to remain in contact with the culture medium for 20 – 30 min. Place in a slanted position so that fluid layers over the colonies.
4. Remove 0.6 mL of the liquid extract from the culture with a sterile capillary pipette and place in a labeled, clean 13 x 75 mm test tube. Stopper the tube.

PROCEDURE

Materials Provided: BBL Taxo TB Niacin Test Strips and/or BBL Taxo TB Niacin Test Control.

Materials Required But Not Provided: Sterile distilled or deionized water or saline (0.85% NaCl), sterile capillary pipette, sterile 1 mL pipette, 13 x 75 mm test tubes, test tube stoppers, forceps, culture media and quality control organisms.

Test Procedure:

A positive control using BBL Taxo TB Niacin Test Control should be run with each series of tests and is prepared as follows:

1. Place 0.6 mL of sterile distilled or deionized water or saline into a 13 x 75 mm test tube.
2. Add one BBL Taxo TB Niacin Control disc, stopper and shake gently three times during 15 min. Allow tube to remain at room temperature.

A negative control is prepared by adding 0.6 mL of the same diluent to a 13 x 75 mm test tube in place of the extract.

1. Process the culture being tested for presence of niacin as stated under “Specimen Collection and Preparation” and place 0.6 mL extract in a 13 x 75 mm test tube as directed.
2. Have on hand a corresponding negative control in a 13 x 75 mm test tube and a positive control prepared as stated above.
3. Using a pair of flamed forceps, drop a BBL Taxo TB Niacin Test Strip with arrow downward into each tube (positive control, negative control and test culture) and stopper immediately.
4. Shake tubes gently to mix the fluid with the reagent on the bottom of the strip, but do not tilt. Repeat this gentle shaking after 5 – 10 min.
5. After 12 – 15 min, but no more than 30 min, compare the color of the extracts.
6. Autoclave the above stoppered tubes after completion of the test.

User Quality Control:

Identity Specifications –

TB Niacin Test Strips - White paper strips which may have slight coloration due to the placement of solutions on the strips. Paper strips approximately 2 1/2" x 1/4" with arrows printed in black.

TB Niacin Test Control – white, round 1/4" paper discs with “N” over “C” printed on both sides.

Cultural Response – Perform the test following the above procedure with the positive and negative control cultures listed. The TB Niacin Test Control should give the appropriate positive reaction when tested with a TB Niacin Test Strip.

Organism	ATCC™	Niacin Production
<i>Mycobacterium tuberculosis</i> H37Ra	25177	+
<i>Mycobacterium intracellulare</i>	13950	-

- + yellow color in the extract
- no color in the extract

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user

refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

A positive test for niacin is indicated by the appearance of a yellow color in the extract of the test culture and positive control, and no color in the negative control.

LIMITATIONS OF THE PROCEDURE

Although a strong positive niacin test on non-chromogenic mycobacteria is strongly indicative of *M. tuberculosis*, it is considered only preliminary identification of this organism. Complete and final identification of *M. tuberculosis* and other clinically significant mycobacteria is based upon results obtained from a battery of information, which includes biochemical tests, drug susceptibility patterns and morphological characteristics. For more information on these tests refer to appropriate texts.^{5-7,9}

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of **BBL Taxo** TB Niacin Test Strips are tested for specific product specifications. Dilutions of 100, 20, 10, 8, 6, 4 and 2 µg/mL of nicotinic acid (niacin) in deionized water are prepared in tubes. To each dilution a TB Test Strip is added with the arrow pointing down. The tubes are gently shaken and allowed to sit for 5 – 10 min after which they are gently shaken again. The color of the extract is compared to a reference control after 12 – 15 min and again after 20 – 25 min.

Additionally, extracts of the following cultures are prepared using sterile saline: *Mycobacterium tuberculosis* H37Ra (ATCC 25177), *M. tuberculosis* (ATCC 27294 and 25618), *M. avium* (ATCC 25291) and *M. intracellulare* (ATCC 13950). Samples of the lot of TB Niacin Test Strips are placed in the culture extract with the arrow on the strip pointed downward. The culture extract tubes are gently shaken and allowed to sit for 5 – 10 min after which they are gently shaken again. The color of the extract is observed after 12 – 15 min and again after 20 – 25 min. A positive (yellow color in the liquid) is observed in all the culture extracts of *M. tuberculosis*. No color change in the liquid (a negative reaction) is observed with *M. intracellulare* or *M. avium*.

All lots of TB Niacin Test Control discs are also tested for specific product specifications. Samples of the lot are placed in tubes containing processed water. The tubes are gently shaken and allowed to sit for 15 min to allow the niacin to extract from the discs. A TB Niacin Test Strip is added to each tube and the tubes are capped. The tubes are gently shaken 3 times in the next 15 min.

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Benex Limited