# Mycobacteria Growth Indicator Tube 7 mL With BD BACTEC™ MGIT™ 960 Supplement Kit

L000180JAA(04)  
2016-11

English: pages 1 – 4      Italiano: pagine 13 – 16  
Français : pages 5 – 8      Português: páginas 17 – 21  
Deutsch: Seiten 9 – 12      Español: páginas 21 – 24

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokupy vám poskytné místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliojotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

## INTENDED USE

The **BD BBL™ MGIT™** Mycobacteria Growth Indicator Tube supplemented with **BD BACTEC™ MGIT™** Growth Supplement and **BD BBL™ MGIT™ PANTA™** antibiotic mixture is intended for the detection and recovery of mycobacteria using the **BD BACTEC™ MGIT™ 960** and **BD BACTEC™ MGIT™ 320** Systems. Acceptable specimen types are digested and decontaminated clinical specimens (except urine), and sterile body fluids (except blood).

## SUMMARY AND EXPLANATION

From 1985 to 1992, the number of reported cases of infection with *Mycobacterium tuberculosis* (MTB) increased 18%. Tuberculosis still kills an estimated 3 million persons a year worldwide, making it the leading infectious disease cause of death.<sup>1</sup> Between 1981 and 1987, AIDS case surveillances indicated that 5.5% of the patients with AIDS had disseminated nontuberculous mycobacterial infections; e.g., MAC. By 1990, the increased cases of disseminated nontuberculous mycobacterial infections had resulted in a cumulative incidence of 7.6%.<sup>2</sup> In addition to the resurgence of MTB, multidrug-resistant MTB (MDR-TB) has become an increasing concern. Laboratory delays in the growth, identification and reporting of these MDR-TB cases contributed at least in part to the spread of the disease.<sup>3</sup>

The U.S. Centers for Disease Control and Prevention (CDC) have recommended that every effort must be made for laboratories to use the most rapid methods available for diagnostic mycobacteria testing. These recommendations include the use of both a liquid and a solid medium for mycobacterial culture.<sup>3,4</sup>

The **MGIT** Mycobacteria Growth Indicator Tube contains 7 mL of modified Middlebrook 7H9 Broth base.<sup>5,6</sup> The complete medium, with OADC enrichment and **PANTA** antibiotic mixture, is one of the most commonly used liquid media for the cultivation of mycobacteria.

All types of clinical specimens, pulmonary as well as extrapulmonary (except blood and urine) can be processed for primary isolation in the **MGIT** tube using conventional methods.<sup>4</sup> The processed specimen is inoculated into a **MGIT** tube, placed into the **BD BACTEC™ MGIT™** instrument for continuous monitoring until positive or the end of the testing protocol.

## PRINCIPLES OF THE PROCEDURE

A fluorescent compound is embedded in silicone on the bottom of 16 x 100 mm round bottom tubes. The fluorescent compound is sensitive to the presence of oxygen dissolved in the broth. Initially, the large amount of dissolved oxygen quenches emissions from the compound and little fluorescence can be detected. Later, actively respiring microorganisms consume the oxygen and allow the fluorescence to be detected.

Tubes entered into the **BD BACTEC MGIT** instrument are continuously incubated at 37 °C and monitored every 60 min for increasing fluorescence. Analysis of the fluorescence is used to determine if the tube is instrument positive; i.e., the test sample contains viable organisms. An instrument positive tube contains approximately 10<sup>5</sup> to 10<sup>6</sup> colony forming units per milliliter (CFU/mL). Culture vials which remain negative for a minimum of 42 days (up to 56 days) and which show no visible signs of positivity are removed from the instrument as negatives and sterilized prior to discarding.

The **BD BACTEC MGIT** Growth Supplement is added to each **MGIT** tube to provide substances essential for the rapid growth of mycobacteria. Oleic acid is utilized by tubercle bacteria and plays an important role in the metabolism of mycobacteria. Albumin acts as a protective agent by binding free fatty acids which may be toxic to *Mycobacterium* species, thereby enhancing their recovery. Dextrose is an energy source. Catalase destroys toxic peroxides that may be present in the medium.

Contamination is reduced when supplementing the **BD BBL MGIT** broth base with **BD BACTEC MGIT** Growth Supplement/ **BD BBL MGIT PANTA** antibiotic mixture prior to inoculation with a clinical specimen.

## REAGENTS

The **BD BBL MGIT** Mycobacteria Growth Indicator Tube contains: 110 µL of fluorescent indicator and 7 mL of broth. The indicator contains Tris 4,7-diphenyl-1, 10-phenanthroline ruthenium chloride pentahydrate in a silicone rubber base. The tubes are flushed with 10% CO<sub>2</sub> and capped with polypropylene caps.

Approximate Formula\* Per L of Purified Water:

Modified Middlebrook 7H9 Broth base .....	5.9 g
Casein peptone .....	1.25 g

**BACTEC MGIT** Growth Supplement contains 15 mL Middlebrook OADC enrichment.

Approximate Formula\* Per L of Purified Water:

Bovine albumin.....	50.0 g	Catalase .....	0.03 g
Dextrose.....	20.0 g	Oleic acid.....	0.1 g
Polyoxyethylene stearate (POES).....	1.1 g		

The **BBL MGIT PANTA** vial contains a lyophilized mixture of antimicrobial agents.

Approximate Formula\* Per Vial Lyophilized **PANTA**:

Polymyxin B .....	6,000 units	Trimethoprim .....	600 µg
Amphotericin B .....	600 µg	Azlocillin .....	600 µg
Nalidixic acid .....	2,400 µg		

\*Adjusted and/or supplemented as required to meet performance criteria.

**Storage of Reagents: BD BBL MGIT** Mycobacteria Growth Indicator Tubes – On receipt, store at 2–25 °C. DO NOT FREEZE. Minimize exposure to light. Broth should appear clear and colorless. Do not use if turbid. **MGIT** tubes stored as labeled prior to use may be inoculated up to the expiration date and incubated for up to eight weeks.

**BD BACTEC MGIT** Growth Supplement – On receipt, store in the dark at 2–8 °C. Avoid freezing or overheating. Do not open until ready to use. Minimize exposure to light.

**BD BBL MGIT PANTA** Antibiotic Mixture – On receipt, store lyophilized vials at 2–8 °C. Once reconstituted, the **PANTA** mixture must be stored at 2–8 °C and used within 5 days.

## WARNINGS AND PRECAUTIONS:

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. “Standard Precautions”<sup>7-10</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Working with *Mycobacterium tuberculosis* grown in culture requires Biosafety Level 3 practices, containment equipment and facilities.<sup>4</sup>

Prior to use, each **MGIT** tube should be examined for evidence of contamination or damage. Discard any tubes if they appear unsuitable.

Dropped tubes should be examined carefully. If damage is seen, the tube should be discarded.

In the event of tube breakage: 1) Close the instrument drawers; 2) Turn off the instrument; 3) Vacate the area immediately; 4) Consult your facility/ CDC guidelines. An inoculated leaking or broken vial may produce an aerosol of mycobacteria; appropriate handling should be observed.

Autoclave all inoculated **MGIT** tubes prior to disposal.

## SPECIMEN COLLECTION AND HANDLING

All specimens should be collected and transported as recommended by the CDC, the *Clinical Microbiology Procedures Handbook* or your laboratory procedure manual.<sup>11</sup>

## DIGESTION, DECONTAMINATION AND CONCENTRATION

Specimens from different body sites should be processed for inoculation of **MGIT** tubes as follows:

**SPUTUM:** Specimens should be processed using the NALC-NaOH method as recommended by the CDC’s *Public Health Mycobacteriology: A Guide for the Level III Laboratory*.<sup>4</sup> Alternatively, use the **BD BBL™ MycoPrep™** kit for processing mycobacterial specimens (see “Availability”).

**GASTRIC ASPIRATES:** Specimens should be decontaminated as for sputum. If the volume of the specimen is more than 10 mL, concentrate by centrifugation. Resuspend the sediment in about 5 mL of sterile water and then decontaminate. Add a small amount of NALC powder (50 to 100 mg) if the specimen is thick or mucoid. After decontamination, concentrate again prior to inoculation into **MGIT** tube.

**BODY FLUIDS:** (CSF, synovial fluid, pleural fluid, etc.): Specimens which are collected aseptically and are expected to contain no other bacteria can be inoculated without decontamination. If the specimen volume is more than 10 mL, concentrate by centrifugation at 3,000 x g for 15 min. Pour off supernatant fluid. Inoculate **MGIT** tube with sediment. Specimens that are expected to contain other bacteria must be decontaminated.

**TISSUE:** Tissue specimens should be processed as recommended by the CDC’s *Public Health Mycobacteriology: A Guide for the Level III Laboratory*.<sup>4</sup>

The routine inoculation of solid media is especially important for optimal recovery of mycobacteria from tissue specimens as these specimen types are particularly susceptible to sporadic organism recovery.

**STOOL:** Suspend 1 g of feces in 5 mL of Middlebrook Broth. Agitate the suspension on a vortex mixer for 5 s. Proceed to the NALC-NaOH procedure as recommended by the CDC’s *Public Health Mycobacteriology: A Guide for the Level III Laboratory*.<sup>4</sup>

**NOTE:** For all specimen processing methods, a phosphate buffer solution (pH 6.8) should be used to QS the sample decontaminant mixture to 50 mL prior to centrifugation. Resuspension of pellet must also be done using a fresh preparation of phosphate buffer solution (pH 6.8).

## PROCEDURE

**Materials Provided:** **BD BBL MGIT** Mycobacteria Growth Indicator Tubes and **BD BACTEC MGIT 960** Supplement Kit, containing **BD BACTEC MGIT** Growth Supplement and **BBL MGIT PANTA** Antibiotic Mixture (see “Availability”).

**Materials Required But Not Provided:** **Falcon™** brand 50 mL centrifuge tubes, 4% sodium hydroxide, 2.9% sodium citrate solution, N-acetyl-L-cysteine powder, phosphate buffer pH 6.8, vortex mixer, 37 °C incubator, 1 mL sterile pipettes, sterile transfer pipettes, **BD BBL** Middlebrook and Cohn 7H10 Agar, **BD BBL MycoPrep** Specimen Digestion / Decontamination Kit, **BD BBL** Middlebrook 7H9 Broth (see “Availability”) or other mycobacterial agars or egg-based media. Tissue homogenizer or sterile swab, **BD BBL** Normal Saline (see “Availability”), microscope and materials for staining slides, adjustable 1,000 µL pipetter, corresponding sterile pipette tips, 5% sheep blood agar plates and tuberculocidal disinfectant.

## INOCULATION OF MGIT TUBES:

**BD BBL MGIT 7 mL** Tubes must be used with a **BD BACTEC MGIT** instrument.

1. Reconstitute a lyophilized vial of **BD BBL MGIT PANTA** Antibiotic Mixture with 15 mL of **BD BACTEC MGIT** Growth Supplement.
2. Label the **MGIT** tube with the specimen number.
3. Unscrew the cap and aseptically add 0.8 mL of Growth Supplement/ **BD BBL MGIT PANTA** Antibiotic Mixture. For best results, the addition of Growth Supplement/ **BD BBL MGIT PANTA** Antibiotic Mixture should be made just prior to specimen inoculation.
4. Add 0.5 mL of the concentrated specimen suspension prepared above. Also add a drop (0.1 mL) of specimen to a 7H10 agar plate or other mycobacterial solid agar or egg-based medium.
5. Tightly recap the tube and mix well.
6. Tubes entered into the instrument will be automatically tested for the duration of the recommended 42 day testing protocol.

For specimens in which mycobacteria with different incubation requirements are suspected, a duplicate **MGIT** tube can be set up and incubated at the appropriate temperature; e.g., 30 or 42 °C.<sup>13</sup> Inoculate and incubate at the required temperature. These tubes must be manually read (refer to the **BD BACTEC MGIT** Instrument *User’s Manual*).

For specimens suspected of containing *Mycobacterium haemophilum*, a source of hemin must be introduced into the tube at the time of inoculation and the tube incubated at 30 °C. These tubes must be manually read (refer to the **BD BACTEC MGIT** Instrument *User’s Manual*).

7. Positive tubes, identified by the **BD BACTEC MGIT** instrument should be subcultured and an acid-fast smear prepared (see “Results”).

**All quality control testing, reprocessing, smear preparations, sub-culturing, etc., of presumptive positive tubes must be performed using bio-safety level (BSL) III practices and containment facilities.**

**Processing a Positive MGIT Tube:** NOTE – All steps should be performed in a biological safety cabinet.

1. Remove the **MGIT** tube from the instrument and transport to an area using BSL III practices and containment facilities.
2. Using a sterile transfer pipet, remove an aliquot from the bottom of the tube (approx. 0.1 mL) for stain preparations (AFB and Gram stains).
3. Inspect smear and preparations. Report preliminary results only after acid-fast smear evaluation.

At the end of six weeks incubation, perform a visual check of all instrument negative tubes. If the tube appears visually positive (i.e., non-homogenous turbidity, small grains or clumps) it should be subcultured, acid-fast stained and treated as a presumptive positive, provided the acid-fast smear result is positive. If the tube shows no signs of positivity, it should be sterilized prior to discarding.

**Reprocessing Contaminated MGIT tubes:** Contaminated **MGIT** tubes may be re-decontaminated and re-concentrated using the procedure in Appendix E - Supplemental Procedures of the **BD BACTEC MGIT Instrument User's Manual**.

**User Quality Control:** Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Quality Control Certificates are provided on the BD website. Quality Control Certificates list test organisms, including ATCC® cultures specified in the CLSI Approved Standard M22-A3, *Quality Control for Commercially Prepared Microbiological Culture Media*.<sup>12</sup>

NOTE: Middlebrook 7H9 Broth (supplemented) is exempt from User QC testing according to CLSI M22-A3.<sup>12</sup>

RESULTS

An instrument-positive sample is determined by the **BD BACTEC MGIT** instrument and confirmed by an acid-fast smear.

REPORTING OF RESULTS

An instrument positive tube must be confirmed by acid-fast smear. A positive AFB smear result indicates the presence of mycobacteria.

**If AFB smear positive, subculture to solid media and report as:** Instrument positive, AFB smear positive, ID pending.

**If microorganisms other than AFB are present report as:** Instrument positive, AFB smear negative. Contaminated.

**If no microorganisms are present:** Reenter the tube into the instrument as an ongoing negative tube within 5 h of removal. Allow tube to complete test protocol. No reportable result.

Perform subculture from the **BD BBL MGIT** tube for identification and drug susceptibility testing.

LIMITATIONS OF THE PROCEDURE

Recovery of mycobacteria in the **MGIT** tube is dependent on the number of organisms present in the specimen, specimen collection methods, patient factors such as presence of symptoms, prior treatment and the methods of processing.

Decontamination with the N-acetyl-L-cysteine Sodium hydroxide (NALC-NaOH) method is recommended. Other decontamination methods have not been tested in conjunction with the **BD BBL MGIT** medium. Digestant/decontaminant solutions may have harmful effects on mycobacteria.

Colony morphology and pigmentation can only be determined on solid media. Mycobacteria may vary in acid-fastness depending on strain, age of culture and other variables. The consistency of microscopic morphology in **BD BBL MGIT** medium has not been established.

An AFB smear-positive **MGIT** tube can be subcultured, to both selective and nonselective mycobacterial media, for isolation to perform identification and susceptibility testing.

**MGIT** tubes which are instrument-positive may contain other non-mycobacterial species. Non-mycobacterial species may overgrow mycobacteria present. Such **MGIT** tubes should be re-decontaminated and re-cultured (refer to the **BD BACTEC MGIT Instrument User's Manual**). Reprocessing is strongly recommended if the original specimen source cannot be easily recollected; e.g. tissue specimen.

**MGIT** tubes which are instrument-positive may contain one or more species of mycobacteria. Faster growing mycobacteria may be detected prior to slower growing mycobacteria; therefore, it is important to subculture positive **MGIT** tubes to ensure proper identification of all mycobacteria present in the sample.

Due to the richness of the **MGIT** broth and to the non-selective nature of the **MGIT** indicator, it is important to follow the stated digestion/decontamination procedure to reduce the possibility of contamination. Adherence to procedural instructions, which includes use of recommended inoculum volume (0.5 mL), is critical for optimum recovery of mycobacteria.

The use of **PANTA** antibiotic mixture, although necessary for all non-sterile specimens, may have inhibitory effects on some mycobacteria.

Seeded culture studies were performed with twenty-four species (ATCC and wild strains) of mycobacteria using inoculum levels ranging from 10<sup>1</sup> to 10<sup>2</sup> CFU/mL. The following species were detected as positive in the **BD BACTEC MGIT 960** System:

<i>M. avium</i> *	<i>M. gordonae</i> *	<i>M. nonchromogenicum</i>	<i>M. terrae</i>
<i>M. abscessus</i>	<i>M. haemophilum</i> †	<i>M. phlei</i>	<i>M. trivale</i>
<i>M. bovis</i>	<i>M. intracellulare</i>	<i>M. simiae</i> *	<i>M. tuberculosis</i> *
<i>M. celatum</i>	<i>M. kansasii</i> *	<i>M. scrofulaceum</i>	<i>M. xenopi</i> *
<i>M. fortuitum</i> *	<i>M. malmoense</i>	<i>M. smegmatis</i>	
<i>M. gastr</i>	<i>M. marinum</i>	<i>M. szulgai</i> *	


\*Species recovered during clinical evaluation of the **BACTEC MGIT 960** System. In addition, *M. mucogenicum* was recovered at one of the clinical sites.

†The *M. haemophilum* was recovered using the addition of a source of hemin to the **MGIT** tube prior to inoculation.

Clinical studies have demonstrated recovery of mycobacteria from respiratory specimens, gastric aspirates, tissue, stool and sterile body fluids except blood; recovery of mycobacteria from other body fluids has not been established for this product.

EXPECTED VALUES

Figure 1 – Frequency distribution of recovery times for clinical trial specimens positive in the **BD BACTEC™ MGIT™ 960** System

PERFORMANCE CHARACTERISTICS

The **BD BACTEC MGIT 960** System was evaluated at six clinical sites including one non-US site, which represented public health laboratories as well as large acute care hospitals in geographically diverse areas. The site population included patients infected with HIV, immunocompromised patients and transplant patients. The **BD BACTEC MGIT 960** System was compared to the **BD BACTEC 460TB** radiometric system and conventional solid growth media for the detection and recovery of mycobacteria from clinical specimens, except blood. A total of 3,330 specimens were tested during the study. A total of 353 specimens were positive which represented 362 isolates recovered during the study. The distribution of positives by specimen type is: respiratory (90%), tissue (7%), body fluids (1%), stool (0.85%) and bone marrow (0.65%). Of the 362 isolates, 289 (80%) were recovered by the **BD BACTEC MGIT 960** System, 271 (75%) were recovered by the **BD BACTEC 460TB** System and 250 (69%) were recovered by conventional solid media. Of the 3,330 specimens tested in the clinical study, 27 (0.8%) **MGIT 960** tubes were determined to be false positive (instrument-positive, smear and/or subculture-negative). Of the 313 **MGIT 960** instrument positive tubes, 27 (8.6%) were determined to be false positive. The false negative rate (instrument-negative, smear and/or subculture-positive) was determined to be 0.5% based on terminal subcultures of 15% of instrument negative vials. The average breakthrough contamination rate for the **BD BACTEC MGIT 960** System was 8.1% with a range of 1.8–14.6%.

Table 1: Detection of Mycobacteria Positive Isolates in Clinical Evaluations

Isolates	Total isolates	Total MGIT 960	MGIT Only	Total BD BACTEC 460TB	BD BACTEC 460TB Only	Total CONV	CONV Only
MTB	132	102	4	119	11	105	3
MAC	172	147	36	123	12	106	3
<i>M. asiaticum</i>	1	0	0	0	0	1	1
<i>M. fortuitum/chelonae</i>	22	18	6	13	1	15	1
<i>M. genavense</i>	1	0	0	1	0	1	0
<i>M. kansasii</i>	5	5	1	4	0	4	0
<i>M. malmoense</i>	1	0	0	1	0	1	0
<i>M. marinum</i>	1	0	0	0	0	1	1
<i>M. mucogenicum</i>	1	1	1	0	0	0	0
<i>M. simiae</i>	1	1	0	1	0	1	0
<i>M. szulgai</i>	2	2	0	2	0	2	0
<i>M. xenopi</i>	2	2	1	1	0	0	0
MOTT	2	1	1	1	1	0	0
<i>Mycobacteria</i> spp.	2	2	1	1	0	1	0
<i>M. goodnae</i>	11	6	3	3	2	6	3
<i>M. nonchromogenicum</i>	6	2	0	1	0	6	4
All MYCO	362	289	54	271	27	250	16

AVAILABILITY

Cat. No.	Description	Cat. No.	Description
245122	<b>BD BBL™ MGIT™</b> Mycobacteria Growth Indicator Tubes, 7 mL, carton of 100 tubes.	240862	<b>BD BBL™ MycoPrep™</b> Specimen Digestion/Decontamination Kit, ten 75 mL bottles of NALC-NaOH solution and 5 packages of phosphate buffer.
245124	<b>BD BACTEC™ MGIT™ 960</b> Supplement Kit, 6 vials, 15 mL, <b>BD BACTEC™ MGIT™</b> Growth Supplement and 6 vials, lyophilized, <b>BD BBL™ MGIT™ PANTA™</b> Antibiotic Mixture. Each Growth Supplement/ <b>BD PANTA™</b> vial sufficient for 15–18 <b>BD MGIT™</b> tubes.	240863	<b>BD BBL™ MycoPrep™</b> Specimen Digestion/Decontamination Kit, ten 150 mL bottles of NALC-NaOH solution and 10 packages of phosphate buffer.
220908	<b>BD BBL™</b> Lowenstein-Jensen Medium Slants, package of 10 (20 x 148 mm tubes with cap).	221174	<b>BD BBL™</b> Middlebrook and Cohn 7H10 Agar, package of 20.
220909	<b>BD BBL™</b> Lowenstein-Jensen Medium Slants, carton of 100 (20 x 148 mm tubes with cap).	221819	<b>BD BBL™</b> Normal Saline, 5 mL, carton of 100.

REFERENCES

1. Bloom, B.R., and C.J.L. Murray. 1992. Tuberculosis: commentary on a reemergent killer. *Science* 257:1055-1064.
2. Horsburg, C.R., Jr., 1991. *Mycobacterium avium* complex infection in the acquired immunodeficiency syndrome. *N. Engl. J. Med.* 324:1332-1338.
3. Tenover, F.C., et al, 1993. The resurgence of tuberculosis: is your laboratory ready? *J. Clin. Microbiol.* 31:767-770.
4. Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. USDHHS, Centers for Disease Control, Atlanta.
5. Cohn, M.L., R.F. Waggoner and J.K. McClatchy. 1968. The 7H11 medium for the cultivation of mycobacteria. *Am. Rev. Respir. Dis.* 98:295-296.
6. Youmans, G.P. 1979. Cultivation of mycobacteria, the morphology and metabolism of mycobacteria, p. 25-35. *Tuberculosis*. W.B. Saunders Co., Philadelphia.
7. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
9. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
11. Isenberg, Henry D. (ed.) 1992. *Clinical microbiology procedures handbook*. vol. 1. American Society for Microbiology, Washington, D.C.
12. Clinical and Laboratory Standards Institute. 2004. Approved Standard M22-A3. Quality control for commercially prepared microbiological culture media, 3rd ed., CLSI, Wayne, Pa.
13. Lindeboom, J. A., et al. 2011. Clinical Manifestations, Diagnosis, and Treatment of *Mycobacterium haemophilum* Infections. *Clinical Microbiology Reviews* 24, 701-717

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or [www.bd.com](http://www.bd.com).

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For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температура / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najniža dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurilimiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite minimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sicaklık alt sınırı / Мінімальна температура / 温度下限

Control / Контроль / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Controllo / Контроль / kontroll / Контроль / 对照

Positive control / Положительный контроль / Pozitívny kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positivne kontroll / Contrôle positif / Pozitívna kontrol / Pozitív kontrol / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Pozitívna kontrol / Positieve controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивный контроль / 阳性对照试剂

Negative control / Отрицательный контроль / Negatívny kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negativne kontroll / Contrôle négatif / Negativna kontrola / Negatív kontroll / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiama kontrolė / Negatívna kontrol / Negatieve controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативный контроль / 阴性对照试剂

Method of sterilization: ethylene oxide / Метод на стерилизация: этиленов оксид / Způsob sterilizace: etylenoxid / Steriliseringsmetode: etylenoxid / Sterilisationsmethode: Etylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Steriliseringsmetode: etüleenoksiid / Méthode de stérilisation: oxyde d'éthylène / Metoda sterilizacije: etilen oksid / Sterilizálás módszere: etilén-oxid / Metodo di sterilizzazione: ossido di etilene / Стерилизация әдісі: этилен тотығы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizācijas metode: etilēnokšids / Gesteriliseerd met behulp van ethyleenoxide: etylenoksid / Metoda sterylizacji: tlenek etylu / Método de esterilização: óxido de etileno / Metodă de sterilizare: oxid de etilenă / Метод стерилизации: этиленоксид / Metóda sterilizácie: etylénoksid / Metoda sterilizacije: etilen oksid / Steriliseringsmetod: etenoxid / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: этиленоксидом / 灭菌方法: 环氧乙烷

Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringsmetode: bestrålning / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseringsmetode: kiirgus / Méthode de stérilisation: irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Стерилизация әдісі: сәуле тудыру / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizācijas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringsmetode: bestrålning / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Метод стерилизации: облучение / Metóda sterilizácie: ožiarenie / Metoda sterilizacije: ozračavanje / Steriliseringsmetod: strålning / Sterilizasyon yöntemi: ırdıyasyon / Метод стерилизації: опромінювання / 灭菌方法: 辐射

Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Biooogilised riskid / Risques biologiques / Biološki rizik / Biológiallag veszélyes / Rischio biologico / Биологиялык төүөкөндөр / 생물학적 위험 / Biologinis pavojus / Biologiskie riski / Biologisch risico / Biologisk risiko / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Биологическая опасность / Biologické riziko / Biološki rizici / Biologisk risk / Biyolojik Riskler / Биологічна небезпека / 生物学风险

Caution, consult accompanying documents / Внимание, направте справка на придружаващите документи / Pozor! Prostudujte si příloženou dokumentaci / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα / Precaución, consultar la documentación adjunta / Ettevaatus! Lugeka kaasnevat dokumentatsiooni / Attention, consulter les documents joints / Uprozejenje, koristí prateću dokumentaciju / Figyelem! Olvassa el a mellékelt tájékoztatót / Attenzione: consultare la documentazione allegata / Абайлаңыз, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dămesio, ziurêkite pridedamus dokumentus / Piesardzība, skatīt pavaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsiktig, se vedlagt dokumentasjon / Należy zapoznać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Výstraha, pozri sprievodné dokumenty / Pažnja! Pogledajte priložena dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlikte verilen belgeleri başvurun / Увага: див. супутню документацію / 小心, 请参阅附带文档。

Upper limit of temperature / Горен лимит на температура / Horní hranice teploty / Øvre temperaturgrænse / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite superior de temperatura / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőmérsékleti határ / Limite superiore di temperatura / Температураның рұқсат етілген жоғарғы шегі / 상한 온도 / Aukščiausia laikymo temperatūra / Augšējā temperatūras robeža / Hoogste temperatuurilimiet / Øvre temperaturgrænse / Gorna granica temperatury / Limite máximo de temperatura / Limită maximă de temperatură / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Øvre temperaturgräns / Sicaklık üst sınırı / Максимальна температура / 温度上限

Keep dry / Πазете сухо / Skladujte v suchém prostředí / Opbevares tørt / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivas / Conserver au sec / Držati na suhom / Száraz helyen tartandó / Tenere all'asciutto / Құрғақ күйінде ұста / 건조 상태 유지 / Laikykite sausi / Uzglabāt sausu / Droog houden / Houdes tørt / Przechowywać w stanie suchym / Manter seco / A se feri de umezeală / Не допускать попадания влаги / Uchovávať v suchu / Držite na suvom mestu / Förvaras tørt / Kuru bir şekilde muhafaza edin / Беретти від вологи / 请保持干燥

Collection time / Време на събиране / Čas odběru / Opsamlingstidspunkt / Entnahmeuhrzeit / Ώρα συλλογής / Hora de recogida / Kogumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жинау уақыты / 수집 시간 / Paėmimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de colheita / Ora colectării / Время сбора / Doba odboru / Vreme prikupljanja / Uppsamlingstid / Toplama zamanı / Час забору / 采集时间

Peel / Обелете / Otefete zde / Abn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Үстіңгі қабатын алып таста / 벗기기 / Plešti čia / Atīmēt / Schillen / Trekk av / Oderwać / Destacar / Se dezlipește / Отклеить / Odtrhnite / Oljuštiti / Dra isår / Ayırma / Відклеїти / 撕下

Perforation / Перфорация / Perforace / Perforering / Διείσδυση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Тесик тесу / 절취선 / Perforacija / Perforācija / Perforatie / Perforaça / Perforaça / Perforação / Perforare / Перфорация / Perforacia / Perforasyon / Перфорация / 穿孔

Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packung nicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Erep paket бүзылган болса, пайдаланба / 패키지가 손상된 경우 사용 금지 / Jei pakuotė pažeista, nenaudoti / Nelietot, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не използвайте при повреждени упаковки / Nepoužívejte, ak je obal poškozený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包装破损, 请勿使用

Keep away from heat / Пазете от топлина / Nevystavujte přilísnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od izvora topline / Óvja a melegtől / Tenere lontano dal calore / Сапқын жерде сақта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no karstuma / Beschermen tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródeł ciepła / Manter ao abrigo do calor / A se feri de căldură / Не награвяте / Uchovávejte mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Берегти від дії тепла / 请远离热源

Cut / Срежете / Odstřihněte / Klip / Schneiden / Κόψτε / Cortar / Lögga / Découper / Reži / Vágja ki / Tagliare / Keciřiz / 잘라내기 / Kirpti / Ougriez / Knippen / Kutt / Odciać / Cortar / Decupați / Отрезать / Odstřihněte / Iseći / Klipp / Kesme / Pozpizati / 剪下

Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskuupäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тізбекүні / 수집 날짜 / Paëmimo data / Savākšanas datums / Verzameldatum / Dato prövetaking / Data pobrania / Data de colheita / Data colectării / Дата сбора / Dátum odberu / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期

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Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA

Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

**Australian Sponsor:**

Becton Dickinson Pty Ltd.  
4 Research Park Drive  
Macquarie University Research Park  
North Ryde, NSW 2113  
Australia

ATCC is a trademark of the American Type Culture Collection.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.