
Revisions

SO 0191-5

Rev from	Rev to	ECO #
1206	2010/06	5348-10

Notes:

1. BD Cat. Number 261187
2. Blank (Sheet) Size : Length: 11" Width: 27"
 Number of Pages: 12 Number of Sheets: 1
 Page Size: Length 11" Width 4.5" Final Folded Size: 4 1/2" x 1 7/8"
3. Style (see illustrations below): # 4

4. See Specification Control Number L001128 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755 (Blue)
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L001128		Category and Description Package Insert, DMACA Indole Reagent Droppers		Sheet: 1 of 13
				Scale: N/A
				A

BD DMACA Indole Reagent Droppers

English: pages 1 – 2 Italiano: pagine 5 – 7
Français : pages 2 – 4 Português: páginas 7 – 8
Deutsch: Seiten 4 – 5 Español: páginas 9 – 10

L001128
2010/06

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Um Anleitungen zu erhalten, wenden Sie sich bitte an Ihren BD-Kundendienst. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Contattare il rappresentante BD di zona per ottenere il foglietto illustrativo. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіңіздің жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

For the detection of indole production as an aid in the identification of aerobic, anaerobic, or facultatively anaerobic organisms.¹

SUMMARY AND EXPLANATION

BBL™ DMACA Indole Reagent Droppers are used in a rapid spot test for determining the ability of bacteria to produce indole by the deamination of tryptophan. The production of indole has been widely used as an aid in differentiating certain genera and species within the same genus. p-Dimethylaminocinnamaldehyde (DMACA) has been shown to be sensitive to the production of indole.² In addition, the color change with DMACA is both rapid and long lasting, making results easier to interpret than with more traditional reagents.³

Each dropper is good for one day's use after the ampule has been broken.

PRINCIPLES OF THE PROCEDURE

The reaction indicates the presence of the enzyme tryptophanase that reacts with tryptophan to produce indole. The indole produced reacts with DMACA to form a blue-green color.

REAGENTS

DMACA Indole Reagent Droppers contain 0.5 mL of aqueous 1% p-dimethylaminocinnamaldehyde in 10% hydrochloric acid.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper laboratory procedures in handling and disposing of infectious materials. DMACA indole reagent is very acidic. Avoid contact with skin, eyes and mucous membranes. Rinse thoroughly with water if spilled.

Storage Instructions: Store at room temperatures 15 – 30°C (59 – 86°F). Protect from light.

Product Deterioration: This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until expiration date. Change in color of the reagent from light yellow to brown indicates improper storage, which may cause weaker reactions.

PROCEDURE

Material Provided: DMACA Indole Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, filter paper, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure: Rapid Spot Method.

To use DMACA Indole Reagent Droppers, hold upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate dropper any further as the plastic may puncture and injury may occur.** Tap bottom of the dropper on the tabletop a few times. When ready for use, invert for convenient drop-by-drop dispensing of reagent.

1. Moisten a piece of Whatman No. 1 filter paper (or equivalent) with a few drops of indole reagent.
2. Remove a well-isolated colony from an 18 – 24 h culture on a Blood Agar Plate with a sterile loop or wooden stick and smear on moistened filter paper.⁴
3. Observe for a blue to blue-green color within 2 min.

User Quality Control:

Examine the reagent for signs of deterioration (See "Product Deterioration").

Positive (*Escherichia coli* ATCC™ 25922) and negative (*Klebsiella pneumoniae* ATCC 33495) controls should be performed along with the organisms to be tested.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Positive: Blue to blue/green. Negative: No color or pink. Consult appropriate references for indole reactions with specific organisms.^{5,6}

LIMITATIONS OF THE PROCEDURE

Colonies to be tested should be isolated on media containing sufficient tryptophan, and no glucose⁷ (e.g., **Trypticase™ Soy Agar**).

Allow at least 2 min for the absence of a reaction to develop before considering the indole reagent test negative.

Colonies should not be tested from media containing dyes (i.e., MacConkey or EMB agars). Dye carryover may give misleading results.

Colonies from mixed cultures should not be tested. Indole positive colonies have been shown to cause nearby (within 5 mm) indole negative colonies to appear weakly positive.⁸

False negative reactions may occur with some organisms, especially *Clostridium* spp. which form indole but break it down more rapidly.⁷

For *Flavobacterium* spp., a 48 h culture should be tested if 24 h culture was negative.⁴

PERFORMANCE CHARACTERISTICS

An internal study was performed comparing BBL DMACA Indole Reagent and another commercially available DMACA Indole Reagent.⁹ As noted in the table, all reactions were identical for the 11 organisms tested using both reagents.

Organism	Expected Results	BBL Reagent	Commercial Reagent
<i>Enterobacter aerogenes</i> ATCC 19615	–	–	–
<i>Acinetobacter anitratus</i> ATCC 43608	–	–	–
<i>Klebsiella pneumoniae</i> ATCC 33495	–	–	–
<i>Flavobacterium breve</i> BD 8130	W+	W+	W+
<i>Providencia stuartii</i> ATCC 33672	+	+	+
<i>Escherichia coli</i> ATCC 25922	+	+	+
<i>Morganella morganii</i> ATCC 25830	+	+	+
<i>Bacteroides thetaiotamicron</i> ATCC 29741	+	+	+
<i>Weeksella virosa</i> BD 5990	W+	W+	W+
<i>Klebsiella oxytoca</i> BD 9214	+	+	+
<i>Vibrio parahaemolyticus</i> ATCC 17802	W+	W+	W+

Key	
–	Negative
+	Positive
W+	Weak Positive

AVAILABILITY

Cat No.	Description
261187	BBL™ DMACA Indole Reagent Droppers, packaged 50 droppers/carton.

REFERENCES

1. Sutter, V.L., and W.T. Carter. 1972. Evaluation of media and reagents for indole-spot tests in anaerobic bacteriology. *Am. J. Clin. Pathol.* 58: 335-338.
2. Lowrance, B.L., P. Reich, and W.H. Traub. 1969. Evaluation of two spot-indole reagents. *Appl. Microbiol.* 17: 923-924.
3. Miller, J.M., and J.W. Wright. 1982. Spot indole test: evaluation of four reagents. *J. Clin. Microbiol.* 15: 589-592.
4. Pickett, M.J. 1989. Methods for identification of *Flavobacteria*. *J. Clin. Microbiol.* 27: 2309-2315.
5. Baron, E.J., and S.M. Finegold. 1990. Baily and Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Company, St. Louis.
6. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
7. MacFaddin, J.F. 2000. Biochemical tests for identification of medical bacteria. 3rd ed. Lippincott, Williams and Wilkins, Baltimore.
8. Vracko, R. and J.C. Sherris. 1963. Indole-spot test in bacteriology. *Am. J. Clin. Pathol.* 39: 429-432.
9. Data on file, BD Diagnostic Systems.

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Benex Limited