
**Revisions**

Rev from	Rev to	ECO #
1206	2010/06	5348-10

**Notes:**

- BD Cat. Number 261181
- Blank (Sheet) Size : Length: 11"      Width: 22.5"  
 Number of Pages: 10      Number of Sheets: 1  
 Page Size: Length 11"      Width 4.5"      Final Folded Size: 4 1/2" x 1 7/8"
- Style (see illustrations below): # 4

- See Specification Control Number L001133 for Material Information
- Ink Colors: Printed two sides  Yes       No  
 No. of Colors: 1      PMS# 2755 Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L001133		Category and Description Package Insert, Oxidase Reagent Droppers	Sheet: 1 of 11 Scale: N/A	<b>A</b>

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## INTENDED USE

Oxidase Reagent Droppers are used in the Kovacs oxidase test as a qualitative reaction in the identification of nonfermenters and miscellaneous gram-negative bacteria.<sup>1,2</sup>

## SUMMARY AND EXPLANATION

The oxidase test is based on the production of an enzyme called indophenol oxidase. This enzyme oxidizes a redox dye (present in the reagent) which results in a color change of yellow to dark purple.<sup>1,2</sup>

## PRINCIPLES OF THE PROCEDURE

Indophenol oxidase, in the presence of atmospheric oxygen, oxidizes the phenylenediamine oxidase reagent to form a dark-purple compound, indophenol.<sup>1,2</sup>

## REAGENTS

Oxidase Reagent Droppers contain 0.5 mL of a 1% aqueous solution of N,N,N',N'-tetramethyl-p-phenylenediamine dihydrochloride which has been formulated with agents to ensure maximum stability.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

Avoid contact with the skin. Rinse thoroughly with water if spilled.

**Storage Instructions:** Store at controlled room temperature 15 – 30°C (59 – 86°F). Protect from light. The reagent droppers need no refrigeration.

**Product Deterioration:** This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until expiration date. Do not use if reagent is other than light yellow.

## PROCEDURE

**Material Provided:** Oxidase Reagent Droppers.

**Materials Required But Not Provided:** Ancillary culture media, filter paper, quality control organisms and laboratory equipment as required for this procedure.

### Test Procedure:

1. Hold reagent dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to crush ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate dropper any further, as the plastic may puncture and injury may occur.**
2. Tap bottom on tabletop a few times. Then invert for convenient drop-by-drop dispensing of reagent.
3. Preparation for testing
  - a. Colonies to be tested must be isolated from other colonies.
  - b. The use of fresh isolates (18 – 24 h cultures) is recommended for routine testing.
  - c. If refrigerated, cultures must be allowed to reach room temperature prior to testing.**Note:** Various conditions can affect the results obtained. Review the conditions listed in "Limitations of the Procedure."
4. Performing the test - Filter Paper Method
  - a. Add a few drops of Oxidase test reagent to a strip of filter paper (Whatman No. 1 or equivalent).
  - b. Streak a loopful of bacteria onto the reagent-saturated paper with a platinum loop or wooden applicator stick. Use of steel or nichrome loops may cause false-positive reactions.<sup>2</sup>

### User Quality Control:

It is recommended that positive (*Pseudomonas aeruginosa* ATCC™ 27853) and negative (*Escherichia coli* ATCC 25922) controls be run simultaneously with the organism to be tested.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidelines and CLIA regulations for appropriate Quality Control practices.

## RESULTS

Positive reactions turn the bacteria violet to purple immediately, or up to 30 s. Negative reactions remain colorless or turn light pink / light purple after 30 s. Delayed reactions should be ignored.

## LIMITATIONS OF THE PROCEDURE

Allow up to 30 s for a positive reaction. Any delayed reaction should be considered negative. Do not add excess reagent, as it may cause the reaction to fade on oxidase-positive organisms. Steel loop, nichrome loop and wire loop containing iron may give a false-positive reaction.<sup>3</sup> A platinum loop or wooden applicator stick is recommended.

Perform the oxidase test on gram-negative bacilli, but **only** on colonies from **nonselective** and/or **non-differential** media to ensure valid results. Selective or differential media can carry over the indicator to the filter paper and cause inaccurate results (false-negative reactions).<sup>2</sup> Colonies grown on media with high concentrations of glucose may inhibit oxidase activity.<sup>3</sup>

False-negative results may occur with mixed cultures containing the two genera *Pseudomonas* and *Neisseria* (not all *Pseudomonas* spp. elaborate oxidase).<sup>2</sup>

Timing is critical for interpretation of results.<sup>2</sup>

Do not test strictly anaerobic organisms.<sup>2</sup>

Viscid colonies may be negative due to poor penetration of the reagent.

Reactions from weak oxidase-positive organisms, e.g., *Pasteurella* species, may be inaccurate. Results inconsistent with other biochemical reactions or with the organisms should be repeated.

#### PERFORMANCE CHARACTERISTICS<sup>4</sup>

The oxidase test of Kovacs was applied to 1,660 strains of various genera, of which 210 strains gave a positive reaction and 64 strains gave a delayed positive reaction. Members of the *Neisseriaceae* and *Pseudomonadaceae* were generally oxidase positive, whereas members of the *Enterobacteriaceae* and, with few exceptions, gram-positive organisms, were oxidase-negative; equivocal results were obtained in the *Brucellaceae* and the genus *Bacillus*.

#### AVAILABILITY

##### Cat. No. Description

261181 BBL™ Oxidase Reagent Droppers, packaged 50 droppers/carton.

#### REFERENCES

1. Kovacs, N. 1956. Identification of *Pseudomonas pyocyanea* by the oxidase reaction. *Nature (London)* 178: 703.
2. York, M.K., M.M. Traylor, J. Hardy, and M. Henry. 2004. Biochemical tests for the identification of aerobic bacteria, p. 3.17.39.1. In H.D. Isenberg (ed.), *Clinical microbiology procedures handbook* 2nd ed., vol. 1, 2 and 3. American Society for Microbiology, Washington, D.C.
3. Chapin, K.C. and T.L. Lauderdale. 2003. Reagents, stains, and media: bacteriology, p. 354–393. In P.R. Murray, E.J. Baron, J.H. Tenover, M.A. Tenover and R.H. Tenover (eds.), *Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.
4. Steel, K.J. 1961. The oxidase reaction as a taxonomic tool. *J. Gen. Microbiol.* 25:297-306.

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Benex Limited