
Revisions

Rev from	Rev to	ECO #
0903	2010/06	5346-10

Notes:

1. BD Cat. Number 261510
2. Blank (Sheet) Size: Length: See VS Width: See VS
 Number of Pages: 12 Number of Sheets: 1
 Page Size: Length See VS Width See VS Final Folded Size: See VS
3. Style (see illustrations below): # 4

4. See Specification Control Number D11743RM for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p> <p style="font-weight: bold; font-size: large;">BD</p> <p>Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA</p>	
Proofer	Date		
Checked By	Date		
Part Number: L001179		Category and Description Package Insert, Bio-Bag Type C Environmental Chamber	Sheet: 1 of 13 <hr/> Scale: N/A

A

BD Bio-Bag™ Type C Environmental Chamber

English: pages 1–2 Italiano: pagine 6–7
Français : pages 3–4 Español: páginas 7–9
Deutsch: Seiten 4–6

 L001179
2010/06

Pokyny vám poskytnie miestni zástupce spoločnosti BD. / Kontakt den lokale BD repräsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Um Anleitungen zu erhalten, wenden Sie sich bitte an Ihren BD-Kundendienst. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Contattare il rappresentante BD di zona per ottenere il foglietto illustrativo. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзгiңiздiң жергiлiктi БД өкiлiне жүгiнiп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

Bio-Bag™ Type C is a disposable environmental system designed to rapidly generate and maintain a carbon dioxide environment for 48 h. Each system contains a sealable translucent **Bio-Bag** and a CO₂ Gas Generator.

SUMMARY AND EXPLANATION

A CO₂-enriched atmosphere is conducive to improved recovery of certain microorganisms; i.e., *Neisseria gonorrhoeae*,^{1,2} *Neisseria meningitidis*,^{2,3} and *Haemophilus* species.^{3,4} Growth and hemolysis on blood agar of some strains of streptococci are enhanced by 5 to 10% CO₂.⁴

Bio-Bag Type C is designed to generate and maintain a CO₂-fortified atmosphere for at least 48 h. This disposable system is intended to be used with appropriate media provided by the user. After incubation, bacterial growth may be macroscopically observed through the **Bio-Bag** at any interval without disrupting the CO₂-enriched atmosphere.

PRINCIPLES OF THE PROCEDURE

A CO₂-enriched atmosphere is provided in each sealed **Bio-Bag Type C** system. A self-contained gas generator consists of an ampule of dilute hydrochloric acid solution and a tablet of sodium bicarbonate. When the ampule is crushed, the tablet is activated producing CO₂. The resulting atmosphere is conducive to the growth of organisms which require a CO₂-enriched environment (e.g., *Neisseria gonorrhoeae*).⁵

REAGENTS

Each **Bio-Bag Type C** system consists of:

- 1 Gas Impermeable Environmental Chamber with a zipper closure.
- 1 Gas Generator consisting of one tablet of 130 mg sodium bicarbonate and one ampule of 0.5 mL of 1.56 N hydrochloric acid.

Warnings and Precautions

For *in vitro* Diagnostic Use.

Do not use if **Bio-Bag** appears damaged.

Bio-Bag must be tightly sealed to maintain proper gaseous environment.

Avoid skin contact with liquid contents of generator ampule.

Do not use generator if package has been previously opened or if generator appears to have been previously activated.

Do not allow generator to come in contact with water prior to use.

Warning: Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁶⁻⁹ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Storage: Store at room temperature 15 – 30°C (59 – 86°F).

PROCEDURE

Materials Provided: **Bio-Bag Type C** System Environmental Chambers and Gas Generators.

Materials Required But Not Provided:

Appropriate medium required for the isolation of organisms requiring a CO₂-enriched atmosphere and incubator.

Quality Control:

A stock strain of *Neisseria gonorrhoeae* or *Haemophilus influenzae* should be tested on a suitable medium in **Bio-Bag** Type C periodically to assure adequate conditions for recovery and characteristic growth.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

Test Procedure

1. Using accepted technique, inoculate specimen onto desired medium.
2. Insert inoculated plate into **Bio-Bag** Type C.
3. Insert CO₂ Gas Generator into **Bio-Bag** making sure pledget (polyester tip) end points toward the open end of the **Bio-Bag**.
4. To properly seal **Bio-Bag**, hold left side in left hand, grip zipper on left side firmly between thumb and fingers of right hand, move along zipper from left to right. Check to be sure continuous lock of zipper has been achieved.
5. Hold **Bio-Bag** so that CO₂ Generator remains upright. Crush lower portion of the ampule in the CO₂ Generator by a quick firm squeeze. Reaction is indicated by bubbling. If no bubbling occurs, flick lower portion of generator gently with finger to start the reaction. The CO₂ Generator must remain upright for at least 30 sec. (Noticeable inflation of the **Bio-Bag** does not occur with the quantity of gas required for this system to function properly.)
6. The **Bio-Bag** system is now ready for appropriate incubation.

RESULTS

CO₂-enriched conditions are achieved rapidly to create an environment suitable for the growth of *Neisseria gonorrhoeae*, *Neisseria meningitidis* and *Haemophilus sp.*

PERFORMANCE CHARACTERISTICS

Bio-Bag Type C was compared with the candle extinction jar, Gono-Pak and JEMBEC™ systems to determine their efficiency and reliability for the isolation of *Neisseria gonorrhoeae*.⁵ A total of 191 anal and 130 urethral specimens were

tested. There were 104 isolates of *N. gonorrhoeae* (24 anal and 80 urethral). The candle jar and **Bio-Bag** systems detected 98 (94%) of the isolates. The Gono-Pak and JEMBEC systems detected 102 (98%) and 100 (96%) of the isolates, respectively. These differences are not statistically different.

AVAILABILITY

Cat. No.	Description
261510	50 Bio-Bag ™ Type C Chambers, 50 sets per carton / 50 generators.

REFERENCES

1. Kellogg, D.S. Jr., Holmes, K., Hill, G.A., 1976. Cumitech 4, Laboratory diagnosis of gonorrhoeae. Coordinating eds., S. Marcus, J.C. Sherris. American Society for Microbiology, Washington, D.C.
2. Morello, J.A., Bohnhoff, M.: *Neisseria* and *Branhamella*, p. 111-130. In E.H. Lennette, A. Balows, W.J. Hausler, J.P. Truant (eds.), Manual of clinical microbiology, 3rd ed. American Society for Microbiology, Washington, D.C., 1980.
3. Bannatyne, R.M., Clausen, C., McCarthy, L.R., 1979. Cumitech 10, Laboratory diagnosis of upper respiratory tract infections. Coordinating ed., I.B.R. Duncan. American Society for Microbiology, Washington, D.C.
4. Kilian, M., *Haemophilus*, p. 330. In E.H. Lennette, A. Balows, W.J. Hausler, J.P. Truant (eds.), Manual of clinical microbiology, 3rd ed., American Society for Microbiology, Washington, D.C., 1980.
5. Carlson, B.L., Haley, M.S., Tisei, N.A., McCormack, W.M., 1980. Evaluation of four methods for isolation of *Neisseria gonorrhoeae*. J. Clin. Microbiol. 12:301-303.
6. National Committee for Clinical Laboratory Standards. 2001. Approved guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
7. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
8. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
9. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official journal L262, 17/10/2000, p. 0021-0045.

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JEMBEC is a trademark of Miles Scientific. BD, BD Logo, and Bio-Bag are trademarks of Becton, Dickinson and Company. ©2010 BD
Benex Limited