
## Revisions

SO 0191-5

Rev from	Rev to	ECO #
0703	2010/06	5347-10

### Notes:


- BD Cat. Number 261194
- Blank (Sheet) Size: Length: 11"      Width: 22.5"  
 Number of Pages: 10      Number of Sheets: 1  
 Page Size: Length 11"      Width 4.5"      Final Folded Size: 4 1/2" x 1 7/8"
- Style (see illustrations below): # 4


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number L001203 for Material Information
- Ink Colors: Printed two sides ☒ Yes      ☐ No  
 No. of Colors: 1      PMS# 2755 Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 <b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L001203		Category and Description  Package Insert, India Ink Reagent Droppers	Sheet: 1 of 11	A
			Scale: N/A	

# BD India Ink Reagent Droppers

English: pages 1 – 2      Italiano: pagine 5 – 6  
Français: pages 2 – 3      Español: páginas 6 – 7  
Deutsch: Seiten 4 – 5

 L001203  
2010/06

Pokyny vám poskytnú miestni zástupcov spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Um Anleitungen zu erhalten, wenden Sie sich bitte an Ihren BD-Kundendienst. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Contattare il rappresentante BD di zona per ottenere il foglietto illustrativo. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіңіздің жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

India Ink Reagent Droppers are used to enhance the microscopic detection of *Cryptococcus* spp. in wet preparations.

## SUMMARY AND EXPLANATION

India ink facilitates the visualization of cryptococcal polysaccharide capsules. The capsular material of cryptococci displaces the colloidal carbon particles of the ink so that the capsule appears as a clear halo around the microorganisms against a black background.<sup>1</sup>

India Ink Reagent Droppers may be used to stain pus, exudate, tissue, sputum and sediment of centrifuged urine and cerebral spinal fluid (CSF) specimens.<sup>2</sup> The India ink smear is particularly useful for CSF specimens because *Cryptococcus neoformans* is the most common cause of fungal central nervous system (CNS) infections throughout the world.<sup>3</sup> The appearance of encapsulated yeast forms in an India ink smear from CSF is diagnostic for cryptococcal meningitis.<sup>4</sup>

## PRINCIPLES OF THE PROCEDURE

India ink is a negative stain resulting in a dark background and unstained organisms. Cryptococci are spherical (sometimes oval) yeast forms which range from 4 µm to 20 µm in diameter. The cryptococcal polysaccharide capsule is present on nearly all *Cryptococcus* spp. and is resistant to coloration by India ink.

## REAGENTS

India Ink Reagent Droppers contain 0.5 mL of an aqueous suspension of India ink with suspending agent and preservative.

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper laboratory procedures in handling and disposing of infectious materials.

**Storage Instructions:** Store at controlled room temperature 15 – 30°C (59 – 86°F). Each dropper is good for one day's use after breaking ampule.

**Product Deterioration:** Reagents are hermetically sealed in an ampule which affords protection of the solution from chemical instability until expiration date. Do not use if extraneous material or organisms are present in stain control slide (see "User Quality Control"). Do not use after the expiration date.

## PROCEDURE

**Material Provided:** India Ink Reagent Droppers.

**Materials Required But Not Provided:** Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

### Test Procedure

Note: CSF, and other body fluids, should be centrifuged at 1,500 x g for 15 min prior to preparation of slide if greater than 1 mL of specimen is available.<sup>2</sup>

1. Hold reagent dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break the ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate dropper any further as the plastic may puncture and injury may occur.**
2. Tap bottom of dropper on tabletop a few times. Then invert for convenient drop-by-drop dispensing of reagent.
3. To a small loopful of sterile saline, water or broth on a precleaned microscope slide, add a small amount of growth from a young agar culture or specimen sediment, or add a drop of CSF directly to the slide.
4. Add one drop of India ink and mix.
5. Cover with a thin, large coverslip. Press gently to allow the fluid to spread as a thin film beneath the coverglass.
6. Scan slide under low power and examine under high power or oil immersion for encapsulated yeast forms.

### User Quality Control

1. Examine the reagent for signs of deterioration (see "Product Deterioration").
2. Positive (a 48-h culture of *Cryptococcus neoformans* ATCC™ 36556), negative (*Candida albicans* ATCC 10231) and stain (1 drop India ink + 1 drop sterile saline) controls should be run simultaneously with the organism or specimen being tested.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

In positive smears, cryptococcal capsules appear as clear, refractile halos either surrounding the yeast cells or positioned between the well-defined cell wall and the surrounding black mass of India ink particles. Capsules may be broad or narrow. The yeast cells may be round, oval, or elongate; buds may be absent, single, or, rarely, multiple. The buds may be detached from the mother cell but enclosed in a common capsule.

## LIMITATIONS OF THE PROCEDURE

1. Care must be taken not to misinterpret fat droplets or lymphocytes as yeast forms. Fat droplets will not have a well-defined cell wall. Lymphocytes have a shaggy border and an eccentric nucleus. The capsules of cryptococci are usually sharply outlined and the well-defined yeast cell wall is centrally located within the capsule.<sup>3</sup>
2. Organisms in AIDS patients may not have capsules.<sup>1</sup>
3. India ink smears performed on CSF sediment are positive in only 60% of all cases of cryptococcal meningitis.<sup>3</sup> All requests for India ink smears should be supplemented with culture or antigenic detection.<sup>3</sup>
4. *Rhodotorula* spp. and rare *Torylospira* spp. will produce capsules.<sup>4</sup>
5. Positive India ink smears may be used as a presumptive diagnosis only. Definitive diagnosis must be accomplished by culture and/or antigenic detection.<sup>1,3</sup>
6. If smear is too dark, dilute the India ink 1:2 with sterile, distilled water and make a new smear. Specimens treated with 20% KOH are not to be used for India ink preparations.

## PERFORMANCE CHARACTERISTICS

*Cryptococcus neoformans* is a causative agent of meningitis. Several studies<sup>5-10</sup> point to the utility of the India Ink stain as one of the tools for diagnosing *C. neoformans*. These articles show that the India Ink stain is useful in identifying *C. neoformans* in cerebrospinal fluid (CSF) samples.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company. © 2010 BD.
Benex Limited