
Revisions

Rev from	Rev to	ECO #
0703	2010/06	5347-10

Notes:

- BD Cat. Number 261204
- Blank (Sheet) Size: Length: 11" Width: 22.5"
 Number of Pages: 10 Number of Sheets: 1
 Page Size: Length 11" Width 4.5" Final Folded Size: 4.5" x 1 7/8"
- Style (see illustrations below): # 4

- See Specification Control Number L001217 for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755 Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L001217		Category and Description Package Insert, Methylene Blue Loeffler Reagent Droppers	Sheet: 1 of 11 Scale: N/A	A

BD Methylene Blue Loeffler Stain Droppers

English: pages 1 – 2 Italiano: pagine 5 – 6
Français: pages 2 – 4 Español: páginas 7 – 8
Deutsch: Seiten 4 – 5

L001217
2010/06

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Um Anleitungen zu erhalten, wenden Sie sich bitte an Ihren BD-Kundendienst. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Contattare il rappresentante BD di zona per ottenere il foglietto illustrativo. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзініздің жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

BBL™ Methylene Blue Loeffler Stain Droppers are used for presumptive identification of *Corynebacterium diphtheriae* and as an adjunct to the Gram stain.

SUMMARY AND EXPLANATION

Methylene Blue is recognized as a simple stain that can be used for a variety of purposes, such as determining bacterial morphology and aiding in the interpretation of the Quellung test.¹ It is the stain of choice for the presumptive identification of *Corynebacterium diphtheriae*. Upon staining, this species is characterized by a beaded or banded appearance of the cell, with deeply stained metachromatic granules contrasted with a lighter blue cytoplasm.^{2,3}

Methylene Blue is used as an adjunct to the Gram stain, especially in the staining of gram-negative bacteria such as *Haemophilus influenzae* and *Neisseria* species.^{4,5} The stain can reveal the morphology of fusiform bacteria and spirochetes (from oral infections such as Vincent's angina) that may be difficult to see with the Gram stain.⁵ The stain is also used to detect the presence of fecal white blood cells that may show the presence of invasive bowel disease.⁵

BBL Methylene Blue Loeffler Stain Droppers can be used as a counterstain with **BBL Flagella Stain Droppers**, which are recommended for staining the flagella of bacteria.

BBL Methylene Blue Loeffler Stain Droppers are supplied in a hermetically sealed ampule for improved stability of the reagent prior to use.

PRINCIPLES OF THE PROCEDURE

Methylene Blue is a cationic dye that imparts a blue color to the negatively charged polyphosphate portions of cells (such as ribonucleic acid).¹ Metachromasia, which is observed with *C. diphtheriae*, is due to a high concentration of polymerized polyphosphates. The deeply blue stained polyphosphate granules, surrounded by a lighter blue stained cytoplasm, are known as metachromatic granules or Babes-Ernst bodies.

REAGENTS

BBL Methylene Blue Loeffler Stain Droppers contain: 0.3 g Methylene Blue, 30 mL Alcohol, 95%, 0.1 mL Potassium Hydroxide, 10%, and 100 mL Distilled Water.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Flammable. Harmful by inhalation and if swallowed. Keep container tightly closed. Avoid contact with skin and eyes. Wear suitable protective clothing.

Follow established laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store Methylene Blue Loeffler Stain Droppers at 15 – 30°C. Methylene Blue Loeffler Stain Droppers are ready to use.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Refer to appropriate references for specific procedures on processing specimens.⁴⁻⁷

The microscopic morphology of *C. diphtheriae* is particularly distinctive when grown on Loeffler Blood Serum.⁸

PROCEDURE

Materials Provided: Methylene Blue Loeffler Stain Droppers.

Materials Required But Not Provided: Glass slides, Bunsen burner, inoculating loop and quality control organisms.

Test Procedure:

1. Prepare a smear of the specimen to be stained on a clean glass slide.
2. Allow the smear to air dry thoroughly.
3. Heat or methanol fix the smear.
4. Hold dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break the inner ampule. **Caution: Break ampule close to its center one time only. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.** Tap bottom on tabletop a few times. Invert for convenient drop by drop dispensing of solution.
5. Flood the slide with stain and allow to stain for 1 – 3 min.
6. Rinse the slide with tap or distilled water and blot dry.
7. Examine the smear with an oil-immersion objective.

User Quality Control:

Identity Specifications – Blue, clear solution, inside a glass ampule contained in a plastic container.

Cultural Response – Test the performance of Methylene Blue Loeffler Stain Droppers using the procedure described above. Observe *C. diphtheriae* for metachromatic granules.

Quality control requirements must be performed in accordance with applicable

local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

Organism	ATCC™	Stain Reaction
<i>Escherichia coli</i>	25922	blue rods
<i>Streptococcus pyogenes</i>	19615	blue cocci
<i>Corynebacterium diphtheriae</i> biotype <i>gravis</i>	8028	deep blue granules

RESULTS

Bacterial cells stain medium blue, while the background stains lighter blue. Stained *C. diphtheriae* appear as pleomorphic beaded rods with swollen ends, producing a club shape, and containing deeply stained metachromatic granules.⁶

LIMITATIONS OF THE PROCEDURE

Over-staining may reduce the contrast between granules and cytoplasm or between bacteria and background.⁷

The Loeffler staining test may be used in the presumptive identification of *Corynebacterium diphtheriae*. Additional biochemical testing is recommended for complete identification because some strains of *Actinomyces*, *Propionibacterium*, and pleomorphic forms of streptococci may demonstrate the characteristic stained appearance of *C. diphtheriae*.⁴

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of BBL Methylene Blue Loeffler Stain Droppers are tested to verify specific product characteristics. One drop of a heavy saline suspension (approximately McFarland standard no. 3) of *Escherichia coli* (ATCC 25922) and *Streptococcus pyogenes* (ATCC 19615) is placed on a microscope slide, allowed to air dry and heat fixed. A smear of *Corynebacterium diphtheriae* biotype *gravis* (ATCC 8028) is placed on a clean microscopic slide, allowed to air dry and heat fixed. The slides are flooded with BBL Methylene Blue Loeffler Stain and allowed to stand for 2 min. They are then rinsed with water and blotted dry. The slides are examined under the oil-immersion lens of a light microscope. Bacterial cells of *E. coli* and *S. pyogenes* stain blue although the color intensity may vary. Cell morphology is readily distinguishable. Metachromatic granules in the cells of *C. diphtheriae* stain a deeper blue than the surrounding cytoplasm.

AVAILABILITY

Cat. No. Description

261204 BBL™ Methylene Blue Loeffler Stain Droppers, 50 x 0.5 mL.

REFERENCES

1. Lillie, R. D. (ed.). 1977. H. J. Conn's biological stains, 9th ed. Williams and Wilkins, Baltimore.
2. Clark, G. (ed.). 1981. Staining procedures, 4th ed. Williams and Wilkins, Baltimore.
3. Saelinger, C. B., and F. F. Bonventre. 1983. *Corynebacterium diphtheriae*. In N. R. Rose and A. L. Baren (ed.), Microbiology: basic principles and clinical applications. Macmillan, New York.
4. Koneman, E. W., S. D. Allen, V. R. Dowell, Jr., and H. M. Sommers. 1983. Color atlas and textbook of diagnostic microbiology, 2nd ed. J. B. Lippincott Company, Philadelphia.
5. Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
6. Isenberg, H. D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
7. Baron, E. J., and S. M. Finegold. 1992. Bailey & Scott's diagnostic microbiology, 8th ed. The C. V. Mosby Company, St. Louis.
8. Lennette, E. H., A. Balows, W. J. Hausler, Jr., and J. P. Truant (ed.). 1980. Manual of clinical microbiology, 3rd ed. American Society for Microbiology, Washington, D.C.

BD Methylene Blue Loeffler Stain Droppers

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Benex Limited