BD BBL™ Flagella Stain Droppers

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INTENDED USE

Flagella Stain Droppers are used for demonstrating bacterial flagella and their arrangement on the cell.

SUMMARY AND EXPLANATION

Due to their narrow diameter, bacterial flagella cannot be seen with the light microscope without the use of special stains. The presence of bacterial flagella, flagellum length, shape, curvature, number per cell and arrangement on the cell are useful in the identification of motile bacteria. The Flagella Stain is based on the tannic acid-crystal violet formula first described by Ryu¹ and later by Kodaka, Armfield, Lombard and Dowell.²

The presence of flagella is especially useful in the identification of nonfermentative, gram-negative bacilli³⁻⁵ and anaerobic bacteria.²

PRINCIPLES OF THE PROCEDURE

The alcoholic solution of crystal violet functions as the primary stain. As alcohol evaporates during the staining process, crystal violet forms a precipitate around the flagella, increasing their apparent size. Tannic acid and aluminum potassium sulfate act as mordants. Phenol retards growth of molds.

REAGENTS

The Flagella Stain Droppers contains: 0.6% Crystal Violet in Ethanol, 2.0% Tannic Acid, 2.5% Phenol and 5.7% Aluminum Potassium Sulfate.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper established laboratory procedures in handling and disposing of infectious materials.

Danger

H315 Causes skin irritation. **H319** Causes serious eye irritation. **H341** Suspected of causing genetic defects. **H350** May cause cancer. **H412** Harmful to aquatic life with long lasting effects. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P302 + P352** IF ON SKIN: Wash with plenty of soap and water. **P305 + P351 + P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store Flagella Stain Droppers at 15 – 30 °C. The Flagella Stain is light sensitive. Protect from light. Flagella Stain Droppers are ready for use.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Collect specimens in sterile containers or with sterile swabs and transport immediately to the laboratory according to recommended guidelines.³⁻⁵

Process each specimen using procedures appropriate for that sample.³⁻⁵

PROCEDURE

Materials Provided: Flagella Stain Droppers and ampule crusher.

Materials Required But Not Provided: Glass microscope slides, inoculating loop, Bunsen burner or incinerator, wax pencil and quality control organisms.

Test Procedure:

Smear Preparation

1. Obtain a pure culture of the organism to be tested. Grow the test organism on blood agar for 18 – 48 h.³
2. Draw a border around the clear portion of a microscope slide with a wax pencil.
3. Place a drop of distilled or deionized water on the slide, approximately 1 cm from the frosted edge.
4. Gently touch a colony of the culture being tested with an inoculating loop and then lightly touch the drop of water without touching the slide. Do not mix.
5. Tilt the slide at a slight angle to allow the drop to flow to the opposite end of the slide.
6. Allow the slide to air-dry at room temperature. Do not fix with heat.

Staining Procedure

1. Holding the Flagella Stain dispenser upright and with the tip pointing away from the user, squeeze gently to crush the glass ampule inside the dispenser.
2. Invert and squeeze slightly to flood the slide with the Flagella Stain.
3. Allow the stain to remain on the slide for approximately 4 min.
4. Carefully rinse off the stain by adding running water to the slide. Do not tip the slide.

- When done rinsing, gently tilt the slide to allow excess water to run off. Allow the slide to air-dry at room temperature.
- Examine the slide microscopically with the oil-immersion objective.

User Quality Control:

Identity Specifications – Dark purple solution with a slight precipitate inside a glass ampule contained in a plastic container.

Cultural Response – Test the performance of Flagella Stain Droppers using the staining procedure described above. Examine microscopically using the oil-immersion lens of a light microscope.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Organism	ATCC™	Flagella Pattern
<i>Klebsiella pneumoniae</i>	13883	no flagella
<i>Proteus vulgaris</i>	13315	peritrichous*
<i>Pseudomonas aeruginosa</i>	27853	monotrichous**

* Flagella arising over the entire cell surface.

** Single polar flagellum.

RESULTS

Bacteria and their flagella should stain purple. Compare morphology against reference descriptions, drawings or pictures.

LIMITATIONS OF THE PROCEDURE

Flagella staining is used to help in the identification of bacteria. Additional biochemical testing is recommended for complete identification.

Many factors, including the age of the culture, the temperature of incubation and the presence of fermentable substances in the growth medium, have been shown to affect flagellation.⁴

Flagella are easily dislodged from a cell. Care is recommended, therefore, in transferring organisms to the droplet of water on the slide.

Staining quality may be affected by the length of staining time. Adjust staining time to achieve optimal results.

PERFORMANCE CHARACTERISTICS²

Kodaka, et al. developed a simple, rapid procedure for staining bacterial flagella using the Flagella stain developed by E. Ryu. In an 18-month period, more than 100 strains of motile and nonmotile bacteria tested using the staining procedure gave satisfactory results. Those included members of the genera *Anaerobiospirillum*, *Bacillus*, *Bacteroides*, *Campylobacter*, *Clostridium*, *Escherichia*, *Fusobacterium*, *Listeria*, *Proteus*, *Pseudomonas*, *Serratia* and *Wolinella*.

AVAILABILITY

Cat. No.	Description
261206	BD BBL™ Flagella Stain Droppers, 50.

REFERENCES

- Ryu, E. 1937. A simple method of staining bacterial flagella. *Kitasato Arch. Exp. Med.* 14:218-219.
- Kodaka, H., A. Y. Armfield, G. L. Lombard, and V. R. Dowell, Jr. 1982. Practical procedure for demonstrating bacterial flagella. *J. Clin. Microbiol.* 16:948-952.
- Gilardi, G. L. 1991. *Pseudomonas* and related genera, p. 429 – 441. In A. Balows, W. J. Hausler, K. L. Herrmann, H. D. Isenberg, and H. J. Shadomy (ed.), *Manual of clinical microbiology*, 5th ed. American Society for Microbiology, Washington, D. C.
- Koneman, E. W., S. D. Allen, V. R. Dowell, Jr., and H. M. Sommers. 1983. *Color atlas and textbook of diagnostic microbiology*, 2nd ed. J. B. Lippincott Company, Philadelphia.
- Baron, E. J., L. R. Peterson, and S. M. Finegold. 1994. *Bailey & Scott's diagnostic microbiology*, 9th ed. Mosby-Year Book, Inc., St. Louis.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Australian Sponsor:
Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company. © 2016 BD.