
Revisions

Rev from	Rev to	ECO #
1204	2010-01	5253-10

Notes:

- BD Cat. Number 220235, 220236
- Blank (Sheet) Size : Length: 23 5/8" Width: 3 15/32"
 Number of Pages: 2 Number of Sheets: 1
 Page Size: Length 23 5/8" Width 3 15/32" Final Folded Size: Special Fold*
- Style (see illustrations below): # 1

- See Specification Control Number OEM - Copan for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 1665C Orange
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

*Peel-off insert / Label. Insert perforated 3/8" from bottom with 3/16" notch for easy peel-off. Top insert size 4 3/16" x 3 15/32" includes 3/8" bottom portion of glued area. Insert is folded under with one horizontal fold and followed by two parallel folds to form the size of 3 1/4" x 3 15/32." Insert to be glued on top 5/16" and bottom 3/8" to roll label, size 3 15/32" (web) x 4 3/8".

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L009332/H284		Category and Description Peel-off Insert / Label, CultureSwab MaxV(+)	Sheet: 1 of 3 <hr/> Scale: N/A	A

Pokyny vám poskytne miestni zástupca spoločnosti BD. / Kontakt den lokale BD representant for at få instruktioner. / Kasutusjuhiste suhtes kontaktteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD igaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie zyskate u miestneho zistupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свяжитесь с местным представителем на BD за инструкциями. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

BD CultureSwab™ MaxV(+) is a sterile ready-to-use system intended for the collection, transport and preservation of clinical specimens for bacteriological examination.

SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of bacterial infections involves the collection and safe transportation of a clinical specimen from the patient to the laboratory. This can be accomplished using the **BD CultureSwab MaxV(+)** collection and transport device. Each **BD CultureSwab MaxV(+)** unit is comprised of a sterile peel pouch containing a rayon-tipped swab applicator used to collect the sample and a tube containing transport medium into which the swab applicator is placed after sampling. The **BD CultureSwab MaxV(+)** is available with Amies Medium Without Charcoal. This transport medium is non-nutritive, buffered with phosphate and provides a reduced environment, due to the formulation containing sodium thioglycollate.¹ Organisms in the sample material are protected from drying by moisture in the transport medium. The medium is designed to maintain the viability of organisms during transit to the laboratory. The **BD CultureSwab MaxV(+)** pouch is made of a plastic film which retards the penetration of atmospheric air into the product.

BD CultureSwab MaxV(+) media are contained in a constricted (venturi) tube. Nitrogen gas is flushed into the transport tube during the media filling and capping process. During final packaging of the swab and tube, air is removed from the pouch by vacuum and nitrogen gas is flushed inside.

REAGENTS

The nominal formula for the medium per L of distilled water is:

Amies Transport Medium Without Charcoal

Sodium Chloride	3.0 g	Monopotassium Phosphate	0.2 g
Potassium Chloride	0.2 g	Disodium Phosphate	1.15 g
Calcium Chloride	0.1 g	Sodium Thioglycollate	1.00 g
Magnesium Chloride	0.1 g	Bacteriological Agar	7.5 g

Precautions: For *in vitro* Diagnostic Use.

It must be assumed that all specimens contain infectious microorganisms; therefore, all specimens should be handled with appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste.

ⓧ **BD CultureSwab MaxV(+)** is for single use only; reuse may cause a risk of infection and/or inaccurate results.

Storage: Store **BD CultureSwab MaxV(+)** at 5–25°C.

Product Deterioration: Contents sterile if unopened or not damaged. Do not use if they show evidence of damage, dehydration or contamination. Do not use past expiration date.

SPECIMEN COLLECTION AND HANDLING

For specific recommendations about collection of specimens for microbiological analysis and primary isolation techniques, consult appropriate references.²⁻⁴

Once a swab sample is collected, it should be placed in the tube of medium, transported to the laboratory as soon as possible and cultured onto appropriate primary isolation media.

PROCEDURE

Materials Required: Fifty (50) units of sterile **BD CultureSwab MaxV(+)** devices contained in each Vi-Pak Pouch. **Materials Required But Not Provided:** Appropriate materials for isolating, differentiating and culturing aerobic and anaerobic bacteria. These materials include culture media plates or tubes and incubation systems, gas jars or anaerobic workstations.

Directions for Use:

The directions for use are printed on each **BD CultureSwab MaxV(+)** unit, along with descriptive diagrams. Directions for use are summarized as follows:

1. Peel open the **BD CultureSwab MaxV(+)** pouch.
2. Remove cap from transport tube.
3. Remove applicator swab and collect specimen.
During specimen collection, the applicator tip should only touch the area where the infection is suspected to minimize potential contamination.
4. Place applicator swab in transport tube.
5. Record patient's name and information on tube label.
6. Send specimen to the laboratory for immediate analysis.

EXPECTED RESULTS

The survival of bacteria in a transport medium depends on many factors. These include the type of bacteria, duration of transport, storage temperature, concentration of bacteria in the sample and formulation of the transport medium. **BD CultureSwab MaxV(+)** will maintain viability of many microorganisms for 24–48 h. For fastidious bacteria such as *Neisseria gonorrhoeae* and *Streptococcus pneumoniae* and anaerobes, swab specimens should be plated directly onto culture medium or transported immediately to the laboratory and cultured within 24 h.

LIMITATIONS OF THE PROCEDURE

BD CultureSwab MaxV(+) Amies Medium Without Charcoal is intended for the collection and transport of bacteriological samples only. Preferred samples for anaerobic investigations are: tissue samples obtained during surgical procedures, biopsies from tissue or bone, fluid, pus, or aspirates collected using a syringe. For detailed information and recommendations for transporting fluid and tissue specimens for anaerobic culture, refer to specific publications.³⁻⁷ Samples containing viruses or chlamydia should be collected and transported using alternative specific transport systems.

Transport media, staining reagents, immersion oil, glass slides and specimens themselves may sometimes contain nonviable organisms visible upon Gram staining. Caution should therefore be used when interpreting Gram stains from either sterile body fluids or specimens from normally sterile body fluids.

PERFORMANCE CHARACTERISTICS

Recovery studies were performed using **BD CultureSwab MaxV(+)** Amies Medium Without Charcoal with a variety of aerobic and anaerobic organisms. Swabs were dosed with inoculum and inserted into the transport tube containing medium. The tubes were stored at room temperature prior to subculturing onto appropriate media.

Aerobic organisms evaluated were *Escherichia coli* (NCTC 9001 and ATCC™ 25922), *Haemophilus influenzae* (ATCC 10211), *Neisseria gonorrhoeae* (ATCC 43069), *Neisseria meningitidis* (NCTC 10025 and ATCC 13090), *Pseudomonas aeruginosa* (NCTC 9332 and ATCC 27853), *Staphylococcus aureus* (NCTC 5532 and ATCC 25923), and *Streptococcus pyogenes* (ATCC 19615). Anaerobic organisms evaluated were *Bacteroides fragilis* (ATCC 25285), *Bacteroides levii* (ATCC 29147), *Bacteroides thetaiotaomicron* (ATCC 29741), *Bacteroides vulgatus* (ATCC 8482), *Clostridium difficile* (ATCC 9689), *Clostridium perfringens* (ATCC 13124), *Clostridium sporogenes* (ATCC 3584), *Clostridium tertium* (ATCC 19405), *Fusobacterium necrophorum* (ATCC 25286), *Fusobacterium nucleatum* (ATCC 25586), *Peptostreptococcus anaerobius* (ATCC 27337), *Peptostreptococcus magnus* (ATCC 29328), *Porphyromonas gingivalis* (ATCC 33277), *Prevotella melaninogenica* (ATCC 25845), *Propionibacterium acnes* (ATCC 6919), and *Streptococcus pneumoniae* (ATCC 6305).

All organisms tested remained viable for more than 24 h when maintained at room temperature.

AVAILABILITY - See end of insert.

REFERENCES

1. Amies CR. 1967. A modified formula for the preparation of Stuart's transport medium. Can. J. Public Health. 58: 296-300.
2. Isenberg H.D., F.D. Schoenkencht, and A. von Graevenitz. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Rubin. American Society for Microbiology, Washington, D.C.
3. Balows, A., W.J. Hausler, Jr, K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
4. Isenberg, H.D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
5. Zavala, M.K., D.M. Citron, E.J.C. Goldstein. 1998. Evaluation of a novel specimen transport system for anaerobic bacteria. Clin. Infect. Dis. 25 (supplement 2): S132-S133.
6. Pery, J.L. 1997. Assessment of swab transport systems for aerobic and anaerobic organism recovery. J. Clin. Microbiol. 35: 1269-1271.
7. Summanen, P., E.J. Baron, D.M. Citron, C.A. Strong, H.M. Wexler, and S.M. Finegold. 1993. Wadsworth anaerobic bacteriology manual, 5th ed. Star Publishing Co., Belmont, Calif.

TECHNICAL INFORMATION: In the United States, telephone Technical Services, toll free (800) 638-8663. ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo and CultureSwab are trademarks of Becton, Dickinson and Company. © 2010 BD.

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