
## Revisions

SO 0191-5

Rev from	Rev to	ECO #
2010-01	2010/12	5525-10

**Notes:**

1. BD Cat. Number 441358
2. Blank (Sheet) Size: Length: See Spec Width: See Spec  
 Number of Pages: 12 Number of Sheets: See Spec  
 Page Size: Length 8.5" Width 5.5" Final Folded Size: See Spec
3. Style (see illustrations below): # See Spec

4. See Specification Control Number 16354RM for Material Information
5. Ink Colors: Printed two sides  Yes  No  
 No. of Colors: 1 PMS# Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	<b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L009851		Category and Description Package Insert, Male Urethral Specimen Collection Kit for the BD ProbeTec CT/GC Q <sup>X</sup> Amplified DNA Assay	Sheet: 1 of 13 <hr/> Scale: N/A	A

# **BD Male Urethral Specimen Collection Kit for the BD ProbeTec™ *Chlamydia trachomatis*/ *Neisseria gonorrhoeae* (CT/GC) Q<sup>x</sup> Amplified DNA Assays**

English: pages 1 – 2 Français : pages 6 – 7  
Dansk: side 2 – 3 Norsk: sider 7 – 8  
Deutsch: Seiten 3 – 4 Svenska: sidan 9 – 10  
Español: páginas 5 – 6

L009851  
2010/12

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontaktteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / **Свяжитесь с местним представител на BD за инструкции.** / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті BD өкіліне жүгініп нұсқа алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

The Male Urethral Specimen Collection Kit for the **BD ProbeTec CT/GC Q<sup>x</sup> Amplified DNA Assays** is used to collect and transport male patient urethral specimens to the laboratory for testing with the **BD ProbeTec CT/GC Q<sup>x</sup> Amplified DNA Assays** on the **BD Viper™ System** in Extracted Mode.

## SUMMARY AND EXPLANATION

The Male Urethral Specimen Collection Kit for the **BD ProbeTec CT/GC Q<sup>x</sup> Amplified DNA Assays** consists of a sterile swab and a tube of Q<sup>x</sup> Swab Diluent for specimen transport. The rayon-tipped collection swab consists of a stainless steel wire in a blue-colored plastic shaft that is scored for breaking. After the specimen is collected, the cap of the diluent tube is removed, the swab is placed into the Q<sup>x</sup> Swab Diluent tube and the shaft is broken at the score mark. The tube is recapped and transported to the testing facility.

## REAGENTS

**Materials Provided:** Each Male Urethral Specimen Collection Kit for the **BD ProbeTec CT/GC Q<sup>x</sup> Amplified DNA Assays** contains one sterile rayon-tipped collection swab and one tube of Q<sup>x</sup> Swab Diluent. The Swab Diluent contains approximately 2 mL of potassium phosphate/potassium hydroxide buffer with DMSO and preservative.

**Materials Required But Not Provided:** Gloves

**Storage Instructions:** Store collection kit at 2 – 33°C. Do not use beyond the expiration date. Swab sterility is guaranteed if primary swab container is intact.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

1. **Do not pre-wet the swab with the Q<sup>x</sup> Swab Diluent prior to collecting specimen.**
2. Pathogenic microorganisms including hepatitis viruses and Human Immunodeficiency Virus may be present in clinical specimens. "Standard Precautions"<sup>1-4</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
3. Optimal performance of the **BD ProbeTec CT/GC Q<sup>x</sup> Amplified DNA Assays** requires proper specimen collection, handling and transport.
4. Use only collection and transport systems labeled for use with the **BD ProbeTec CT/GC Q<sup>x</sup> Amplified DNA Assays**.
5. Proper labeling should accompany each specimen to the laboratory.
6. The Q<sup>x</sup> Swab Diluent tube contains dimethyl sulfoxide (DMSO). DMSO is harmful by inhalation, in contact with skin and if swallowed. Keep container in a well-ventilated place. Avoid contact with eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. Wear suitable protective clothing and gloves.
7. When breaking the shaft of the swab, take care to avoid splashing of contents.
8. Dispose of unused Q<sup>x</sup> Swab Diluent tubes according to the Material Safety Data Sheet.
9. Take care to avoid cross-contamination during the specimen handling steps. Specimens can contain high levels of organisms. Ensure that specimen containers do not contact one another, and discard used materials without passing over open containers. If gloves come in contact with specimen, change gloves to avoid cross-contamination.
10. Specimens must be collected and tested before the expiration date of the Q<sup>x</sup> Swab Diluent tube.
11. **Ⓢ Male Urethral Specimen Collection Kit for the **BD ProbeTec *Chlamydia trachomatis*/ *Neisseria gonorrhoeae* (CT/GC) Q<sup>x</sup> Amplified DNA Assays** is for single use only; reuse may cause a risk of infection and/or inaccurate results.**

## SPECIMEN COLLECTION AND TRANSPORT

1. Remove the swab from the packaging.
2. Insert the swab 2 – 4 cm into the urethra and rotate for 3 – 5 s.
3. Withdraw the swab.
4. Uncap the Q<sup>x</sup> Swab Diluent tube.
5. Fully insert the collection swab into the Q<sup>x</sup> Swab Diluent tube.
6. Break the shaft of the swab at the score mark. Use care to avoid splashing of contents.
7. **Tightly** recap the tube.
8. Label the tube with patient information and date/time collected.
9. Transport to the laboratory.

### Swab Storage and Transport:

The Q<sup>x</sup> Swab Diluent tube with swab must be stored and transported to the laboratory and/or test site at 2 – 30°C within 30 days of collection. Alternatively, the Swab Diluent tube may be stored and transported frozen (-20°C) to the laboratory and/or test site within 180 days of collection.

For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

## REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guide for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
3. U. S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5<sup>th</sup> ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/31/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

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Becton Dickinson France S.A.S.