BD GeneOhm™ StrepB

REF 441240

Table of Contents / Tables des Matières / Inhaltstabelle / Índice / Indice

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English

Intended use

BD GeneOhm™ StrepB assay is a qualitative *in vitro* diagnostic test for the rapid detection of Group B streptococcus (GBS) DNA in vaginal/rectal specimens from prepartum or intrapartum women. The test performed on the SmartCycler® automated analyzer utilizes polymerase chain reaction (PCR) for the amplification of a *cfb* gene sequence of GBS recovered from clinical samples and fluorogenic target-specific hybridization for the detection of the amplified DNA.

BD GeneOhm™ StrepB assay can be used to establish GBS colonization status of prepartum and intrapartum women.

Summary and explanation of the test

A vaginal/rectal specimen is collected and transported to the laboratory using a recommended swab with Liquid Stuart Medium (refer to Materials required but not provided). The swab is eluted in sample buffer; an aliquot of the specimen is then lysed and added to the PCR reagents which contain the GBS-specific primers used to amplify the genetic target of the assay (*cfb* gene), if present. The assay also includes an internal control (IC) to detect PCR inhibitory specimens and to confirm the integrity of assay reagents. Amplified targets (*cfb* gene and IC) are detected with hybridization probes labelled with quenched fluorophores (molecular beacons). The amplification, detection and interpretation of the signals are done automatically by the Cepheid SmartCycler® software. The whole procedure takes about 45 minutes.

Recently, the incidence of perinatal group B streptococcal disease in the United States has been decreasing because of antibiotic prophylaxis used for prevention of GBS diseases^{1,2}. A national survey of ACOG members in 2000 found that 73.5% of respondents used a screening-based strategy, showing that antepartum screening highly contributed to the 70% decrease of early-onset GBS disease since the 1996 Centers for Disease Control and Prevention (CDC) guidelines were implemented³. Moreover, in cases of short labor, the availability of antepartum screening results allows for the administration of intravenous antibiotics as soon as possible after admission (< 5 hours). Based on current diagnostic tools, the US CDC recommends universal prenatal culture-based screening for vaginal and rectal GBS colonization of all women at 35-37 weeks gestation².

Although this recommendation is based on documentation of a strong protective effect of the culture-based screening strategy relative to the risk-based strategy, antepartum screening has limitations^{2,4}. The results of the screening are sometimes not known to the physician at the time of delivery, colonization status may change between testing and delivery, and some women do not undergo antepartum screening^{2,4}. In each case, these women may not be treated appropriately for GBS colonization, thereby increasing the risk of infecting the newborn. A sensitive rapid-screening test for GBS that can provide intrapartum screening results in time to administer antibiotic prophylaxis would be useful to several groups of women including those who have had inadequate prenatal care, those whose results of antepartum screening are unknown at the time of delivery, and those facing preterm delivery.

Warning

- The use of BD GeneOhm™ StrepB for intrapartum screening should not preclude the use of other strategies (e.g., antepartum testing). Intrapartum BD GeneOhm™ StrepB results are useful to identify candidates for intrapartum antibiotic prophylaxis when results can be available in time to administer intravenous antibiotics at least 4 hours before delivery.
- The BD GeneOhm™ StrepB assay does not provide susceptibility results. Additional time is needed to culture and perform susceptibility testing that would be recommended for penicillin-allergic women.

Principle of the procedure

The genetic target of the BD GeneOhm™ StrepB assay is the *cfb* gene. This gene encodes for the CAMP factor, a diffusible extra-cellular protein which is present in virtually all GBS isolates. Detection of the CAMP factor is used for the presumptive identification of GBS by biochemical methods. Phenotypic and molecular characterizations showed that the *cfb* gene is well conserved within this species^{5,6}. Primers complementary to conserved and specific regions of the gene have been developed to amplify a 154 bp fragment of the *cfb* gene.

A vaginal/rectal swab is collected and transported to the laboratory using the recommended swab with Liquid Stuart Medium (refer to Materials required but not provided). The swab is placed in sample buffer

to elute its content and an aliquot is transferred to the lysis tube. Lysis occurs through a combination of chemical and physical actions and takes less than 15 minutes. A sample of the lysate is added directly to the PCR reagents contained into the SmartCycler® reaction tube, and run with the SmartCycler® instrument, a random access real-time fluorescence thermal cycler.

In specimens containing GBS, the 154 bp region of the *cfb* gene will be amplified and detected. Amplification of the IC, a 180-bp DNA fragment consisting of a 134-bp sequence not found in GBS flanked by the sequence of each of the two GBS-specific primers, will also take place.

The amplified DNA targets are detected with molecular beacons, a hairpin-forming single-stranded oligonucleotides labelled at one end with a quencher and at the other end with a fluorescent reporter dye (fluorophore). In the absence of target, the fluorescence is quenched. In the presence of target, the hairpin structure opens upon beacon/target hybridization, resulting in emission of fluorescence. For the detection of GBS amplicons, the molecular beacon contains the fluorophore FAM at the 5' end and the non-fluorescent quencher moiety DABCYL at the opposite end of the oligonucleotide. For the detection of the IC amplicons, the molecular beacon contains the fluorophore TET at the 5' end and the quencher DABCYL at the 3' end. Each beacon-target hybrid fluoresces at a wavelength characteristic of the fluorophore used in the particular molecular beacon. The amount of fluorescence at any given cycle, or following cycling, depends on the amount of specific amplicons present at that time. The SmartCycler® simultaneously monitors the fluorescence emitted by each beacon, interprets all data and at the end of the cycling program provides a final result (see Interpretation of Results).

Reagents

BD GeneOhm™ StrepB	50 tests
Sample Buffer	60 x 1 mL
Tris-EDTA buffer	
Lysis tube	50 tubes
Glass beads	
Master Mix	50 tubes
< 0.001% DNA polymerase complex	
< 0.001% Internal Control - non-infectious DNA containing GBS primers binding sequences and a unique sequence for probe hybridization	
< 0.001% primers	
< 0.002% molecular probes	
< 0.05% dATP, dCTP, dGTP, dTTP	
Bovine serum albumin	
Carbohydrate	
< 0.005% non-infectious genomic DNA from <i>Streptococcus pneumoniae</i>	
Positive Control	50 tubes
< 0.001% DNA polymerase complex	
< 0.001% Internal Control - non-infectious DNA containing GBS primers binding sequences and a unique sequence for probe hybridization	
< 0.001% primers	
< 0.002 % molecular probes	
< 0.05% dATP, dCTP, dGTP, dTTP	
Bovine serum albumin	
Carbohydrate	
< 0.001% non-infectious genomic GBS DNA with the <i>cfb</i> gene	
Negative Control	50 tubes
< 0.001% DNA polymerase complex	

- < 0.001% Internal Control - non-infectious DNA containing GBS primers binding sequences and a unique sequence for probe hybridization
- < 0.001% primers
- < 0.002% molecular probes
- < 0.05% dATP, dCTP, dGTP, dTTP
- Bovine serum albumin
- Carbohydrate
- < 0.005% non-infectious genomic DNA from *Streptococcus pneumoniae*

Diluent**50 x 1 mL**

- Tris-HCl buffer
- MgCl₂
- (NH₄)₂SO₄

Precautions

This test is for *in vitro* diagnostic use only.

- Do not use the kit if the outer carton safety seal is broken.
- Do not use reagents if the protective pouches are open or torn upon arrival.
- Close protective pouches of Master Mix and Controls quickly with the zip seal after each use.
- Do not use reagents if desiccant is not present inside Master Mix and Control pouches.
- Do not remove desiccant from Master Mix and Control pouches.
- Reagents are not interchangeable between lots.

Never pool reagents from different tubes even if they are from the same lot.

- Do not use the reagents after their expiration date.
- Do not interchange caps among reagents as contamination may occur and compromise test results.
- Avoid microbial and deoxyribonuclease (DNAse) contamination of reagents when removing aliquots from tubes. The use of sterile DNAse-free disposable filter-blocked or positive displacement pipettor tips is recommended.
- To avoid contamination of the environment with GBS amplicons, do not open the reaction tubes post-amplification.
- Use a new tip for each specimen or reagent.
- Performing the assay outside the time ranges recommended can produce invalid results. Assays not falling within the specified time ranges should be repeated.
- Additional controls may be tested according to guidelines or requirements of local, state, provincial and/or federal regulations or accrediting organizations.
- In cases where open-tube PCR tests are also conducted by the laboratory, separated and segregated working areas should be used for specimen preparation and amplification / detection activities. Supplies and equipment should be dedicated to each area and should not be moved from one area to another. Gloves must always be worn and must be changed before changing from one area to another or before manipulating lyophilized reagents.
- Always handle specimens as if they are infectious and in accordance with safe laboratory procedures such as those described in *Biosafety in Microbiological and Biomedical Laboratories*⁷ and in the CLSI Document M29⁸.
- Wear protective clothing and disposable gloves while handling kit reagents. Wash hands thoroughly after performing the test.
- Do not pipet by mouth.
- Do not smoke, drink, or eat in areas where specimens or kit reagents are being handled.
- Dispose of unused reagents and waste in accordance with country, federal, provincial, state and local regulations.

Materials provided

- Sample buffer
- Lysis tube
- Master mix
- Positive Control (PC)
- Negative Control (NC)
- Diluent
- Specimen identification labels

Storage, handling and stability

Collected specimen

Specimens should be kept between 2°C and 30°C during transport. Protect against freezing or exposure to excessive heat.

Specimens that can be tested within 24 hours can be kept at room temperature; if not, it is recommended that they be refrigerated. Specimens stored between 2°C and 8°C are stable for up to 6 days.

Reagents

Note: Storage conditions must follow the specifications written on each pouch. Tubes outside of their protective bag and unused within specified time limit should be discarded.

Kit Component		Master mix and Controls (white, red or green labels)	Lysis tube (yellow cap)	Sample buffer, and Diluent (blue cap, and black strip label respectively)
Sealed pouch	Temperature	2-25°C	2-25 °C	2-25 °C
	Stability	Expiration date	Expiration date	Expiration date
Opened pouch ¹	Temperature	2-8 °C	2-25 °C ²	2-25 °C ²
	Stability	1 month ³	Expiration date	2 months ³

¹ Once the original seal on the pouch is broken, carefully close the pouch with the zip seal after each use and store at 2-8°C.

² Although these reagents can be stored at room temperature they should be kept with their accompanying reagents of the same lot at 2-8°C.

³ Provided the bag is properly closed with the zip seal after each use.

Kit Component outside of their protective pouch		Master mix and Controls (white, red or green labels)
Not reconstituted tubes	Temperature	15-25 °C
	Stability	2 hours

Materials required but not provided

- **BBL™ CultureSwab™ Liquid Stuart** (Becton Dickinson catalog no. 220099), **Copan Transystem™ Liquid Stuart** (Copan Italia International catalog no. 141C.USE), **Copan Venturi Transystem™ Liquid Stuart** (Copan Diagnostics Inc. catalog no. 141C.US), **HealthLink TransPorter™ single Liquid Stuart** (HealthLink catalog no. 4432)
- **Vortex** Genie 2 (Fisher) with microtube holder or equivalent; for processing multiple samples, adapters with multiple holding sites can be used

- **Micropipettors** (accurate range between 1-50 µL)
- **Sterile DNase**-free filter-blocked or positive displacement pipettor tips
- **Scissors**
- **Gauzes**
- Disposable **gloves**, powderless
- **Microcentrifuge** for low speed centrifugation
- **Dry heating block** for 1.5 mL tubes or water bath
- Ice or **cooling block** for 1.5 mL tubes
- Stopwatch or **timer**
- **SmartCycler® starter system** with **Dx Software** (processing block, user manual, accessory kit, and desktop computer) (Cepheid, Sunnyvale, CA, USA)

Instructions for use

Specimen collection

In order to obtain an adequate specimen, the procedure for specimen collection must be followed closely.

Using the recommended swab with Liquid Stuart Medium (refer to Materials required but not provided), vaginal-rectal specimens are collected according to the following procedure:

1. **Wipe away excessive amount of secretion or discharge from the vaginal area;**
2. Carefully insert the swab into the **lower one-third part of the vagina**, and sample secretions from the mucosa;
3. Carefully insert the same swab, approximately **2.5 cm beyond the anal sphincter**, and gently rotate to sample anal crypts;
4. Replace the swab in its container;
5. Label the container;
6. Ship the swabs to the laboratory according to hospital standard operating procedures. For specimens sent to a laboratory off-site, the specimen should be protected against exposure to excessive heat.
7. At the laboratory, keep all specimens at room temperature until testing (if within 24 hours from collection);
8. Refer to the section entitled "Storage, handling and stability - Collected specimens" for storage and handling.

Specimen preparation

Note: One Sample Buffer tube (**blue cap**) and one Lysis tube (**yellow cap**) are required **for each specimen** to be tested. Remove the required number of tubes from their protective pouches, **remove the excess air and close the pouches quickly with the zip seal.**

1. **Place the collection device (swab) in a Sample Buffer tube (blue cap).**
Identify the Sample Buffer on the cap and/or the tube label.
2. **Break the swab stem and close the tube tightly.**
Hold the swab by the stem near the rim of the tube (use gauze to minimize risks of contamination). Lift the swab a few millimeters (mm) from the bottom of the tube and bend the stem against the edge of the tube to break it. **Alternative method: use clean scissors to cut the stem.** Make sure the cap will close tightly.
3. **Let stand for 5 minutes.**
4. **Vortex at high speed for 15 seconds.**
For processing multiple samples, adapters with multiple holding sites can be used.
5. **Transfer 50 µL of cell suspension to the Lysis tube (yellow cap); close tightly.**
Use a new micropipettor tip for each specimen.
6. **Vortex at high speed for 5 minutes.**
For processing multiple samples, adapters with multiple holding sites can be used.
7. **Centrifuge the Lysis tube briefly (quick spin).**

At low speed for 2 to 5 seconds; to bring the content to the bottom of the tube.

8. Heat at 95 ± 2°C for 2 minutes.

Use a dry heating block for 1.5 mL tubes or a water bath.

9. Put the Lysis tube on ice or on a cooling block.

BD GeneOhm™ StrepB assay procedure

Note: One Master Mix tube (**SmartCycler® tube with a white label**) is required for **each specimen** to be tested. One Positive Control tube (**red label**) and one Negative Control tube (**green label**) are required per assay run. One Diluent tube (**label with black strip**) is required for the preparation of up to 40 PCR reactions. Remove the required number of tubes from their protective pouch, **remove the excess air, and close the pouch quickly with the zip seal.**

Prepare only enough SmartCycler tubes to fill available I-CORE® modules on the SmartCycler® instrument.

1. Place Master Mix, Positive and Negative Control tubes to be tested on the SmartCycler® cooling block.

Identify the Master Mix tube(s) (white label) on the cap with the specimen identification labels provided with the kit.

2. Open the Master Mix, Positive and Negative Control tubes.

Hold the tube firmly by the ribbed upper portion, and lift the cap with an upward vertical movement. Do not lift horizontally. **Avoid touching the optical detection windows at the bottom edges of the tube and the lower diamond-shaped area.**

3. Add 25 µL of Diluent (black label) to all tubes; partially close the Master Mix tubes.

Deliver the Diluent into the reservoir (upper part) of each tube.

4. Add 1.5 µL of each lysate to a different Master Mix tube; close the tubes tightly.

Take care not to aspirate beads when pipetting into the Lysis tube. After the addition of the specimen, pipette up and down 2-3 times in the reservoir to ensure transfer of the complete volume. Close the Lysis tube and the Master Mix tube. Use a new micropipettor tip for each specimen.

5. Close the Positive and Negative Control tubes.

The Negative Control should be the last sample prepared for a given assay (i.e. all samples run simultaneously on the SmartCycler®).

6. Centrifuge all reaction tubes for 5-10 seconds.

Use the specially adapted microcentrifuge provided with the SmartCycler®.

7. Keep the reaction tubes at 2-8 °C on the SmartCycler® cooling block before loading on the instrument.

The remaining lysates should be frozen at -20 ± 5 °C for later use, if necessary.

8. Just before loading the tubes on the instrument, take the COOLING BLOCK WITH THE TUBES STILL IN PLACE and VORTEX UPSIDE DOWN for 5-10 seconds.

9. Create a run with the BD GeneOhm™ StrepB assay protocol.

Refer to the SmartCycler® Dx Software Operator Manual if needed. You should enter the identification parameters for the specimens before starting the run.

10. Insert each reaction tube in an I-CORE® module of the SmartCycler® and close the I-CORE® lid.

Place the Positive and Negative Controls at their appropriate position (see the section entitled "Quality control"). Press down all the tubes firmly into place.

11. Start the run.

Quality Control

Positive and negative Controls

Quality control procedures are designed to monitor assay performance. The Positive Control is intended to monitor for substantial reagent failure. The Negative Control is used to detect reagent or environmental contamination (or carry-over) by either *S. agalactiae* or GBS amplicons. Positive and Negative Controls are assay controls (run controls). An invalid control invalidates the run. Finally, an internal control incorporated into each reaction mixture is intended to monitor PCR inhibition in each specimen.

One Positive Control and one Negative Control must be processed for each assay run on the SmartCycler®. The software automatically assigns the position of the controls on the instrument (refer to the SmartCycler® Dx Software Operator Manual).

Specimen processing controls

Control strains may be tested according to guidelines or requirements of local, state and/or federal regulations or accrediting organizations. A culture of group B *Streptococcus* (e.g., *S. agalactiae*, American Type Culture Collection, ATCC 12973) or a well characterized clinical isolate of *S. agalactiae* may be used as a specimen processing control while a culture of *Streptococcus bovis* (e.g. ATCC 33317) or of any other non-group B *Streptococcus* may be used as an external negative control.

Transfer 3 colonies of 1-2 mm in size from a fresh sheep blood agar plate into 3 mL of Tryptic Soy Broth and grow to an OD of 0.6 at 600 nm. Prepare serial dilutions in saline to obtain a bacterial suspension of approximately 10⁶ CFU/mL. Dip the recommended swab with Liquid Stuart Medium (refer to Materials required but not provided) into the bacterial suspension, replace the swab into its container (to allow contact with the transport medium), let stand at room temperature for 5 min and then process and test as a clinical specimen (Refer to sections entitled Specimen preparation and BD GeneOhm™ StrepB assay procedure). All specimens and controls should yield valid results (no invalid positive or negative control and no failed internal control).

This procedure may also serve as a quality control procedure for specimen collection devices. In such cases, testing should be done in triplicate according to the package insert instructions (Specimen preparation), including controls. All specimens and controls should yield valid results (no invalid positive or negative control and no failed internal control).

For general QC guidance, the user may wish to refer to CLSI MM3⁹ and C24¹⁰.

Antimicrobial susceptibility testing for penicillin-allergic patients

Note: The Sample Buffer should not affect organism viability. However, use of the buffer in culturing has not been evaluated. A second swab may be obtained from women in order to conduct antimicrobial susceptibility testing. Laboratories may choose to validate other approaches. Procedures for culturing and antimicrobial susceptibility testing are available from the CDC².

Interpretation of results

The decisional algorithm for the BD GeneOhm™ StrepB assay is embedded in the SmartCycler® software. Interpretation of assay results is done according to the following criteria:

Assay result reported	IC result reported	Interpretation of result
NEG	PASS	No GBS DNA detected, presumed negative for GBS or number of organisms may be below the limit of detection of the assay
POS	NA	GBS DNA detected, presumed positive for GBS
Unresolved	FAIL	Unresolved—inhibitory specimen or reagent failure
ND	ND	Not determined due to I-Core® Module failure (with Warning or Error Codes ^A)

IC = Internal Control; NA = not applicable; ND = not determined.

^A Refer to the SmartCycler® Dx Software Operator Manual for the interpretation of warning and error codes.

An invalid positive or negative control invalidates the assay run. In such cases, assay results obtained in that run are invalid and must not be reported. Invalid assay run or instrument error codes or warnings are flagged on-screen and on reports. Before reporting GBS results, always verify that the assay run is valid. Refer to the SmartCycler® Dx Software Operator Manual for printing results.

Invalid assay Run

Using frozen lysate(s), prepare new reaction tubes for all clinical specimens within that assay run along with new control tubes.

Unresolved specimens

Repeat testing with the corresponding specimen frozen lysate. The effect of the freeze-thaw cycle has been shown to reduce PCR inhibitory substances.

Specimen not determined due to I-CORE® module failure

Repeat testing with the corresponding specimen frozen lysate. For the interpretation of warning or error code messages, refer to the SmartCycler® Dx Software Operator Manual.

Limitations of the procedure

- Performance of this test has been established with the SmartCycler® instrument, with vaginal-rectal specimen from antepartum and intrapartum patients collected with Copan Venturi Transystem® with Liquid Stuart Medium. Therefore, this product can only be used with the SmartCycler®; moreover, the use of specimen collection and transport system other than those listed in the Material required but not provided section is not recommended. Other clinical sources have not been assessed and performance characteristics of this test are unknown on other specimen types.
- Negative test results may occur from improper specimen collection, handling or storage, presence of inhibitor, technical error, sample mix-up or because the number of organisms in the specimen is below the analytical sensitivity of the test. Careful compliance to the instructions given in this insert and in the SmartCycler® Dx Software Operator Manual is necessary to avoid erroneous results. Use of this test should be limited to personnel trained on the procedure and on the use of the SmartCycler®.
- Since detection of group B *Streptococcus* is dependent on the number of organisms present in the sample, reliable results are dependent on proper specimen collection, handling and storage.
- BD GeneOhm™ StrepB assay results may sometimes be unresolved or invalidated due to an invalid control, and require retesting that can result in a delay for obtaining the results.
- False negative results may also occur due to the presence of polymerase inhibitory substances; the internal control present in the Master Mix allows the detection of such substances.
- Mutations in primer or probe binding regions may affect detection, especially when organisms are present at less than 10⁴ organisms. There are no published reports of GBS strains or isolates lacking the *cfb* gene. If such a case occurred, BD GeneOhm™ StrepB assay would yield a false negative result. No phenotypically CAMP-negative GBS isolates have been evaluated with BD GeneOhm™ StrepB assay.
- A positive test result does not necessarily indicate the presence of viable organism. It is however presumptive for the presence of group B *Streptococcus*.
- Methods for culturing from the BD GeneOhm™ StrepB Sample Buffer have not been evaluated. Laboratories must validate their own culturing procedures or collect a second specimen to hold in reserve.
- Results from BD GeneOhm™ StrepB assay should be used as an adjunct to clinical observations and other information available to the physician. The test is not intended to differentiate carriers of group B *Streptococcus* from those with streptococcal infection. Test results may also be affected by concurrent antibiotic therapy, therefore therapeutic success or failure cannot be assessed using this test since DNA may persist following antimicrobial therapy.
- Although there is no need for reagent preparation and the main technical operations include pipetting, good laboratory technique is essential to the proper performance of this assay. Due to the high analytical sensitivity of this test, extreme care should be taken to preserve the purity of all reagents, especially in cases where multiple aliquots are taken from a single tube.
- Good laboratory practices and the use of gloves are recommended to avoid contamination of specimens or reagents.

Interfering substances

Potentially interfering substances include, but are not limited to the following: amniotic fluid, meconium, blood, feces, lubricant, urine, and vaginal mucosa. The presence of excessive vaginal secretions (antepartum or intrapartum) may inhibit PCR and give unresolved results. Therefore, it is recommended that the vaginal area be wiped away prior to sampling. The presence of excessive blood may also lead to PCR inhibition and an unresolved result.

In an investigational study involving 803 vaginal/rectal specimens from intrapartum maternity patients, potentially interfering substances mentioned above were reported for 63% of collected specimens. Only 1% (10/803) of the specimens yielded unresolved results. Of those, there was no potentially interfering

substance reported for 6 specimens, amniotic fluid for 1 and a combination of substances for 3 of them. All but one were resolved following a freeze-thaw cycle of the specimen. For the specimen that could not be resolved, no potentially interfering substance was observed.

Expected values

Approximately 10% to 30% of pregnant women are colonized with GBS in the vagina or rectum². GBS colonization can be transient, chronic or intermittent. Culture screening of both the vagina and rectum for GBS late in gestation during prenatal care can detect women who are likely to be colonized with GBS at the time of delivery. In the investigational study for BD GeneOhm™ StrepB, the overall colonization rate at the time of delivery determined with the culture technique was 18.6% with a range from 9.1% to 28.7% at the different clinical sites. With BD GeneOhm™ StrepB, overall, 20.1% of intrapartum women were positive. In addition to intrapartum screening, a subset of women (674) also had preterm culture screening. Of these, 11.7% were positive during antepartum testing.

Overall, 10.2% of these women had antepartum culture results different than intrapartum culture results and 11.6% had antepartum culture results different than BD GeneOhm™ StrepB intrapartum results. In other studies, sensitivities of 87%¹¹ (83-92% CI) and 69%¹² (57-79% CI) and specificities of 96%¹¹ (95-98% CI) and 92%¹² (89-94% CI) have been reported for late-prenatal cultures for identifying colonization status at delivery.

Performance characteristics

Clinical performances

Performance characteristics of BD GeneOhm™ StrepB assay were determined in a multi-site prospective investigational study: five institutions with maternity services, 2 in Canada and 3 in the U.S. Each institution had in place a culture-based screening program. Testing was done in clinical laboratories affiliated with each institution. To be enrolled in the study, women had to provide written consent, be in labor and have no contraindication to vaginal examination (e.g. bleeding). There was also no evidence of placenta previa, no urgent indication to proceed to delivery and no antibiotic used in the week prior to admission. The study did not account for times to delivery for the patients.

The method of reference used was the culture technique recommended by the Centers for Disease Control and prevention, that is, microbiological culture in selective broth medium (Todd-Hewitt broth medium supplemented with 15 µg/mL of nalidixic acid, and 10 µg/mL of colistin, or with 8 µg/mL of gentamicin and 15 µg/mL of nalidixic acid: other commercially available media include SBM broth or Lim Broth), followed by overnight incubation and subculture onto solid blood agar medium. Specific identification of colonies suggestive of GBS was done with slide agglutination tests.

Eight hundred and eighty one (881) women consented to participate in the study. Of those, 78 were excluded because the woman had either changed her mind or was not in labor or because of deviations to the protocol. Therefore, the performance characteristics of BD GeneOhm™ StrepB were determined from results of 803 intrapartum maternity patients.

Two vaginal/rectal specimens from each intrapartum maternity patient were collected with the recommended swab with Liquid Stuart Medium (refer to Materials required but not provided) using the CDC recommended procedure². One swab was tested with the CDC recommended culture technique and one swab with the BD GeneOhm™ StrepB assay. The results of the study are presented in Tables 1 to 4.

Table 1. Results obtained with BD GeneOhm™ StrepB assay in reference to the culture technique.

		BD GeneOhm™ StrepB assay		
		Positive	Negative	Total
Culture technique	Positive	140 ^A	9	149
	Negative	27 ^B	626 ^C	653
Total		167	635	802 ^D

^A Fourteen (14) specimens were initially culture negative but upon investigation were found to be culture positive; 1 of the 14 had initially tested unresolved but upon re-testing gave a positive result; 3 specimens that were initially BD GeneOhm™ StrepB assay positive were retested because of invalid controls (positive and negative) and tested positive.

^B One (1) specimen that was initially positive with BD GeneOhm™ StrepB assay was retested because of an invalid control and retested positive.

^C Twelve (12) specimens that were initially negative with BD GeneOhm™ StrepB assay were retested because of an invalid control and all tested negative; 8 specimens that initially gave an unresolved result gave a negative result upon retesting.

^D One (1) specimen that gave an initially unresolved result remained unresolved upon retesting and was not included in the table above.

Investigation of specimens that yielded culture positive/BD GeneOhm™ StrepB negative results (n=9) revealed that all specimens had a bacterial load below the detection limit of BD GeneOhm™ StrepB. It was also shown that all isolates recovered from the culture positive plates yielded a positive result with BD GeneOhm™ StrepB.

Investigation of specimens that yielded culture negative/ BD GeneOhm™ StrepB positive results (n=27) revealed that all specimens did contain the proper target region (including primer binding regions) of the *cfb* gene and that the amplified product was that expected with BD GeneOhm™ StrepB.

The overall prevalence of the study population was 18.6%, with a range from 9.1% to 28.7% at the different clinical sites. The negative predictive value was 98.6% (95% CI, 97.3% - 99.3%) and the positive predictive value was 83.8% (95% CI, 77.4% - 89.1%).

Table 2. Clinical performance of BD GeneOhm™ StrepB assay in reference to the CDC recommended culture technique.

	Clinical Sensitivity	Clinical Specificity	No. of unresolved specimens	Invalid/ total no. of runs
Site 1	93% (n=15) (68%-100%) ^A	93% (n=150) (87%-96%) ^A	5	2/55
Site 2	88% (n=25) (69%-97%)	100% (n=62) (97%-100%)	0	0/38
Site 3	99% (n=77) (93-100%)	97% (n=350) (94%-98%)	1	1/56
Site 4	85% (n=13) (54%-98%)	100% (n=35) (98%-100%)	1	1/21
Site 5	89% (n=19) (67%-99%)	93% (n=56) (82%-98%)	3	1/22
Total :	94% (n=149) (89%-97%)	96% (n=653) (94% - 97%)	10 ^B	5/192

^A Binomial 95% confidence intervals.

^B All specimens were unresolved due to failed internal controls indicative of inhibition or reagent failure. Nine (9) of the 10 were resolved upon re-testing.

Table 3. Stratification of results according to the time elapsed between specimen collection and testing with BD GeneOhm™ StrepB assay.

Elapsed time (hours)	Culture positive		Culture negative		% agreement (95% CI)
	BD GeneOhm™ positive	BD GeneOhm™ negative	BD GeneOhm™ negative	BD GeneOhm™ positive	
0 to 4	11	3	88	4	93% (87-97%)
4 to 8	32	3	220	10	95% (92-97%)
8 to 12	9	0	29	2	95% (83-99%)
12 to 24	41	2	143	4	97% (93-99%)
24 to 48	35	0	127	7	96% (92-98%)
> 48	12	1	19	0	97% (84-100%)
Total	140	9 ^A	626	27 ^B	96% (94-97%)

^A All specimens had a bacterial load below the detection limit of BD GeneOhm™ StrepB.

^B All specimens did contain the proper target region (including primer binding regions) of the *cfb* gene and that the amplified product was that expected with BD GeneOhm™ StrepB.

Membrane rupture did not affect performances of BD GeneOhm™ StrepB assay. Using the presence or absence of amniotic fluid and/or meconium in collected specimens as an indication of membrane rupture, the sensitivity and specificity of BD GeneOhm™ StrepB assay are the same used antepartum or intrapartum (Table 4).

Table 4. Performance of BD GeneOhm™ StrepB assay antepartum and intrapartum.

	Sensitivity	Specificity	% agreement
Antepartum	95% (53/56) (85-99%) ^A	96% (301/315) (93-98%)	95% (354/371) (93-97%)
Intrapartum	94% (87/93) (86-98%)	96% (325/338) (94-98%)	96% (412/431) (93-97%)

^A Binomial 95% confidence intervals.

Analytical specificity

Genomic DNA from 99 ATCC strains representing 27 species of streptococci, other species phylogenetically related to *S. agalactiae*, other bacteria and yeasts commonly found in vaginal and rectal flora and human DNA were tested. Among those were 9 organisms that had not been reported to carry the *cfb* gene but elicit CAMP-like activity. For microbial DNA, 1.5 ng (2×10^5 equivalent genome copies per PCR reaction or 10^8 equivalent genome copies/mL) were used. For human DNA, 75 to 233 ng (up to 1.4×10^5 genome copies per reaction or 10^8 copies/mL) were used. The specificity was 100%.

Analytical sensitivity

The analytical sensitivity (Limit of Detection or LOD) of BD GeneOhm™ StrepB assay was determined with 12 strains of *S. agalactiae* representing 11 known serotypes or variants. Quantitated culture and purified genomic DNA diluted in the Sample Buffer of BD GeneOhm™ StrepB assay were tested in 5 replicates. The LOD is defined as the smallest concentration at which all replicates test positive.

The LOD of BD GeneOhm™ StrepB assay for *S. agalactiae*, serotypes and variants Ia, Ib/c, II, IIc, IIR, III, IIIc, IIIR, V and VR, ranges between 10 and 50 genome copies per reaction with a median of 25 genome copies per reaction. The LOD in CFU is 3 to 9 CFU/reaction. Taking into account the dilution factor due to specimen processing, this translates into 10^3 to 10^4 CFU/swab. The following table details the lowest concentration of each subtype giving a positive result in all 5 replicates.

Table 5. LOD obtained for each serotype tested.

Serotype	Genome copies /reaction	CFU /reaction
ATCC 12973 (II)	10	3
IIIR	25	6
Ia/c	50	4
VR	50	3
III	25	3
Ia	25	3
V	50	8
Ib/c	25	3
IIR	25	4
II	25	9
IIIc	25	3
IIc ^A	50	4

^A For serotype IIc, a mismatch exists in one of the primer binding region.

Reproducibility

A panel of 10 simulated specimens (R1 to R10) with varying concentrations of GBS and the two controls (positive and negative) supplied with BD GeneOhm™ StrepB assay were tested in triplicate on three different days at each of 3 sites (10 specimens plus 2 controls tested X 3 X 3 days X 3 sites). One lot of reagent was used for the study.

Table 6. Summary of reproducibility results.

Specimen ID ^A	Site 1	Site 2	Site 3	Total agreement	Total % agreement
Negative	9/9	9/9	9/9	27/27	100%
Negative	9/9	9/9	9/9	27/27	100%
Weak positive	3/9	7/9	8/9	18/27	66%
Weak positive	8/9	9/9	9/9	26/27	96%
Positive	8/9	9/9 ^B	9/9	26/27	96%
Positive	9/9 ^A	9/9 ^B	9/9	27/27	100%
Positive	9/9 ^A	9/9 ^B	9/9	27/27	100%
Strong positive	9/9	9/9	9/9	27/27	100%
Strong positive	9/9	9/9	9/9	27/27	100%
Strong positive	9/9	9/9	9/9	27/27	100%
Pos Control	3/3 ^C	9/9	9/9	21/21	100%
Neg Control	3/3 ^C	9/9	9/9	21/21	100%
Total agreement	88/96	106/108	107/108	301/312	96.5%
% agreement	91.6 %	98.1 %	99.1 %	96.5%	

^A All positive specimens have varying amounts of GBS.

^B Specimen initially unresolved but upon retesting gave the expected result.

^C All replicates we tested as a single run instead of three separate runs of 10 specimens and 2 controls.

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References / Références / Referenzen / Referencias / Riferimenti

- ¹ Schrag SJ, Zywicki S, Farley MM, Reingold AL, Harrison LH, Lefkowitz LB, Hadler JL, Danila R, Cieslak PR, Schuchat A. Group B streptococcal disease in the era of intrapartum antibiotic prophylaxis. *N Engl J Med* 2000; 342:15-20.
- ² Centers for Disease Control and Prevention. Prevention of perinatal group B streptococcal disease: Revised guidelines from CDC. *MMWR* 2002; 51 (RR-11):1-22.
- ³ Watt JP, Schuchat A, Erickson K, Honig JE, Gibbs R, Schulkin J. Group B streptococcal disease prevention practices of obstetrician-gynecologists. *Obstet Gynecol* 2001; 98:7-13.
- ⁴ Hager WD, Schuchat A, Gibbs R, Sweet R, Mead P and Larsen WL. Clinical commentary. Prevention of perinatal group B streptococcal infection: Current controversies. *Obstet Gynecol* 2000; 96:141-45.
- ⁵ Wilkinson HW. CAMP-disk test for presumptive identification of group B streptococci. *J Clin Microbiol* 1977; 6:42-5.
- ⁶ Podbielski A, Blankenstein O, Luticken R. Molecular characterization of the *cfb* gene encoding group B streptococcal CAMP-factor. *Med Microbiol Immunol* 1994; 183:239-56.
- ⁷ Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories. Richmond JY and McKinney RW (eds) (1993). HHS Publication number (CDC) 93-8395.
- ⁸ Clinical and Laboratory Standards Institute. Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Document M29 (Refer to latest edition).
- ⁹ Clinical and Laboratory Standards Institute. Molecular Diagnostic Methods for Infectious Diseases; Approved Guideline. CLSI document MM3-A2 (Refer to the latest edition).
- ¹⁰ Clinical and Laboratory Standards Institute. Statistical quality control for quantitative measurements: Principles and definitions; Approved Guideline. Document C24 (Refer to latest edition).
- ¹¹ Yancey MK, Schuchat A, Brown LK, Ventura VL, Markenson GR. The accuracy of late antenatal screening cultures in predicting genital group B streptococcal colonization at delivery. *Obstet Gynecol* 1996; 88:811-15.
- ¹² Melin P, Schmitz M, Tsobo C, Hayette MP, De Mol P. Rapid intrapartum test (StrepB OIA®) and prenatal cultures for identification of Group B streptococcal (GBS) carriers at delivery: a prospective study. Interscience Conference on Antimicrobial Agents and Chemotherapy. Toronto, 2000. Poster # 357.

This kit is sold under license from the Public Health Research Institute of the City of New York, Inc. and may be used under the PHRI patent rights only for human *in vitro* clinical diagnostics.

The purchase of this product allows the purchaser to use it for amplification and detection of nucleic acid sequences for providing human *in vitro* diagnostics. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby.

With sale of any kits (label license in package insert): "The purchase of this product includes a limited, non-transferable license under U.S. Patent No. 6,787,338, owned by the University of Utah Research Foundation and licensed to Idaho Technology, Inc., to use only the enclosed amount of product according to the specified protocols. No right is conveyed, expressly, by implication, or by estoppel, to use any instrument or system under any claim of U.S. Patent No. 6,787,338, other than for the amount of product contained herein."

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BD and BD logo are trademarks of Becton, Dickinson and Company. © 2014 BD
GeneOhm Sciences Canada, Inc.
Benex Limited
Becton Dickinson Pty Ltd.