
Revisions

Rev from	Rev to	ECO #
G	H	3538-05

Notes:

- BD Cat. Number 442018, 442121
- Blank (Sheet) Size : Length: N/A Width: N/A
 Number of Pages: N/A Number of Sheets: N/A
 Page Size: Length N/A Width N/A Final Folded Size: N/A
- Style (see illustrations below): N/A

- See Specification Control Number VS0250 for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS # Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: PP043JAA		Category and Description Package Insert, BACTEC 13A Mycobacteria Culture Vials	Sheet: 1 of 20 Scale: 1 : 1	A

BD BACTEC™ 13A Mycobacteria Culture Vials

Middlebrook 7H13

English: pages 1 – 4 Italiano: pagine 9 – 12
 Français : pages 4 – 6 Portuguais: páginas 12 – 14
 Deutsch: Seiten 6 – 9 Español: páginas 15 – 17

 PP043JAA(H)
2005/05

See symbol glossary at end of insert. / Viz popis symbolů na konci příbalového letáku. / Se symbolglossaret i slutningen af indlægssedlen. / Zie lijst met symbolen aan het einde van de bijsluiters. / Vaadake sūbolite seletust infolehe lõpus. / Katsõ pakkausselosteen lopussa olevaa kuvamerkkien sanastoa. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / A jelmagyarázat a használati utasítás végén található. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Zr. informacino lapelo pabaigoje patelkiamą simbolių glosarijų. / Se i symbolforklaringen på slutten av produktvedlegget. / Zobacz objaśnienie symboli na końcu ulotki. / Consulte o glossário de símbolos no fim do folheto informativo. / Pozri slovník symbolov na konci letáka. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Pokyny vám poskytné místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőtől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Instrukcie ziskate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar.

INTENDED USE

BACTEC™ 13A Mycobacteria Medium is for mycobacteria cultures. Principal use is with the **BACTEC 460TB System**.

A qualitative procedure for the culture and recovery of microorganisms (mainly mycobacteria) from blood. The primary objective of this medium is to detect mycobacteria in patients with depressed immune systems (e.g., Acquired Immunodeficiency Syndrome [AIDS]).

SUMMARY AND EXPLANATION

The **BACTEC Radiometric Technique** has been widely used for the rapid recovery of mycobacteria from sputum and other clinical specimens.¹⁻⁴ In this technique, ¹⁴C₂ produced by mycobacteria from the metabolism of a ¹⁴C-labeled substrate, present in **BACTEC 7H12** brand medium, is quantitatively measured. Several investigators have successfully used the **BACTEC 460TB** system for isolation of mycobacteria from blood specimens from AIDS patients.⁵⁻⁸ However, processing of blood is required to lyse the cells and concentrate the specimen because only 1 mL of specimen can be inoculated into **BACTEC 7H12** medium. The use of **BACTEC 13A Mycobacteria Medium** eliminates the time-consuming and potentially hazardous processing step by allowing a larger volume (up to 5 mL) of blood to be inoculated directly into the 30 mL vial of **BACTEC 13A Mycobacteria Medium**. The medium is supplemented with SPS and other additives to prevent clotting and enhance the growth of mycobacteria, especially *Mycobacterium avium* complex.

The sample to be tested is inoculated into one or more vials and incubated. The culture vial is periodically inserted into the **BACTEC 460TB** instrument for testing, which consists of aspiration of the supernatant gas and assay of its radioactive content. A positive reading indicates the presence of viable microorganisms in the vial. It is limited to detection of microorganisms that will grow in a particular type of culture vial.

PRINCIPLES OF THE PROCEDURE

BACTEC 13A Mycobacteria Medium is a modified **BACTEC 7H12 Medium**, which when supplemented with **BACTEC Enrichment** supports growth of mycobacteria. If microorganisms are present in the test sample inoculated into the **BACTEC** vial, ¹⁴C₂ will be produced and liberated into the vial head space when the organisms metabolize the substrates present in the vial. The instrument analyzes the vial head space gas for radioactivity and, if a threshold level is exceeded, indicates that the vial is positive; i.e., that the sample contained viable organisms.

REAGENTS

The **BACTEC** culture vials contain the following reactive ingredients prior to processing:

BACTEC 13A Mycobacteria Medium, 30 mL each

Amount added per Vial of Processed Water

7H9 Broth Base	0.47% w/v
Casein hydrolysate	0.1% w/v
Sodium polyanetholesulfonate (SPS)	0.025% w/v
Polysorbate 80	0.02% w/v
Catalase	1440 units
¹⁴ C Substrate	5 µCi (185 kBq)

BACTEC Enrichment, 10 mL each

Bovine Serum Albumin	15.0% w/v
----------------------	-----------

All **BACTEC** media are dispensed with added CO₂. Composition may have been adjusted to meet specific performance requirements.

Warning and Precautions:

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

Pathogenic microorganisms including, hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁹⁻¹² and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Prior to use, each vial should be examined for evidence of contamination such as cloudiness, bulging or depressed stopper, or leakage. DO NOT USE any vial showing evidence of contamination. A contaminated vial could contain positive pressure. Vial contamination may not be readily apparent.

Prior to use, the user should examine the vials for evidence of damage or deterioration. Vials displaying turbidity, contamination, or discoloration (darkening) should not be used. On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions a vial may not be sealed sufficiently. In both cases, the contents of the vials may leak or spill, especially if the vial is inverted. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: Before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling should be performed in a biological cabinet if possible, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or securely fastened Luer-Lok™ brand tips.

Storage Instructions

BACTEC 13A Mycobacteria vials and the BACTEC 13A Enrichment should be stored in a cool, dry place at 2° – 25°C.

The medium and the Enrichment are stable at the recommended storage conditions until the labeled expiration date.

Specimen Collection

The blood specimen must be collected using sterile techniques to reduce the chances of contamination. The recommended blood volume is 5 mL. It is recommended that the specimen be inoculated into the **BACTEC** vials at bedside. Most commonly, a 10cc or 20cc syringe with a Luer-Lok brand tip is used to draw the sample. **The inoculated BACTEC vial should be transported to the laboratory as quickly as possible.**

PROCEDURE

Material Provided:

BACTEC 13A Mycobacteria Medium

Materials Required But Not Provided:

BACTEC 12B Medium

BACTEC 460TB Instrument

10-20cc syringe

Tuberculin syringe with permanently attached needle

Disinfectant solution

Alcohol

Remove the flip-off cap from the **BACTEC** vial top and inspect the vial for cracks, contamination, excessive cloudiness, and bulging or indented septum. DO NOT USE if any defect is noted. Before inoculating, swab the septum with alcohol (*iodine* is not recommended). Aseptically inject 5 mL of specimen per vial.

BACTEC Enrichment is provided in separate vials and must be added to each **BACTEC 13A Mycobacteria Medium** vial in order to achieve optimum growth of mycobacteria. Aseptically add 0.5 mL of the Enrichment to each 13A medium vial. This addition may be done prior to specimen inoculation, as the supplemented medium is stable at room temperature for several months. If specimen inoculation is done outside the laboratory or at bedside, the Enrichment must be added after receiving the inoculated vials in the laboratory using a separate syringe for each inoculation. A uniform policy of Enrichment addition should be adopted to ensure the addition.

Inoculated vials should be incubated at 37° ± 1°C without shaking and should be tested on a **BACTEC 460TB** instrument using 5% – 10% CO₂ in air. Test the vials every two to three days for the first two weeks and weekly thereafter, for a total of six weeks.

A Growth Index (GI) reading of 20 or more is considered positive. Smears should be made from all positive vials. Stain smears with the Ziehl-Neelsen or fluorochrome method. If the smear is negative, continue testing the vial daily and repeat smears. In addition, subculture to blood or chocolate agar plates may be performed to investigate the presence of non-mycobacterial growth.

If a smear from a **BACTEC 13A Mycobacteria Medium** vial is positive for AFB, report as culture positive and discontinue testing on the **BACTEC 460TB** instrument. Subculture to **BACTEC 12B (7H12) Medium** and solid medium. Incubate subcultures at 37° ± 1°C. Test **BACTEC 12B** vials on the **BACTEC 460TB** instrument daily.

TB Differentiation: When a subculture 12B vial GI reaches about 50, perform the **BACTEC NAP** test (Cat. #442103) following the recommended procedure. This will differentiate TB from mycobacteria other than tuberculosis (MOTT). The NAP test cannot be performed directly from the 13A vial because of the presence of blood in the medium.

Speciation of mycobacteria should be carried out by the conventional techniques as recommended by Centers for Disease Control and Prevention.¹³

Subculturing: Prior to sampling or subculturing or staining, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, insert a sterile needle with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

Drug Susceptibility Testing: Any drug susceptibility test for *M. tuberculosis* may be set up at GI 300 + 1 additional day of incubation or GI ≥ 500 (see **BACTEC 460TB** system Product and Procedure Manual for details).

Quality Control

The user should examine the culture vials for evidence of deterioration prior to use. Vials displaying turbidity, contamination, or discoloration (darkening) should not be used.

Quality Control Certificates are provided with each carton of media. Quality control certificates show test organisms, including ATCC™ cultures specified in the NCCLS standard, *Quality Assurance for Commercially Prepared Microbiological Culture Media*.

As part of the daily maintenance schedule, Performance Test Vials (BACTEC Performance Test Kit, Catalog No. 440534) should be tested on the BACTEC instrument. If low readings are obtained, yet Performance Test Vials were tested properly, possible instrument problems may be indicated (see BACTEC PTK package insert PP-046).

Do not use vials past their expiration date.

Interpretation of Results

The instrument reads out directly on a Growth Index scale of 0 – 999. The Growth Index is actually a measure of metabolic activity of organisms in the medium. When the GI reaches 20 or more, the culture is considered presumptively positive. The AFB smear confirms the presence of mycobacteria and can be made as soon as the vial shows a positive GI value. When AFB are observed on the smear made from a positive vial, the specimen may be reported as culture positive for AFB (ID pending). Negative results should be reported after six weeks.

GI numbers below 20 are representative of background readings and indicate negative cultures. However, when a GI \geq 20 is reached, the number of AFB in the medium is high enough to visualize the organism on smear. If the smear is negative, positive vials should be tested daily and smears repeated until AFB are seen. Subculture to a solid medium from a positive BACTEC 13A Mycobacteria vial will provide colonies for identification tests.

LIMITATIONS OF THE PROCEDURE

Contamination

Care must be taken to prevent contamination of the sample during collection and inoculation into the BACTEC vial. A contaminated sample will give a positive reading but will not indicate a relevant clinical sample. Such a determination must be made by the user based on such factors as type of organism recovered, occurrence of the same organism in multiple cultures, patient history, appropriate equipment checks, etc. Contamination may also result from improper maintenance of BACTEC needles and tubing, poor septum cleansing or malfunction of the BACTEC needle heating unit (see Instrument Operation and Maintenance Manual).

The GI level reflects the rate and amount of growth occurring in the medium but cannot be as accurately quantified as by colony counts on solid media. Colonial morphology and pigmentation characteristics can only be determined on solid media.

PERFORMANCE CHARACTERISTICS

In a study evaluating the recovery and detection time of mycobacteria conducted by Strand et al., 13A Mycobacteria Medium was compared to a manual lysis-centrifugation method utilizing 7H11 medium. The recovery rate for 13A Mycobacteria Medium was 96.9% and for the manual method the recovery rate was 98.5%. The mean times to recovery were essentially the same (23.6 and 23.8 days, respectively). Thus, the direct-inoculation technique using 13A Mycobacteria Medium showed performance equivalence to the lysis-centrifugation method.¹⁴

AVAILABILITY

Cat. No.	Description
442018	BACTEC™ 13A Mycobacteria Medium 50 vials Medium, 30 mL each and 4 vials Enrichment, 10 mL each
442121	BACTEC™ 13A Mycobacteria Medium 12 vials Medium, 30 mL each and 2 vials Enrichment, 10 mL each
442004	BACTEC™ 12B Medium, 4 mL, 100 per carton
442103	BACTEC™ NAP Differentiation Kit, 5 µg NAP per vial, 10 vials per carton
442102	BACTEC™ SIRE Drug Kit, 2 vials lyophilized Streptomycin, 2 vials lyophilized Isoniazid, 2 vials lyophilized Rifampin, 2 vials lyophilized Ethambutol
442104	BACTEC™ Diluting Fluid, 9.9 mL, 10 vials per carton
442146	BACTEC™ Isoniazid, 1 vials lyophilized Isoniazid

REFERENCES

- Middlebrook, G. et al. Automated radiometric detection of growth of *Mycobacterium tuberculosis* in selective media. *Am. Rev. Respir. Dis.* 115: 1066–1069, 1977.
- Roberts, G.D. et al. Evaluation of the BACTEC radiometric method for recovery of mycobacteria and drug susceptibility testing of *Mycobacterium tuberculosis* from acid-fast smear-positive specimens. *J. Clin. Microbiol.* 18: 689–696, 1983.
- Fadda, G. and S. Roe. Recovery and susceptibility testing of *Mycobacterium tuberculosis* from extrapulmonary specimens by the BACTEC radiometric method. *J. Clin. Microbiol.* 19: 720–721, 1984.
- Kirihara, J. M. et al. Improved detection times for *Mycobacterium avium* complex and *Mycobacterium tuberculosis* with the BACTEC Mycobacterium radiometric system. *J. Clin. Microbiol.* 22: 841–845, 1985.
- Macher, A. M. et al. Bacteremia due to *Mycobacterium avium-intracellulare* in the Acquired Immunodeficiency Syndrome. *Annal. Int. Med.* 99: 782–785, 1983.
- Gill, V.J. and F. G. Witebsky. Blood cultures for mycobacteria. *J. Med. Technology*, 1:8: 611–614, 1984.
- Gill, V.J. et al. Use of lysis-centrifugation (isolator) and radiometric (BACTEC) blood culture systems for the detection of mycobacteremia. *J. Clin. Microbiol.* 22: 543–546, 1985.
- Kiehn, T.E. et al. Use of BACTEC radiometric 12A broth for the recovery of clinically significant *Mycobacterium avium* complex (MAC) from patients with AIDS. Abstract C-47. Annual Meeting, ASM, Washington, D.C. 1986.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 4th ed. U.S. Government Printing Office, Washington, D.C.

12. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
13. Kent, P. J., et al. Public Health Mycobacteriology: A Guide for the Level III Laboratory. Centers for Disease Control, Division of Laboratory Training and Consultation, Atlanta, GA, 1985.
14. Strand, C.L., C. Epstein, S. Verzosa, E. Effatt, P. Hormozi and S. Siddiqi, 1989. Evaluation of a new blood culture medium for Mycobacteria. J. Clin. Pathol. 97:316-318.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BACTEC, and Luer-Lok are trademarks of Becton, Dickinson and Company. ©2005 BD.
BENEX Limited