
Revisions

Rev from	Rev to	ECO #
0607	2010/07	5394-10

Notes:

- BD Cat. Number 440534
- Blank (Sheet) Size : Length: N/A Width: N/A
 Number of Pages: N/A Number of Sheets: N/A
 Page Size: Length N/A Width N/A Final Folded Size: N/A
- Style (see illustrations below): N/A

- See Specification Control Number VS-0250 for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 199 Red
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: PP046JAA		Category and Description Package Insert, BACTEC PTK Performance Test Kit	Sheet: 1 of 15 <hr/> Scale: 1:1	A

BD BACTEC™ PTK Performance Test Kit

English:	pages	1 – 2	Italiano:	pagine	7 – 8
Français :	pages	3 – 4	Português:	páginas	9 – 10
Deutsch:	Seiten	5 – 6	Español:	páginas	11 – 12

PP046JAA
2010/07

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзініздің жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

To provide an overall performance and calibration test for the BACTEC™ 460TB instrument, and to set forth appropriate radiation safety precautions and to suggest instructions for the handling, use, storage and disposal of the radioactive material.

SUMMARY AND EXPLANATION

The Performance Test Kit provides reagents to evaluate the overall performance and calibration of the BACTEC 460TB instrument. Successful performance of this test assures that the electronics and displays are correctly calibrated and operating properly. Performance of this test verifies that the instrument needles are not clogged and there are no air leaks. The test should be performed daily and the result recorded as part of the Quality Assurance program for the BACTEC 460TB system.

PRINCIPLES OF THE PROCEDURE

A standard amount of ¹⁴C-bicarbonate solution is injected into a vial containing acid. This liberates a predetermined quantity of ¹⁴CO₂ which, when measured on a BACTEC brand 460TB instrument, will give a standard reading if the instrument is operating correctly. Instrument malfunctions such as air leaks, clogged needles or electronic faults will result in a non-standard reading.

REAGENTS

The following are included in the BACTEC PTK Performance Test Kit:

11A ¹⁴ C Solution	1 vial (20 mL)	1 μCi (37 kBq) of ¹⁴ C radioactive material in bicarbonate solution
11B Test acid	3 vial (10 mL)	2N hydrochloric acid solution

Precautions

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

Prior to use, each vial should be examined for evidence of contamination such as cloudiness, bulging or depressed stopper, or leakage. **DO NOT USE** any vial showing evidence of contamination.

Prior to use, the user should examine the vials for evidence of damage or deterioration. On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions, a vial may not be sealed sufficiently. In both cases the contents of the vials may leak or spill, especially if the vial is inverted.

The Test Acid vials contain strong acid. Flush with water immediately and seek medical assistance if it splashes into mouth or eyes.

To minimize the potential of leakage during inoculation into vials, use syringes with permanently attached needles or Luer-Lok™ brand tips.

HANDLING

Exercise care in handling the glass vials to prevent breakage and also damage to the seal. Do not drop the vials; they are rugged but will break.

STORAGE

These vials should be stored in a cool, dry place (2 – 25°C). Note the "Use by" date on the outside of each package and use the oldest stock first. It is good practice to segregate radioactive materials by storage in a cabinet accessible only to authorized personnel. The cabinet should be labeled with the words, "Caution: Radioactive Materials."

BROKEN VIALS

If a vial is accidentally broken, flood with disinfectant if necessary and then pick up any large pieces of glass and the metal seal with plastic gloved hands using tongs or tweezers to avoid being cut by broken glass. With

disposable plastic gloves on your hands, use a pad of paper towels to mop up the liquid and to pick up the small pieces of glass. When all the free liquid and glass fragments have been removed, use more paper towels and a solution of detergent in hot water to clean the affected area. Rinse with clean water and dry. Put all pieces of glass, used paper towels and gloves into the same receptacle used for the disposal of used **BACTEC** vials.

PROCEDURE

Materials Provided

- 1 vial of 11A ¹⁴C Solution
- 3 vials of 11B Test Acid

Materials Required But Not Provided

BACTEC 460TB Instrument

- 1 mL syringe with permanently attached needle

Inject 0.2 mL of the ¹⁴C solution into one Test Acid vial. Shake vigorously for at least 10 s. Test on the **BACTEC** 460TB instrument. Tested vial should give a reading of approximately 45 – 65 Growth Index units. On the very first test of a new Test Acid vial, the reading will be 40 – 50. Each Test Acid vial can be used for approximately 30 tests before discarding. When the ¹⁴C Solution is exhausted, discard the kit and order a replacement.

EXPECTED VALUES AND PERFORMANCE CHARACTERISTICS

If the reading obtained is out of range, repeat the test being careful to accurately measure the ¹⁴C solution. The most probable causes of low readings are blocked needles and loose fittings or hose connections. Replace the needles, tighten any loose connections and repeat the test. A persistently low reading may be due to an electronic malfunction.

The most probable cause of a high reading is excess ¹⁴C Solution in the Test Acid vial. Repeat the test. If still high, it may be due to an electronic malfunction.

- CAUTION: RADIOACTIVE MATERIAL – NOT FOR INTERNAL OR EXTERNAL USE IN HUMANS OR ANIMALS.
- THE CONTENTS OF THESE **BACTEC** TEST VIALS ARE EXEMPT FROM NRC OR AGREEMENT STATE LICENSING REQUIREMENTS.
- INTRODUCTION INTO FOODS, BEVERAGES, COSMETICS, DRUGS, MEDICINALS OR INTO PRODUCTS MANUFACTURED FOR COMMERCIAL DISTRIBUTION IS PROHIBITED.
- EXEMPT QUANTITIES SHOULD NOT BE COMBINED.
- THE CONTENTS OF THESE **BACTEC** TEST VIALS DO NOT EXCEED 10 MICROCURIES OF CARBON-14.

NOTICE TO GENERAL LICENSEES

The radioactive material may be received, acquired, processed and used only by physicians, clinical laboratories or hospitals and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and a general license of the U.S. Nuclear Regulatory Commission or of the State with which the Commission has entered into an agreement for the exercise of regulatory authority.

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Benex Limited