
Revisions

SO 0191-5

Rev from	Rev to	ECO #
0607	2010/07	5394-10

Notes:

1. BD Cat. Number 442143
2. Blank (Sheet) Size: Length: N/A Width: N/A
 Number of Pages: N/A Number of Sheets: N/A
 Page Size: Length N/A Width N/A Final Folded Size: N/A
3. Style (see illustrations below): N/A

4. See Specification Control Number VS-0250 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 363
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p> <p style="font-weight: bold; font-size: large;">BD</p> <p>Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA</p>	
Proofer	Date		
Checked By	Date		
Part Number: PP072JAA		Category and Description Package Insert, BACTEC Pyrazinamide (PZA) Drug Kit	Sheet: 1 of 13 <hr/> Scale: N/A

A

BD BACTEC™ Pyrazinamide (PZA) Drug Kit

PP072JAA
2010/07

English: pages 1 – 3 Italiano: pagine 6 – 7
 Français : pages 3 – 4 Português: páginas 8 – 9
 Deutsch: Seiten 4 – 6 Español: páginas 10 – 11

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свѣржете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкции свяжитесь с местным представителем компании BD. / Өзініздің жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

BACTEC™ Pyrazinamide (PZA) Drug Kit is for pyrazinamide susceptibility testing. Principal use is with BACTEC PZA Test Medium in conjunction with the BACTEC 460TB instrument.

Pyrazinamide is provided in lyophilized form along with a special Reconstituting Fluid (RF) for use with BACTEC PZA Test Medium for susceptibility testing of *Mycobacterium tuberculosis*.

SUMMARY AND EXPLANATION

For conventional PZA drug susceptibility testing of *M. tuberculosis* a special medium with a pH of 5.5 is required.¹ For the BACTEC PZA susceptibility test, BACTEC PZA Test Medium with a pH of 6.0 should be used.² Lyophilized PZA along with a Reconstituting Fluid is provided to make the BACTEC procedure more convenient. The drug should be rehydrated with Reconstituting Fluid provided with the kit. RF contains a growth promotion substance, polyoxyethylene stearate (POES). When the reconstituted PZA is added to the medium, POES is also introduced into the medium. RF should also be added to BACTEC 12B Medium used to establish the culture prior to testing. RF is also added to the PZA Test Medium control vial.

PRINCIPLES OF THE PROCEDURE

For BACTEC PZA susceptibility test, BACTEC PZA Test Medium at pH 6.0 should be used.² At this pH, PZA activity against mycobacteria can be determined without inhibiting the growth of most *M. tuberculosis* isolates. To compensate for the increase in the pH, a higher concentration of PZA (100 µg/mL) is used in the BACTEC test as compared to the conventional test (25 – 50 µg/mL).

The principle of the BACTEC PZA susceptibility test procedure is the same as that in the BACTEC susceptibility test of the antituberculosis drugs.³⁻⁵ Once the lyophilized PZA is reconstituted with 5 mL of Reconstituting Fluid, the desired PZA concentration in the PZA Test Medium is achieved by adding 0.1 mL of the prepared drug stock to PZA Test Medium.

REAGENTS

Lyophilized Drug
 Pyrazinamide 20,000 µg/vial
 Reconstituting Fluid (RF)
 Polyoxyethylene stearate 4,000 µg/mL
 Processed water 10 mL

Warnings and Precautions:

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁶⁻⁹ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Contaminated material should be sterilized by autoclaving before discarding. Reconstitute the drug vial only with Reconstituting Fluid (refer to the BACTEC PZA Test Medium package insert, PP-071JAA).

Storage and Stability

Prior to reconstitution, the PZA Drug Kit should be stored at 2–8°C. Once reconstituted, the drug may be stored frozen in small aliquots at -5 – -70°C for up to six months. Avoid repeated thawing and freezing. The expiration date applies to the unopened lyophilized bottle stored as recommended.

Disposal

All used bottles should be sterilized by autoclaving before discarding.

SPECIMEN PREPARATION

Refer to the **BACTEC PZA Test Medium** package insert, PP-071JAA.

PROCEDURE

Materials Provided:

Lyophilized PZA drug
Reconstituting Fluid

Materials Required But Not Provided:

BACTEC PZA Test Medium
BACTEC 12B Medium
BACTEC 460TB Instrument
Tuberculin Syringe with permanently attached needle
5 mL Syringe
70% Isopropyl Alcohol

} See "Availability"

Preparation of PZA Solution:

Clean the tops of vial with an alcohol swab. Using a disposable syringe, aseptically transfer 5 mL of RF into the PZA drug vial. Mix thoroughly until dissolved. This stock solution contains 4,000 µg PZA/mL.

Precautions: Do not use the lyophilized PZA vial if it is cracked or displays any moisture. Do not prepare PZA with water or any diluent other than RF. The volume of RF in the vial may vary. Measure exactly 5 mL from the supplied vial when reconstituting the PZA powder. Addition of the correct amount of the test drug into the medium is critical to the performance of the susceptibility test. Add exactly 0.1 mL of stock solution using a 1 mL tuberculin syringe.

Refer to the PZA Test Medium package insert (PP-071JAA) for further instructions.

QUALITY CONTROL

Refer to the PZA Test Medium package insert (PP-071JAA).

LIMITATIONS OF THE PROCEDURE

The PZA Drug concentration has been designed for use with the **BACTEC PZA** susceptibility testing of *M. tuberculosis*. A single concentration is recommended for the **BACTEC** procedure; however, if additional concentrations are desired (other than the recommended concentration) the amount of Reconstituting Fluid used to reconstitute the drug vial may be manipulated. All other desired dilutions of the drug should be made with RF.

EXPECTED VALUES

When PZA is added to the PZA Test Medium along with a test culture of *M. tuberculosis*, there will be a decrease or increase in the daily GI output depending upon the susceptibility of the culture. The PZA Test Medium serving as a growth control (drug-free) will demonstrate a gradual daily increase in the GI. (For details see the PZA Test Medium package insert, PP-071JAA.)

PERFORMANCE CHARACTERISTICS

See the **BACTEC PZA Test Medium** Package insert (PP-071JAA) for details.

AVAILABILITY

Cat. No.	Description
442139	BACTEC™ PZA Test Medium Culture Vials, 4 mL, ctn. of 10
442143	BACTEC™ PZA Drug/Reconstituting Fluid Kit (50 tests)
442004	BACTEC™ 12B Medium
445203	BACTEC™ 460TB System Instrument

REFERENCES

1. Butler, W.R. and J.O. Kilburn. 1982. "Improved method for testing susceptibility of *Mycobacterium tuberculosis* in pyrazinamide." *J. Clin. Microbiol.* 16:1106-1109.
2. Salfinger, M. et al. "Rapid radiometric method for pyrazinamide susceptibility testing of *Mycobacterium tuberculosis*." Abstract 1210, 20th ICAAC Meeting, Los Angeles, CA, 1988.
3. Middlebrook, G. et al. 1977. "Automated radiometric detection of growth of *Mycobacterium tuberculosis* in selective medium." *Am. Rev. Resp. Dis.* 115: 1066-1069.
4. Roberts, G.D. et al. 1983. "Evaluation of the **BACTEC** radiometric method for recovery of mycobacteria and drug susceptibility testing of *M. tuberculosis* from acid-fast smear positive specimens." *J. Clin. Microbiol.* 18:689-696.
5. Siddiqi, S.H. et al. 1981. "Evaluation of a rapid radiometric method for drug susceptibility testing of *Mycobacterium tuberculosis*." *J. Clin. Microbiol.*, 13:908-912.
6. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, PA.

- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

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Benex Limited