BD BACTEC™ Peds Plus™/F Culture Vials Soybean-Casein Digest Broth with Resins

English:	pages	1 – 3	Español:	páginas	12 – 14
Français :	pages	4 – 6	Dansk:	side	14 – 17
Deutsch:	Seiten	6 – 9	Português:	páginas	17 – 20
Italiano:	pagina	9 – 11	Svenska:	sidan	20 – 22

PP091JAA(02)
2015-05

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD BACTEC Peds Plus™/F culture vials (enriched Soybean-Casein Digest broth with CO₂) are for aerobic blood cultures. Principal use is with the BD BACTEC fluorescent series instruments for the qualitative culture and recovery of aerobic microorganisms (mainly bacteria and yeast) from pediatric and other blood specimens which are generally less than 3 mL in volume.

SUMMARY AND EXPLANATION

The sample to be tested is inoculated into the vial which is then inserted into the BD BACTEC fluorescent series instrument for incubation and periodic reading. Each vial contains a sensor which can detect increases in CO₂ produced by the growth of microorganisms. The sensor is monitored by the instrument every ten minutes for an increase in its fluorescence, which is proportional to the amount of CO₂ present. A positive reading indicates the presumptive presence of viable microorganisms in the vial. Detection is limited to microorganisms that will grow in a particular type of medium.

Resins have been described for the treatment of blood specimens both prior to and after their inoculation into culture media. Resins have been incorporated into BD BACTEC culture media to enhance recovery of organisms without a need for special processing.¹⁻³

PRINCIPLES OF THE PROCEDURE

If microorganisms are present in the test sample inoculated into the BD BACTEC vial, CO₂ will be produced when the organisms metabolize the substrates present in the vial. Increases in the fluorescence of the vial sensor caused by the higher amount of CO₂ are monitored by the BD BACTEC fluorescent series instrument. Analysis of the rate and amount of CO₂ increase enables the BD BACTEC fluorescent series instrument to determine if the vial is positive; i.e., that the test sample contains viable organisms.

REAGENTS

The BD BACTEC Peds Plus/F culture vials contain the following active ingredients prior to processing:

List of Ingredients

Processed Water	40 mL
Soybean-Casein Digest Broth	2.75% w/v
Yeast Extract	0.25% w/v
Animal Tissue Digest	0.10% w/v
Sodium Pyruvate	0.10% w/v
Dextrose	0.06% w/v
Sucrose	0.08% w/v
Hemin	0.0005% w/v
Menadione	0.00005% w/v
Sodium Polyanetholsulfonate (SPS)	0.020% w/v
Pyridoxal HCl (Vitamin B6)	0.001% w/v
Nonionic Adsorbing Resin	10.0% w/v
Cationic Exchange Resin	0.6% w/v

All BD BACTEC media are dispensed with added CO₂.

Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. Standard Precautions⁴⁻⁷ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Prior to use, each vial should be examined for evidence of damage, contamination or deterioration. Vials displaying evidence of damage or contamination such as leakage, cloudiness, discoloration (darkening), bulging or depressed septum should not be used.

A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, contaminated culture media could be refluxed into the patient's vein. Vial contamination may not be readily apparent. When using direct draw procedures, monitor the process closely to avoid refluxing materials into patient.

On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions, a vial may not be sealed sufficiently. In both cases the contents of the vial may leak or spill. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: Before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling should be performed in a biological safety cabinet if possible, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or Luer-Lok™ brand tips.

Storage Instructions

The BD BACTEC vials are ready for use as received and require no reconstitution or dilution. Store at 2° – 25 °C, in a dry location **out of direct light**.

SPECIMEN COLLECTION

The specimen must be collected using sterile techniques to reduce the chance of contamination. The range of blood volume which can be cultured is 0.5 to 5.0 mL. Optimum results are obtained with 1.0 to 3.0 mL. If the volume of blood cultured is less than 0.5 mL, recovery of some fastidious organisms, such as *Haemophilus* species, may require the use of an appropriate supplement, as described later in this Package Insert. It is recommended that the specimen be inoculated into the **BD BACTEC** vials at bedside. Most commonly, a syringe with a **Luer-Lok™** brand tip is used to draw the sample. If appropriate, a **Vacutainer™** brand Needle Holder and a **Vacutainer** brand Blood Collection Set, **Vacutainer Safety-Lok™** Blood Collection Set or other tubing "butterfly" set may be used. If using a needle and tubing set (direct draw), carefully observe the direction of blood flow when starting sample collection. The vacuum in the vial will usually exceed 5 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. When the desired 1 – 3 mL has been drawn, the flow should be stopped by crimping the tubing and removing the tubing set from the **BD BACTEC** vial. **The inoculated BD BACTEC vial should be transported as quickly as possible to the laboratory.**

PROCEDURE

Remove the flip-off cap from **BD BACTEC** vial top and inspect the vial for cracks, contamination, excessive cloudiness, and bulging or indented stoppers. **DO NOT USE** if any defect is noted. Before inoculating, swab the septum with alcohol (iodine is **not** recommended). Aseptically inject or draw directly a maximum of 5 mL of specimen per vial (see the section on Limitations of the Procedure). **Inoculated aerobic vials should be placed in the BD BACTEC fluorescent series instrument as soon as possible** for incubation and monitoring. If placement of an inoculated vial into the instrument has been delayed and visible growth is apparent, it should not be tested in the **BD BACTEC** fluorescent series instrument, but rather it should be subcultured, Gram-stained and treated as a presumptively positive bottle.

Vials entered into the instrument will be automatically tested every ten minutes for the duration of the testing protocol period. Positive vials will be determined by the **BD BACTEC** fluorescent series instrument and identified as such (see the appropriate **BD BACTEC** fluorescent series instrument User's Manual). The sensor inside the bottle will not appear visibly different in positive and negative vials, however the **BD BACTEC** fluorescent series instrument can determine a difference in fluorescence.

If, at the end of the testing period, a negative vial appears visually positive (i.e., chocolateized blood, bulging septum, lysed and/or very darkened blood), it should be subcultured and Gram-stained and treated as a presumptive positive.

Positive vials should be subcultured and a Gram-stained slide prepared. In a great majority of cases, organisms will be seen and a preliminary report can be made to the physician. Subcultures to selective media and a preliminary direct antimicrobial susceptibility test may be prepared from fluid in the **BD BACTEC** vials.

Subculturing: Prior to subculturing, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, insert a sterile needle with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

For maximum yield of isolates, negative cultures may be checked by stain and/or subcultured at some point prior to discarding as negative.

QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Do not use culture vials past their expiration date.

Do not use culture vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each carton of media. Quality Control Certificates list test organisms, including ATCC™ cultures specified in the CLSI Standard, *Quality Control for Commercially Prepared Microbiological Culture Media*.⁸ The range of time-to-detection in hours was ≤ 72 hours for each of the organisms listed on the Quality Control Certificate for this medium:

Organism
<i>Streptococcus pyogenes</i> ATCC 19615
<i>Escherichia coli</i> ATCC 25922
<i>Streptococcus pneumoniae</i> * ATCC 6305
<i>Pseudomonas aeruginosa</i> ATCC 27853
<i>Candida albicans</i> ATCC 18804
<i>Neisseria meningitidis</i> ATCC 13090
<i>Alcaligenes faecalis</i> ATCC 8750
<i>Haemophilus influenzae</i> ATCC 19418
<i>Staphylococcus aureus</i> ATCC 25923

*CLSI Strain

For information on Quality Control for the **BD BACTEC** fluorescent series instrument, refer to the appropriate **BD BACTEC** fluorescent series instrument User's Manual.

LIMITATIONS OF THE PROCEDURE

Contamination

Care must be taken to prevent contamination of the sample during collection and inoculation into the **BD BACTEC** vial. A contaminated sample will give a positive reading, but will not indicate a relevant clinical result. Such a determination must be made by the user based on such factors as type of organisms recovered, occurrence of the same organism in multiple cultures, patient history, etc.

Recovery of SPS Sensitive and Fastidious Organisms from Blood Samples

Because blood can neutralize the toxicity of SPS toward organisms sensitive to SPS (such as some *Neisseria* species), the presence of optimum volumes of blood (1 – 3 mL) is of benefit in the recovery of these organisms.

Some fastidious organisms, such as certain *Haemophilus* species, require growth factors, such as NAD, or factor V, which are provided by the blood specimen. If the blood specimen volume is very small (0.5 mL or less), an appropriate supplement may be required for recovery of these organisms. **BD BACTEC** brand **FOS™** Fastidious Organism Supplement may be used as a nutritional supplement.

Nonviable Organisms

A Gram-stained smear from a culture medium may contain small numbers of nonviable organisms derived from medium constituents, staining reagents, immersion oil, glass slides, and specimens used for inoculation. In addition, the patient specimen may contain organisms that will not grow in the culture medium or in media used for subculture. Such specimens should be subcultured to special media as appropriate.⁹

Antimicrobial Activity

Neutralization of the antimicrobial activity by resins varies depending on dosage level and timing of specimen collection.

Studies have demonstrated that the resins present in this medium do not adequately neutralize imipenem-cilastatin antimicrobial preparations.

Recovery of *Streptococcus pneumoniae*

In aerobic media, *S. pneumoniae* will typically be visually and instrument positive, but in some cases no organism will be seen on Gram stain or recovered on routine subculture. If an anaerobic vial was also inoculated, the organism can usually be recovered by performing an aerobic subculture of the anaerobic vial, since this organism has been reported to grow well under anaerobic conditions.¹⁰

General Considerations

Optimum recovery of isolates will be achieved by adding 1 – 3 mL of blood. Use of lower or higher volumes may adversely affect recovery and/or detection times. Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms. False negative readings may result when certain organisms are present which do not produce enough CO₂ to be detected by the system or if significant growth has occurred before placing the vial into the system. False positivity may occur when the white blood cell count is high.

EXPECTED RESULTS

Seeded culture studies were performed using inoculum levels targeted at 10 to 50 CFU per culture vial of both ATCC™ and wild microbial strains. The following is a list of the organisms which were detected as positive in **BD BACTEC Peds Plus™/F Culture Medium** within a five (5) day period.

<i>Acinetobacter baumannii</i>	<i>Haemophilus influenzae</i>	<i>Serratia marcescens</i>
<i>Acinetobacter Iwoffii</i>	<i>Haemophilus parainfluenzae</i>	<i>Staphylococcus aureus</i>
<i>Alcaligenes faecalis</i>	<i>Klebsiella pneumoniae</i>	<i>Staphylococcus epidermidis</i>
<i>Candida albicans</i>	<i>Neisseria gonorrhoeae</i>	<i>Stenotrophomonas (Xanthomonas) maltophilia</i>
<i>Candida (Torulopsis) glabrata</i>	<i>Neisseria meningitidis</i>	<i>Streptococcus pneumoniae</i>
<i>Corynebacterium J-K</i>	<i>Proteus mirabilis</i>	Group A <i>Streptococcus</i>
<i>Cryptococcus neoformans</i>	<i>Providencia stuartii</i>	Group D <i>Streptococcus</i>
<i>Enterobacter cloacae</i>	<i>Pseudomonas aeruginosa</i>	
<i>Escherichia coli</i>	<i>Salmonella typhimurium</i>	

Internal studies have demonstrated that antimicrobials are effectively neutralized by the resins used in **BD BACTEC** resin media. In these tests, antimicrobials were added in clinically relevant concentrations directly to resin media prior to inoculation with susceptible strains. These tests were performed in parallel using non resin media as controls. Antimicrobials representative of the following categories were found to be neutralized by the resins: penicillins, cephalosporins (1st, 2nd and 3rd generation), macrolides, aminoglycosides, fluoroquinolones, tetracycline and chloramphenicol.

PERFORMANCE CHARACTERISTICS

In an external 2-site clinical study comparing the performance of **BD BACTEC Peds Plus/F Culture Medium** to the **BD BACTEC NR Peds Plus Culture Medium**, a total of 3,249 compliant paired specimens were tested. A total of 430 organisms were recovered. Of these, 342 (80%) were deemed clinically significant and 88 (20%) were considered not clinically significant. Of the clinically significant isolates, 219 (64%) were recovered in both media, 63 (18%) were recovered in the **BD BACTEC Peds Plus/F Culture Medium** only, and 60 (18%) were recovered into the **BD BACTEC NR Peds Plus Culture Medium** only. Table 1 shows the isolates recovered by media type. Four false positive culture vials and no false negative culture vials were identified in the **BD BACTEC Peds Plus/F Culture Medium** during this study.

The average time-to-detection for the **BD BACTEC Peds Plus/F Culture Medium** was 25 hours, while the average time-to-detection for **BD BACTEC NR Peds Plus Culture Medium** was 33 hours.

Table 1: Clinical Study Isolate Recovery – Media Type

Organism	Recovered in Peds Plus/F Only	Recovered in NR Peds Plus Only	Recovered in both
Gram Negatives	9	6	50
Gram Positives	51	43	146
Yeast	3	11	23

List of Organisms Recovered in BD BACTEC Peds Plus/F during Clinical Trials:

<i>Actinobacillus lignieresii</i>	<i>Escherichia coli</i>	<i>Streptococcus constellatus</i>
<i>Burkholderia (Pseudomonas) cepacia</i>	<i>Klebsiella pneumoniae</i>	<i>Streptococcus pneumoniae</i>
<i>Candida albicans</i>	<i>Neisseria meningitidis</i>	<i>Streptococcus salivarius</i>
<i>Candida tropicalis</i>	<i>Pseudomonas fluorescens</i>	<i>Streptococcus intermedius</i>
<i>Citrobacter freundii</i>	<i>Salmonella</i> spp.	
<i>Enterobacter aerogenes</i>	<i>Serratia liquefaciens</i>	
<i>Enterobacter cloacae</i>	<i>Staphylococcus aureus</i>	
<i>Enterococcus</i> spp.	<i>Staphylococcus</i> spp. coag. neg.	
	<i>Streptococcus</i> Group B	

AVAILABILITY

Cat. No.	Description
442194	BD BACTEC™ Peds Plus™/F Culture Vials , case of 50 vials

REFERENCES:

- Wallis, C. et al. Rapid isolation of bacteria from septicemic patients by use of an antimicrobial removal device. *J. Clin. Microbiol.* 1980, 11:462-464.
- Applebaum, P.C. et al. Enhanced detection of bacteremia with a new **BD BACTEC** resin blood culture medium. *J. Clin. Microbiol.* 1983, 17:48-51.
- Jungkind, D.L. et al. Evidence for a second mechanism of action of resin in **BD BACTEC** NR16A aerobic blood culture medium. Abstracts of the Annual Meeting of Amer. Soc. for Microbiol. 1989.
- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
- U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
- Clinical and Laboratory Standards Institute. 2004. Approved Standard M22-A3. Quality control for commercially prepared microbiological culture media, 3rd ed., CLSI, Wayne, Pa.
- Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller. 2007. Manual of clinical microbiology, 9th ed., American Society for Microbiology, Washington, D.C.
- Howden, R.J. *J. Clin. Path.* 1976, 29:50-53.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Becton Dickinson Pty Ltd.
Benex Limited